**Human/Rat  $\beta$ Amyloid (42) ELISA Kit Wako, High-Sensitive****TABLE OF CONTENTS**

1. Introduction .....	2
2. Kit Contents .....	2
3. Required Apparatuses .....	2
4. Principle .....	2
5. Assay Method	
5-1. Preparation of reagents to be used .....	2
5-2. Test sample preparation .....	3
5-3. Assay procedure .....	3
5-4. Standard curve .....	4
6. Test performance	
6-1. Sensitivity .....	4
6-2. Reproducibility .....	4
6-3. Specificity .....	4
6-4. Linearity of dilution .....	5
6-5. Spike recovery .....	5
7. Troubleshooting .....	5
8. Method of sample preparation .....	6
9. References .....	8
10. Assay Summary .....	9

## 1. Introduction

Alzheimer's Disease (AD) is characterized by the presence of extracellular senile plaques (SPs) and intracellular neurofibrillary tangles (NFT) in the brain. The major protein component of SPs is  $\beta$ Amyloid peptide (A  $\beta$ ). The peptide is 40 or 42 amino acid peptide cleaved from the Amyloid Precursor Protein (APP) by  $\beta$ -Secretase and  $\gamma$ -Secretase. A  $\beta$ 42 is more prone to aggregate than A  $\beta$ 40. Therefore the initial A  $\beta$  deposition begins with A  $\beta$ 42 but not with A  $\beta$ 40. A  $\beta$ 42 positive and A  $\beta$ 40-negative plaques may represent early-stage diffuse type SPs, and A  $\beta$ 40-positive plaque appears in the advanced stage, especially more often in cored portion of the mature plaque.

In this kit, we use the monoclonal antibody BNT77, which epitope is A  $\beta$ (11-28) and the monoclonal antibody BC05, which specifically detects the C-terminal portion of A  $\beta$ 42. Therefore this kit is designed to be used for the quantitative determination not only Human or Rat A  $\beta$ (1-42) but also A  $\beta$ (x-42) with a truncated or modified N-terminus in samples such as tissue culture medium, tissue homogenate, CSF and plasma. Additionally, the sensitivity of the kit increased about 10 times compared with Human/Rat (42)  $\beta$ Amyloid ELISA Kit Wako (290-62601).

## 2. Kit Contents

① Antibody (BNT77)-coated Microtiter Plate (Anti Human A $\beta$ (11-28) MoAb (Clone No.BNT77))	1 plate
② Standard Solution (Human $\beta$ Amyloid (1-42), 20 pmol/L)	2 mL $\times$ 2 vials
③ Standard Diluent	30 mL $\times$ 1 vial
④ Wash Solution (20 $\times$ )	50 mL $\times$ 1 vial
⑤ HRP-conjugated Antibody (BC05) Solution (Anti Human A $\beta$ (35-43) MoAb (Clone No.BC05)Fab'-HRP)	12 mL $\times$ 1 vial
⑥ TMB Solution	12 mL $\times$ 1 vial
⑦ Stop Solution	12 mL $\times$ 1 vial
⑧ Plate Seal	3 sheets

**Storage** : Keep at 2–10℃.

**Expiration date** : Printed on the package

**Package** : 96 tests

## 3. Required Apparatuses

1. A Microplate Reader
2. A Micropipette
3. A Microplate Washer

## 4. Principle

This kit is constructed as a sandwich ELISA format with two kinds of antibodies. The monoclonal antibody BNT77, which epitope is A  $\beta$ (11-28) is coated on 96 well surfaces of a separable microplate and acts as a capture antibody for A  $\beta$ (x-42). Captured A  $\beta$ (x-42) is recognized by another antibody, BC05 (Fab' fragment), which is specifically detects the C-terminal portion of A  $\beta$ (x-42), labeled with HRP. After addition of TMB solution, positive samples will develop a blue color. The reaction is terminated by the addition of stop solution, which produces a yellow color. The absorbance is then measured at 450 nm.

## 5. Assay Method

### 5-1. Preparation of reagents to be used

Number	Reagent Name	Preparation
①	Antibody (BNT77)-coated Microplate	Ready to use
②	Standard Solution (Human $\beta$ Amyloid (1-42), 20 pmol/L)	Dilute the Standard with Standard Diluent. (e.g. 20, 10, 5, 2, 1, 0.5, 0.1 (pmol/L))
③	Standard Diluent	Ready to use
④	Wash Solution (20 $\times$ )	Working solution is prepared by dilution of 50 mL of the Wash Solution(20 $\times$ ) with distilled deionized water to 1000 mL.
⑤	HRP-conjugated Antibody (BC05)	Ready to use
⑥	TMB Solution	Ready to use
⑦	Stop Solution	Ready to use
⑧	Plate Seal	Ready to use

### 5-2. Test sample preparation

Test samples, containing A  $\beta$ (x-42) at more than 20 pmol/L, should be diluted with the Standard Diluent. If the concentration of A  $\beta$ (x-42) in samples cannot be estimated in advance, a pre-assay with several different dilutions is recommended to determine the proper dilution of samples. Additionally, in the case of plasma samples, they should be diluted 2-4 times with the Standard Diluent to disrupt interaction of A  $\beta$ (x-42) with masking proteins. In the case of culture medium with FCS, A  $\beta$ (x-42) like substances in FCS may be measured. Therefore negative control should be taken.

### 5-3. Assay procedure

- ① Take the aluminium package, containing the antibody-coated microplate, from a refrigerator which is maintained at 2-10 °C, and leave it until it reaches room temperature. Separate the required wells from the microplate, and seal the remaining wells in the aluminium package and store them in the refrigerator.
- ② Dispense 100  $\mu$ L standard diluent into zero blank well (*See Note 1*).
- ③ In the assays of standard solution and test samples (e.g. tissue culture medium, tissue homogenate, cerebrospinal fluid (CSF) and plasma) dispense 100  $\mu$ L of solution into each well (*See Note 2,3*).
- ④ Seal the wells with the plate seal, and leave refrigerated overnight.
- ⑤ After removal of the plate seal from the plate, discard the solutions from the wells by microplate washer or aspirator, and then wash 5 times with wash solution (*See Note 4,5*).
- ⑥ Dispense 100  $\mu$ L HRP-conjugated antibody solution into the wells, seal them with the plate seal, and let the microplate stand in the refrigeration or cold room for 1 hour.
- ⑦ After removal of the plate seal, remove the antibody solution from the wells and wash them 5 times with wash solution (*See Note 4,5*).
- ⑧ Add 100  $\mu$ L TMB solution to the wells within a short interval, thus starting the HRP reaction at room temperature in the dark (*See Note 1,6*).
- ⑨ Add 100  $\mu$ L stop solution to the wells, 30 minutes after addition of the TMB solution, in order to terminate the reaction (*See Note 7*).
- ⑩ Read the absorbance of each well at 450 nm having blank the plate reader against a TMB solution blank. Read the plate within 30 minutes after adding the stop solution.
- ⑪ Read the Human or Rat  $\beta$ Amyloid (x-42) concentrations for unknown samples and controls from the standard curve.

### Notes


1. All reagents except for the wash solution and HRP-conjugated Antibody Solution should be used only after returning to room temperature.
2. At least more than duplicate assays are recommended for each sample.
3. The standard curve should be made at each study.
4. Try to complete washing within 2 or 3 minutes in total. The temperature of the wash solution should be equal to or lower than room temperature. And differences in soaking time (i.e., time taken for pouring the wash solution into wells) between the wells in the washing process may cause inconsistent results. Try to perform the washing process in such a way as to ensure that the soaking time for each well is as similar as possible (e.g., start washing from the right row of the plate for the first wash and from the left row in the second wash).
5. If using an aspirator for removal of the solutions from microplate wells, pay attention to not touch and damage the surfaces of the wells.
6. The vessels and apparatus should be well cleaned, prior to use. Contamination results in color development of the solutions.
7. When handling many samples, systematic dispensing of the solutions should be required in order to attain a constant reaction time in each sample assay.
8. The kit is constructed with well-adjusted combinations of the antibody preparations and biological reagent solutions for each lot. Combining the biological reagent solutions and the plate from different lots may cause unexpected results.
9. The quantitation may indicate abnormally high measurements caused by the cross reaction of mouse IgG1 antibodies to unidentified components in human plasma. In that case, the plasma sample should be absorbed with a column which fixes refined mouse IgG. Then you can minimize the phenomenon. (*see Reference 8,9*)
10. This kit is for research use only, and not for use in diagnostic procedures.

**Assay procedure table**

	Test Sample	Standard	Test Sample Blank
Reagents	Test Sample 100 $\mu$ L	Diluted standard 100 $\mu$ L	Standard Diluent 100 $\mu$ L
Incubate for overnight in the refrigerator with plate seal			
Washing 5 times			
HRP-conjugated Antibody	100 $\mu$ L	100 $\mu$ L	100 $\mu$ L
Incubate for 1 hour in the refrigerator with plate seal			
Washing 5 times			
TMB Solution	100 $\mu$ L	100 $\mu$ L	100 $\mu$ L
Incubate for 30 minutes at room temperature in the dark with plate seal			
Stop Solution	100 $\mu$ L	100 $\mu$ L	100 $\mu$ L
Read the absorbance of each well at 450 nm with a microplate reader.			

**5-4. Standard curve**

Standard (pmol/L)	Mean (OD at 450 nm)	CV (%)
0.0	0.046	1.26
0.1	0.056	2.74
0.5	0.097	2.39
1.0	0.154	0.99
2.0	0.264	1.31
5.0	0.582	1.07
10.0	1.099	0.48
20.0	2.092	1.01

**6. Test performance****6-1. Sensitivity**

- Dynamic range.....0.1-20.0 (pmol/L)
- Sensitivity.....0.024 (pmol/L)  
 ※The sensitivity was determined by adding two standard deviations to the mean absorbance obtained when the zero standard was assayed 24 times.

**6-2. Reproducibility**

## 1) Intra-assay (n=24)

	Sample 1	Sample 2	Sample 3
Mean (pmol/L)	3.16	7.40	16.41
SD	0.02	0.06	0.15
CV(%)	0.77	0.82	0.96

SD = Standard Deviation, CV = Coefficient of Variation

## 2) Inter-assay (n=5)


	Sample 1	Sample 2	Sample 3
Mean (pmol/L)	3.17	7.23	16.38
SD	0.05	0.13	0.18
CV (%)	1.50	1.86	1.12

**6-3. Specificity (n=4)**

	Cross-reaction (%)
Human $\beta$ Amyloid (1-40)	$\leq 0.1$
Human $\beta$ Amyloid (1-42)	100.0
Human $\beta$ Amyloid (1-43)	12.7
Rat (mouse) $\beta$ Amyloid (1-40)	$\leq 0.1$
Rat (mouse) $\beta$ Amyloid (1-42)	156.0

#### 6-4. Linearity of dilution (n = 3)

Sample	Dilution	Measured (pmol/L)
Plasma A	1/2	3.193
	1/4	1.566
	1/8	0.736
	1/16	0.377
Plasma C	1/2	2.743
	1/4	1.361
	1/8	0.686
	1/16	0.335
Plasma E	1/2	3.033
	1/4	1.627
	1/8	0.780
	1/16	0.409

#### 6-5. Spike recovery (n = 3)

Sample (4 fold dilution)	Spiked (pmol/L)	Measured (pmol/L)	CV (%)	Expected (pmol/L)	Recovery (%)
Plasma A	0.00	1.50	1.35		
	1.25	2.79	4.03	2.75	101.57
	2.50	4.06	0.96	4.00	101.50
	5.00	6.62	1.67	6.50	101.90
Plasma B	0.00	2.24	2.14		
	1.25	3.54	0.16	3.49	101.23
	2.50	4.74	1.13	4.74	99.90
	5.00	7.08	0.90	7.24	97.71
Plasma C	0.00	1.29	1.31		
	1.25	2.47	0.99	2.54	97.35
	2.50	3.80	0.90	3.79	100.22
	5.00	6.18	0.83	6.29	98.30
Plasma D	0.00	1.91	1.06		
	1.25	3.05	0.80	3.16	96.43
	2.50	4.41	1.17	4.41	99.96
	5.00	6.64	0.75	6.91	96.01
Plasma E	0.00	1.55	2.21		
	1.25	2.75	0.20	2.80	98.17
	2.50	3.93	0.74	4.05	97.07
	5.00	6.09	0.49	6.55	92.96
Plasma F	0.00	1.80	1.90		
	1.25	2.86	0.71	3.05	93.85
	2.50	4.01	0.24	4.30	93.24
	5.00	6.07	0.70	6.80	89.21

## 7. Troubleshooting

Q1 : The OD in the standard solution is too low to draw a calibration curve. Why does this happen ?

A1 : It is possible that more time than necessary was taken for washing. Try to complete washing within 2 or 3 minutes in total. The temperature of the wash solution should be equal to or lower than room temperature. If you do not use a microplate washer or aspirator, discard the solution in the wells first by decanting (be careful not to cause inter-well contamination at this time), fill the wells with the wash solution using a bottle washer, etc. and discard the solution again. Repeat the procedure 5 times. The total washing time in this case should also be no more than 2 to 3 minutes.

Q2 : The absorbance of 0 pmol/L is higher than standards with the low concentration range, so the calibration curve cannot be drawn. Why does this happen ?

A2 : There is a possibility of insufficient washing. Samples, standard solutions, or HRP conjugated antibody solutions can attach to the upper part of the sides of the wells, and may not be removed by the washing operation. Be careful to avoid touching the sides of the wells when adding the solution. During the movement of the plate be careful not to drop or tilt it. During the reaction time, the plate should remain still.

In case of manual washing operations, make sure at least 300  $\mu$ L of the wash solution is dispensed into each well. When dispensing the first two times, do not overflow the wash solution from the wells in order to avoid well-to-well contamination. The third time, thoroughly wash by overflowing of the wash solution. When using a plate washer, we recommend washing with the maximum amount that can be set.

Q3 : Does A  $\beta$  peptide in the standard solution aggregate ?

A3 : The standard solution is designed not to aggregate. It is stable when stored in a refrigerator.

Q4 : The measurement values are not consistent. What are the possible causes of this ?

A4 : Differences in soaking time (i.e., time taken for pouring the wash solution into wells) between the wells in the washing process may cause inconsistent results. Manual washing using a pipet is particularly prone to variations in timing. Try to perform the washing process in such a way as to ensure that the soaking time for each well is as similar as possible (e.g., start washing from the right row of the plate for the first wash and from the left row in the second wash).

Q5 : Storage method

A5 : Human specimens such as cerebral tissues and cerebrospinal fluid can be stored in a freezer, although measurement values for some specimens may decrease even after such storage. If the concentration of A  $\beta$  is low, then freeze-thawing of the specimen should be avoided. For a culture supernatant, mix it with 0.2% cow serum albumin and 0.075% CHAPS (final concentrations for both) to minimize losses by adhesion to the vial wall, etc. before storage in a freezer.

Q6 : Is it possible to measure a serum sample ?

A6 : Yes. However, it is not recommended because the value in serum is often lower than that in plasma that is collected simultaneously. The serum sample also tends to deteriorate to a greater degree on refrigeration.

Q7 : Which anticoagulant can I use for plasma samples ?

A7 : We recommend EDTA2K. EDTA2Na and heparin have also been used.

Q8 : Can I use concentrations for dilution of the standard solution other than those illustrated ?

A8 : No problem at all. You can prepare a series of dilutions suitable for your experimental system.

Q9 : Is a calibration curve of a type other than log-log acceptable ?

A9 : No problem at all.

Q10 : While it is recommended that the assays be performed at least in duplicate, is it also necessary to perform duplicate assays for standard samples ?

A10 : We recommend that the assays are at least duplicated for standard samples, but you may decide the number of standard samples, such as  $n=1$ , to suit to each experimental system.

Q11 : How should we store the reagent after opening ?

A11 : The reagent is stable in a refrigerator provided that the microplates are returned to the aluminum packs properly with a desiccating agent.

## 8. Method of sample preparation

Following methods are examples of how to prepare the sample.

### ◆ Brain tissue sample

#### 1. Aggregated A $\beta$ (brain tissue in patients with Alzheimer's disease) (see Reference 6)

① Mince about 1 g of brain tissue, add 5 volumes of TBS\* and homogenize (10 strokes) using a Teflon homogenizer. Centrifuge this homogenate at  $500,000 \times g$ , 4  $^{\circ}$ C for 20 minutes.

\* (50 mmol/L Tris-HCl, pH 7.6, 150 mmol/L NaCl, protease inhibitors [0.1 mmol/L diisopropyl fluorophosphate, 0.5 mmol/L phenylmethylsulfonyl fluoride, 1 g/mL N $\alpha$ -P-tosyl-L-lysine chloromethyl ketone, 1 g/mL antipain, 0.1 g/mL leupeptin])

- ② Suspend the precipitate in 5 volumes of TBS/protease inhibitors (containing 1.0 mol/L sucrose) and centrifuge at  $500,000 \times g$ , 4°C for 20 minutes.
- ③ Add 3 volumes of 1% Triton X-100/TBS/protease inhibitors to the precipitate and homogenize (10 strokes). Incubate this homogenate at 37°C for 15 minutes and centrifuge at  $500,000 \times g$  for 20 minutes.
- ④ Add 3 volumes of 2% sodium dodecyl sulfate (SDS)/TBS/protease inhibitors to the precipitate and homogenize. Incubate the homogenate at 37°C for 15 minutes and centrifuge at  $500,000 \times g$ , 25°C.
- ⑤ Add 1 mL of 70% formic acid to the precipitate, ultrasonicate and centrifuge at  $500,000 \times g$ , 4°C.
- ⑥ Collect the supernatant and dry by Speed Vac. Add 100  $\mu$ L of DMSO and perform ultrasonic disintegration for a short time for suspension. Store this DMSO-dissolved, formic acid-extracted A $\beta$  at -80°C until use.
- ⑦ For measurement, dissolve the sample in the Standard Diluent in the kit as necessary and follow the instructions in the package insert. The recommended dilution is 1,000 times for A $\beta$  40 and 100 times for A $\beta$  42.

## 2. Soluble A $\beta$ (normal brain tissue)

We recommend the following methods of 2-1 and 2-2 in the case of soluble A $\beta$  in the normal brain tissues. Because we can detect the difference between soluble A $\beta$  and aggregated A $\beta$  clearly compared with the methods of 2-3.

### 2-1. The method using RIPA buffer

- ① Homogenize (25 strokes) about 100 mg the brain tissue in 1 mL of RIPA buffer (250 mM NaCl, 1% NP-40, 0.5% Na-deoxycholate, 0.1% SDS, 50 mM Tris HCl, pH=8.0) with protease inhibitors. Centrifuge this homogenate at  $100,000 \times g$ , 4°C, for 20min.
- ② Dilute the supernatant with the Standard Diluent in the kit, if necessary, and measure as directed in the package insert.

### 2-2. The method using Tris buffer

- ① Homogenize (25 strokes) about 100 mg the brain tissue in 1 mL of Tris buffer (150 mM NaCl, 50 mM Tris HCl, pH=7.6) with protease inhibitors. Centrifuge this homogenate at  $100,000 \times g$ , 4°C, for 20min.
- ② Dilute the supernatant with the Standard Diluent in the kit, if necessary, and measure as directed in the package insert.

### 2-3. The method using 70% formic acid (see Reference 7)

- ① Homogenize (6 strokes) about 150 mg of tissue in 1 mL of 70% formic acid. Centrifuge this homogenate at  $100,000 \times g$  for 1 hour.
- ② Dilute the supernatant 20-fold with 1 mol/L Tris Base to neutralize.
- ③ After neutralization, dilute the sample with the Standard Diluent in the kit, if necessary, and measure as directed in the package insert.

## ◆ Culture supernatant sample

- ① To measure, dilute the sample with the Standard Diluent in the kit as necessary and follow the instructions in the package insert. The recommended dilution is 5-fold for A $\beta$  40 and no dilution for A $\beta$  42.

## ◆ Cerebrospinal fluid sample

- ① To measure, dilute the sample with the Standard Diluent in the kit as necessary and follow the instructions in the package insert. The recommended dilution is 50-fold.

## ◆ Plasma sample (see Reference 8,9)


- ① Centrifuge the blood collected in a vacuum blood collection tube with EDTA2K at  $5,000 \times g$ , 4°C for 15 minutes to separate the plasma, which is then stored at -80°C until use.
- ② To measure, dilute the plasma sample 4-fold with the Standard Diluent in the kit to avoid the effect of interfering substances in the plasma and follow the instructions in the package insert.

## 9. References :


- 1) Suzuki, N., Cheung, TT., Cai, XD., Odaka, A., Otvos, L. Jr., Eckman, C., Golde, TE. and Younkin, SG. : *Science*, **264**, 1336 (1994)
- 2) Iwatsubo, T., Odaka, A., Suzuki, N., Mizusawa, N. and Ihara, Y. : *Neuron*, **13**, 45 (1994)
- 3) Asami-Odaka, A., Ishibashi, Y., Kikuchi, T., Kitada, C. and Suzuki, N. : *Biochemistry*, **34**, 10272 (1995)
- 4) Scheuner, D., Eckman, C., Jensen, M., Song, X., Citron, M., Suzuki, N., Birb, TD., Hardv, J., Hutton, M., Kukull, W., Larson, E., Levv-Lahad, E., Viitanen, M., Peskind, E., Poorkai, P., Schellenberg, G., Tanzi, R., Wasco, W., Lannfelt, L., Selkoe, D. and Younkin, S. : *Nature Med.* **2**, 864 (1996)
- 5) Kosaka, T., Imagawa, M., Seki, K., Arai, H., Sasaki, H., Tsuji, S., Asami-Odaka, A., Fukushima, T., Imai, K. and Iwatsubo, T. : *Neurology* **48**, 741 (1997)
- 6) Hosoda, R., Saido, TC., Otvos, L., Jr., Arai, T., Mann, DMA., Lee, VM-Y, Trojanowski, JQ. and Iwatsubo, T. : *J. Neuropathol. Exp. Neurol.*, **57**, 1089 (1998)
- 7) Borchelt, D., R., *et al.* : *Neuron*, **17**, 1005 (1996)
- 8) Scheuner, D., Eckman, C., Jensen, M., Song, X., Citron, M., Suzuki, N., Bird, TD., Hardy, J., Hutton, M., Kukull, W., Larson, E., Levy-Lahad, E., Viitanen, M., Peskind, E., Poorkaj, P., Schellenberg, G., Tanzi, R., Wasco, W., Lannfelt, L., Selkoe, D. and Younkin, S. : *Nat Med.*, **2**, 864 (1996)
- 9) Kosaka, T., Imagawa, M., Seki, K., Arai, H., Sasaki, H., Tsuji, S., Asami-Odaka, A., Fukushima, T., Imai, K. and Iwatsubo, T. : *Neurology*, **48**, 741 (1997)

## 10. Assay Summary


Add 100  $\mu$ L of Standards and Samples and incubate in the refrigerator with plate seal for overnight.


Wash 5 times.

Add 100  $\mu$ L of HRP-conjugated Antibody Solution and incubate in the refrigerator for 1 hour.

Wash 5 times.

Add 100  $\mu$ L of TMB Solution and incubate for 30 minutes at room temperature in the dark with plate seal.

Add 100  $\mu$ L of Stop Solution and read the absorbance of each well at 450 nm with a microplate reader.

### FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

#### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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