FUJIFILM

Wako

〈For Research Use Only〉

Code No. 292-67801 (50 reactions)

for Genomic Research
**8-OHdG Assay Preparation
Reagent Set**

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The 8-OHdG Assay Preparation Reagent Set contains nuclease, buffer and reagents which are used for the preparation of the assay of 8-OHdG (8-Hydroxy-2'-deoxyguanosine), a marker for oxidative stress. This set is able to reduce the variation of the enzyme reaction by addition of each reagent to the sample DNA in accordance with the protocol. Therefore, you are able to get stable assay results for a marker of oxidative stress.

【Features】

- 1) All reagents for the preparation of 8-OHdG assay are included.
- 2) Get stable assay results for a marker of oxidative stress with reduced variation of the reaction.

【Kit contents】

- | | |
|----------------------------------|------------------------|
| 1) Acetic Acid Buffer | 950 μ L \times 1 |
| 2) Nuclease P ₁ | 500 units \times 1 |
| 3) Tris Buffer | 1 mL \times 1 |
| 4) Alkaline Phosphatase Solution | 50 μ L \times 1 |

【Materials needed】

Reagent :

- 1) Deionized distilled water (DDW) : Code No. 316-90101

Equipments :

- 1) Centrifugal ultrafiltration device (10,000 M.W.) (e.g. : Sartorius VIVASPIN 500 code. VS0102)
- 2) High speed refrigerated microcentrifuge
- 3) Incubator (98°C, 37°C)

【Storage】

Store at -20°C.

【Package size】

50 reactions

【Method】

(Preparation)

Add 500 μ L of Nuclease P₁ (lyophilized form) to 500 μ L of DDW in order to prepare 1 unit/ μ L Nuclease P₁ Solution. Refrigerate the Nuclease P₁ Solution. The Nuclease P₁ Solution may be used for approx. 6 months after preparation.

(Protocol) Standard conditions

It is possible to change the concentration of DNA sample for assay of 8-OHdG by HPLC.

- (1) Heat the DNA sample which is prepared 200 μ g/150 μ L (= 1.33 mg/mL) at approx. 98°C for 2 minutes.
- (2) Chill in ice water and leave for 5-10 minutes.
- (3) Flash centrifuge.
- (4) Add 19 μ L of Acetic Acid Buffer and 10 μ L of Nuclease P₁ Solution.
- (5) Incubate at 37°C for 30 minutes.
- (6) Flash centrifuge.
- (7) Add 20 μ L of Tris Buffer and 1 μ L of Alkaline Phosphatase Solution. (Final concentration of DNA : 1 mg/mL)
- (8) Incubate at 37°C for 30 minutes.
- (9) Flash centrifuge.
- (10) Transfer to a filter cup of a Centrifugal ultrafiltration device (10,000 M.W.) (e.g. : Sartorius VIVASPIN 500 code. VS0102).
- (11) Centrifuge at 15,000 \times G at 4°C for approx. 20 minutes.
- (12) Collect the filtrate. Proceed to the assay of 8-OHdG.

The Oxidation of DNA advances early in the state of the filtrate after the treatment of Nuclease P₁ Solution. We recommend you to proceed to the assay of 8-OHdG as soon as possible, on the day if possible.

Flow chart

Standard conditions

200 μ g DNA/150 μ L DDW

|
98°C, 2 min.

|
Chill on ice, 5-10 min.

|
← 19 μ L **Acetic Acid Buffer**
← 10 μ L **Nuclease P₁ Solution**

37°C, 30 min.

|
← 20 μ L **Tris Buffer**
← 1 μ L **Alkaline Phosphatase Solution**

37°C, 30 min.

|
Centrifugal ultrafiltration (10,000 M.W.), 15K \times G, 4°C, 20 min.

|
Filtrate (1 mg/mL DNA)

|
8-OHdG assay (ELISA or HPLC)

It is possible to change the concentration of DNA sample for assay of 8-OHdG by HPLC.

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Code No.	Description	Package Size
307-07921	High Sensitive 8-OHdG Check	96 tests

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://www.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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