Code No. 296-69901 (200 cm²)
 292-69903 (1000 cm²)
 290-69904 (2000 cm²)

ImmunoStar[®] LD

ImmunoStar[®] LD is a chemiluminescence reagent which includes a unique luminol derivative L-012 as the substrate. It allows detection of targets with a high signal to noise ratio, minimizing background noise. The substrate produces a signal of long enough duration to be captured by image sensors.

[Kit Components]


	200 cm ²	1,000 cm ²	2,000 cm ²
Chemiluminescence Solution A	10 mL	50 mL	100 mL
Chemiluminescence Solution B	10 mL	50 mL	100 mL

[Storage]


2-10°C in the dark.

〈Protocol Summary〉

- ① Just before use, prepare Luminescence Working Solution by mixing Chemiluminescence Solution A and Chemiluminescence Solution B at a ratio of 1:1.


- ② Add the Luminescence Working Solution onto a plastic wrap by pipette. The volume needed for 100 μL/cm².


- ③ Place a blot protein side down with the membrane on the wrap.

- ④ Incubate the membrane for a few seconds to 5 minutes.


- ⑤ Drain excess the Luminescence Working Solution.

— 1/8 —

- ⑥ Cover the membrane with a clean wrap and remove any air bubbles.

- ⑦ Expose the wrapped membrane to X-ray film or image analyzer.

[Protocol]

Membrane Preparation

1. Prepare a protein blotted membrane by electrophoretic transfer or dot blotting.
2. Introduce the membrane into 50 mL of blocking buffer, and incubate for 1 hour at room temperature (RT) or incubate overnight at a cool place.
3. Wash the membrane with 50 mL of wash Solution (e.g. PBS-T or TBS-T) with agitation for 10 minutes, followed by 5 minutes washing, two or three times.
4. Prepare primary antibody with antibody diluent in an appropriate concentration and incubate the membrane with the antibody solution for more than 1 hour at RT or overnight at a cool place.
5. Wash the membrane with 50 mL of wash Solution with agitation for 10 minutes, followed by 5 minutes washing, two or three times.
6. Prepare secondary antibody with antibody diluent in an appropriate concentration and incubate the membrane with the antibody solution with agitation for 1 hour at RT.
7. Wash the membrane with 50 mL of wash Solution with agitation for 10 minutes, followed by 5 minutes washing, two or three times.

Detection Procedure

1. Just before use, prepare Luminescence Working Solution by mixing Chemiluminescence Solution A and Chemiluminescence Solution B at a ratio of 1:1.
2. Add the Luminescence Working Solution onto a plastic wrap by a pipette. The volume needed for 100 μL/cm².
3. Place a blot membrane side down with the membrane on the wrap. Be sure to touch the entire protein side with the Luminescence Working Solution on the wrap.
4. Incubate the membrane in the Luminescence Working Solution for a few seconds to 5 minutes.
5. Drain excess the Luminescence Working Solution from the membrane. Use absorbent tissue to remove the excess liquid.
6. Cover the membrane with a clean wrap or polyethylene sheet and remove any air bubbles.
7. a) CCD imager protocol
 - ① Place the membrane on an imager tray with protein side up.
 - ② Place the tray in the CCD imager.
 - ③ Operate the imager as recommended by manufacturer.
- b) X-ray film protocol
 - ① Place the membrane in a film cassette with protein side up.
 - ② Operate the X-ray film processor as recommended by manufacturer in the dark room using red safe lights. A recommended first

— 2/8 —

exposure time is 15 or 60 seconds as a test to determine optimal exposure time.

[Notes]

- Wash solution : Select wash solution A or B depending on your purpose.
Wash Solution A : PBS-0.1% Tween 20
Wash Solution B : TBS-0.1% Tween 20
- The optimal antibody concentration can be varied depending on the antibody titer to be used. Perform a systematic dot blot analysis to determine the appropriate concentration.
- Immediately before use, prepare blocking buffer, antibody solution and Luminescence Working Solution.
- The signal cannot be sufficiently detected in streptavidin-biotin system according to circumstances.

[Trouble Shooting]

Problems	Possible causes	Solutions
Weak or no signal	Insufficient protein was loaded to the gel.	Increase the amount of protein loaded to the gel.
	Transfer efficiency to a membrane was poor.	Stick the gel, membrane and blotting papers each other. Use fresh transfer buffer. Check protein transfer by staining the gel and membrane.
	The concentration of antibody is too low.	Increase the amount of antibody.
	Incubation time of antigen-antibody reaction was too short.	Increase incubation time of antigen-antibody reaction.
	The Luminescence Working Solution activity is reduced.	Immediately before use, prepare the Luminescence Working Solution.
	Exposure time was too short.	Increase exposure time.
Excessive signal	Too much protein of interest was loaded to the gel.	Decrease the amount of protein loaded to the gel.
	The concentration of antibody was too high.	Decrease the amount of antibody.
	The incubation time of antigen-antibody reaction was too long.	Shorten incubation time of antigen-antibody reaction.
Spots in the protein bands	Too much protein was loaded on the gel and/or the too high concentration of antibody was used.	Decrease the amount of protein loaded on the gel and/or lower the concentration of antibody. It is recommended to perform a systematic dot blot analysis.
High background	The exposure time is too long.	Shorten the exposure time.
	The concentration of antibody was too high.	Decrease the amount of antibody.
	Incubation time of antigen-antibody reaction was too long.	Shorten incubation time of antigen-antibody reaction.
	The volume of antibody solution was inadequate.	Increase the volume of antibody solution to cover the whole membrane and agitate the membrane in antibody solution sufficiently.
	Blocking was insufficient.	Increase blocking duration, numbers or volume of washes. Use a different blocking reagent.
	Washing was insufficient.	Add additional washing steps or increase washing duration.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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