## CTGF (Full+N-terminal region) ELISA Kit Wako

### [1. Introduction]

CTGF is a secretory protein of approximately 38 kDa produced by umbilical vein and vascular endothelial cells. CTGF is composed of 4 domains, designated Modules 1 through 4. Module 1 is the IGF-binding domain, Module 2 the VWC domain, Module 3 the TSP-1 domain, and Module 4 the CT domain, and different factors bind to each domain. CTGF is reportedly involved in the proliferation and differentiation of chondrocytes and in cell adhesion, as well as being a major factor in tissue fibrosis, and has been proposed as a candidate biomarker for various fibrotic disorders.

Meanwhile, it has been reported that CTGF in blood contains an N-terminal region in which the C-terminal region has been cleaved, giving rise to the problem that the amount of N-terminal CTGF in blood cannot be measured because it is mixed with platelet-derived full-length CTGF at the time of blood sample collection<sup>1)</sup>. N-terminal CTGF has been proposed as a potential biomarker for pulmonary fibrosis<sup>2)</sup>.

This product is an ELISA kit for detection of CTGF (Full) + CTGF (N-terminal region) using monoclonal antibodies that recognize CTGF Module 1 and Module 2. The amount of N-terminal CTGF in blood can be measured when this product is used in combination with a kit for detecting CTGF (Full) (Code No. 290-84701 CTGF (Full) ELISA Kit Wako)<sup>1)</sup>.

1) Miyazaki, O. *et al.* : *Ann. Clin. Biochem.*, **47**, 205 (2010).

2) Kono, M. *et al.* : *Clin. Chim. Acta*, **412**, 2211 (2011).

### [2. Kit performance]

Calibration curve range	7.81 to 500 pM
Analyte	CTGF (Full+N-terminal region)
Sample	Human serum, human plasma (EDTA)
Sample amount required	Human serum, 5 $\mu$ L Human plasma (EDTA), 10 $\mu$ L
Assay time	2 hours 50 minutes
Detection method	Colorimetry

### [3. Materials supplied]

Components	State	Amount
Antibody-coated Plate	Use after washing	96 wells (8 $\times$ 12)/plate
CTGF Standard	Lyophilized product, use after reconstitution	1 vial
Buffer	Ready-to-use	20 mL/1 vial
Biotin-conjugated Antibody Solution	Use after preparation	100 $\mu$ L/1 vial

Peroxidase-conjugated Streptavidin Solution	Use after preparation	100 $\mu$ L/1 vial
TMB Solution	Ready-to-use	12 mL/1 vial
Stop Solution	Ready-to-use	12 mL/1 vial
Wash Solution (10 $\times$ )	Use after preparation	100 mL/1 vial
Standard Sample Buffer	Ready-to-use	100 mL/1 vial
Plate Seal	Ready-to-use	4 sheets
Instruction Manual	—	1 copy

### [4. Assay principle]

Anti-CTGF Module 1 antibody is coated in each well of the assay plate. The standard solution or sample is added to these wells and allowed to react with the antibody. Next, biotin-conjugated anti-CTGF Module 2 antibody is added to the wells, followed by peroxidase-conjugated streptavidin. Finally, the concentration of CTGF (Full+N-terminal region) in the sample can be determined by measuring the peroxidase activity in the wells.

### [5. Apparatus, equipment, and materials required]

- ☐ Purified water (distilled water)
- ☐ Test tubes for dilution of standard solution and samples
- ☐ Glass utensils for dilution of Wash Solution (graduated cylinder, beaker)
- ☐ Pipettes with disposable tips (one capable of pipetting 10  $\mu$ L of liquid accurately and one capable of pipetting 200 to 500  $\mu$ L)
- ☐ Repeater pipette, capable of repeatedly dispensing 100  $\mu$ L
- ☐ Water-absorbent material such as a paper towel (to remove any solution remaining on the plate after washing)
- ☐ Vortex-type mixer
- ☐ Microplate shaker (range approx. 500 to 800 rpm)
- ☐ Automatic washer for 96-well plate (if available) or washing bottle
- ☐ Microplate reader capable of measuring at  $450 \pm 10$  nm, with the correction wavelength set at 600 to 650 nm
- ☐ Software for data analysis

### [6. Reagent preparation]

Equilibrate the reagents in this kit to room temperature (20°C to 25°C) before use (over about 2 hours).

#### 6-1. Preparation of Standard Solution

Reconstitute the CTGF Standard with a separately specified\* volume of purified water to prepare the standard stock solution (5,000 pM). Next, add the room-temperature Standard Sample Buffer included in this kit to the standard stock solution to prepare the standard solutions at the concentrations below.

\*As the volume of purified water to be added will differ depending on the lot, see the volume specified in the accompanying sheet. An example follows below.

Concentration (pM)	Volume of Standard Solution	Standard Sample Buffer
500	Standard Stock Solution : 30 $\mu$ L	270 $\mu$ L
250	500 pM solution : 150 $\mu$ L	150 $\mu$ L
125	250 pM solution : 150 $\mu$ L	150 $\mu$ L
62.5	125 pM solution : 150 $\mu$ L	150 $\mu$ L
31.3	62.5 pM solution : 150 $\mu$ L	150 $\mu$ L
15.6	31.3 pM solution : 150 $\mu$ L	150 $\mu$ L
7.81	15.6 pM solution : 150 $\mu$ L	150 $\mu$ L
0 (Blank)	—	150 $\mu$ L

#### 6-2. Biotin-conjugated Antibody Solution

Dilute 100-fold with buffer.

#### 6-3. Peroxidase-conjugated Streptavidin Solution

Dilute 100-fold with buffer.

#### 6-4. Wash Solution (10 $\times$ )

Dilute 10-fold with purified water (distilled water).

Example for entire 96-well microplate : Add 100 mL of Wash Solution (10 $\times$ ) to 900 mL of purified water (distilled water).

○ Other reagents are ready to use.

### [7. Stability and storage of each reagent]

#### 7-1. Antibody-coated Plate

Unused antibody-coated strips should be returned into the zip-lock bag provided in the kit and stored at 2°C to 10°C. These are stable until the expiration date.

#### 7-2. CTGF Standard

Store reconstituted standard solution, which is standard stock solution 5,000 pM, at 2°C to 10°C and use up within 1 week. Use the standard solution diluted to each concentration immediately after preparation and do not store any remaining solution.

#### 7-3. Buffer

If only part of the Buffer is to be used, transfer slightly more than the volume needed to another container, immediately close the cap tightly without equilibrating the remaining Buffer to room temperature, and store at 2°C to 10°C. It is stable until the expiration date.

#### 7-4. Biotin-conjugated Antibody Solution

If the kit is divided for multiple assays, prepare these solutions by dilution after removal from the refrigerator, immediately close the cap tightly without equilibrating the remaining stock solution to room temperature, and store at 2°C to 10°C. It is stable until the expiration date. Discard any remaining diluted solution after use.

#### 7-5. Peroxidase-conjugated Streptavidin Solution

If the kit is divided for multiple assays, prepare these solutions by dilution after removal from the refrigerator, immediately close the cap tightly without equilibrating the remaining stock solution to room temperature, and store at 2°C to 10°C. It is stable

until the expiration date. Discard any remaining diluted solution after use.

#### 7-6. TMB Solution

If only part of the TMB Solution is to be used, transfer slightly more than the volume needed to another container, immediately close the cap tightly without bringing the remaining TMB Solution to room temperature, and store at 2°C to 10°C. It is stable until the expiration date.

#### 7-7. Stop Solution

Store any unused Stop Solution at 2°C to 10°C with the cap tightly fastened. It is stable until the expiration date.

#### 7-8. Wash Solution (10 $\times$ )

Store Wash Solution (10 $\times$ ) with the cap tightly closed at 2°C to 10°C. It is stable until the expiration date. Discard any remaining diluted Wash Solution.

#### 7-9. Standard Sample Buffer

If only part of the Standard Sample Buffer is to be used, transfer slightly more than the volume needed to another container, immediately close the cap tightly without equilibrating the remaining Standard Sample Buffer to room temperature, and store at 2°C to 10°C. It is stable until the expiration date.

### [8. Sample preparation]

Serum : Dilute 31-fold with the Standard Sample Buffer contained in the kit and assay.

Plasma (EDTA) : Dilute 11-fold with the Standard Sample Buffer contained in the kit and assay.

### [9. Assay procedure]

- (1) Remove the protective solution and fill each well with Wash Solution prepared beforehand, and wash 3 times. Next, invert the plate and gently blot it with a paper towel or similar product to remove any liquid remaining in the wells.
- (2) Add 50  $\mu$ L of standard solution at each concentration to the wells for assay of the standard.
- (3) Add 50  $\mu$ L of sample solution diluted with the Standard Sample Buffer to the wells for assay of the sample.
- (4) Agitate the plate on a microplate shaker or similar device.
- (5) Cover with a plate seal, and allow to stand for 1 hour at room temperature (20°C to 25°C).
- (6) After the reaction is complete, discard the reaction solution and fill each well with wash solution and wash 3 times. Next, invert the plate and gently blot it with a paper towel or similar product to remove any liquid remaining in the wells.
- (7) Add 50  $\mu$ L of Biotin-conjugated Antibody Solution to each well. Agitate the plate on a microplate shaker or similar device.
- (8) Cover with a plate seal, and allow to stand for 1 hour at room temperature (20°C to 25°C).
- (9) After the reaction is complete, repeat the washing procedure

in Step 6.

- (10) Add 50  $\mu\text{L}$  of Peroxidase-conjugated Streptavidin Solution to each well.  
Agitate the plate on a microplate shaker or similar device.
- (11) Cover with a plate seal, and allow to stand for 30 minutes at room temperature (20°C to 25°C).
- (12) After the reaction is complete, repeat the washing procedure in Step 6.
- (13) Add 100  $\mu\text{L}$  of TMB Solution to each well.  
Agitate the plate on a microplate shaker or similar device.
- (14) Cover with a plate seal, and allow to stand for 20 minutes at room temperature (20°C to 25°C).
- (15) Add 50  $\mu\text{L}$  of Stop Solution to each well to stop the color development reaction.
- (16) Agitate the plate, then measure the absorbance at 450 nm (sub-wavelength : 620 nm) using a microplate spectrophotometer. A range of 600-650 nm can also be used for the sub-wavelength.

#### [10. Calculation method]

Prepare a standard curve by plotting the absorbance readings on the Y-axis against the concentrations of the standard solution (pM) on the X-axis\*\*.

Read the unknown sample concentration (pM) corresponding to the absorbance of the diluted sample from the standard curve.


Multiply this concentration by the dilution factor to determine the assay value.

To determine the CTGF (N-terminal region) concentration, obtain the measured value with the CTGF (Full + N-terminal region) Kit *Wako* and subtract the measured value of CTGF (Full) from the measured value of CTGF (Full + N-terminal region).  
Calculation formula : CTGF (N-terminal region) concentration (pM) = CTGF (Full + N-terminal region) concentration (pM) - CTGF (Full) concentration (pM)

\*\*For calculations using computer software, we recommend the use of a cubic polynomial, and 4 or 5 parameters.

#### [11. Typical standard curve]

Standard curve (example)

— 5/16 —

#### [12. PRECAUTIONS]

1. This kit should only be used by persons who have completed training in ELISA techniques, or under the direction of an instructor.
2. When performing the assay using manual procedures, ensure that the reproducibility of the pipetting is stable.
3. Wear protective gloves, eyewear, and clothing during preparation and when using this kit.
4. Avoid contact of skin with kit reagents. If any reagents come in contact with the eyes, mouth, wounds, skin, etc., immediately wash with copious amounts of tap water and seek treatment from a physician if necessary.
5. Do not drink, eat, or smoke in any location where this kit is used.
6. Handle the sample with thoroughgoing care, being aware that the sample may present an infection risk. This kit contains animal-derived ingredients.
7. Used samples and other used consumable items should be immersed in 1% formalin, 2% glutaraldehyde or more than 0.1% sodium hypochlorite solution for more than 1 hour or be sterilized by autoclaving before disposal. Dispose of used consumable items and unused reagents in accordance with the applicable laboratory rules and local regulations.
8. Do not use reagents with different lot numbers in the same assay procedure. When allowing the plate to stand in each step, ensure that a plate seal is applied to prevent the wells drying out, contamination with foreign matter, deviations in temperature, and evaporation of the dispensed reagents.
9. The ELISA is affected by measurement environment. Strictly adhere to the measurement operations and maintain the room temperature at the place of static reaction (20 to 25°C as the temperature on the laboratory table or in the incubator). Also, avoid taking measurements in a windy, low-humidity environment (including air conditioner wind). Do not place a heat source (e.g., personal computer, incubator) near the reaction plate.

#### [13. Assay procedure summary]

Ensure that you read the instruction manual and check the sample conditions, assay conditions, and assay procedure before beginning the assay.

- |                          |   |
|--------------------------|---|
| <input type="checkbox"/> | Thoroughly equilibrate the plates and reagents to room temperature (20°C to 25°C) before use (approx. 2 hours).   |
| <input type="checkbox"/> | Dilution of Wash Solution (10×) : Dilute 10-fold with purified water brought to room temperature.   |
| <input type="checkbox"/> | Preparation of standard solution (example) : Reconstitute CTGF Standard with the amount of purified water specified in the accompanying sheet*** to prepare the Standard Stock Solution (5,000 pM).<br>Next, prepare the standard solution with the room-temperature Standard Sample Buffer supplied with the kit. Refer to the accompanying sheet for the amount of purified water to be added.<br>***The amount of purified water to be added varies depending on the lot, so first check the amount specified on the accompanying sheet. An example follows below. |

— 6/16 —

Dilution example	Concentration (pM)	500	250	125	62.5	31.3	15.6	7.81	0
	Standard Solution (μL)	30	150*	150*	150*	150*	150*	150	—
	Standard Sample Buffer (μL)	270	150	150	150	150	150	150	150

\* : Standard solution at 1 step higher concentration

<input type="checkbox"/>	<b>Antibody-coated Plate</b>
<input type="checkbox"/>	↓ Wash 3 times (*1)
<input type="checkbox"/>	<b>Sample or Standard Solution</b> 50 μL/well
<input type="checkbox"/>	↓ Agitate (*2), then allow to stand for 1 hour at room temperature (20°C to 25°C) (*3).
<input type="checkbox"/>	*Preparation of Biotin-conjugated Antibody Solution (Dilute 100-fold with room-temperature buffer.)
<input type="checkbox"/>	↓ Wash 3 times (*1)
<input type="checkbox"/>	<b>Biotin-conjugated Antibody Solution</b> 50 μL/well
<input type="checkbox"/>	↓ Agitate (*2), then allow to stand for 1 hour at room temperature (20°C to 25°C) (*3).
<input type="checkbox"/>	*Preparation of Peroxidase-conjugated Streptavidin Solution (Dilute 100-fold with room-temperature buffer)
<input type="checkbox"/>	↓ Wash 3 times (*1)
<input type="checkbox"/>	<b>Peroxidase-conjugated Streptavidin Solution</b> 50 μL/well
<input type="checkbox"/>	↓ Agitate (*2), then allow to stand for 30 minutes at room temperature (20°C to 25°C) (*3).
<input type="checkbox"/>	↓ Wash 3 times (*1)
<input type="checkbox"/>	<b>TMB Solution</b> 100 μL/well
<input type="checkbox"/>	↓ Agitate (*2), then allow to stand for 20 minutes at room temperature (20°C to 25°C) (*3).
<input type="checkbox"/>	<b>Stop Solution</b> 50 μL/well
<input type="checkbox"/>	↓ Agitate (*2)
<input type="checkbox"/>	<b>Measure the absorbance (main wavelength, 450 nm ; sub-wavelength, 620 nm ; 600 to 650 nm)</b>

(\*1) After dispensing the washing solution into the wells for each wash, shake the plate gently in the palm of your hand for 10 seconds to discard. After 3 consecutive washes, blot the plate upside down on a paper towel to completely remove the wash solution. Following removal of the wash solution, immediately dispense the next solution, ensuring that the wells do not dry out. The recommended volume of washing solution to be added when using a pipette is 300 μL/well.

(\*2) Three 10-second cycles of agitation at 600 to 800 rpm are recommended.

(\*3) After agitation, cover with a plate seal and allow to stand. Peel off the protective paper and affix the plate seal with the adhesive side facing the plate. Do not reuse used plate seals.

#### [Storage]

Store at 2°C to 10°C

#### [Expiration date]

Indicated on the label

#### [Packaging]

For 96 assays

### FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
 Telephone : + 81-6-6203-3741  
 Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

#### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
 Richmond, VA 23237  
 U.S.A.  
 Telephone : + 1-804-271-7677  
 Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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