Human Chromogranin A ELISA Kit *Wako* «For Research Use Only»

[1. Introduction]

Chromogranin A (CgA) is an acidic glycoprotein originally isolated from the chromaffin granules of the adrenal medulla. Human CgA consists of 439 amino acids. CgA is widely distributed in the endocrine and nervous systems, with particularly high levels in the adrenal medulla and pituitary gland. CgA coexists and is co-released with catecholamines and correlates with the secretion of catecholamines in the blood. Thus, CgA can be used as a biomarker of the activity of the sympathetic-adrenal system, making immunodetection of CgA in the blood extremely important. Blood CgA levels in patients with neuroendocrine tumors, especially pheochromocytomas and pituitary tumors, are markedly higher than in healthy individuals. CgA is also present in the submandibular duct and is released into saliva upon stimulation of the autonomic nervous system, making salivary CgA a new biomarker of mental stress. Thus, plasma CgA can be used as a biomarker of tumor, and salivary CgA as a biomarker of mental stress.

This ELISA kit is designed for the measurement of CgA in human plasma and saliva.

[2. Kit Performance]

Calibration curve range	0.14-33.33 pmol/mL
Sample	Human plasma, human saliva
Sample amount required	50 μ L (n=2)
Assay time	O/N + 2.5 hours
Intra-assay CV (%)	CV < 14%
Inter-assay CV (%)	CV < 16%

[3. Materials supplied]

Components	State	Amount
Antibody-coated Plate	Use after washing	96 wells (8 \times 12)/plate
Chromogranin A Standard	Lyophilized product, use after reconstitution	1 vial
Biotin-conjugated Chromogranin A	Lyophilized product, use after reconstitution	1 vial
Anti-Human Chromogranin A Antibody	Lyophilized product, use after reconstitution	1 vial
Peroxidase-conjugated Streptavidin Solution	Ready-to-use	12 mL \times 1 vial
Substrate Buffer	Use after preparation	25 mL \times 1 vial
OPD tablet	Use after preparation	2 tablets
Stop Solution	Ready-to-use	12 mL \times 1 vial
Buffer (5 \times)	Use after preparation	12 mL \times 1 vial

Components	State	Amount
Wash Solution (20 \times)	Use after preparation	50 mL \times 1 vial
Plate Seal	Ready-to-use	3 sheets
Instruction Manual	—	1 copy

[4. Assay Principle]

Goat anti-rabbit IgG antibody is immobilized in each well of the assay plate. The standard solutions or samples, labeled antigen (Biotin-conjugated Chromogranin A), and Anti-Human Chromogranin A Antibody are sequentially added in the wells to allow for competitive binding. Next, Peroxidase-conjugated Streptavidin Solution is added. Finally, the concentrations of CgA in the samples can be determined by measuring the peroxidase activity in the wells.

Human CgA* synthetic peptide is used as the antigen.

*Human CgA (344-374) :

E-E-E-D-N-R-D-S-S-M-K-L-S-F-R-A-R-A-Y-G-F-R-G-P-G-P-Q-L-R-R

[5. Supplies and Equipment]

- Micropipettes and disposable tips (25 to 1,000 μ L) (8- or 12-channel pipettes are recommended.)
- Microplate reader which can read extinction 2.5 at 490 nm
- Microplate shaker or conventional shaker
- Test tubes for dilution of standard solutions and samples (glass or polypropylene)
- Microplate washer (For manual operation, a continuous dispenser, needle dispenser, and an aspirator, or vacuum pump are recommended)
- Graduated cylinder (1,000 mL)
- Distilled or deionized water

[6. Reagent preparation]

The reagents in the kit must be brought to room temperature (20-30°C) before use (about 2 hours).

6.1 Preparation of Standard Solution

Add 1 mL of the Buffer to the vial containing the standard and dissolve.

Take 0.1 mL of the standard stock solution and dilute with 0.2 mL of the Buffer to make a 33.33 pmol/mL standard solution. Repeat the same dilution procedure to prepare standard solutions of 11.11, 3.70, 1.23, 0.41, and 0.14 pmol/mL. For 0 pmol/mL standard solution, use the Buffer.

Concentration (pmol/mL)	Volume of standard solution	Buffer
33.33	Standard stock solution : 0.1 mL	0.2 mL
11.11	33.33 pmol/mL standard solution : 0.1 mL	0.2 mL
3.70	11.11 pmol/mL standard solution : 0.1 mL	0.2 mL
1.23	3.70 pmol/mL standard solution : 0.1 mL	0.2 mL
0.41	1.23 pmol/mL standard solution : 0.1 mL	0.2 mL
0.14	0.41 pmol/mL standard solution : 0.1 mL	0.2 mL
0 (Blank)	—	0.2 mL

6.2 Biotin-conjugated Chromogranin A

Add 6 mL of distilled water to the vial containing Biotin-conjugated Chromogranin A and dissolve.

6.3 Anti-Human Chromogranin A Antibody

Add 12 mL of distilled water to the vial containing Anti-Human Chromogranin A Antibody and dissolve.

6.4 Substrate Solution

Add one OPD tablet to 12 mL of the Substrate Buffer and dissolve.

Note : Prepare immediately before the end of the 2-hour incubation after the addition of the Peroxidase-conjugated Streptavidin Solution.

6.5 Buffer (5×)

Dilute 10 mL of the Buffer (5×) with 40 mL of distilled water.

6.6 Wash Solution (20×)

Dilute 50 mL (the whole amount) of the Wash Solution (20×) with 950 mL of distilled water.

○ Other reagents are ready to use.

[7. Assay Procedure]


1. Bring the kit components to room temperature (20-30°C). Prepare Standard Solutions, Biotin-conjugated Chromogranin A, Anti-Human Chromogranin A Antibody, Buffer (5×) and Wash Solution (20×) as described above.
2. Add 350 μ L of the Wash Solution into each well and aspirate the solution. Alternatively, invert the plate, discard the buffer, and then remove the remaining buffer by blotting on a paper towel. Repeat this washing procedure further twice (total 3 times).
3. Add 50 μ L of the Buffer into each well, and add 25 μ L of standard solutions (0, 0.14, 0.41, 1.23, 3.70, 11.11, 33.33 pmol/mL) or samples. Add 50 μ L of Biotin-conjugated Chromogranin A and 100 μ L of Anti-Human Chromogranin A Antibody.
4. Cover with the plate seal and incubate overnight (16-20 hours) at room temperature with shaking (approximately 100 rpm).
5. After the reaction is complete, discard the solution from each well and wash as in step 2, 3 times in total.
6. Add 100 μ L of Peroxidase-conjugated Streptavidin Solution into each well.
7. Cover with the plate seal and incubate for 2 hours at room temperature with shaking (approximately 100 rpm).
8. Immediately before the end of the incubation in step 7, dissolve an OPD tablet in the Substrate Buffer to prepare the Substrate Solution.
9. Discard the solution from each well and wash as in step 2, 4 times in total.
10. Add 100 μ L of the Substrate Solution into each well and allow to stand in dark place for 30 minutes at room temperature.
11. Add 100 μ L of the Stop Solution to each well.

12. Measure the absorbance at 490 nm using a microplate reader.

13. Using commercially available software, create a standard curve from the measured values of each CgA standard solution using a 4 (or 5)-parameter or Log-Logit regression equation and determine CgA concentrations of the samples.

When using semi-log graph paper, plot the concentration of the standard solution on the x axis (log scale) and the absorbance of the standard solution at each concentration on the y axis (linear scale) to create a standard curve. Determine the CgA concentration that corresponds to the absorbance of each sample using the standard curve.

[8. Typical Standard Curve]

[9. Precautions]

1. Plasma and saliva samples should be assayed immediately. If samples cannot be assayed immediately after collection, they should be divided into small portions and stored frozen at -30°C or below. Avoid repeated freezing and thawing of samples.
2. In general, reagents should be freshly prepared. In particular, the standards, Biotin-conjugated Chromogranin A, and Anti-Human Chromogranin A Antibody should be used immediately after preparation. If the kit is used in portions, the prepared standards, Biotin-conjugated Chromogranin A, and Anti-Human Chromogranin A Antibody should be divided into small portions and stored frozen at -30°C or below.
3. Precipitates may form in the Wash Solution (20×) during storage. The precipitates will dissolve with dilution, and the quality will not be affected.
4. When diluting standard solution, be sure to use a new tip for each dilution step. Dispensing to each well should be performed accurately, as it affects the accuracy of the assay. When adding samples into wells, use a new tip for each sample and take care to avoid cross-contamination.
5. If the concentration of a sample exceeds 33.33 pmol/mL, dilute the sample with the Buffer in the kit and use for the assay.
6. Since the reaction efficiency of each reaction step will decrease if the assay plate is left standing, be sure to use a microplate shaker to shake the plate during the incubation (except for the chromogenic reaction). Shaking should be done slowly (about 100 rpm) to prevent splashing of the reaction solution on the plate seal.

7. All measurements should be performed in duplicate.
8. Absorbance should be measured immediately after the reaction is stopped.
9. The coloration level of the substrate may be slightly affected by the reaction temperature, time, degree of shaking of the assay plate, etc. Therefore, a standard curve must be prepared for each assay.
10. Avoid exposure to strong light during storage or use of each reagent.
11. Do not combine kits from different lots.

[10. Measurement Examples]

Spiked Recovery Test

	Spiked CgA (pmol/mL)	Assay value (pmol/mL)	Recovery (%)
Plasma	0	0.54	—
	0.25	0.86	108.86
	1	1.93	125.32
	4	6.58	144.93

	Spiked CgA (pmol/mL)	Assay value (pmol/mL)	Recovery (%)
Saliva	0	0.47	—
	0.25	0.81	112.50
	1	1.40	95.24
	4	3.66	81.88

[11. Assay Procedure Summary]

Ensure that you read the instruction manual and check the sample conditions, assay conditions, and assay procedure before beginning the assay.

<input type="checkbox"/>	Thoroughly equilibrate the plates and reagents to room temperature (20℃ to 30℃) before use (approx. 2 hours)																								
<input type="checkbox"/>	Preparation of the reagents : Biotin-conjugated Chromogranin A : Add 6 mL of distilled water to dissolve the content. Anti-Human Chromogranin A Antibody : Add 12 mL of distilled water to dissolve the content. Buffer (5×) : Dilute 10 mL of the Buffer (5×) with 40 mL of distilled water. Wash Solution (20×) : Dilute 50 mL (the whole amount) of the Wash Solution (20×) with 950 mL of distilled water.																								
<input type="checkbox"/>	Preparation of standard solutions (example) : Add 1 mL of the Buffer to the vial containing CgA standard to dissolve the content. Dilute with the Buffer as below to prepare standard solutions of 33.33, 11.11, 3.70, 1.23, 0.41, and 0.14 pmol/mL.																								
Dilution example	<table><tr><td>Concentration (pmol/mL)</td><td>33.33</td><td>11.11</td><td>3.70</td><td>1.23</td><td>0.41</td><td>0.14</td><td>0</td></tr><tr><td>Standard solution (mL) : Original</td><td>0.1</td><td>0.1*</td><td>0.1*</td><td>0.1*</td><td>0.1*</td><td>0.1*</td><td>—</td></tr><tr><td>Buffer (mL)</td><td>0.2</td><td>0.2</td><td>0.2</td><td>0.2</td><td>0.2</td><td>0.2</td><td>0.2</td></tr></table>	Concentration (pmol/mL)	33.33	11.11	3.70	1.23	0.41	0.14	0	Standard solution (mL) : Original	0.1	0.1*	0.1*	0.1*	0.1*	0.1*	—	Buffer (mL)	0.2	0.2	0.2	0.2	0.2	0.2	0.2
Concentration (pmol/mL)	33.33	11.11	3.70	1.23	0.41	0.14	0																		
Standard solution (mL) : Original	0.1	0.1*	0.1*	0.1*	0.1*	0.1*	—																		
Buffer (mL)	0.2	0.2	0.2	0.2	0.2	0.2	0.2																		
	* : Standard solution from preceding step																								

<input type="checkbox"/>	Assay plate with immobilized antibody
<input type="checkbox"/>	↓ Wash 3 times.
<input type="checkbox"/>	Buffer 50 μL
<input type="checkbox"/>	Sample or Standard Solutions 25 μL
<input type="checkbox"/>	Biotin-conjugated Chromogranin A Solution 50 μL
<input type="checkbox"/>	Anti-Human Chromogranin A Antibody Solution 100 μL
<input type="checkbox"/>	↓ Cover with the plate seal and incubate overnight with shaking (16-20 hours, approximately 100 rpm) at room temperature.
<input type="checkbox"/>	↓ Wash 3 times.
<input type="checkbox"/>	Peroxidase-conjugated Streptavidin Solution 100 μL
<input type="checkbox"/>	↓ Cover with the plate seal and incubate for 2 hours at room temperature with shaking (approximately 100 rpm). Immediately before the end of incubation, prepare Substrate Solution (dissolve one OPD tablet in 12 mL of the solvent for substrate).
<input type="checkbox"/>	↓ Wash 4 times.
<input type="checkbox"/>	Substrate Solution 100 μL
<input type="checkbox"/>	↓ Allowed to stand for 30 minutes at room temperature.
<input type="checkbox"/>	Stop Solution 100 μL
<input type="checkbox"/>	Absorbance measurement (primary wavelength at 490 nm)

[Storage]

Store at 2-10°C

[Expiration date]

Indicated on the label

[Package]

96 reactions

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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E-E-E-E-D-N-R-D-S-S-M-K-L-S-F-R-A-R-A-Y-G-F-R-G-P-G-P-Q-L-R-R

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