Cytotoxic Fluoro test *Wako*

This test provides a simple protocol and the necessary reagents to measure viability of cells treated with various cytotoxic substrates. It is based on fluorometric assay on cell lysis and screening (FACLS). Proliferation and aging of cells can also be estimated by this procedure. The test, which consisted of a unique fluorescent dye, a cell lysis reagent, utilizes properties of the fluorescent dye that differentiates between intact cells and damaged cells, and reacts with DNA to become fluorescent. The dye permeates only damaged or dead cell membrane, and reacts with DNA in the cells. The intensity of fluorescence emissions of the DNA-dye complex before and after addition of the cell lysis reagent to the cell sample represents quantities of damaged cells and total cells applied to the assay, respectively. The procedure is applicable to adhesive cells as well as floating cells.

The followings are fields to be applicable;

- Cytotoxicity assays
- Apoptosis and necrosis, determination of ratio of dead cells and living cells
- Cell proliferation and aging


[Characteristics]

1. Determination of viability of cells by percentage of intact and / or living cells per total cells due to unique properties of the dye
2. A simple 2-step procedure, performed in about 5 minutes, applicable to simultaneous screening assay for many samples
3. Cell pre-treatments such as washing of cells, exchange of media, and cell separation by centrifugation are not necessary. Difficult preparation of dye solution like MTT and XTT are avoided.
4. Sensitive and accurate measurement based on cell membrane permeability and the DNA-affinity of the unique fluorescent dye, avoiding interference with changes of pH by addition of drug, temperature, reaction time and serum containing in culture medium, all of which are inevitably associated with use of enzyme reaction in the other assays.
5. No hazardous wastes are generated from use, disuse of RI, toxic organic solvent, or dye.

[Principle of the test]

The fluorescence reagent has unique properties. The dye emits intense fluorescence at 100-1,000 times of the original intensity only when the dye binds to DNA. In addition, the dye enters into cells with damaged cell membranes or dead cells, but not intact or living cells. Thus, fluorescence emission of damaged and dead cells can be measured by staining their DNA with the dye. In the assay system, therefore, one can estimate.

- 1) Damaged cell number by measuring fluorescence intensity [F_D] generated after addition of the dye to the cell sample.
- 2) The number of total cells in the assay by measuring fluorescence intensity [F_T] after subsequent addition of the cell lysis reagent to the sample to disrupt membrane permeability of intact cells.
- 3) The number of intact or living cells by the fluorescence intensity resulting from subtracting [F_D] from [F_T].
- 4) Viability of the cell sample expressed by the percentage of intact cells per total cells.

[Contents]

Fluorescence Reagent	500 μ L	1 vial
Lysis Solution	10 mL	1 vial

※ It is highly recommended to use a known quantity of the reagent.
(Example : Fluorescence Reagent 50 μ L/100 assays)

[Storage] Store at -20°C

[Materials and apparatuses to be used]

1. Materials
 - Culture medium
 - Cells used to be Distilled assayed
 - Distilled water
2. Apparatuses
 - Multiplate reader, Infinite Series (Tecan, Austria)
 - Microdispensors and microtipes
 - Blood cell counter slide glasses
 - Cell culture dishes

[Procedures]

<Setting of Fluorescence/Absorbance multiplate reader in the case of Infinite Series, Tecan Austria>

- (1) Turn on the reader and personal computer connected to the reader. Open the software i-Control regulating the reader system.
- (2) Set assay conditions on the board as follows.

Road Mode	: fluorescence
Excitation	: 420 nm
Emission	: 460 nm
Read	: Top
Shake	: linear, high, 120 sec
Gain	: optimal
Delay Time	: 10
Lag Time	: 0
Integration Time	: 40
Flash Time	: 100

<Preparation of working fluorescence solution>

Add 50 μ L of Fluorescence Reagent to a vial containing 950 μ L diluent with distilled water, and vortex well. Storing on ice until use.

<Assay>

- (1) Dispense 200 μ L of the culture medium to be assayed after culture or treatment of drug to each of well and add 10 μ L of fluorescence working solution.
- (2) Set the plate to fluorescence microplate reader and measure the fluorescence intensity.
- (3) Draw the plate from the reader and add 10 μ L of Lysis Solution to each of the well assayed, gently mixing the suspension by finger tapping.
- (4) Measure the fluorescence intensity.

[Making the standard curve]

<Preparation of cell suspension>

- (1) Collect cells from the culture dish at late log phase of growth and transfer them to a 50mL-centrifuge tube.

- (2) Centrifuge the tube at 1,500 r.p.m. for 5 minutes at room temperature to precipitate the cells, and then carefully remove the supernate from the tube.
- (3) Suspend the cell pellet in 10 mL culture medium by repeats of gentle pipetting several times.
- (4) Repeat (2) followed by (3) with the cell suspension.
- (5) Count cell number of 10 μ L of the cell suspension using blood cell counter slide glass and calculate cell density (cells/mL) of the suspension.
- (6) Adjust cell density of the suspension by diluting the suspension with the culture medium to the desired cell density, such as 5×10^4 to 4×10^4 cells/mL. Usually, 3 mL of cell suspension is needed per plate.

<Dispensing dilution series of cell>

- (1) Flash airduster to the microplate for fluorescence measurement to remove the dust.
- (2) While gently agitating the cell suspension described in the above to distribute cells uniformly, dispense 400 μ L of the cell suspension to 6 wells coded as B11 to G11 in the plate, as shown in figure 1.

	1	2	3	4	5	6	7	8	9	10	11	12
A											400 μ L	
B												
C												
D												
E												
F												
G												
H												

Fig. 1

- (3) Dispense 200 μ L of the culture medium to each of the wells coded as B2-G10. for dilution, as shown in figure 2

	1	2	3	4	5	6	7	8	9	10	11	12
A		200 μ L/well										
B												
C												
D												
E												
F												
G												
H												

Fig. 2


- (4) Transfer 200 μ L of the cell suspension in the wells to adjacent wells coded as B10-G10.
- (5) Continue the dilution as described in (4) until the wells coded as B4-G4., and then remove 200 μ L of the last diluted suspensions in B4-G4 wells, which contain 400 μ L. The dilution series contains 8 density of cells, and the culture medium in the residual wells are for blanks, as shown in figure 3.

	1	2	3	4	5	6	7	8	9	10	11	12
A	cells	0	0	78.13	156.3	312.5	625	1250	2500	5000	10000	
B												
C												
D												
E												
F												
G												
H												

Fig. 3

<Measurement>

- (1) Add 10 μ L of fluorescence working solution to each of the wells to be assayed, and lightly tap the side of the microplate to distribute the dye uniformly.
- (2) Set the plate on the microplate reader and measure the fluorescence intensity.
- (3) Remove the plate from the reader and add 10 μ L of Lysis Solution to each of the wells to be assayed, gently mixing the suspension by tapping the side of the microplate.
- (4) Measure the fluorescence intensity.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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