Negative Gel Stain MS Kit

It is known that protein bands separated by SDS/polyacrylamide gel electrophoresis (SDS-PAGE) can be visualized as transparent bands against the background of milky white gel stained by the negative gel stain with a Zn/imidazol reagent. We have improved the method using a new imidazol derivative reagent, which allows a clear and stable image of protein bands on the gel as sensitive as that by the silver staining, with such the short times as 10 minutes. The staining technique is useful to obtain the clear and sensitive resolution pattern of the gel before immunoblotting as well as to excise and purify the band of interest from the gel without significant deterioration of amino acid residues for the subsequent studies of protein such as sequencing and mass analysis of peptide.

Features

1. Visualize protein bands on SDS-PAGE gel much more sensitively and simply with a few steps than by CBB R-250 staining.
2. Optimize the staining intensity by either de-staining the excess stain or re-staining following de-stain of the gel.

Contents

- | | | |
|-------------------------|--------|----------|
| 1. Staining Solution A | 500 mL | 1 bottle |
| 2. Staining Solution B | 500 mL | 1 bottle |
| 3. De-staining Solution | 500 mL | 1 bottle |

Materials and apparatuses to prepare

A transparent 50-100 mL-vessel for stain of the gel

A rotator shaker

A sheet of black paper

Procedure

1. Immerse a SDS-PAGE gel into 25mL of Staining Solution A in a transparent vessel for staining, and gently shake it by a rotator shaker for 5-10 minutes.
2. Drain off the staining solution from the vessel, and then rinse the gel for 1 minute three times with 50mL of distilled deionized water.
3. Add 25mL of Staining Solution B to the vessel, and shake it for 10-50 seconds, while checking staining intensity of the gel against a sheet of black paper.
4. Stop the staining by drain of the staining solution followed by wash of the gel with distilled deionized water for 1 minute three times as described above in 2. In case of excess staining of the gel, the gel could be de-stained with 2% EDTA3Na solution prepared until a desirable density, and washed with distilled deionized water as described above.

The image of stained gel rapped in a vinyl bag is stable for about 3 weeks at 2-10 °C.

When used in subsequent experiments such as Western Blotting and de-staining, gel is de-stained by the follow treatment.

5. Immerse the gel into 25mL of De-staining Solution provided followed by gentle stirring the vessel until the staining image reaches a desirable intensity. It takes about 15 minutes for complete removal of the staining, which can be used for the subsequent experiments such as Western Blotting as well as re-staining after wash of the gel as described above.

Notes

1. Be careful not to mix the solutions provided to each other, otherwise, the solutions are inactive or insensitive.
2. Wear the gloves on staining.

Reference

1. Fernandez-patron, *et al.*: *Anal. Biochem.*, **224**, 263 (1995).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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