<For research Use Only>

Code No. 293-71501 (40 reactions)

For apoptosis research

Apoptosis *in situ* Detection Kit Wako**1. The principle of Apoptosis *in situ* Detection Kit** Wako

Apoptosis *in situ* Detection Kit is designed for detection of apoptotic cells appearing in paraffin-embedded sections, frozen sections of tissue and neutralized formalin-fixed culture cells, enabling an easy and speedy assay.

The kit is based on TUNEL [Terminal deoxynucleotidyl Transferase(TdT)-mediated dUTP nick end labeling] procedure, that is the addition of fluorescein -dUTP to 3'-terminals of apoptotically fragmented DNA with TdT followed by immunochemical detection using anti-fluorescein antibody conjugated with horseradish peroxidase (POD) and DAB as a substrate.

2. Features

- Rapid detection can be performed.
The whole process from the de-paraffinizing step to the microscopic examination can be completed in about 2 hours.
- Complicated preparations of various reagent are not needed.
The kit contains the essential reagents required for detection of apoptosis.
- The kit shows a clear positive image with low background.

3. Applicable samples

- paraffin-embedded tissue sections
- frozen tissue sections
- neutralized formalin-fixed culture cells

The specimen should be mounted on a silane-coated glass slide, preferably pre-treated with a reagent for protection of detachment of the section from the glass slide.

4. The contents of the kit (approx. 1cm² × 40 reactions)

- | | |
|-----------------------------------|------------|
| (1) Protein Digestion Enzyme | 1 mL |
| (2) TdT | 40 μ L |
| (3) TdT Substrate Solution | 4.4 mL |
| (4) 100 × POD-Conjugated Antibody | 44 μ L |
| (5) DNase I | 4 μ L |
| (6) 10 × DNase I Reaction Buffer | 40 μ L |

5. Storage

Store the kit at -20 °C
Store the remaining components at indicated temperatures below.

Reagents (2), (3), (5), (6) : -20 °C
Reagents (1), (4) : 2~10 °C
For long term storage, Reagent (3) should be aliquoted and stored at -80 °C.


6. Apparatuses and reagents to be prepared

[Apparatuses]

- (1) 12 to 19 Coplin-staining jars (for protein digestion, counter staining, washing steps, etc.)
- (2) A moist chamber, shown below, which is pre-warmed with an incubator maintained at 37 °C before starting the experiment.

- (3) A water bath and an incubator maintained at 37°C
- (4) Vortex mixer
- (5) Microscope furnished with a camera.

Moist Chamber

[Reagents]

<For approximately 10 specimens using 100 mL Coplin-staining jar>

	Paraffin-embedded section	Cultured cell (or frozen section)
(1) Lemosol [®] , Lemosol [®] A or xylene	400 mL	200 mL
(2) Ethanol (for hydration and dehydration)	620 mL	260 mL
(3) 30% hydrogen peroxide solution	10 mL	10 mL
(4) PBS(-) at pH 7.2~7.4	1,300 mL	1,300 mL
(5) Counter staining solution (methyl green solution etc.)	100 mL	100 mL
(6) 99% 1-butanol (for a substitute of dehydration by ethanol)	200 mL	200 mL
(7) Permanent mounting medium	several mL	several mL
(8) 10% Formalin neutral buffer solution for fixation	—	necessary volume for sample numbers
(9) Chromogenic substrate (DAB Solution)	1 mL	1 mL

We recommend DAB Immunohistochemistry Substrate.

(Genway Biotech Inc. Code No. 517-88601)

This product includes Reagent A (3 mL) and Reagent B (100 mL).

7. Preparation of working solutions

(1) Preparation of the reagents for pre-treatment of the specimen

(a) Protease solution (for paraffin-embedded section)

Just before dehydration of specimen, dilute Protein Digestion Enzyme with distilled deionized water (DDW) in a Coplin-staining jar, and warm the jar in a water bath at 37°C.

<Working solution>

Protein Digestion Enzyme 50 μ L
DDW 100 mL

(b) Permeabilizing solution (for frozen section and formalin-fixed culture cells)

Prepare a permeabilizing solution of 0.1% sodium citrate and 0.1% Triton X-100 in a Coplin-staining jar, and cool it on ice.

<Working solution>

Sodium citrate, anhydrate 0.1 g
Triton X-100 0.1 mL
DDW 100 mL

(2) Preparation of TdT reaction solution

Prepare TdT reaction solution just before use, and store it on ice until use.

(a) For the standard procedure

<Working solution for a single assay>

TdT	1 μ L
TdT Substrate Solution	99 μ L

(b) For the enhanced procedure

<Working solution for a single assay>

TdT	1 μ L
TdT Substrate Solution	49 μ L

(3) Preparation of 3% hydrogen peroxide solution

Just before use, prepare 3% hydrogen peroxide solution in a Coplin-staining jar.

<Working solution>

30% hydrogen peroxide solution	10 mL
PBS(-)	90 mL

(4) Preparation of POD-conjugated antibody solution

Dilute 100 $\times$ POD-Conjugated Antibody with PBS(-) just before use, and store it on ice until use.

<Working solution for a single assay>

100 $\times$ POD-Conjugated Antibody	1 μ L
PBS(-)	99 μ L

(5) Preparation of Chromogenic Solution

Prepare DAB Solution just before use, and store it at room temperature until use.

When use DAB Immunohistochemistry Substrate (Genway Biotech Inc. Code No. 517-88601)

<Working solution for a single assay>

Reagent A	3 μ L
Reagent B	100 μ L

(6) Preparation of DNase I Solution for the positive control

Prepare DNase I solution just before use, and store it on ice until use.

<Working solution for a single assay>

10 $\times$ DNase I Reaction Buffer	10 μ L
DNase I	1 μ L
DDW	89 μ L

8. Procedures

For frozen section

- (1) Immerse the specimen on the glass slide into a 4% formalin neutral buffer solution (made by dilution of 10% formalin neutral buffer solution with PBS, as an example) in a Coplin-staining jar, and leave it for 5 minutes at room temperature.
- (2) Immerse the glass slide into PBS, and leave it for 5 minutes as a washing step, which is repeated once.
- (3) Proceed to the step (4) for cultured cells.

For cultured cells

- (1) Suspend cultured cells in 4% formalin neutral buffer solution (made by dilution of 10% formalin neutral buffer solution with PBS, as an example) at a concentration of about 2×10^6 cells/mL, and leave them for 10 min at room temperature.
- (2) Drop 50 μ L of the cell suspension onto a slide glass, and dry it well.

Note : The suspension should be mounted on a silane-coated glass slide, preferably, pre-coated with a reagent for protection of detachment of cells from the glass slide.

- (3) Wash the glass slide twice by immersing it into PBS for 5 min.
- (4) Immerse the glass slide into the cold Permeabilizing solution, and leave it for 2 minutes on ice.
- (5) Repeat a brief wash (for about 10 seconds) of the glass slide by PBS(–) twice, then wipe the residual solution from around the specimen and the back of the glass slide using a paper tissue (e.g. Kimwipe).
- (6) Proceed to the step (4) for paraffin-embedded section described below.

For paraffin-embedded section

- (1) Deparaffinize and hydrate the specimen with distilled water by transfer of the glass slide among a series of the solvent/distilled water in Coplin-staining jars as follows.

Lemosol[®] or Lemosol[®] A : 5 min. $\times$ 2 times
ethanol (99.5% or more) : 5 min. $\times$ 2 times
90% ethanol : 3 min.
70% ethanol : 3 min.
DDW : 3 min.

- (2) Immerse the glass slide into the pre-warmed Protease solution in a Coplin-staining jar, and digest protein of the specimen by incubation at 37°C for 5 minutes.

Note : The factors of the reaction time and temperature will greatly affect the result, so the reliability of the method is dependent on optimizing the reaction time in your specimen.

- (3) Prepare 4 Coplin-staining jars, each containing an appropriate volume of PBS, and wash the glass slide with PBS every 2 minutes at room temperature by transferring the slide from one Coplin-staining jar to the other. After washing, wipe the residual solution around the specimen and the back of the slide with a paper tissue (e.g. Kimwipe).
- (4) Drop 100 μ L (or 50 μ L for the enhanced procedure, mentioned above) of TdT reaction solution onto the specimen on the glass slide, and label 3'-terminals of DNA by incubation for 10 minutes in a moist chamber at 37°C.

Note : This reaction time and temperature are critical, so optimize the reaction time in your specimen.

- (5) Prepare 3 Coplin-staining jars, each containing an appropriate volume of PBS. After removing TdT reaction solution on the specimen, as much as possible, wash the glass slide with PBS every 5 minutes at room temperature by transferring the slide glass from one Coplin-staining jar to the other.

Note : Collect the TdT reaction solution used as the buffer solution containing potassium cacodylate, which is hazardous, and discard properly.

- (6) Eliminate intrinsic peroxidase activity of the specimen by incubating the glass slide in a Coplin-staining jar containing 3% hydrogen peroxide solution for 5 minutes at room temperature.

- (7) Prepare 2 Coplin-staining jars, each containing an appropriate volume of PBS, and wash the glass slide with PBS as described in step (5). After washing, wipe the residual solution around the specimen and the back of the slide with a paper tissue.
- (8) Drop 100 μ L of POD-conjugated antibody solution onto the specimen on the glass slide, and incubate it for 10 minutes in a moist chamber at 37 °C.
- (9) After removing the antibody solution, wash the glass slide with PBS using 3 Coplin-staining jars and wipe the residual solution as described in step (7).

Note : Prepare DAB Solution during washing in step (9), and store it at room temperature until use.

- (10) Drop 100 μ L DAB Solution onto the specimen, and incubate it for 5 minutes at room temperature.
- (11) After removing the substrate solution, wash the specimen by moving the slide up and down 10 times in DDW, and repeat twice in fresh DDW.

Note : If the background arises, adjust the incubation time, observing the specimen under a microscope.

- (12) Counter-stain the specimen by immersing the glass slide in an appropriate staining solution, such as methyl green solution, in a Coplin-staining jar for 5 minutes at room temperature.
- (13) Dehydrate the specimen by either method according to the following:
 - a) Immerse the glass slide in the order of 70%, 90%, and anhydrous (99.5% or more) ethanol every 10 seconds.
 - b) Immerse the glass slide in 1-butanol (99% or more) for about 10 seconds until the specimen changes from a deep blue to light green. Otherwise, repeat in the same manner with fresh solvent until the specimen becomes light green.
- (14) Immerse the specimen into either fresh Lemosol[®], Lemosol[®] A or xylene for 10 seconds, and repeat.
- (15) Drop an appropriate permanent mounting medium, such as Softmount, onto the specimen, cover it with a cover glass, and inspect it under a microscope after the mounting medium hardens.

Note : Make sure the specimen is not observed by fluorescent microscope because fluorescein is used as a tag for the antibody but not for signaling.

For Controls

Negative and positive controls are essential for an accurate judgment.

- 1) Positive control
 1. Drop 100 μ L of DNase I solution onto the specimen after step (3) for paraffin-embedded section or step (5) for cultured cells, and incubate for 15 minutes in a moist chamber at 37 °C.
 2. After removal of the enzyme solution, wash the glass slide with PBS and wipe off the residual solution as described in step (7) for paraffin-embedded section.
 3. Proceed to step (4) for paraffin-embedded sections.
- 2) Negative control
 1. Incubate the specimen with TdT Substrate Solution without TdT at step (4) for paraffin-embedded section.
 2. Proceed to the next step.


9. Reference

- (1) Gavrieli, Y., Sherman, Y. and Ben-Sasson, S. A. : *J. Cell Biology*, **119**, 493 (1992).
- (2) Funaki, N., Sasano, H., Shizawa, S., Nio, M., Iwami, D., Ohi, R. and Nagura, H. : *J. pathol.*, **186**, 429 (1998).
- (3) Li, H. C., Yashiki, S., Sonoda, J., Lou, H., Ghosh, S. K., Byrnes, J. J., Lema, C., Fujiyoshi, T., Karasuyama, M. and Sonoda, S. : *Jpn. J. Cancer Res.*, **91**, 34 (2000).
- (4) Shibata, MA., Morimoto, J. and Otsuki, Y. : *Cancer Gene Therapy.*, **9**, 16 (2002).
- (5) Doi, H., Shibata, MA., Shibata, E., Morimoto, J., Akao, Y., Inuma, M., Tanigawa, N. and Otsuki, Y. : *Anticancer Res.*, **29**, 2485 (2009).

10. Appendix(Flow Chart of the Experiments)

Paraffin-embedded tissue section

Neutralized 4% formalin-fixed culture cells (Neutralized 4% formalin-fixed frozen section)

Allied products

Wako catalog No.	Description	Package Size
122-03991	Lemosol [®]	1 L
128-03993	Lemosol [®]	3 L
120-04411	Lemosol [®] A	1 L
126-04413	Lemosol [®] A	3 L
192-16301	Softmount	250 mL
162-19321	Phosphate Buffered Saline Powder (0.01 mol/L, pH 7.2–7.4)	1 L × 20
138-12701	Methyl Green Solution (0.5 v/w%)	100 mL
293-71403	Apoptosis Ladder Detection Kit <i>Wako</i>	96 reactions
297-71401		24 reactions
297-55901	Annexin V-Fluorescein Staining Kit	50 reactions
517-88601	DAB Immunohistochemistry Substrate	1 set

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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