LabAssay™ Cholesterol

Please, read this instruction carefully before use.

1. Introduction

Cholesterol concentration in serum is known to be closely related with production, absorption and dissimulation in liver and intestinal tract, and lipoprotein metabolism in blood.

LabAssay™ Cholesterol is able to measure total cholesterol in mouse and human serum by an enzyme reaction using *N*-Ethyl-*N*-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline sodium salt (DAOS). It is a simultaneous multi-sample assay format using a microplate, but it measurements can also be made using a test tube.

2. Kit Contents

	Components	State	Amount
(1)	Buffer Solution	liquid	150 mL/1 bottle
(2)	Chromogen Substrate	Freeze dry	For 150 mL/1 bottle
(3)	Standard Solution	liquid	5 mL/1 bottle

※ Assay in a microplate ; 500 tests. Assay in a test tube ; 50 tests.

3. Storage and Expiration

When the complete kit is stored at 2°C - 10°C, the kit is stable until the expiration date shown on the label on the container.

4. Assay Principle

Cholesterol esters in the sample are decomposed into free cholesterol and fatty acid by cholesterol esterase, when the sample reacts with chromogen reagent. The cholesterol is oxidized with existing free cholesterol by cholesterol oxidase, and simultaneously hydrogen peroxide is produced. The produced hydrogen peroxide let DAOS and 4-Aminoantipyrin oxidize and condensate quantitatively by peroxidase (HRP), which produces a blue pigment. Quantitation of total cholesterol in the sample can be made by measurement of the absorbance.

5. Materials and Apparatuses to be prepared

- 96 wells microplate (transparent type) ·Micropipette ·Incubator maintaining at 37°C *
- Plate mixer *
- Microplate reader with 600 nm wavelength filter (* if the microplate reader is not equipped.)

(For Test Tube method)

- Test tube ·Pipette ·Incubator maintaining at 37°C
- Spectrophotometer or colorimeter with 600 nm wavelength filter

6. Preparation of Reagents to be used

① Chromogen Reagent :

Chromogen reagent is prepared by dissolving 1 vial of Chromogen Substrate with 1 vial of Buffer Solution.
After adjustment, it may be stored at 2°C - 10°C for up to 3 weeks.

② Standard Solution

Standard solution is prepared by dilution of the provided Standard Solution.

No.	Provided Standard	Distilled or deionized water	Sample volume	Cholesterol concentration
1	50 μ L	150 μ L	2 μ L	50.0 mg/dL
2	100 μ L	100 μ L	2 μ L	100.0 mg/dL
3	undiluted	—	2 μ L	200.0 mg/dL
4	undiluted	—	4 μ L	397.4 mg/dL ^{*1}
5	undiluted	—	6 μ L	592.2 mg/dL ^{*1}

^{*1} The test sample volume is usually 2 μ L, but 4 or 6 μ L are taken in this case. The cholesterol concentration must be corrected accordingly as indicated in the table above because of volume increasing.

7. Preparation of Samples

- Serum sample should be used immediately after collecting blood.
- Hemolysis causes a slightly increase of the cholesterol value.

- Ascorbic acid and bilirubin may not significantly affect the assay.

8. Procedure

(1) Assay in a Microplate

Perform the assay in the wells according to the following table scheme.


	Test	Standard	Blank
Sample	Serum 2 μ L	Standard solution 2 μ L	—
Chromogen reagent	300 μ L	300 μ L	300 μ L
Mix well and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Main wavelength 600 nm / Sub wavelength 700 nm			

(2) Assay in a Test Tube

Perform the assay in a test tube according to the following table scheme.

	Test	Standard	Blank
Sample	Serum 20 μ L	Standard solution 20 μ L	—
Chromogen reagent	3 mL	3 mL	3 mL
Mix well and incubate at 37°C for 5 min. Measure the absorbance of the test sample and standard solution with the blank solution as the control. Colorimeter with 600 nm wavelength filter Spectrophotometer : Main wavelength 600 nm (subwavelength 700 nm)			

9. Standard Curve

10. Performance

(1) Sensitivity

- The absorbance is below 0.11, when measuring purified water as a sample.
- The absorbance is 0.13 to 0.65, when measuring 200 mg/dL cholesterol as a sample.

(2) Specificity

The cholesterol concentration is less than $\pm 10\%$, when measuring the known concentration of control serum as a sample.

11. Usage Notes

(1) Notes on the assay

- Measure the total cholesterol according to this insert sheet.
- Use measurement apparatuses according to the use manual.
- Keep the provided reagents under the indicated conditions before the expiration. Avoid using frozen reagents. Each of

the reagents should be used immediately after opening, or store them under the indicated conditions.

- Vials should be opened carefully, as the inner is decompressed.
- The chromogen reaction is completed in ca. 2 min. Incubation up to 20 min. may not affect the assay. The reaction is stable and may not significantly affect the absorbance after 2 hours of the reaction.
- This kit is for research use only, but not for diagnostic use.

(2) Safety Precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water.
- A pipette should be used with a safety pipette filler.
- Vials should be opened carefully.

(3) Waste

- The waste should be processed appropriately according to law.

[Reference]

1. Allain, C. C., Poon, L. S., Chan, C. S. G., Richmond, W. and Fu, P. C.: *Clin. Chem.*, 20, 470 (1974).
2. Richmond, W. : *Clin. Chem.*, 19, 1350 (1973).
3. Richmond, W. : *Scand. J. Clin. Lab. Invest.*, 29 (Supl.), 126 (1972)

LabAssayTM Cholesterol

[Storage]	Store at 2°C - 10°C
[Expiration date]	Indicated on the label
[Package]	For 500 tests
[Cat #]	293-93601

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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