## For Tissue Optical Clearing Reagent AbSca/e-G Solution Set

### [Introduction]

Dr. Atsushi Miyawaki *et al.* developed a sorbitol-based optical clearing method, Sca/e S.

AbSca/e-G Solution Set is a reagent used for the AbSca/e-G, an improved version of AbSca/e, that is Sca/e S's three-dimensional immunohistochemical staining method.

Compared to AbSca/e, it has improved antibody staining ability and enables observation using more antibodies and fluorescent dyes. It can be used immediately by mixing the reagents in the contents.

It can also be used for non-tissue transparency samples such as regular tissue sections.

### [Set Contents]

- (1) Solution A 50 mL × 1 bottle
- (2) Solution B 2 mL × 1 bottle

### [Storage]

Store in the dark at 2-10°C.

### [Additional Required Materials]

Reagents :

Code No.	Description	Size
193-18455	SCALEVIEW-A2	500 mL
196-18521	SCALEVIEW-S0	250 mL
194-18561	SCALEVIEW-S4	250 mL
041-34425	deSca/e Solution	500 mL
160-16061	Paraformaldehyde	100 g
164-25511	1 × PBS (-)	5 L
215-00611	Urea	100 g

Equipments :

- 1) 5 mL/15 mL/30 mL/50 mL conical centrifuge tube (for each sample)
- 2) 2 mL microcentrifuge tube (for each sample)
- 3) Culture dish (100 mm, 60 mm, and 35 mm diameter) (for each sample)
- 4) Seesaw shaker
- 5) Fluorescence microscope and Recommend Objective lenses for tissue clearing samples

### [Procedure] (1-2 mm brain samples)

(Fixation)

- 1) Transcardially perfuse an anesthetised mouse with 4% paraformaldehyde (PFA)/PBS (pH 7.6-7.8).

- 2) Remove the whole brain and subject it to post-fixation in 4% PFA/PBS at 4°C for 72 hrs.
- 3) Wash the samples in PBS.
- 4) Prepare vibratome sections.

(Pretreatment)

- 1) Transfer to 15 mL tube\*<sup>1</sup>, add 5 mL of deSca/e Solution, and rinse it at room temperature × 2 times.
- 2) Replaced with 6 mL of SCALEVIEW-S0 and incubate it with gentle shaking (10 rpm) at 37°C for 6 hrs.
- 3) Replaced with 6 mL of SCALEVIEW-A2 and incubate it with gentle shaking (10 rpm) at 37°C for 6-8 hrs.
- 4) Replaced with 6 mL of 8M Urea Solution and incubate it with gentle shaking (10 rpm) at 37°C for 6-8 hrs.
- 5) Replaced with 6 mL of SCALEVIEW-A2 and incubate it with gentle shaking (10 rpm) at 37°C for 3-4 hrs.
- 6) Replaced with 10 mL of deSca/e Solution and incubate it with gentle shaking (10 rpm) at Room temperature for 6 hrs.

\*1...Use tubes that match the size of the sample.

(Antibody staining, Re-fixation, Clearing)

- 1) Prepare antibody-diluted AbSca/e-G Solution.
  - Add 75  $\mu$ L of Solution B to 1.5 mL of Solution A.
  - Add antibody and mix well.
- 2) Transfer the sample to 2 mL tube, add 1), and incubate it with gentle shaking (10 rpm) at 37°C for 3 days.
- 3) Transfer the sample to 15 mL tube\*<sup>2</sup>, add 5 mL of AbSca/e-O Solution\*<sup>3</sup>, and incubate it with gentle shaking (10 rpm) at room temperature for 5 min × 3 times.
- 4) Replaced with 5 mL of 4% PFA/PBS (pH 7.6-7.8) and incubate it with gentle shaking (10 rpm) at room temperature for 30 min.
- 5) Replaced with 8 mL of SCALEVIEW-S4 and incubate it with gentle shaking (10 rpm) at 37°C for 6 hrs.
- 6) Replaced with flesh SCALEVIEW-S4 and store at room temperature in the dark.

\*2...Use tubes that match the size of the sample.

\*3...AbSca/e-O Solution : 0.33 M Urea and 0.1% (wt/vol) Triton X-100 in PBS (-) Solution

(Imaging)

- 1) Mount onto an imaging chamber. Use appropriate amount of SCALEVIEW-S4 solution to mount the sample. Observe the SCALEVIEW-S cleared samples using fluorescence microscopy.

(Note)

- For thicker brain slices (>2 mm) or hemisphere mouse brain whole mouse brain, incubation time in each step can become longer and SCALEVIEW-S solutions in each step can increase amount (e.g., hemisphere brain 25 mL, whole brain 40 mL each).
- SCALEVIEW-S4 mounted samples will not be solidified. If

necessary for each sample, 1.5-3.0% melting agarose can be added to minimize movement artifacts during imaging.

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