For Tissue Optical Clearing Reagent SeeDB2 Trial Kit

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Drs. Ke and Imai developed a new type of tissue clearing agent, SeeDB2. SeeDB2 is ideal for high-resolution three-dimensional imaging of fluorescent proteins. SeeDB2G and S are designed match the refractive indices of glycerol (1.46) and oil (1.52), minimizing spherical aberrations for high-NA glycerol- and oil- immersion objective lenses, respectively. SeeDB2 is particularly powerful in confocal and super-resolution microscopy using high-NA objective lenses. It is noteworthy that many of fluorescent proteins are highly preserved in SeeDB2, much better than in PBS or other commercialized mounting media optimized for fluorescent dyes. Thus, SeeDB2 is also an ideal mounting medium for samples labeled with fluorescent proteins. SeeDB2 is useful not only for thick brain tissues, but also for thin samples for cell biology or tissue sections.

【Kit contents】

- | | |
|----------------------|------------------|
| (1) SeeDB2G solution | 50 mL × 1 bottle |
| (2) SeeDB2S solution | 50 mL × 1 bottle |
| (3) PBS (-) | 50 mL × 1 bottle |
| (4) Saponin | 5 g × 1 bottle |

【Storage】

Store below 25°C in the dark.

【Additional Required Materials】

Reagents :

- 1) Paraformaldehyde (Code No. 160-16061)
- 2) 1 × PBS (-) (Code No. 164-25511)
- 3) Sodium Azide (Code No. 195-11092)
- 4) 5 mL Eppendorf tube
- 5) Immersion solution :
Immersol 518F (Carl Zeiss) or Immersion oil type F (Leica, Olympus)
for oil-immersion lens ; glycerol (Leica) for glycerol-immersion lens

Equipments :

- 1) Ring forceps (or paint brush for thin slices) or paint brush
- 2) Overhead tube rotator
- 3) Sterilization filter (optional ; for example, Steradisc 0.2 μ m from Kurabo and Millex 0.22 μ m from Merck-Millipore)
- 4) Vibratome (optional ; for example Dosaka MicroSlicer)
- 5) SeeThrough Chamber
- 6) Glass slides and glass coverslips (No. 1.5H, for example, Precision cover glasses from Marienfeld)

Microscope :

Cofocal Microscope with high-NA objective lenses (NA>1.2).
Super-resolution Microscope
(e.g., TCS SP8X STED, LSM880-Airyscan, TCS SP8 Hyvolution2, etc.)
Light-sheet microscope (e.g., TCS SP8 DLS with clearing objective lens)

〈Note〉

- It is highly recommended to use oil-immersion lens for SeeDB2S and glycerol-immersion lens for SeeDB2G. Coverslips should be No. 1.5H (0.170 μ m thick). In light-sheet microscopy, it is recommended to use a specialized objective lens for cleared samples (index optimized).
- White floating matter maybe observed in SeeDB2S Solution, but there is no problem in quality. We recommend that make sure to mix SeeDB2S Solution lightly before you use.

【Procedure】

(Solution preparation)

Clearing reagent preparation (for 3 mL reaction in 5 mL tube)

1. 20% saponin stock solution
Dissolve 0.6 g Saponin in 3 mL distilled water. It is recommended to filter-sterilize the solution for longer storage.
2. Permeabilization solution
300 μ L 20% saponin
2700 μ L PBS
3. Clearing solution 1 (2 : 1 mixture)
300 μ L 20% saponin
1700 μ L distilled water
1000 μ L SeeDB2G solution
4. Clearing solution 2 (1 : 1 mixture)
300 μ L 20% saponin
1200 μ L distilled water
1500 μ L SeeDB2G solution
5. Clearing solution 3 (for SeeDB2G)
Dissolve 60mg saponin in 3 mL SeeDB2G solution
6. Clearing solution 4 (for SeeDB2S, optional) :
Dissolve 60 mg saponin in 3 mL SeeDB2S solution

〈Note〉

- These solutions are stable up to 1 week at room temperature and 4 weeks at 4°C. For longer storage, it is recommended to filter-sterilize the 20% saponin solution and/or add sodium azide to each solution at 0.01% to prevent bacterial contamination. SeeDB2S solution includes sodium azide at 0.01%.
- 3 mL each is sufficient to clear up to 2 mm-thick mouse brain slices.

Procedure (brain slices)

(Sample preparation)

- 1) Fix the dissected sample in 4% PFA in PBS (4°C, o/n).
- 2) Wash the samples in PBS.
- 3) Prepare vibratome sections (optional, 0.2-2 mm thick).

(Clearing)

- 4) Transfer the sample into 3 mL of permeabilization solution and incubate it with gentle rotation (\sim 4 rpm) at room temperature for 12-16 hours.
- 5) Transfer the sample to a new tube filled with 3 mL of clearing solution
1. Place the tube on an overhead rotator at room temperature for 6-10 hours.
- 6) Transfer the sample to a new tube filled with 3 mL of clearing solution
2. Place the tube on an overhead rotator at room temperature for

6-10 hours.

- 7) Transfer the sample to a new tube filled with 3 mL of clearing solution 3. Place the tube on an overhead rotator at room temperature for at least 12 hours.
- 8) (optional for SeeDB2S) Transfer the sample to a new tube filled 3 mL of clearing solution 4. Place the tube on an overhead rotator at room temperature for at least 12 hours.
- 9) Transfer the sample into 3 mL of SeeDB2G solution (from step 7) or SeeDB2S solution (from step 8).

(Imaging)

- 10) Pick the cleared sample with ring forceps or paint brush and mount onto an imaging chamber (e.g., SeeThrough Chamber from Wako) or a glass slide. Use appropriate amount of SeeDB2G/S solution to mount the sample and seal the sample with a glass coverslip. Excess amount of SeeDB2G/S at the edge of the coverslip will be air dried within minutes.
- 11) Observe the SeeDB2G/S-cleared samples using fluorescence microscopy. Use glycerol immersion to image SeeDB2G samples with glycerol-immersion objective lens ; and oil to image SeeDB2S samples with oil-immersion objective lens.

(Note)

- The cleared sample should not be exposed to the air for long term, because the surface may become sticky due to water evaporation. SeeDB2G/S solution cannot be used as an immersion solution.
- For the preparation of imaging chambers, see SeeDB Resources (<https://sites.google.com/site/seedbresources/>).
- For very thin samples and tissue sections, some steps before mounting can be omitted. However, direct soaking in SeeDB2G/S solution may result in shrinkage and deformation of samples. Serial incubation is recommended to best maintain the morphology.
- For thinner brain slices (<1 mm) or smaller samples (e.g., Drosophila brain), incubation time in each step can be shortened (e.g., 1h or 3h each).
- SeeDB2G/S-mounted samples will not be solidified. If necessary for small samples, 0.1% low-melting agarose (e.g., UltraPure Low Melting Point Agarose from Thermo Fisher) can be added to minimize movement artifacts during imaging.
- While SeeDB2 solution protect fluorescent proteins from photobleaching, some of chemical dyes (e.g., Alexa Fluor dyes) are sensitive to photobleaching in SeeDB2 solution.


Figure 1. Clearing and mounting procedure

[Example images]


Figure 2. Mouse brain slice (adult, 1.5 mm-thick) before and after clearing with SeeDB2.


Figure 3. Confocal images of brain slices from a Thy1-YFP (line H) mouse. An oil-immersion objective lens (NA 1.4) was used. Note that fluorescence level did not decrease even though laser power was constant throughout depth. Scale bars are 2 μ m.

[Related products]

Code No.	Description	Size
294-35631	See Through Chamber, 0.3 mm thick	10 set
291-35641	See Through Chamber, 0.5 mm thick	10 set
295-35661	See Through Chamber, 1.0 mm thick	10 set
292-35671	See Through Chamber, 2.0 mm thick	10 set
299-35681	See Through Chamber, 3.0 mm thick	10 set
160-16061	Paraformaldehyde	100 g
164-25511	1 \times PBS (-)	5 L

[References and Notes]

Ke et al., : *Cell Reports* **14**, 2718 (2016)

SeeDB Resources (<https://sites.google.com/site/seedbresources/>) : updated information and technical TIPS from the authors.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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