Code No. 294-83501 (10 tests)

Human ES/iPS Cell Staining Kit-BF

This kit is the product which can check the undifferentiated state of human ES/iPS cells under bright-field using the recombinant lectin rBC2LCN, which has a high affinity for undifferentiated human ES/iPS cells.

rBC2LCN is the recombinant lectin of the N-terminal domain of BC2L-C, which derived from *Burkholderia cenocepacia*. rBC2LCN has a high affinity for H type 3 (Fuc α 1-2Gal β 1-3GalNAc), which is the mucin-like O-type sugar chain on podocalyxin present in the cell membrane of human ES/iPS cells, so rBC2LCN is reported as an undifferentiated marker of human ES/iPS cells.

Peroxidase-labeled rBC2LCN staining procedures that produce a colored human ES/iPS cells can be visualized under bright-field microscopy. Because peroxidase-labeled rBC2LCN is used, undifferentiated human ES/iPS cells can be stained with a very simple procedure, and the undifferentiated state of human ES/iPS cells can be confirmed under bright-field without using a fluorescence microscope.

This product is for laboratory use only.

[Kit Components]

	Code No.	Components	Volume/Quantity
1	291-83511	rBC2LCN-POD	100 μ L/1 tube
2	298-83521	DAB Buffer	15 mL/1 bottle
3	295-83531	DAB Reagent	600 μ L/1 tube
4	292-83541	Washing Buffer (10 \times)	10 mL/1 bottle

[Reagents Required, not Supplied]

- Purified Water (Distilled Water)
- D-PBS (-)
- Hydrogen Peroxide Solution (if necessary)

[Reagents Preparation]

Please prepare each reagent according to the following procedures before use. Also, please prepare only the required volume at the time of use.

Preparation of rBC2LCN-POD Reaction Solution

Add 1/100 volume of rBC2LCN-POD to D-PBS (-) and mix well. After mixing, store in the dark and use within 2 hours.

Preparation of DAB Reaction Solution

Add 40 μ L of DAB Reagent to 1 mL of DAB Buffer, and mix well. After mixing, store in the dark and use within the day.

Preparation of Washing Buffer (1 \times)

Dilute the Washing Buffer (10 \times)^{*1} 10-fold with purified water (distilled water) and mix well.

[Procedures]

Please prepare the reagents according to [Reagents Preparation].

- (1). Prepare the vessels containing cells to be stained.^{*2, 3}
- (2). Discard the solution in the vessels^{*4}, and add the rBC2LCN-POD Reaction Solution to the vessels by referring to the Table 1.

Table 1

Vessel	rBC2LCN-POD Reaction Solution (mL)
6 well plate	1
12 well plate	0.5
24 well plate	0.25
96 well plate	0.1
35 mm dish	1
60 mm dish	2
100 mm dish	5

- (3). Incubate for 1 hour at room temperature in the dark.
- (4). Wash 3 times with Washing Buffer (1 \times).
- (5). Discard the Washing Buffer (1 \times) from the vessels, and add DAB Reaction Solution to the vessels by referring to the Table 2.

Table 2

Vessel	DAB Reaction Solution (mL)
6 well plate	1
12 well plate	0.5
24 well plate	0.25
96 well plate	0.1
35 mm dish	1
60 mm dish	2
100 mm dish	5

- (6). Incubate vessels at room temperature until staining develops. 2-10 minutes generally provides an acceptable staining intensity.
- (7). Wash the samples in Purified Water.

※1 The components of the Washing Buffer (10 \times) may precipitate when stored in a refrigerator. Therefore, after returning to room temperature, check that there are no precipitates before using. If there are precipitates, please use it after dissolution.

※2 If the cells are in culture, please fix the cells with 4% Para-formaldehyde.

※3 When cells with endogenous peroxidase are included, a non-specific reaction occurs. Please inactivate endogenous peroxidase with 3% Hydrogen Peroxidase Solution according to following procedures.

〈Procedures to inactivate peroxidase〉

- 1). Dilute the Hydrogen Peroxide Solution with D-PBS (-) to prepare a 3% Hydrogen Peroxidase Solution.
- 2). After discarding the solution in the vessels containing cells, add 3% Hydrogen Peroxide Solution.
- 3). Incubate at temperature for 5 minutes.
- 4). Wash 3 times with D-PBS (-).

※4 If the solution in which the cells are stored is not buffer solution such as D-PBS (-), please wash them before adding rBC2LCN-POD Reaction Solution.

【Package】

10 tests

【Storage】

Store at 2-10℃ in the dark (Don't freeze)

Please use opened reagents as soon as possible since they may be affected by storage conditions. The remained solution should be stored at 2-10℃, in the dark, with the cap tightened.

【References】

1. Watanabe, T., Yasuda, S., Kusakawa, S., Kuroda, T., Futamura, M., Ogawa, M., Mochizuki, H., Kikkawa, E., Furukawa, H., Nagaoka, M. and Sato, Y. : *Cytotherapy* (2020). DOI:https://doi.org/10.1016/j.jcyt.2020.07.009
2. Onuma, Y., Tateno, H., Hirabayashi, J., Ito, Y. and Asashima, M. : *Biochem. Biophys. Res. Commun.*, **431**, 524 (2013).
3. Tateno, H., Matsushima, A., Hiemori, K., Onuma, Y., Ito, Y., Hasehira, K., Nishimura, K., Ohtaka, M., Takayasu, S., Nakanishi, M., Ikehara, Y., Nakanishi, M., Ohnuma, K., Chan, T., Toyoda, M., Akutsu, H., Umezawa, A., Asashima, M. and Hirabayashi, J. : *Stem Cells Transl. Med.*, **2**, 265 (2013).
4. Tateno, H., Toyota, M., Saito, S., Onuma, Y., Ito, Y., Hiemori, K., Fukumura, M., Matsushima, A., Nakanishi, M., Ohnuma, K., Akutsu, H., Umezawa, A., Horimoto, K., Hirabayashi, J. and Asashima, M. : *J. Biol. Chem.*, **286**, 20345 (2011).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : +1-804-271-7677
Facsimile : +1-804-271-7791
<http://www.wakousa.com>

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