LBIS™ Rat LH ELISA Kit

Please, read this instruction carefully before use.

1. Intended use

LBIS™ Rat LH ELISA Kit is a sandwich ELISA system for quantitative measurement of rat LH. This is intended for research use only.

2. Introduction

LH is produced and stored in basophilic cells called gonadotrophs in the anterior pituitary gland. LH is found in all vertebrates from fishes to mammals. Some reports indicated the presence of LH-like substance in the testis.

LH has a heterodimer structure consisted of α -subunit and β -subunit with a total molecular weight of approximately 29000. α -Subunit is a glycoprotein which is common to other pituitary hormones follicle-stimulating hormone (FSH) and thyroid-stimulating hormone (TSH), while β -subunit is specific to LH, and is also a glycoprotein.

The receptor of LH is a transmembrane receptor of G-protein coupled type, which relates to PKA system and penetrates the cell membrane 7 times.

In females, LH acts on matured granulosa cells in ovarian follicle, where it causes, in cooperation with FSH, follicular maturation, estrogen production, and ovulation, and after ovulation, acts on corpora lutea and promotes progesterone production and secretion.

In males, LH acts on the interstitial cells (Leydig cells) in the testis, causing production and secretion of androgens, and secondarily promotes spermatogenesis via androgen.

Deficiency of LH leads to the lower secretion of sex steroids, atrophy of interstitial cells, and failure of ovulation and luteinization, while excessive LH causes hyperplasia of testicular interstitial cells followed by atrophy, increased secretion of estrogen or androgen, super-ovulation, and accelerated sexual maturation.

GnRH (LHRH) secreted from hypothalamus act directly on gonadotrophs and stimulate LH secretion. Physiologically, positive feedback of estrogen at preovulatory period causes LH surge via GnRH and leads to ovulation. LH secretion is promoted in postmenopausal period women and in aging in men.

LH secretion was minimized by increased blood levels of sex steroids, opioid peptides especially endorphin, in childhood, and pregnancy.

3. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Rat LH ELISA Kit, standards or samples are incubated in monoclonal anti-LH β antibody coated wells to capture LH. After 2 hours' incubation and washing, biotinylated anti-LH α antibody is added and incubated further for 1 hour to bind with captured LH. After washing, Peroxidase-conjugated Streptavidin Solution is added, and incubated for 30 minutes. After washing, HRP-complex remaining in wells are reacted with a chromogen (TMB) for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is proportional to LH concentration. The standard curve is prepared by plotting absorbance against standard LH concentrations. LH concentrations in unknown samples are determined using this standard curve.

4. Performance Characteristics

- Assay range

The assay range of the kit is 0.313 ng/mL - 10 ng/mL (for 5 \times sample dilution, 1.565 ng/mL - 50 ng/mL)

- Specificity

The kit uses monoclonal antibodies specific to LH.

- Precision of assay

Within assay variation (2 samples, 8 replicates assay) Mean CV was within 10%.

- Reproducibility

Between assay variation (3 samples, 4 days, 4 replicates assay) Mean CV was within 10%.

- Recovery test

Standard LH was added in 3 concentrations to 2 serum samples and was assayed in duplicates.

The recoveries were 95.1% - 106%.

- Dilution test

Two serum samples were serially diluted by 3 steps.

The dilution curves showed excellent linearity. ($R^2 = 0.999$).

5. Reference Assay Data

Rat LH assay data : Mean 1.93 ng/mL, SD 1.22 ng/mL (Highest : 4.33 ng/mL/Lowest : 0.925 ng/mL)

Subspecies : CD, female, n=9, 8 weeks-old, Serum, fed ad libitum

These are reference data. Blood LH levels may change due to breeding, sampling, and sample storage conditions.

6. References

- 1) Neonatal immune challenge alters reproductive development in the female rat
L Sominsky *et al.*
Horm Behav. 2012 Feb 14.
- 2) Gonadotropin-Inhibitory Hormone Inhibits GnRH-Induced Gonadotropin Subunit Gene Transcriptions by Inhibiting AC/cAMP/PKA-Dependent ERK Pathway in L β T2 Cells.
Son Y. L. *et al.* : *Endocrinology*, **153** (5), 2332 (2012).

7. Precautions

- For professional use only, beginners are advised to use this kit under the guidance of experienced person.
- Wear gloves and laboratory coats when handling assay materials.
- Do not drink, eat or smoke in the areas where assays are carried out.
- In treating assay samples of animal origin, be careful for possible biohazards.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the Stop Solution because it is sulfuric acid. The Stop Solution and the substrate solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Avoid contact with the acidic Stop Solution and TMB Solution, which contains hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye-glasses and protective clothing when handling these reagents.
- The materials must not be pipetted by mouth.
- Unused samples and used tips should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- Use clean laboratory glassware.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.

8. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated Plate	Use after washing	96 wells/1 plate
(B) Rat LH Standard	Concentrated. Use after dilution	200 μ L/1 bottle
(C) Buffer	Ready for use.	60 mL/1 bottle
(D) Biotin-conjugated Antibody Solution	Concentrated. Use after dilution.	100 μ L/1 bottle*
(E) Peroxidase-conjugated Streptavidin Solution	Concentrated. Use after dilution.	100 μ L/1 bottle*
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(G) Buffer for Sample	Ready for use.	12 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 \times)	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate Seal	—	4 sheets

* You can take out 100 μ L from vials.

【Storage and Stability】

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C.

[(B) Rat LH Standard Solution]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

[(C) Buffer], [(F) TMB Solution] and [(G) Buffer for Sample]

If not opened, store at 2°C - 10°C. Once opened, we recommend using them as soon as possible to avoid influence by environmental condition.

[(D) Biotin-conjugated Antibody Solution] & [(E) Peroxidase-conjugated Streptavidin Solution]

Unused working solution (already diluted) should be disposed.

[(H) Stop Solution]

Close the stopper tightly and store at 2°C - 10°C.

[(I) Wash Solution (10×)]

The rest of undiluted buffer : store at 2°C - 10°C. Dispose any unused diluted buffer.

9. Equipments or supplies required but not supplied ☐ Use as a check box

- ☐ Deionized water (or Distilled water) ☐ Test tubes for preparation of standard solution series.
- ☐ Glassware for dilution of Wash Solution (10×) (a graduated cylinder, a bottle)
- ☐ Pipettes (disposable tip type). One should be able to deliver 10 μ L precisely, and another for 10 μ L - 100 μ L and 100 μ L - 500 μ L. ☐ Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 50 μ L. ☐ Paper towel to remove washing buffer remaining in wells. ☐ A vortex-type mixer. ☐ A shaker for 96 well-plate (600 rpm - 1200 rpm)
- ☐ An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. ☐ A 96 well-plate reader (450 nm \pm 10 nm, 620 nm : 600 nm - 650 nm) ☐ Software for data analysis.

10. Preparation of Samples

This kit is intended to measure rat serum (do not use serum-separation-accelerant such as Serum Gel for fear of low assay value) and plasma.

The necessary sample volume for the standard procedure is 10 μ L.

- **Dilute rat serum (plasma) to 2× using (G) Buffer for Sample. Mix and let it stand still for 10 min at 20 - 25°C, then further dilute it to 2.5× using (C) Buffer (final dilution rate is 5×),** and add it to wells for assay. If sample is in good condition, the limit of final dilution rate would be 2.5×. In this case, dilute sample to 2× using (G) Buffer for Sample, mix and let it stand still for 10 min and dilute it further 1.25× using (C) Buffer.
 - * **To avoid influence of blood (high lipid or hemolysis, etc.), 5× at the final dilution is recommended for assay. Even if the final dilution is 5×, abnormal value might be obtained when original sample had heavy chyle or hemolysis.**
- When plasma samples are used, we recommend EDTA-2Na (at a final concentration of 1 mg/mL) to keep samples' pH and avoid interference of Ca²⁺. Heparin is not suitable.
- Anesthesia while sampling may influence the assay system. We do not recommend ether anesthesia.
- Don't use sample tubes for human when collecting blood. We have not checked all kinds of sample tubes. Please contact us if necessary.
- Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.
- If presence of interfering substance is suspected, examine by dilution test at more than 2 points.

Storage and stability

Samples should be immediately assayed or stored below -35°C. Defrosted samples should be mixed thoroughly for best results. LH in samples will be inactivated if stored at 2°C - 10°C. If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35°C. Avoid repeated freeze-thaw cycles. Sample dilution should be made just before your assay.

11. Preparation of Reagents

- ◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.
- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Rat LH Standard]

Make a serial dilution of master standard solutions using (C) Buffer to prepare each standard solution.

Volume of standard solution	(C) Buffer	Concentration (ng/mL)
Original solution : 50 μ L	450 μ L	10.0
10 ng/mL solution : 200 μ L	200 μ L	5.0
5 ng/mL solution : 200 μ L	200 μ L	2.5
2.5 ng/mL solution : 200 μ L	200 μ L	1.25
1.25 ng/mL solution : 200 μ L	200 μ L	0.625
0.625 ng/mL solution : 200 μ L	200 μ L	0.313
0 (Blank)	200 μ L	0

[(D) Biotin-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with (C) Buffer to **1 : 100**.

[(E) Peroxidase-conjugated Streptavidin Solution]

Prepare working solution by dilution of (E) with (C) Buffer solution to **1 : 100**.

[(I) Wash Solution (10 \times)]

Dilute 1 volume of the concentrated Wash Solution (10 \times) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated Wash Solution (10 \times) and 900 mL of deionized water (or distilled water).

12. Assay Procedure

Remove the cover sheet of the Antibody-coated Plate after bringing up to 20°C - 25°C.

- (1) Wash the (A) Antibody-coated Plate by filling the well with 300 μ L of washing buffer and discard 4 times (*①), then strike the plate upside-down onto several layers of paper towels to remove residual buffer in the wells (*⑤).
- (2) Pipette 50 μ L of properly diluted sample to the designated sample wells. (Limit of dilution is 5 \times in standard assay procedure)
- (3) Pipette 50 μ L of standard solution to the wells designated for standards.
- (4) Shake the plate gently on a plate shaker (*②).
- (5) Stick a Plate Seal (*③) on the plate and incubate for 2 hours at 20°C - 25°C.
- (6) Discard the reaction mixture and rinse wells as step (1).
- (7) Pipette 50 μ L of Biotin-conjugated Antibody Solution to all wells, and shake as step (4).
- (8) Stick a Plate Seal (*③) on the plate and incubate the plate for 1 hour at 20°C - 25°C.
- (9) Discard the reaction mixture and rinse wells as step (1).
- (10) Pipette 50 μ L of Peroxidase-conjugated Streptavidin Solution to all wells, and shake as step (4).
- (11) Stick a Plate Seal (*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
- (12) Discard the reaction mixture and rinse wells as step (1).
- (13) Pipette 50 μ L of (F) TMB Solution to wells, and shake as step (4).
- (14) Stick a Plate Seal (*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
- (15) Add 50 μ L of the (H) Stop Solution to all wells and shake as step (4).
- (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (*④)) immediately using a plate reader.

*Refer to 16. Summary of Assay Procedure for notes of *①, *②, *③, *④, and *⑤.


13. Technical Tips

- For manual operation, proficiency in pipetting technique is recommended.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- The TMB Solution should be almost clear pale yellow before use. It turns blue during reaction, and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the Antibody-coated Plate is module type of 8 wells \times 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.

14. Calculations

- (1) Prepare a standard curve using semi-logarithmic or two-way logarithmic section paper by plotting absorbance* (Y-axis) against LH concentration (ng/mL) on X-axis.
- (2) Using the standard curve, read the LH concentration of a sample at its absorbance*, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, repeat the assay after proper dilution of samples with (C) Buffer.

* We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.

Rat LH assay standard curve (an example)

Absorbance may change due to assay situation.

* Absorbance at 450 nm minus absorbance at 620 nm.

15. Trouble Shooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not be added.
- 2) Reagents necessary for coloration such as Biotin-conjugated Antibody Solution, Peroxidase-conjugated Streptavidin Solution, or TMB Solution might not be added.
- 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotin-conjugated Antibody Solution or Peroxidase-conjugated Streptavidin Solution.
- 4) Contamination of enzyme inhibitor (s).
- 5) Influence of the temperature under which the kits had been stored.
- 6) Excessive hard washing of the well plate.
- 7) Addition of TMB Solution soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- Intense coloration in all wells including blank

Possible explanations :

- 1) Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.)
- 2) Overdeveloping. Incubation time with TMB Solution should be decreased before addition of Stop Solution.
- 3) Too high incubation temperature. Adjust the temperature to 20°C - 25°C.

- High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon

- Q-2 : I found there contains liquid in 96 well-plate when I opened the box. What is it?


A-2 : When we manufacture 96 well-plate, we insert preservation stabilizer in wells.

16. Summary of Assay Procedure ☐ : Use as a check box

* First, read this instruction manual carefully and start your assay after confirmation of details.

- ☐ Bring the (A) Antibody-coated Plate and all reagents to 20°C - 25°C for 2 hours.

- ☐ Wash Solution (10×) concentrate must be diluted to **10 times** by deionized water (or distilled water) that returned to 20°C - 25°C.
- ☐ Standard LH solution dilution example : Use (C) Buffer returned to 20°C - 25°C.

- ☐ Sample dilution-1 : Dilute rat sample to 2× using (G) Buffer for Sample, mix and let it stand still for 10 minutes at 20°C - 25°C.
- ☐ Sample dilution-2 : Dilute the above prepared sample to 2.5× using (C) Buffer returned to 20°C - 25°C.

- ☐ Antibody-coated Plate
- ☐ ↓ Washing 4 times (*①), (*⑤)
- ☐ Properly diluted Samples /Standards 50 µL
- ☐ ↓ Shaking (*②), Incubation for 2 hours at 20°C - 25°C. (Standing (*③))
- ☐ Dilute (D) Biotin-conjugated Antibody Solution to 1 : 100 with (C) Buffer returned to 20°C - 25°C. Dilute reagents during the first reaction.
- ☐ ↓ Washing 4 times (*①), (*⑤)
- ☐ Biotin-conjugated Antibody Solution 50 µL
- ☐ ↓ Shaking (*②), Incubation for 1 hour at 20°C - 25°C. (Standing (*③))
- ☐ Dilute (E) Peroxidase-conjugated Streptavidin Solution to 1 : 100 with (C) Buffer returned to 20°C - 25°C. Dilute reagents during the second reaction.
- ☐ ↓ Washing 4 times (*①), (*⑤)
- ☐ Peroxidase-conjugated Streptavidin Solution 50 µL
- ☐ ↓ Shaking (*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (*③))
- ☐ ↓ Washing 4 times (*①), (*⑤)
- ☐ TMB Solution 50 µL
- ☐ After dispense, the color turns to blue depending on the concentration
- ☐ ↓ Shaking (*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (*③))
- ☐ Stop Solution 50 µL
- ☐ After dispense, the color turns to yellow depending on the concentration.
- ☐ ↓ Shaking (*②), Immediately shake.
- ☐ Measurement of absorbance (450 nm, Ref 620 nm (*④)) immediately.
- ☐ Ref. wave cancels the dirt in the back of plate.

*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume : 300 µL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.


Standard of plate-washing pressure : 5 mL/min - 25 mL/min. (Adjust it depending on the nozzle's diameter.)

- *② Guideline of shaking : 600 rpm - 1200 rpm 3 times for 10 seconds.
- *③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.
- *④ 600 nm - 650 nm can be used as reference wavelength.
- *⑤ After removal of wash buffer, immediately dispense the next reagent.

Worksheet Example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	10 ng/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
B	5 ng/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
C	2.5 ng/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
D	1.25 ng/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
E	0.625 ng/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
F	0.313 ng/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
G	0 (Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40
H	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33	Sample 41

Assay Worksheet

17. Storage and Expiration

The complete kit is stored at 2°C - 10°C (Do not freeze). Opened reagents should be used as soon as possible to avoid loss in optimal assay performance caused by storage environment.

LBIS™ Rat LH ELISA Kit

[Storage] Store at 2 - 10°C (Do not freeze).
 [Expiration date] Indicated on the label
 [Package] For 96 tests
 [Cat #] 294-88501

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://flwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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