## Limulus ES-II Single Test *Wako* 0.015

### Package insert

---

**[Summary and General Information]**

Endotoxin, lipopolysaccharide (LPS) is a component of the outer membrane of Gram-negative bacteria. Since endotoxin causes serious symptoms such as fever and shock in the human body, the bacterial endotoxin test for pharmaceuticals, especially parenteral drugs, is important.

In 1956, F. B. Bang reported clotting of the horseshoe crab body fluid due to Gram negative bacteria<sup>1)</sup>. In 1964, J. Levin and F. B. Bang discovered that horseshoe crab amebocytes (Limulus amebocyte lysate, or LAL) contained all the components needed to clot due to endotoxin<sup>2)</sup>. Subsequently, bacterial endotoxin detection by the Limulus test has been widely used. This method known as the Bacterial Endotoxin Test was recognized in 1980 by United States Pharmacopoeia and, in 1988 by Japanese Pharmacopoeia<sup>3), 4)</sup>.

Kakinuma et al reported that LAL reacted with not only endotoxin but also  $\beta$ -1,3-glucan<sup>5)</sup> and Iwanaga et al reported that the cascade system in LAL activated by  $\beta$ -1,3-glucan was different from the one activated by endotoxin<sup>6), 7)</sup>. It was also reported that a LAL-reactive material was eluted from cellulose membranes, which activated LAL<sup>8)</sup>. The specificity of the LAL test for the detection of endotoxin became a problem.

This kit contains LAL ES-II reagent and Control Standard Endotoxin (CSE) and can specifically and rapidly assay and detect endotoxin with high sensitivity. LAL ES-II reagent is supplied in single use tubes able to assay just 0.2 mL sample with the LAL reagent including the buffer ingredient. In addition, the labelled gel-clot sensitivity was determined using the Japanese Pharmacopoeia Reference Standard Endotoxin (JP-RSE), and this reagent is suitable for an analysis by the gel-clot technique and the kinetic-turbidimetric assay using Toxinometer. CSE contains 500 ng LPS (Lipopolysaccharide, LPS) which has been extracted and purified from *E. coli* UKT-B by phenol method<sup>9), 10)</sup>, and combined with mannitol and glycine as additives, and is tested in comparison to the Japanese Pharmacopoeia Reference Standard Endotoxin and the potency is printed on the label as a reference.

**[Features]**

1. Highly sensitive for specific detection of endotoxin without the interference of  $\beta$ -1,3-glucan in samples.
2. Gel forms strongly and easy to determine.
3. The reagent can assay by the gel-clot technique and the kinetic-turbidimetric assay using Toxinometer. Toxinometer can detect lower endotoxin than the gel-clot technique.
4. Store at 2 ~ 10°C.
5. Single test vials mean there is no reagent wasted.


**[Principle]**

In the gelation of LAL by endotoxin, serine proteases are activated in the cascade reactions. In the last reaction, the gel-forming protein precursor (coagulogen) is hydrolyzed to coagulin, which makes an insoluble gel. (Fig. 1)  $\beta$ -1,3-glucan also activates LAL. However, activation by  $\beta$ -1,3-glucan is dependent on the concentration of  $\beta$ -1,3-glucan itself. In a very high

concentration.  $\beta$ -1,3-glucan does not activate LAL. In such a high concentration  $\beta$ -1,3-glucan does not interfere with the activation of LAL by endotoxin. The LAL ES-II contains a large amount of  $\beta$ -1,3-glucan so that activation by  $\beta$ -1,3-glucan is completely inhibited and therefore, endotoxin can be specifically assayed.

The gel-clot technique is that LAL reagent is mixed with a sample. The mixture is incubated at 37°C for one hour. It is assayed by turning the tube upside down to determine if a gel has formed.

The kinetic-turbidimetric assay is when a sample containing endotoxin is mixed with the LAL reagent. During gelation, turbidity in the reaction mixture increases (i.e., transmittance decreases). The transmittance decrease is measured with the Toxinometer, which is programmed with a threshold value for the percent transmittance. The time it takes for the transmittance to reach the threshold value is also measured by the Toxinometer. This time is defined as gelation time (Tg). Endotoxin is calculated by Tg and endotoxin concentration.

**Fig. 1 the cascade system of Limulus Amebocyte Lysate**

#### **[Kit Contents]**

1. LAL ES-II, lyophilized 25 tubes containing lyophilized extract of *Limulus polyphemus* amebocytes containing Tris buffer and  $\beta$ -1,3-glucan derivatives.  
..... 25 vials (for 0.2 mL)

Gel-clot sensitivity was determined using the Japanese Pharmacopoeia Reference Standard Endotoxin

Store at 2 to 10°C

2. Control Standard Endotoxin, lyophilized 1 vial  $\times$  500 ng (as purified LPS)  
Endotoxin purified from *E. coli* UKT-B, containing mannitol and glycine as additives.  
Store at 2 to 10°C

#### **[Usage]**

##### **A. The gel-clot technique**

- I. Materials required but not supplied
  1. Pipettes (for 0.2 mL, 0.5 mL, 2.0 mL and 5.0 mL)
  2. Test tubes with aluminum caps for dilution of standard solutions
  3. Block heater or an incubator able to maintain  $37 \pm 1^\circ\text{C}$ .
  4. Endotoxin-free water

NOTE : A test tube and the pipette should be dehydrogenated by 250°C and more than 30 min. Water should be Endotoxin-free water. The plastic products such as tips for micropipets should be confirmed that they are no endotoxin contamination, and adsorption or interference with the measurement.

##### **II. Procedure**

1. Preparation of Control Standard Endotoxin (CSE) Solution  
Remove the aluminum cap on the vial. The inside of the vial is a

vacuum. Release the vacuum by pulling slowly up on the rubber stopper and remove the stopper from the vial. The EU value (EU/vial) is printed on the label.

Determine the volume of water to be added into the vial to make a 1,000 EU/mL solution. Add the determined volume of endotoxin-free water into the vial. Put the rubber stopper on the vial. After turning the vial upside down several times, agitate the vial vigorously with a Vortex mixer for approximately 2 minutes. The Control Standard Endotoxin Solution can be stored at 2 ~ 10°C for 1 month. When the stored solution is used, agitate the solution vigorously with Vortex mixer for more than 1 minute before use.

## 2. Dilution of standard endotoxin solutions, the pH adjustment and dilution of samples

A dilute standard solution should be agitated with a Vortex mixer for more than 30 seconds before use in the next dilution procedure. One dilution should not be more than 10-fold.

When the pH of the LAL ES- II solution reconstituted with 0.2 mL of sample solution is out of the range between 6.0 and 8.0, adjust the pH range of sample solution with a dilute sodium hydroxide solution or dilute hydrochloric acid.




## 3. Dispensing of samples

Take LAL ES- II tubes out of the package and put the tubes in a tube rack. When powder of LAL is attached to the rubber stopper, knock on several times of bottoms of the vial. A little LAL attaches to a rubber stopper and does not influence the result of the examination. Each rubber stopper is loosened by releasing the vacuum and slowly pulling up the rubber stopper slightly. (Fig. 2)

Dispense 0.2 mL of a sample into each LAL ES- II tube and put the rubber stopper on the tube. Mix the single test vial with a vortex mixer for 5 seconds to dissolve. Check to dissolved. Incubate the reaction mixture at  $37 \pm 1^\circ\text{C}$  for  $60 \pm 2$  min without vibration. (Fig. 3)

## 4. Judgment

After 60 minute incubation, slowly turn each tube upside down to determine if a gel is formed. A positive result is defined as the formation of a firm gel capable of maintaining its integrity. A negative result is defined as the absence of gel or the formation of a viscous mass, which does not hold. (Fig. 4)

### III. Interfering Testing

The inhibition or enhancement of the LAL test should be determined for each sample formulation before the LAL test is used to assess the endotoxin content of any sample. When the pH of the LAL ES-II solution reconstituted with 0.2 mL of sample solution is out of the range between 6.0 and 8.0, adjust the pH range of sample solution with a dilute sodium hydroxide solution, or dilute hydrochloric acid.

#### (Method)

1. Prepare a series of two-fold dilutions of endotoxin with endotoxin-free water to give concentrations of  $2\lambda$ ,  $\lambda$ ,  $1/2\lambda$ , and  $1/4\lambda$ , where  $\lambda$  is the labeled sensitivity of the LAL reagent. Add 50  $\mu$ L of 200  $\lambda$  endotoxin solution to 4.95 mL of sample solution, which makes a final endotoxin concentration of  $2\lambda$ . Prepare another series of two-fold dilutions of the  $2\lambda$  endotoxin-spiked sample solution diluted with the sample solution to give final endotoxin concentrations of  $2\lambda$ ,  $\lambda$ ,  $1/2\lambda$ , and  $1/4\lambda$ .
2. Perform the test on the above 2 series of dilutions, an unspiked sample solution and a negative control (endotoxin-free water). Please follow "Inhibition and Enhancement Testing" of the station "bacterial endotoxin test" in the Japanese Pharmacopoeia about the number of the repetition.
3. The unspiked sample solutions and the negative controls must be negative. The geometric mean endpoint concentration of the first series of dilutions (diluted with endotoxin-free water) must be greater than or equal to  $1/2\lambda$  and less than or equal to  $2\lambda$ .  
If the negative controls are negative and the unspiked sample solutions are positive, remove endogenous endotoxin by ultrafiltration/appropriate way or dilute the sample solution. And then retest.  
If all the conditions in the above 3. are met and the geometric mean endpoint concentration of the second series of dilutions (diluted with sample solution) is greater than or equal to  $1/2\lambda$  and less than or equal to  $2\lambda$ , it is confirmed that the sample concentration in the sample solution does not interfere with the LAL test.

### IV. Preliminary Test for Dilution Factor

Prepare a series of dilutions (determine an appropriate dilution factor according to the degree of interference) of the sample solution with endotoxin-free water. Add endotoxin solution to each of the dilutions making the final endotoxin concentration  $2\lambda$ . Perform the test on 2 series of dilutions (a series of unspiked sample solutions and a series of  $2\lambda$  endotoxin-spiked sample solutions). If a dilution of the spiked sample solution is positive and the same dilution of the unspiked sample solution is negative, it is confirmed that no interference is observed the sample concentration of that dilution.

### V. Routine Testing (Limit Test)

Perform the test on the following specimens at least in duplicate ( $n=2$ ).

- Negative control (Endotoxin-free water)
- A dilution of endotoxin ( $2\lambda$ ) (Positive Control : PC)
- Sample solution
- Positive product control ( $2\lambda$  endotoxin-spiked sample solution) (PPC)

If all the following conditions are met, the endotoxin concentration of the sample solution is passed.

- Negative control is negative
- PC and PPC are positive

- Sample is negative

If the sample solution is positive, perform the test on a series of 2-fold dilutions of sample solution diluted with endotoxin-free water.

The endotoxin concentration in the sample solution is determined by :  

$$[\lambda \text{ (EU/mL)}] \times [\text{Maximum dilution factor in the series of dilutions of positive results}]$$

## **B. The kinetic-turbidimetric assay using Toxinometer<sup>®</sup>**<sup>12), 13), 14)</sup>

### **1. Preparation of Reagents**

Prepare reagents as stated in the "Preparation of Reagents" section.

### **2. Toxinometer<sup>®</sup> Parameters**

Set the Toxinometer<sup>®</sup> parameters as shown :

Target temperature : 37°C  
 Threshold value (Th) : 94.9%  
 Count : 3  
 Wait time<sup>\*1</sup> : 5 minutes  
 Measurement time<sup>\*2</sup> : 60 minutes

\*1 : When measuring at high concentrations, it is recommended to shorten the Wait time and change the setting to "Use the first count as Reaction time."

\*2 : Please note that the higher the sensitivity, the more time may be required for the endotoxin to be measured by the Toxinometer<sup>®</sup>.

### **3. Assay procedure**

- 1) Remove the appropriate number of LAL ES-II reagent out of the kit box and place each in the test tube rack.
- 2) Remove the cap by using the cap opener provided in the kit box.
- 3) If the lyophilized product is attached to the rubber stopper, gently tap the bottom of the vial on a firm surface before opening the vial. Small amounts of lyophilized product attached to the rubber stopper will not affect testing.
- 4) Use cap opener to pull up rubber stopper releasing the vacuum. If necessary, remove the rubber stopper and cover vial with an endotoxin-free aluminum cap.
- 5) Dispense 0.2 mL of water for BET (negative control), standard solutions, and/or samples into each of the lyophilized LAL ES-II reagent.
- 6) Return the rubber stopper or the aluminum cap on to the single test vial.
- 7) Mix the single test vial with a vortex mixer for 5 seconds to dissolve. Check to dissolved.
- 8) Place each vial into the Toxinometer<sup>®</sup> to be measured by the Toximaster<sup>®</sup> software.

### **4. Standard curve and Data analysis**

The endotoxin concentration of the sample is based on the sample's gel time using the generated standard curve.

A standard curve will be generated by taking the logarithm of the endotoxin concentration on the x axis and the logarithm (log-log) or twice the logarithm (log-loglog) of the gel time on the y axis. Linear or quadratic regression can be used for regression calculation. If the range of the standard curve is wide, improvement of the standard curve fitting can be expected by using twice the logarithm of gel time or quadratic regression.


Figure 2 : An example of standard curve of endotoxin measured with LAL ES-II reagent using a Toxinometer®. This figure shows that a high correlation standard curve can be obtained over a wide concentration range.

#### **【Precautions】**

1. Be very careful not to cause endotoxin contamination. Use endotoxin-free equipment and water
2. Do not use the reagent with a large quantity of insoluble material or when it changes color.
3. For research used only. For in vitro use only. Not for use in diagnostic procedures. Not for internal use in human or animals.
4. Exercise caution when handling LAL, because its toxicity is not known.
5. Even a trace amount of  $\beta$ -1,3-glucan derivatives, which is included in the LAL ES-II, activates normal LAL reagent that is sensitive to  $\beta$ -glucan. Be careful not to contaminate normal LAL with LAL ES-II.

#### **【Procedure to Determine Potency of CSE (Reference)】**

1. Prepare a series of dilutions of Japanese Pharmacopoeia Reference Standard Endotoxin (JP-RSE) according to Table 1. Perform the test on the dilutions of 0.0039 ~ 0.0313 EU/mL in duplicate using the Toxinometer.  
A standard curve,  $\log (Tg \text{ (min)}) = A \times \log (\text{JP-RSE conc. (EU/mL)}) + B$ , will be obtained by the least squares method. The coefficient of correlation ( $r$ ) of the standard curve must be less than  $-0.980$ . The valid range of the standard curve is between mean Tg of the maximum concentration and the minimum concentration.
2. Prepare a series of dilutions of each of 4 vials of CSE according to Table 2. Perform the test on all the dilutions in duplicate by the Toxinometer. Calculate mean Tg (Tg [mean]) for each concentration. Tg's of 3 or more concentrations of CSE solution must be within the valid Tg range of the JP-RSE standard curve.
3. Calculate mean Tg [mean] of 4 vials for each CSE concentration (ng/mL). Convert ng/mL to EU/mL. Calculate EU/ng. See Tables 3 ~ 5 for a calculation example.  
The potency which calculated with the lot of LAL ES-II and CSE is printed on the label as a reference.  
Note : If a LAL lot is changed, the CSE potency test must be repeated because EU/ng of the CSE may be changed.

Table 1. Endotoxin dilution series of JP-RSE

JP-RSE soln. (mL)	H <sub>2</sub> O (mL)	final JP-RSE conc. (EU/mL)
0.5 (10000 EU/mL)	4.5	1000
0.5 (1000 EU/mL)	4.5	100
0.5 (100 EU/mL)	4.5	10
0.5 (10 EU/mL)	4.5	1.0
1.0 (1.0 EU/mL)	3.0	0.25
1.0 (0.25 EU/mL)	3.0	0.0625
2.0 (0.0625 EU/mL)	2.0	0.0313
2.0 (0.0313 EU/mL)	2.0	0.0156
2.0 (0.0156 EU/mL)	2.0	0.0078
2.0 (0.0078 EU/mL)	2.0	0.0039

Table 2. Endotoxin dilution series of CSE

CSE soln. (mL)	H <sub>2</sub> O (mL)	final CSE conc. (ng/mL)
0.5 (100 ng/mL)	4.5	10
0.5 (10 ng/mL)	4.5	1.0
0.5 (1.0 ng/mL)	4.5	0.1
1.0 (0.1 ng/mL)	3.0	0.025
2.0 (0.025 ng/mL)	2.0	0.0125
2.0 (0.0125 ng/mL)	2.0	0.00625
2.0 (0.00625 ng/mL)	2.0	0.00313
2.0 (0.00313 ng/mL)	2.0	0.00156

CSE calculation example

Table 3. Endotoxin conc. and Tg of JP-RSE

JP-RSE (EU/mL)	Tg (min)
0.0039	53.0 52.6
0.0078	44.4 44.0
0.0156	37.4 37.0
0.0313	29.6 29.4

$$\log (\text{Tg (min)}) = -0.276 \log (\text{JP-RSE conc. (EU/mL)}) + 1.061$$

$$r = -0.9973$$

Table 4. CSE conc. (ng/mL) and Tg [mean] (min)

CSE (ng/mL)	Tg [mean] (min)			
	vial 1	vial 2	vial 3	vial 4
0.00156	49.5	52.6	49.4	49.3
0.00313	39.4	45.3	40.0	41.1
0.00625	32.0	34.2	33.0	31.1
0.0125	25.8	27.3	27.3	26.2
0.025	21.2	22.4	21.2	21.0

Table 5. CSE Potency

CSE (ng/mL)	Tg [mean] (min)	Potency (EU/mL)	Potency (EU/ng)
0.00156	50.2	0.0049	3.14
0.00313	41.5	0.0097	3.10
0.00625	32.6	0.0231	3.71
0.0125	26.7	do not use	do not use
0.025	21.5	do not use	do not use

Calculated EU [mean] 3.32 (EU/ng)

$500 \text{ (ng/vial)} \times 3.32 \text{ (EU/ng)} = 1660 \text{ (EU/vial)}$

Round to 3 significant digits

Therefore CSE potency is 1700 (EU/vial).

#### [References]

1. Bang, F.B. : *Bull. Johns Hopkins Hosp.*, **98**, 325 (1956).
2. Levin, J. and Bang, F. B. : *Bull. Johns Hopkins Hosp.*, **115**, 265 (1964).
3. *The United States Pharmacopeia 41st Edition*, U.S. Pharmacopeial Convention Inc.
4. The Japanese Pharmacopoeia, Eighteenth Edition (2021).
5. Kakinuma, J., Asano, T., Torii, H. and Sugino, Y. : *Biochem. Biophys. Res. Commun.*, **101**, 434 (1981).
6. Nakamura, T., Morita, T., Hiranaga, M., Miyata, T., and Iwanaga, S. : *Nihon Saikingaku Zasshi.*, **38**, 781 (1983). Japanese.
7. Iwanaga, S., Morita, T., Miyata, T., Nakamura, T., Hiranaga, M. and Ohtusbo, S. : *Bacterial Endotoxin*, Eds. Homma, J.Y., Kanegasaki, S., Luderitz, O., Shiba, T., and Westpha, O., p.365 Verlag Chemie (1984).
8. Pearson, F.C., Weary, M. and Bohon, J. : *Endotoxin and Their Detection with the Limulus Amebocyte Lysate Test*, Eds. Watson, S.W., Levin, J. and Novitsky, T.J., p.247, Alan R. Liss, Inc. (1982).
9. Westphal, O., Luderiz, O., Eichenberger, E. and Keiderling, W. : *Z. Naturforsch.*, **7b**, 536 (1952).
10. Akama, K., Kuratsuka, K., Homma, R., Kanoh, S., Niwa, M., Iwanaga, S. and Nakahara, C. : *Bacterial Endotoxin*, Eds. Homma, J.Y. *et al.*, p.395, Verlag Chemie (1984).
11. Tsuchiya, M., Takaoka, A., Tokioka, N., and Matsuura, S. : *Nihon Saikingaku Zasshi.*, **45**, 903 (1990). Japanese.
12. Oishi, H., Hatayama, Y., Shiraishi, H., Yanagisawa, K., and Sakata, Y. : *Yakugakuzasshi.*, **105**, 300 (1985). Japanese.
13. Oishi, H., Takaoka, A., Hatayama, Y., Matsuo, T. and Sakata, Y. : *J. Parenter. Sci. Technol.*, **39**, 194 (1985).
14. Oishi, H., Fusamoto, M., Hatayama, Y., Tsuchiya, M., Takaoka, A. and Sakata, Y. : *Chem. Pharm. Bull.*, **36**, 3012 (1988).

---

#### Manufactured by FUJIFILM Healthcare Manufacturing Corporation

2613-2, Ogohara, Komono-Cho, Mie-Gun, Mie 510-1222, Japan

#### Distributed by FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://ffwk.fujifilm.co.jp>

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