FUJIFILM

< For Research Use Only >

Wako

Code No. 295-71701

For Genetic Research
microRNA Extractor® SP Kit

[Introduction]

microRNA Extractor® SP Kit is designed for the extraction of microRNA from serum and plasma of human and rodents (rat etc.). The extraction method of this kit is based on the procedure using chaotropic reagent, which allows for a safe extraction without hazardous phenol/chloroform, and spin column. This kit is useful microRNA extraction kit for circulating microRNA in serum and plasma.

[Features]

- microRNA extraction from serum and plasma of human and rodents
- Safe procedure without hazardous phenol/chloroform
- Easy protocol by using protease, chaotropic reagent and spin column
- Preparation of RNA fraction for real-time qPCR (TaqMan® MicroRNA Assay)

[Kit contents] (50 reactions)

This kit includes a total of 9 components.

(1) Lysis Solution	19 mL	× 1 bottle
(2) Reducing Agent	50 µL	× 1 vial
(3) Protease Solution	500 µL	× 1 vial
(4) Enhancer	500 µL	× 1 vial
(5) Wash Solution I	12 mL	× 1 bottle
(6) Wash Solution II	21 mL	× 1 bottle
(7) Elution Solution	5 mL	× 1 bottle
(8) Spin Column		50 columns
(9) Collection Tube		50 tubes

[Additional materials required]

Reagents :

- 1) 2-Propanol (Code No. 168-21675)
- 2) Ethanol (Code No. 054-07225)
- 3) 1-Butanol (Code No. 022-16035)

Equipment :

- 1) 1.5mL Micro centrifugation tube
- 2) Micro pipette
- 3) Pipette tip
- 4) High speed micro-centrifuge (Max > 14,500 rpm, 20,000 × G)
- 5) Vortex mixer

[Important information before experiment]

1. Equipment

PCR tube, micro centrifugation tube, pipette tip must be sterilized by autoclave. There is no problem with using commercially available RNase free products. We recommend that you must use gloves and a mask to avoid contamination of RNase.

2. Reagents

Please use RNase free reagents for RNA extraction procedure.

3. Handling of infectious wastes

Please correctly handle human serum, plasma, extraction equipment after the procedure as an infectious waste.

[Procedure]

Preparation of serum and plasma sample

- 1) Add enough volume of serum or plasma^{※1} into 1.5 mL microtube.
- 2) Centrifuge at 20,000 × G for 5 min at 4 °C.
- 3) Collect 200 µL of supernatant after centrifugation. ([**SAMPLE**])
The 200 µL is appropriate volume for an extraction.

※1 : We recommend you use EDTA or citric acid as an anticoagulant. We can NOT recommend to use heparin, because heparin will inhibit PCR after RNA extraction. To avoid RNA degradation by RNases, serum and plasma sample has been handled on ice until addition of Lysis solution.

Preparation of reagents

- [**Wash Solution I (+ 2-Pro)**]
Add 18 mL of 2-Propanol into the bottle of Wash Solution I and mix gently.
 - [**Wash Solution II (+ Etha)**]
Add 63 mL of Ethanol into the bottle of Wash Solution II and mix gently.
 - [**Lysis Solution (+ Reducing Agent)**]
Add 2.6 µL of Reducing Agent into 1 mL of Lysis Solution^{※2} and mix by vortex mixer^{※3}. The required volume of Lysis Solution will be changed by the number of samples.
Therefore, please calculate the required volume of Lysis Solution and transfer this into a centrifugation tubes before mix with Reducing Agent.
- ※2 : The precipitants will be appeared in Lysis Solution during the storage.
In this case, the precipitants have to be completely dissolved on 37 °C by water bath etc.
- ※3 : [**Lysis Solution (+ Reducing Agent)**] is stable in 2~10 °C for 1 week.

microRNA extraction procedures

- 1) Add 10 µL of Protease Solution into a new 1.5 mL micro centrifugation tube.
- 2) Add 200 µL of [**SAMPLE**] into a tube of step 1).
- 3) Add 200 µL of [**Lysis Solution (+ Reducing Agent)**]^{※4} and mix by vortex mixer and spin down by mini centrifuge.
- 4) Incubate for 30 minutes on heat block at 37 °C and spin down by mini centrifuge.
- 5) Add 200 µL of 1-Butanol^{※5} and mix by vortex mixer.
- 6) Centrifuge at 12,000 × G for 10 minutes at 4 °C^{※6}.
- 7) Add 10 µL of Enhancer into a new 1.5 mL micro centrifugation tube and transfer the all of supernatant of step 6) into this tube by decantation.
- 8) Add 300 µL of 2-Propanol and mix by vortex mixer.
- 9) Spin down for 2-3 seconds by mini centrifuge^{※7}.
- 10) Mix by pipetting^{※8} and transfer a half volume of the sample solution into Spin Column.
- 11) Centrifuge at 8,000 × G for 1 minute at room temperature. Remove the filtrated solution from Collection Tube by decantation.
- 12) Mix by pipetting^{※8} and add the remaining sample solution into Spin Column.

- 13) Repeat Step 11).
- 14) Add 500 μ L of **Wash Solution I (+2-Pro)** into Spin Column.
- 15) Repeat Step 11).
- 16) Add 700 μ L of **Wash Solution II (+Ethanol)** into Spin Column.
- 17) Repeat Step 11).
- 18) Repeat Step 16) and Step 11).
- 19) Set the Spin Column into Collection Tube^{※9} and centrifuge at $8,000 \times G$ for 3 minutes at room temperature. Remove the Wash Solution completely.
- 20) Set the Spin Column into a new 1.5 mL micro centrifugation tube.
- 21) Add 50 μ L^{※10} of Elution Solution into the center of Spin Column membrane and incubate for 5 minutes at room temperature.
- 22) Centrifuge at $8,000 \times G$ for 1 minute at room temperature. Remove the Spin Column from 1.5 mL micro centrifugation tube. This sample is final RNA fraction^{※11}.

※4 : Do NOT add the **Lysis Solution (+Reducing Agent)** into Protease Solution directly. We recommend you add these reagents by the following.

Protease Solution $\rightarrow$ **SAMPLE** $\rightarrow$ **Lysis Solution (+Reducing Agent)**

※5 : If you could not prepare 1-Butanol, you can use 2-Propanol as an alternative reagent. In the case of using 2-Propanol, the quantity of contaminating protein will be increased in final RNA fraction. We recommend you use 1-Butanol as much as possible.

※6 : The centrifuge has to be set up from 4°C to room temperature for the following procedure.

※7 : Avoid the spin down for a long time.


※8 : This pipetting procedure is necessary.

※9 : To avoid the contamination of alcohol in final RNA fraction, we recommend you use a new 1.5 mL micro centrifugation tube.

※10 : In the case of using 100 μ L of Elution Solution, recovery efficiency of microRNA will be increased 10~30% more than standard procedure. However, the concentration of microRNA will be decreased. Please optimize the adding volume of Elution Solution to your following experiment.

※11 : Store the final RNA fraction at -80°C until the analysis by real-time qPCR.

〔Extraction flow chart〕

〔Storage〕

Store at $2 \sim 10^{\circ}\text{C}$.

〔Products〕

Code No.	Description	Package Size
295-71701	microRNA Extractor [®] SP kit	50 extractions
168-21675	2-Propanol	500 mL
054-07225	Ethanol(99.5)	500 mL
022-16035	1-Butanol	500 mL

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : +1-804-271-7677
 Facsimile : +1-804-271-7791
<http://www.wakousa.com>

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