Code No. 295-73901(100 assays)

LysoPure™ Nuclear and Cytoplasmic Extractor Kit

LysoPure™ Nuclear and Cytoplasmic Extractor Kit enable stepwise extraction of cytoplasmic and nuclear fractions (soluble and insoluble) from mammalian cells or tissue.

A sample is separated into cytoplasmic and nuclear fraction with Nuclear Fractionation Buffer (1). Soluble intact nuclear fraction is extracted from the nuclear fraction with Nuclear Extraction Buffer (2). Insoluble nuclear fraction is solubilized and extracted as a denatured protein fraction from the remaining nuclear fraction with SDS Lysis Buffer (3). The denatured proteins are allowed to apply western blot assay.

[Kit Components]

- ① Nuclear Fractionation Buffer 100 mL × 1 bottle
- ② Nuclear Extraction Buffer 75 mL × 1 bottle
- ③ SDS Lysis Buffer 20 mL × 1 bottle


[Storage]

2-10°C

[Package]

100 assay

[Protocol Flowchart]

[Additional Materials Required]

-Instruments/labware

- Refrigerated microcentrifuge
- Vortex mixer
- Desktop centrifuge
- Ultrasonic sonicator
- 1.5 mL microcentrifuge tubes
- Micropipette
- Pipette tips

- Cell scraper (for cells sample)
- Ultrasonic sonicator
- Dounce or potter homogenizer (for tissue sample)

-Reagents

- PBS (-)
- Trypsin-EDTA (for cells sample)

[Please Read Before Starting]

-Instruments/labware

- Use sterilized microcentrifuge tubes and pipette tips.
- To avoid protease contamination, be sure to wear a mask and plastic gloves during experiment.
- Microcentrifuge tubes must be pre-chilled on ice.

-Reagents

- Be careful not to contaminate reagents with protease.
- Keep reagents in this kit on ice during experiment.

-Treatment of tissue sample from human and laboratory waste

- Please deal with a human tissue sample as a susceptibility infectious sample. This includes waste liquid and equipments such as microcentrifuge tubes, pipette tips and gloves after experiment. Dispose of these biohazardous wastes according to the guidelines of the institution which belong to.

[Protocol]

Sample Preparation

○ For adherent cells

- 1) Culture cells to be $2 \times 10^6 - 1 \times 10^7$ *1
- 2) Remove culture medium from cell culture dishes or plates and wash the cells with PBS (-) twice.
- 3) Harvest cells with trypsin-EDTA or a cell scraper.
- 4) Centrifuge at $200 \times g$ at 4°C for 5 minutes and discard the supernatant.
- 5) Suspend the cells pellet with 1 mL of PBS (-) and transfer it into pre-chilled 1.5 mL microcentrifuge tube.
- 6) Centrifuge at $200 \times g$ at 4°C for 5 minutes and discard the supernatant.
- 7) Suspend the cell pellet with 1 mL of PBS (-).
- 8) Centrifuge at $200 \times g$ at 4°C for 5 minutes and discard the supernatant.

○ For suspension cells

- 1) Culture cells to be $2 \times 10^6 - 1 \times 10^7$ *1
- 2) Centrifuge at $200 \times g$ for 5 minutes at 4°C and discard the supernatant.
- 3) Suspend the cells pellet with 1 mL of PBS (-) and transfer it into a pre-chilled 1.5 mL microcentrifuge tube.
- 4) Centrifuge at $200 \times g$ at 4°C for 5 minutes and discard the supernatant.
- 5) Suspend the cells pellet with 1 mL of PBS (-).
- 6) Centrifuge at $200 \times g$ at 4°C for 5 minutes and discard the supernatant.

*1 It is recommended to determine an appropriate culture medium and cell number depending on cell lines of interest. The maximum number of cell is approximately 2×10^7 /assay.

○ For tissue

- 1) Cut 20-100 mg of tissue into small pieces and transfer it into a pre-chilled 1.5 mL microcentrifuge tube.
- 2) Suspend the tissue with 1mL of PBS (-).
- 3) Centrifuge at $500 \times g$ at 4°C for 5 minutes and discard the supernatant.
- 4) Repeat steps 2 and 3.

Nuclear and Cytoplasmic Protein Extraction

Table.1 Reagent volumes for different cell numbers and tissue weight

Cell number	① Nuclear Fractionation buffer (μ L)	② Nuclear Extraction buffer (μ L)	③ SDS Lysis buffer (μ L)
2×10^6	100	50	50
5×10^6	250	125	100
1×10^7	500	250	200

Tissue weight (mg)	① Nuclear Fractionation buffer (μ L)	② Nuclear Extraction buffer (μ L)	③ SDS Lysis buffer (μ L)
20	100	50	100
50	500	250	200
100	1000	500	400

○ For Cells

Cytoplasmic Fraction

- 1) Add the appropriate volume of “① Nuclear Fractionation Buffer” depending on the cell number (See Table. 1) to the cells pellet and vortex the tube.
- 2) Incubate on ice for 10 minutes.
- 3) Mix by pipetting or vortexing.
- 4) Centrifuge at $500 \times g$ at 4°C for 10 minutes.
Please use the pellet as a sample for nuclear protein extraction.
- 5) Transfer the supernatant into a pre-chilled 1.5 mL microcentrifuge tube and centrifuge at $20,000 \times g$ at 4°C for 10 minutes.
- 6) Transfer the supernatant into a pre-chilled 1.5 mL microcentrifuge tube (Cytoplasmic Fraction).^{*2}
Proceed the next step to extract nuclear fraction if needed.

Soluble Nuclear Fraction

- 7) Add $500 \mu\text{L}$ of “① Nuclear Fractionation Buffer” to the pellet prepared in the above-mentioned Step 4) and vortex the tube. Be careful not to use “② Nuclear Extraction Buffer” in this step.
- 8) Centrifuge at $500 \times g$ at 4°C for 10 minutes and discard the supernatant carefully.
- 9) Add the appropriate volume of “② Nuclear Extraction Buffer” (See Table. 1) to the pellet and vortex sufficiently.
- 10) Incubate on ice for 30 minutes, as vortex every 10 minutes.
- 11) Centrifuge at $20,000 \times g$ at 4°C for 10 minutes.
- 12) Transfer the supernatant into a pre-chilled 1.5 mL microcentrifuge tube (Soluble Nuclear Fraction).^{*2, *3}

Insoluble Nuclear Fraction

- 13) Add $500 \mu\text{L}$ of “② Nuclear Extraction Buffer” to the pellet produced in the above-mentioned Step 11) and vortex the tube. Be careful not to use “③ SDS Lysis Buffer” in this step.
- 14) Centrifuge at $20,000 \times g$ at 4°C for 10 minutes and discard the supernatant carefully.
- 15) Add the appropriate volume of “③ SDS Lysis Buffer” (See Table. 1) to the pellet and vortex the tube.
- 16) Sonicate the tube with ultrasonic sonicator until the pellet is disappear.
- 17) Centrifuge at $20,000 \times g$ at 4°C for 10 minutes.
- 18) Transfer the supernatant into a pre-chilled 1.5 mL microcentrifuge tube (Insoluble Nuclear Fraction).^{*2}

○ For Tissue

- 1) Add the appropriate volume of “① Nuclear Fractionation Buffer” depending on the tissue weight (See Table. 1) to the tissue and homogenize it with a dounce or potter homogenizer.^{*4}
- 2) Transfer the tissue suspension into a pre-chilled 1.5 mL microcentrifuge

fuge tube.

- 3) Incubate on ice for 10 minutes.
- 4) Mix by pipetting or vortexing.
- 5) Centrifuge at $500 \times g$ at 4°C for 10 minutes.
Please use the pellet as a sample for nuclear protein extraction.
- 6) Transfer the supernatant into a pre-chilled 1.5 mL microcentrifuge tube and centrifuge at $20,000 \times g$ at 4°C for 10 minutes.

〈Optional protocol〉

After centrifugation, if cloudiness appears in the supernatant, it is mostly removable by a centrifugal filter (Millipore, Ultrafree-MC $0.45 \mu\text{m}$, UFC30HV00).

- 7) Transfer the supernatant into a pre-chilled 1.5 mL microcentrifuge tube (Cytoplasmic Fraction).^{*2}
Proceed the next step to extract nuclear fraction if needed.
- 8) Add $500 \mu\text{L}$ of “① Nuclear Fractionation Buffer” to the pellet produced in the previous step 5) and vortex. Be careful not to use “② Nuclear Extraction Buffer” in this step.
- 9) Proceed step 8) in Soluble Nuclear Fraction for cells.

- *2 Store extracts -80°C until use. Use protease inhibitors as needed.

*3 Soluble nuclear fraction is extracted under high salt concentration. Dilute it or exchange buffers as needed.

*4 The maximum weight of tissue is approximately 100 mg/assay.

〈Troubleshooting〉

Problem	Possible Cause	Solution
Low cytoplasmic protein	Cells or tissue was not sufficiently lysed.	Increase amount of "Nuclear Fractionation Buffer". Vortex thoroughly and increasing vortexing time.
Low soluble nuclear protein	Reagents were confused.	Check the protocol and reagents' discription.
	Pellet was not sufficiently lysed.	Increase amount of "Nuclear Extraction Buffer". Vortex thoroughly and increasing vortexing time.
Low insoluble nuclear protein	The reagents were confused.	Check the protocol and reagents' discription.
	Pellet was not disappear.	Increase amount of "SDS Lysis Buffer". Increasing sonicating time.
No or low protein activity detected	Proteases were activated.	Keep sample and reagents on ice during experiment. Immediately, store extracts at -80°C . Use protease inhibitor as needed.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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