LBIS™ Mouse OVA-IgE ELISA Kit

Please, read this instruction carefully before use.

1. Intended use

LBIS™ Mouse OVA-IgE ELISA Kit is an ELISA system for quantitative measurement of mouse anti-OVA-IgE antibody titer. This is intended for research use only.

2. Introduction

IgE (immunoglobulin E) is the 5th immunoglobulin found and is composed of 2 heavy chains which contain 5 domains (VH, CH ϵ 1-4) and 2 light chains. It is a glycoprotein with molecular weight of 190000 and in electrophoresis it moves to γ 1 region. Metabolic half life of IgE is about 3 days in man, and normal serum level in human is about 300 ng/mL, however, it is markedly increased in parasite infection and hay fever. Allergy-related IgE is called reagin. The Fc region of reagin increased after sensitization with allergens will bind Fc ϵ R1 receptor of basophilic granulocytes and mast cells in the skin, respiratory, and digestive organs, and sensitizes those cells. Those IgE-sensitized cells will be degranulated when the second allergens bind the surface IgE, and release histamine, serotonin, protease, heparin, chemotactic factors, prostaglandins, leukotrienes, etc., causing type I allergy reactions. FUJIFILM Wako Pure Chemical Corporation's OVA-IgE ELISA Kit is a useful tool for studying mouse immune system by measuring specific anti-OVA IgE after immunization of mice with OVA (ovalbumin).

3. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Mouse OVA-IgE ELISA Kit, biotinylated anti-mouse IgE antibody, standards or samples are incubated in OVA-coated wells to capture OVA-IgE. After 1 hour incubation and washing, HRP (horse radish peroxidase)-conjugated avidin is added, and incubated for 30 minutes together with captured anti-mouse OVA-IgE antibody. After washing, HRP-complex remaining in wells is reacted with a chromogen (TMB) for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to anti-mouse OVA-IgE antibody titer. The standard curve is prepared by plotting absorbance against standard OVA-IgE concentrations. The concentrations in unknown samples are determined using this standard curve.

4. Performance Characteristics

- Assay range
The assay range of the kit is 1.88 U/mL - 120 U/mL.
- Specificity
The Biotin-conjugated Antibody Solution of this kit is specific to mouse IgE.
- Precision of assay
Within assay variation (2 samples, 5 replicates assay) Mean CV was within 10%.
- Reproducibility
Between assay variation (3 samples, 4 days, duplicate assay) Mean CV was within 10%.
- Recovery test
Anti-mouse OVA-IgE was added in 3 concentrations to 2 serum samples and was assayed.
The recoveries were 95.8% - 106%
- Dilution test
2 serum samples were serially diluted by 3 steps.
The dilution curves showed excellent linearity. ($R^2=0.9987$ - 0.9999)

5. Reference Assay Data

Mouse OVA-IgE antibody titer's mean assay value : 139 U/mL, SD : 22.5 U/mL

Strain : BALB/c, 3 males, 8 week-old

OVA administration

: Equal volumes of alum (20 mg/mL) and OVA (50 μ g/mL) were mixed, and 0.2 mL/head of the mixture was intraperitoneally injected twice with 1 week interval, then blood sampling was made 3 weeks later. OVA was first solubilized with 0.1 M carbonate buffer pH 8.5 at a concentration of 1 mg/mL, then diluted with saline to make 50 μ g/mL.

These data should be considered as guidance only. Each laboratory should establish its own normal and pathological reference ranges for OVA-IgE levels independently.

6. Precautions

- For professional use only, beginners are advised to use this kit under the guidance of experienced person.
- Wear gloves and laboratory coats when handling assay materials.
- Do not drink, eat or smoke in the areas where assays are carried out.
- In treating assay samples of animal origin, be careful for possible biohazards.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the Stop Solution because it is sulfuric acid. The Stop Solution and the TMB Solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Avoid contact with the acidic Stop Solution and TMB Solution, containing hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye and clothing protection when handling these reagents.
- The materials must not be pipetted by mouth.
- Residual samples and used tips should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- Use clean laboratory glassware.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.

7. Reagents supplied

Components	Use Status	Amount
(A) Antigen-coated Plate	Use after washing	96 wells/1 plate
(B) Anti OVA- IgE Standard	Concentrated. Use after dilution.	100 μ L/1 bottle
(C) Buffer	Ready for use.	60 mL/1 bottle
(D) Biotin-conjugated Antibody Solution	Concentrated. Use after dilution.	200 μ L/1 bottle
(E) Peroxidase-conjugated Avidin Solution	Concentrated. Use after dilution.	200 μ L/1 bottle
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 $\times$)	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate Seal	—	3 sheets

【Storage and stability】

[(A) Antigen-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C.

[(B) Anti OVA- IgE Standard (1200 U/mL)]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

Unused working solution (already diluted) should be disposed.

[(C) Buffer] and [(F) TMB Solution]

If not opened, store at 2°C - 10°C. It maintains stability until expiration date. Once opened, we recommend using them as soon as possible to avoid influence by environmental condition.

[(D) Biotin-conjugated Antibody Solution] & [(E) Peroxidase-conjugated Avidin Solution]

Unused working solution (already diluted) should be disposed.

[(H) Stop Solution]

Close the stopper tightly and store at 2°C - 10°C.

[(I) Wash Solution (10 $\times$)]

The rest of undiluted buffer : store at 2°C - 10°C. Dispose any unused diluted buffer.

8. Equipments required but not supplied ☐ Use as a check box

☐ Deionized water (or Distilled water) ☐ Test tubes for preparation of standard solution series.

☐ Glassware for dilution of Wash Solution (10×) (a graduated cylinder, a bottle)

☐ Pipettes (disposable tip type). One should be able to deliver 5 μ L - 10 μ L precisely, and another for 100 μ L. ☐ Syringe-type repeating dispenser like Eppendorf multipipette plus which can dispense 50 μ L and 100 μ L. ☐ Paper towel to remove washing buffer remaining in wells. ☐ A vortex-type mixer. ☐ A shaker for 96 well-plate (600 rpm - 1200 rpm) ☐ An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. ☐ A 96 well-plate reader (450 nm $\pm$ 10 nm, 620 nm : 600 nm - 650 nm) ☐ Software for data analysis.

9. Preparation of Samples

This kit is intended to measure anti-mouse OVA-IgE antibody titer in mouse serum or plasma.

Samples should be immediately assayed or stored below -35°C for several days. Defrosted samples should be mixed thoroughly for best results.

Hemolytic and hyperlipemic samples are not suitable.

***To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis, do not use them for assay. Abnormal value might be obtained with hemolysis above 160 mg/dL with this kit.**

If presence of interfering substance is suspected, examine by dilution test at more than 2 points. Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.

Make sure to dilute samples more than 10 times. Recommended is 10 - 50 times depending on the antibody titer. Dilution should be carried out with the Buffer of the kit using small test tubes before assay.

Example of dilution :	Rate	10 times	20 times	50 times	
	Sample (μ L)	5	25*	20*	
	Buffer (μ L)	45	25	30	*One rank lower diluted sample

Storage and Stability

Sample is stable at 2°C - 10°C within a week. If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35°C. Avoid repeated freezing and thawing cycles.

10. Preparation of Reagents

◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.

◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated Reagents]

[(B) Anti OVA- IgE Standard (1200 U/mL)]

Make a serial dilution of original standard solution to prepare each standard solution. Example is shown below.

Volume of standard solution	(C) Buffer	Concentration (U/mL)
Original solution : 10 μ L	90 μ L	120
120 U/mL solution : 50 μ L	50 μ L	60
60 U/mL solution : 50 μ L	50 μ L	30
30 U/mL solution : 50 μ L	50 μ L	15
15 U/mL solution : 50 μ L	50 μ L	7.5
7.5 U/mL solution : 50 μ L	50 μ L	3.75
3.75 U/mL solution : 50 μ L	50 μ L	1.88
0 (Blank)	50 μ L	0

*In this kit, 1 U/mL is prescribed to antigen-binding constant (K_a) 2.0×10^8 M⁻¹ antibody 1.3 ng/mL.

[(D) Biotin-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the (C) Buffer to **1 : 100**.

[(E) Peroxidase-conjugated Avidin Solution]

Prepare working solution by dilution of (E) with the (C) Buffer to **1 : 100**.

[(I) Wash Solution (10×)]

Dilute 1 volume of the Wash Solution (10×) to 10 volume with deionized water (or distilled water) to prepare working solution. Example : 100 mL of Wash Solution (10×) and 900 mL of deionized water (or distilled water).

11. Assay Procedure

Remove the cover sheet of the 96 well-plate after bringing up to room temperature.

- (1) Wash the Antigen-coated Plate (A) by filling the well with 300 μ L of washing buffer and discard 3 times (*①), then strike the plate upside-down onto several layers of paper towels to remove residual buffer in the wells.
- (2) Pipette 50 μ L of Biotin-conjugated Antibody Solution (D) to the designated wells. Shake the plate gently on a plate shaker (*②).

- (3) Pipette 10 μ L of diluted samples to the wells designated for samples.
 - (4) Pipette 10 μ L of standards to the wells designated for standards.
 - (5) Shake the plate on a plate shaker (*②).
 - (6) Stick a Plate Seal (*③) on the plate and incubate for 1 hour at 20°C - 25°C.
 - (7) Discard the reaction mixture and rinse wells as step (1).
 - (8) Pipette **100 μ L** of Peroxidase-conjugated Avidin Solution (E) to all wells, and shake as step (5).
 - (9) Stick a Plate Seal (*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
 - (10) Discard the reaction mixture and rinse wells as step (1).
 - (11) Pipette **100 μ L** of TMB Solution (F) to wells, and shake as step (5).
 - (12) Stick a Plate Seal (*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
 - (13) Add **100 μ L** of the Stop Solution (H) to all wells and shake as step (5).
 - (14) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm*) immediately using a plate reader.
- *Refer to 15. Summary of Assay Procedure for notes of *①, *②, and *③.


12. Technical Tips

- In manual operation, proficiency in pipetting technique is recommended.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2°C - 10°C.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- The TMB Solution should be almost clear pale yellow before use. It turns blue during reaction, and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the Antigen-coated Plate is module type of 8 wells $\times$ 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.
- The standard of this kit is anti-OVA-IgE monoclonal antibody. Therefore, it is possible to compare the assay results even if this kit is used in a different laboratory. Principally it is not possible to compare the assay results with other assay kits' results because the standards don't always have the affinity to the same OVA.

13. Calculations

- (1) Prepare a standard curve using semi-logarithmic or two-way logarithmic section paper by plotting absorbance* (Y-axis) against anti-OVA-IgE concentration (U/mL) on X-axis. Physiological or pathological situation of animals should be judged comprehensively taking other examination results into consideration.
- (2) Using the standard curve, read the anti-OVA-IgE concentration of a sample at its absorbance*, and multiply the assay value by dilution factor. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the buffer solution.

*We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.

Mouse anti-OVA-IgE assay standard curve (an example)
Absorbance may change due to assay environment.

14. Trouble Shooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not be added.
- 2) Reagents necessary for coloration such as Biotin-conjugated Antibody Solution, Peroxidase-conjugated Avidin Solution, or TMB Solution might not be added.
- 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotin-conjugated Antibody Solution or Peroxidase-conjugated avidin Solution.
- 4) Contamination of enzyme inhibitor (s).
- 5) Influence of the temperature under which the kits had been stored.
- 6) Excessive hard washing of the well plate.
- 7) Addition of TMB Solution soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- Blank OD was higher than that of the lowest standard concentration (1.88 U/mL).

Possible explanations :

Improper or inadequate washing. (Change washing frequency from 3 times to 4 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated avidin Solution.)

- High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.

- Q-2 : I found 96 well-plate is empty when I opened the box.

A-2 : As this kit is dried type, not preservation stabilizer is added.


15. Summary of Assay Procedure ☐ : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

☐ Bring the well-plate and all reagents to **20°C - 25°C for 2 hours**.

☐ Wash Solution (10×) concentrate must be diluted to **10 times** by deionized water (or distilled water) that returned to 20°C - 25°C.

☐ Standard solution dilution example :

- ☐ Dilute Biotin-conjugated Antibody Solution to **100×** with Buffer (C) returned to 20°C - 25°C.
- ☐ OVA-coated Plate (Dried-plate)
- ☐ ↓ Washing 3 times (*①), (*⑤)
- ☐ Biotin-conjugated Antibody Solution 50 µL
- ☐ ↓ Shaking (*②)
- ☐ Diluted Samples / Standards 10 µL
- ☐ ↓ Shaking (*②), Incubation for 1 hour at room temp. (Standing (*③))
- ☐ Dilute Peroxidase-conjugated avidin Solution (E) to **1 : 100** with Buffer (C) returned to 20°C - 25°C.
- ☐ ↓ Washing 3 times (*①), (*⑤)
- ☐ Peroxidase-conjugated avidin Solution 100 µL
- ☐ ↓ Shaking (*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (*③))
- ☐ ↓ Washing 3 times (*①), (*⑤)
- ☐ TMB Solution 100 µL
(After dispense, the color turns to blue depending on the concentration.)
- ☐ ↓ Shaking (*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (*③))
- ☐ Stop Solution 100 µL
(After dispense, the color turns to yellow depending on the concentration.)
- ☐ ↓ Shaking (*②) (Immediately shake.)
- ☐ Measurement of absorbance (450 nm, Ref 620 nm(*④)) immediately using a plate reader.
(Ref. wave cancels the dirt in the back of plate.)

*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer.
Guideline of washing volume : 300 µL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 3 times to 4 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated avidin Solution.

Standard of plate-washing pressure : 5 mL/min. - 25 mL/min. (Adjust it depending on the nozzle's diameter.)

*② Guideline of shaking : 600 rpm - 1200 rpm 3 times for 10 seconds.


*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the Plate Seal used once.

*④ 600 nm - 650 nm can be used as reference wavelength.

*⑤ After removal of wash buffer, immediately dispense the next reagent.

Worksheet Example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	120 U/mL	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
B	60 U/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
C	30 U/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
D	15 U/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
E	7.5 U/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
F	3.75 U/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
G	1.88 U/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
H	0 (Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40

Assay Worksheet**16. Storage and Expiration**

When the intact kit is stored at 2°C - 10°C. Reagents, once opened, should be used as soon as possible to avoid losing its optimal assay performance by storage environment.

LBIS™ Mouse OVA-IgE ELISA Kit

[Storage] Store the kit at 2 - 10°C (Do not freeze).

[Expiration date] Indicated on the label.

[Package] For 96 tests

[Cat #] 295-91601

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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