LBIS™ Mouse Urinary Albumin Assay Kit

Please, read this instruction carefully before use.

1. Intended use

LBIS™ Mouse Urinary Albumin Assay Kit is a research reagent for an autoanalyzer to measure albumin in mouse urine or serum using anti-mouse albumin antibody. This is intended for research use only.

2. Introduction

Albumin is mostly a simple hydrophilic protein present in cells and body fluids. Albumin is synthesized in the liver, and serum albumin (M.W. 69,000, pI 4.9) occupies 56% - 60% of total serum proteins. Because of its large population, albumin is very important in maintaining plasma osmotic pressure. Albumin can bind hydrophobic physiological substances e.g. fatty acids, bilirubin, and thyroxine and contributes the transfer of these substances. Concentration of albumin in serum is lowered in liver cirrhosis, malnutrition, and pyrexial diseases due to decreased biosynthesis or increased consumption of albumin in, and also by secretion into urine in renal damage. In normal human subjects, excretion of albumin into urine is very little, about 30 mg/day, but increased in glomerulonephritis, nephritic syndrome, and diabetic nephropathy. Urinary albumin sometimes increases in pyrexia, hypertension, congestive heart failure (CHF), and urinary tract infection (UTI). Even in healthy human, a transient upraise of urinary albumin is observed after hard exercise, muscle work, bathing in high temperature water, mental excitement, stress, intake of much protein, and before menstruation. These are called physiological or functional or sports proteinuria (albuminuria). Orthostatic proteinuria (albuminuria) is observed in teenagers. Rare hereditary analbuminemia in human is really albumin deficiency. A model analbuminemia rat has been established by Dr. Sumi Nagase which has been derived from Sprague-Dawley strain, and is called Nagase analbuminemia rat (NAR). Routine measurement of serum albumin is conveniently made with FUJIFILM Wako Pure Chemical's TIA assay kit for automatic analyzers. FUJIFILM Wako Pure Chemical's Albumin ELISA Kit can measure urinary albumin with high sensitivity, and is also applied to *in vitro* albumin biosynthesis system, checking of contamination with albumin in biological active substance preparations obtained by culture system, and liver transplantation experiments using NAR.

3. Assay Principle

In FUJIFILM Wako Pure Chemical's LBIS™ Mouse Urinary Albumin Assay Kit, albumin in samples and antibody of Reagent 2 bring about antigen-antibody reaction, and it makes turbidity. Since this turbidity depends on its antigen concentration, it makes possible to detect albumin concentration by measuring this turbidity at 340 nm/700 nm.

4. Performance Characteristics

- Assay range
The assay range of the kit is 6.17 $\mu\text{g/mL}$ - 500 $\mu\text{g/mL}$.
- Precision of assay
Within assay variation (3 samples, 5 replicates assay), Mean CV was within 7.5%
- Reproducibility
Between assay variation (3 samples, 4 days, 3 replicates assay), Mean CV was within 7.5%
- Dilution linearity test
Sample was serially diluted by 4 steps and measured.
The dilution curve showed linearity with $R^2 = 0.999$.

5. Reference Data

Strain	Week old	Gender	Number	Sample	Mean	Standard deviation
BALB/c	5 w	Male	4	Urine	16.8 $\mu\text{g/mL}$	1.5 $\mu\text{g/mL}$

*Urine was not diluted.

*These data should be considered as guidance only. Each laboratory should establish its own normal and pathological reference ranges for albumin levels independently.

6. Precautions

- Beginners are advised to use this kit under the guidance of experienced person for autoanalyzer.
- Do not use this product for other purpose or in other ways except stated in this manual.
- Note that operation of autoanalyzer is different in each machine. Follow the protocol of your machine.

- Do not drink, eat or smoke in the area where assays are carried out.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Wear gloves and eye and clothing protection when handling these reagents.
- The materials must not be pipetted by mouth.
- Use clean laboratory glassware.
- Do not use with different lot numbers. Performance or character is different in each lot number. Also avoid extension of reagents within same lot as it may influence assay results.
- Store all reagents in a refrigerator at 2°C - 10°C strictly. Do not freeze reagents.

7. Reagents supplied

Components	Use Status	Amount	
		60 tests	250 tests
1) Reagent 1	Ready for use.	18.5 mL / 1 bottle	57.8 mL / 1 bottle
2) Reagent 2	Ready for use.	6.6 mL / 1 bottle	16.5 mL / 1 bottle
3) Mouse Urinary Albumin Standard	Concentrated. Use after dilution.	2 mL / 1 bottle	1.5 mL / 2 bottles
4) Buffer for Standard	Ready for use.	4 mL / 1 bottle	5 mL / 2 bottles

• Storage and Stability

[1) Reagent 1]

Store at 2°C - 10°C.

[2) Reagent 2]

Store at 2°C - 10°C.

[3) Mouse Urinary Albumin Standard]

Store the unused reagents at 2°C - 10°C.

[4) Buffer for Standard]

Store at 2°C - 10°C.

8. Equipments or supplies required but not supplied ☐ Use as a check box

- ☐ An autoanalyzer
- ☐ Purified water (distilled water)
- ☐ Physiological saline
- ☐ Test tubes for preparation of standard solution series
- ☐ Pipettes (disposable tip type)
- ☐ A vortex-type mixer.

9. Preparation of Samples

This kit is intended to measure albumin in mouse urine or serum.

Urinary samples

Run your assay soon after collection of urinary samples. If you would like to store urinary samples for a long time, store them below -35°C. Turbid samples or those containing insoluble materials should be centrifuged or filtered before testing to remove any particulate matter.

Serum samples

Serum samples should be immediately assayed after collection or stored below -35°C for a long-time storage. Dilute samples to the proper rate for assay. Hemolytic and hyperlipemic serum samples are not suitable.

- * If presence of interfering substance is suspected, examine by dilution test at more than 2 points.
- * Dilute samples with purified water for an autoanalyzer, or dilution of samples should be made in a test tube using saline solution. Obtained value should be raised its dilution rate to be assay value.
- * Before starting assay, shake thawed samples sufficiently. Do not repeat freeze-and-thaw cycles.
- * Adjust samples' pH within 6.5 - 7.5.
- * **To avoid influence of blood (high lipid or hemolysis, etc.) , if your original samples have heavy chyle or hemolysis, do not use them for assay. Abnormal value might be obtained with hemolysis above 200 mg/dL with this kit.**

10. Preparation of Reagents

- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

◆ Do not use the reagents after expiration date indicated on the label of the box.

• **Reagents ready for use**

- 1) Reagent 1
- 2) Reagent 2
- 4) Buffer for Standard

• **Concentrated reagents**

- 3) Mouse Urinary Albumin Standard

Dilute 3) Mouse Urinary Albumin Standard serially 3× using 4) Buffer for Standard.

To make 6 points of standard curve specimen, use Buffer for Standard for the 0 μg/mL.

The example is shown below. (* One rank higher standard.)


Concentration (μg/mL)	500	167	55.6	18.5	6.17	0
Standard solution (μL)	orig.sol.: 500	250*	250*	250*	250*	0
Buffer solution (μL)	0	500	500	500	500	500

Diluted standard solutions can be stable within a week if stored 2°C - 10°C.

11. Assay Procedure

For Hitachi 7180 autoanalyzer

* **Note that operation of autoanalyzer is different depending on the type of machines**

12. Technical Tips

- Be careful to avoid any contamination of assay samples and reagents. We recommend using disposal pipette tips.
- Be careful not to make bubbles in samples and every reagent. Check if there are bubbles before setting into an autoanalyzer. If so, remove them using filter paper before starting assay.
- Calibrate the analyzer before assay.
- Used samples and consumables should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal. Also, dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

13. Storage and Expiration

The complete kit is stored at 2°C - 10°C. Opened reagents should be used as soon as possible to avoid loss in optimal assay performance caused by storage environment.

LBIS™ Mouse Urinary Albumin Assay Kit

[Storage]	Store at 2 - 10°C (Do not freeze)
[Expiration date]	Indicated on the label
[Package]	For 60 tests / 250 tests
[Cat #]	295-92701 (60 tests) 291-92703 (250 tests)

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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