FM Fix (s)

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Tissue fixation is an important procedure to protect biological specimens from degradation due to autolysis and decay, ensuring stable preservation of proteins and other components within the tissues. Formalin fixation is widely used in research and pathology, including microscopy specimen preparation, however, it can affect the steric structure and antigenicity of some proteins, leading to deterioration in staining quality and other properties.

FM Fix is a fixative that has superior tissue staining properties while retaining tissue morphology and three-dimensional structure. It helps to maintain the three-dimensional structure of proteins even after tissue fixation, enabling the stable detection of antigens such as complements. In addition, it is characterized by minimal reduction of fluorescence protein signals such as GFP, making it suitable for tissue transparency techniques.

FM Fix (s) is a fixative optimized for small tissues such as needle biopsies and tissues with good permeability, while maintaining the fixative characteristics of FM Fix (L).

【Features】

- Optimized fixative for needle biopsy tissue
- Superior tissue morphology preservation and staining
- Tissue cryopreserved after fixation with this fixative can be thawed and stained with PAS, PAM, or MT

【Kit contents】

This kit consists of two components.

Solution A (s) : 0.5 mL × 10 bottles

Solution B (s) : 7.5 mL × 10 bottles

【Storage】

Store at -20°C

【Procedure】

Note : Fixing B (s) solution may turn slightly yellow, but this does not affect performance.

1. Reagent Preparation

Prepare the reagent just before use, and use the mixed solution on the same day.

- (1) Thaw each solution immediately before use.
- (2) After thawing, shake or stir the solution lightly to make sure there are no undissolved residues.
- (3) Add all of solution A (0.5 mL × 1 bottle) to solution B (7.5 mL × 1 bottle).
- (4) Close the container lid and shake gently to mix Solution A and Solution B.

2. Fixation Procedure

- (1) Fix a kidney tissue fragment collected by needle biopsy or a wedge-shaped excised kidney tissue fragment by immersing it in the prepared fixative solution for several hours or overnight.
 - * Adjust the fixation time, etc. according to the size and type of tissue specimen.
- (2) After confirming the presence or absence of lesions, separate the blocks for each application and proceed with each step.

<Preparation of paraffin-embedded sections for optical microscopy>

Cut the fixed tissue into specimens of appropriate size, and follow the usual steps for paraffin-embedded sections using an automated processing device.

* If immunostaining is performed on paraffin-embedded sections, antigen activation may be required. In such cases, perform heat treatment with citric acid or EDTA to activate the antigen.

<Preparation of frozen sections for immunostaining>

Cut the fixed tissue into specimens of appropriate size and freeze them after the usual embedding steps using sucrose replacement, OCT compound, etc.

【Precautions】

When handling, wear a gas mask or hose mask, protective gloves, and glasses as necessary.

(For details, please refer to the handling precautions on the label.)

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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