## LabAssay™ Total Protein

Please, read this instruction carefully before use.

### 1. Intended use

This product is a research-use reagent designed for the quantification of total protein using the Biuret method.

Total protein is mainly composed of albumin and globulin. These proteins play important roles in maintaining osmotic pressure and transporting substances within the body. Abnormal levels of total protein can indicate pathological conditions such as liver disease, kidney disease, and chronic inflammation.

When used in conjunction with albumin concentration data, the albumin-to-globulin (A/G) ratio can also be calculated. The A/G ratio is regarded as a useful diagnostic indicator for liver and kidney disorders, nutritional deficiencies, multiple myeloma, and chronic inflammatory diseases. Therefore, total protein and albumin are often measured together.

### 2. Storage and expiration date

Store at 2°C - 10°C and do not freeze.

Store the reconstituted standard at -35°C or below.

The expiration date is indicated on the label on the outer box of the kit.

### 3. Kit Component Reagents

Components		Use Status	Amount
(1)	Chromogen Reagent	Solution	65 mL/1 bottle
(2)	Standard	Freeze-dried	1 bottle

\* For 250 tests

### 4. Principle of the method

When the Chromogen Reagent is added to the sample, proteins form a coordination complex with copper ions and develop a blue-violet color. The total protein concentration in the sample is determined by measuring the absorbance of this blue-violet color (Biuret method).

### 5. Equipment or supplies required but not provided in the kit

- 96-well microplate (transparent type)
- Micropipette
- Microtube
- Pipette
- Plate mixer
- Microplate reader with 540 nm wavelength filter

### 6. Preparation of reagents

① Chromogen Reagent : Ready to use. After opening the bottle, store at 2°C - 10°C.

② Standard : Reconstitute ALT Standard with purified water in the amount specified in "Reconstitution of standard" \* to prepare the original standard solution (8.0 g/dL). Then prepare a dilution series using water that has been allowed to warm up to room temperature.

\* Use purified water equilibrated to room temperature to dissolve the standard.

\* Add the specified volume of purified water and allow the solution to stand at room temperature for 10 minutes. Then mix using a roller mixer or an equivalent device until completely dissolved.

\* Standard solutions diluted to each concentration should be used immediately and must not be stored.

Example of preparation of dilution series of standard solution

Standard Concentration (g/dL)	11.9	8.0	4.0	2.0	1.0	0.0
Standard Solution (μL)	Original standard solution : 50	Original standard solution : 50	Original standard solution : 50	50 *	50 *	—
Water (μL)	—	—	50	50	50	50
Well Volume (μL)	7.5**	5	5	5	5	5

\*One rank higher standard solution.

\*\* : The standard solution is typically dispensed at 5 μL per well ; however, 7.5 μL is used for the 11.9 g/dL standard. Accordingly, the displayed value has been corrected to account for the increased volume.

## 7. Preparation of specimen

### Serum/Plasma

- Analyze specimens immediately after collection.
- The anticoagulant heparin, EDTA · 2Na and citrate does not significantly influence the assay when used in normal amounts.
- Serum samples with a high degree of hemolysis may affect the measurement results. In such cases, collecting the blood sample again is recommended or the results should be used for reference only.
- Chyle and ascorbic acid affect the measured values.
- High levels of bilirubin in the specimen may affect the measured values.
- If the measured value exceeds the measurable range, dilute the sample with saline solution before measurement. Multiply the result by the dilution factor.

## 8. Assay procedure

Ensure that the reagents are brought up to room temperature (20°C - 25°C) before use.


- (1) Dispense 250  $\mu$ L of Substrate-Enzyme Solution into each well of a 96-well microplate.
  - (2) Dispense 5  $\mu$ L of each concentration of standard solution into the well in which the standard is to be measured.  
(Dispense the 11.9 g/dL standard at the volume specified in the standard dilution example.)
  - (3) Dispense 5  $\mu$ L of sample into the well in which the sample is to be measured.
  - (4) Shake the 96-well microplate on a plate mixer (\*①).
  - (5) Allow to react at room temperature at 20°C - 25°C for 30 minutes.
  - (6) Immediately after reaction, measure the absorbance of each well at 540 nm using a spectrophotometer for microplates.
- (\*①) Guideline for shaking : 600 rpm - 700 rpm for 10 seconds, repeated 3 times.

## 9. Calculation

- (1) Create a calibration curve for each measurement, with the standard concentration (g/dL) on the X-axis and absorbance on the Y-axis.
  - (2) From the calibration curve, read the total protein (g/dL) corresponding to the absorbance of the sample. If a diluted sample is used, the concentration reading is multiplied by the sample dilution ratio to obtain the measurement value.
- \*If the absorbance of a sample is outside the standard curve, dilute the sample with saline to an appropriate dilution ratio and repeat the assay.
- \*In computer software calculations, we recommend using 4 parameters or 5 parameters.

## 10. Standard Curve

- Measurement range : 1.0 - 11.9 g/dL

## 11. Performance

### ● Accuracy

- When control serum of known concentration is measured, the measured value is within  $\pm 15.0\%$  of the known concentration.

### ● Repeatability

- When the same sample is measured simultaneously 5 times, the CV (%) of measured values is 15.0% or less.

## Notes

- Avoid using this product at temperatures or reaction times other than those specified.
- The color development is stable, and the absorbance shows almost no change even after standing at room temperature for 1 hour.

- If any reagent accidentally comes into contact with the mouth, eyes, or skin, flush immediately with a large amount of water. Consult a physician if necessary.
- When using a pipette, do not aspirate by mouth; always use a safety pipette filler or equivalent device.
- Store reagents under the specified conditions and do not use expired reagents.
- Do not use reagents that have been accidentally frozen because accurate results may not be obtained.
- After opening, use the reagent as soon as possible. If it is stored, close the lid and store under the specified conditions.
- Do not use the container or accessories included with this product for anything other than their intended purposes.
- Handle samples being cognizant of the risk of infection (e.g., viruses).
- Wear disposable gloves to avoid the risk of infection when performing a test.
- Be careful regarding contamination when collecting reagents with a pipette.
- Do not mix and use reagents with different lot numbers.
- When disposing of the reagents, dispose of them according to Waste Management and Public Cleansing Law (Waste Disposal and Cleaning Act) and wastewater standards. The chromogenic solution contains copper (II) sulfate pentahydrate at 3 g/dL (equivalent to 760 mg/L as copper).
- Any reagent or reagent bottle that comes into contact with a specimen should be considered potentially infectious and handled accordingly.
- Any scattering/splashes of samples or waste fluids should be removed with a disinfectant such as sodium hypochlorite (effective chlorine concentration 1,000 ppm).

LabAssay™ Total Protein

[Storage] Store at 2°C - 10°C (Do not freeze)

[Term of validity] Indicated on the label.

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