
Presep<sup>®</sup> PolyChelate (250 mg/15 mL)

This product is a solid-phase extraction column packed with a chelating resin containing carboxymethylated polyethyleneimine as a functional group. Compared to the commonly used chelating resin containing IDA (Iminodiacetic acid) as a functional group, it can capture and recover many elements (metals) with high efficiency over a wide pH range. It can also capture and recover metal oxoacids (Mo, V, W) at high rates under acidic to neutral conditions, without being affected by alkali metals and alkaline earth metals.

## 【Example of solid-phase extraction procedure\*】

\* These conditions are only an example and may require further study depending on the sample.

## 【Instruction for use and handling】

- Products must be conditioned prior to use. Drying the packing agent during conditioning may result in loss of reproducibility.
- Add a buffer solution that does not interact with elements (metals) such as ammonium acetate before pH adjustment of the sample solution.
- When passing the sample solution, reduce or increase the pressure to maintain a constant flow rate.
- To avoid adsorption and elution of elements (metals), please refrain from using metal utensils. Use glass or plastic utensils that have been thoroughly rinsed with nitric acid.

## 【References】

Analysis example are posted on our website.  
 TOP>Analytical Chemistry>Solid-Phase Extraction  
 Columns>SPE Column for Metal Scavenger  
<https://labchem-wako.fujifilm.com/us/category/00995.html>

## FUJIFILM Wako Laboratory Chemicals site

<https://labchem-wako.fujifilm.com>

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
 Telephone : + 81-6-6203-3741  
 Facsimile : + 81-6-6201-5964  
[ftwk-cs@fujifilm.com](mailto:ftwk-cs@fujifilm.com)

## FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
 Richmond, VA 23237  
 U.S.A.  
 Telephone : + 1-804-271-7677  
 Facsimile : + 1-804-271-7791  
[wkuslabchem@fujifilm.com](mailto:wkuslabchem@fujifilm.com)

## FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12  
 D-41468 Neuss  
 Germany  
 Telephone : + 49-2131-311-0  
 Facsimile : + 49-2131-311100  
[labchem\\_wkeu@fujifilm.com](mailto:labchem_wkeu@fujifilm.com)

## FUJIFILM Wako Chemicals (Hong Kong) Limited

Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road, Tsuen Wan, N.T., Hong Kong  
 Telephone : + 852-2799-9019 Facsimile : + 852-2799-9808  
[wkhk.info@fujifilm.com](mailto:wkhk.info@fujifilm.com)

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## FUJIFILM Wako (Guangzhou) Trading Corporation