for genetic research
DNA Extractor SP Kit
for Serum and Plasma

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“DNA Extraction SP Kit” is a novel extraction kit developed for rapid preparation of circulating DNA in human serum and plasma.

Recently, tumor-specific genes have been amplified and detected in the serum and plasma of patients with various diseases such as lung, breast and colon cancer. Many of these reference articles are currently being published. The research of gene diagnosis using plasma DNA is proceeding to early diagnosis, prognostic value and postoperative surveillance of cancer.

DNA Extractor SP Kit is a powerful pretreatment reagent for the detection and analysis of target DNA because of its high quality and yield. This procedure, using Sodium Iodide (NaI) as a chaotropic agent, can be performed in a simple and safe extraction without hazardous phenol/chloroform. Consequently, the entire DNA extraction procedure is done in a single centrifuge tube.

【Advantages】

- High DNA yield from small amount of serum and plasma

The DNA extraction procedure is based on a simple centrifugal separation using Sodium Iodide and alcohol, that results in less loss of DNA than in solid-phase methods using a silica matrix.

- Safety of operation - Hazardous phenol and chloroform are not employed as reagents for extraction.

- Minimizing the cross and extraneous contamination

The entire extraction of DNA is done in a single centrifuge tube.

- Complete removal of lipids derived from blood using a specialized alcohol solution.

- Extraction of high quality DNA without PCR inhibition

【Kit Contents】

1) Enzyme Reaction Solution	10 mL × 1
2) Protein Digestion Solution	250 μ L × 1
3) Sodium Iodide Solution	15 mL × 1
4) Alcohol Solution	30 mL × 1
5) Washing Solution (A)	50 mL × 1
6) Washing Solution (B)	50 mL × 1

[Materials Needed]**Reagent :**

TE (pH 8.0) or sterile distilled water (D.W.)

Equipments :

- 1) microcentrifuge tubes (1.5 mL or 2 mL)
- 2) vortex mixer
- 3) incubator
- 4) high speed microcentrifuge
- 5) block heater

[Storage and Stability]

Store at 2 ~ 10°C.

The expiry date is printed on the kit label.

[Reactions]

50 reactions (100 μ L of start sample)

25 ~ 30 reactions (200 μ L of start sample)

[Method]**1. Important note before starting**

- If salt precipitation appears in the Sodium Iodide Solution, incubate at 50°C ~ 60°C until the precipitation has dissolved.
- The Protein Digestion Solution should be placed on ice, during operation. Other reagents can be handled at room temperature.
- Samples of serum and plasma should be placed on ice for the avoidance of activation of the endogenous DNase. We recommend that Enzyme Reaction Solution should be added as soon as possible.

2. Standard protocol : 100 μ L sample

1. Dispense 100 μ L of serum or plasma to a plastic 1.5 mL microcentrifuge tube.
2. Add 200 μ L of Enzyme Reaction Solution and mix briefly.
3. Add 5 μ L of Protein Digestion Solution and vortex.
4. Incubate at 56°C for 10 min.

It is possible to store the solution of this step at 4°C or - 20°C for up to one week.
5. Add 300 μ L of Sodium Iodide Solution and mix briefly.
6. Add 600 μ L of Alcohol Solution and vortex.
7. Incubate at room temperature for 10 min.
8. Centrifuge at 12,000 ~ 20,000 \times g for 10 min. at room temperature.
9. Decant and remove the supernatant as much as possible from the tube that is inverted and placed on a paper towel.
10. Add 1mL of Washing Solution (A) to the pellet and vortex.

It is possible to store the solution at - 20°C for long term.
11. Centrifuge at 12,000 ~ 20,000 \times g for 5 min. at room temperature.
12. Repeat Step 9 and remove as much supernatant as possible.
13. Add 1 mL of Washing Solution (B) to the pellet and vortex.

14. Centrifuge at $12,000 \sim 20,000 \times g$ for 5 min. at room temperature.
15. Repeat Step 9 and remove as much supernatant as possible.
16. Dry the pellet well for about 5 min. on the block heater at about 65°C .
17. Add adequate volume ($10 \sim 20 \mu\text{L}$) of TE (pH 8.0) or D.W. to the pellet.
18. Dissolve the pellet completely with vortex mixer and incubation at about 65°C for 3 ~ 5 min.
19. Store at -20°C for long-term storage.

3. Alternative protocol : 200 μL of sample

If you start the experiment from 200 μL of sample, refer to the below procedure.
This protocol may decrease the DNA yield and purity compared to the Standard protocol.

1. Dispense 200 μL of serum or plasma to a plastic 2 mL microcentrifuge tube.
2. Add 300 μL of Enzyme Reaction Solution and mix briefly.
3. Add 8 ~ 10 μL of Protein Digestion Solution cooled on ice and vortex.
4. Incubate at 56°C for 30 min.

It is possible to store the solution of this step at 4°C or -20°C for up to one week.

5. Add 400 μL of Sodium Iodide Solution and mix briefly.
6. Add 900 μL of Alcohol Solution and vortex.
7. Incubate at room temperature for 10 min.
8. Centrifuge at $12,000 \sim 20,000 \times g$ for 10 min. at room temperature.
9. Decant and remove the supernatant as much as possible from the tube that is inverted and placed on a paper towel.
10. Add 1.5 mL of Washing Solution (A) to the pellet and vortex.

It is possible to store the solution at -20°C for long term.

11. Centrifuge at $12,000 \sim 20,000 \times g$ for 5 min. at room temperature.
12. Repeat Step 9 and remove as much supernatant as possible.
13. Add 1.5 mL of Washing Solution (B) to the pellet and vortex.
14. Centrifuge at $12,000 \sim 20,000 \times g$ for 5 min. at room temperature.
15. Repeat Step 9 and remove as much supernatant as possible.
16. Dry the pellet well for about 5 min. on the block heater at about 65°C .
17. Add adequate volume ($10 \sim 20 \mu\text{L}$) of TE (pH 8.0) or D.W. to the pellet.
18. Dissolve the pellet completely with vortex mixer and incubation at about 65°C for 3 ~ 5 min.
19. Store at -20°C for a long-term storage.

[Application data] Result of amplification of p53-Exon5 region:

By using DNA Extractor SP Kit, DNA was extracted from 10 human sera and one plasma sample. The extracted DNA was finally resolved in 20 μL of TE (pH 8.0), and 5 μL of that was added to the PCR. Commercial available kit based on a glass binding method (spin column method) was used as a comparative method.

Result of amplification of p53-Exon5 region


All samples extracted by DNA Extractor SP Kit was able to be amplified in the region of p53-Exon 5. In this experiment, we found that the amount of fragment amplified from our kit was larger than that of a competitor's kit. In the competitor's kit, the target fragment could not be observed in the sample of serum 10 and plasma.

PCR condition : 10 μ L reaction

DNA sample	5 μ L			
Exon5 forward-primer	4 pmol	96°C	30 sec.	
Exon5 reverse-primer	4 pmol	94°C	30 sec.	40 cycles
dNTP	200 μ mol/L	65°C	30 sec.	
10 \times Gene <i>Taq</i> Universal Buffer	1 μ L	72°C	1 min.	
Gene <i>Taq</i> NT	0.5 units	72°C	5 min.	
<hr/>				
D.W.	fill up to 10 μ L			

Flow chart

- ※1) Grouping symbol means the protocol from 200 μ L of sample.
- ※2) Samples of serum and plasma should be placed on ice for the avoidance of activation of the endogenous DNase. We recommend Enzyme Reaction Solution to be added as soon as possible.
- ※3) Protein Digestion Solution should be placed on ice, during operation.
- ※4) Decant and remove the supernatant as much as possible from the tube that is in verted and placed on paper towel.
- ※5) It is possible to store the solution at -20°C for long term.
- ※6) Mix vigorously until the pellet is removed from wall of tube.
- ※7) Dissolve the pellet completely with vortex mixer and incubation.

[References]

- 1) Ishizawa, M., Kobayashi, Y., Miyamura, T. and Matsuura, S. : *Nucleic Acids Res.*, **22**, 1774 (1994)
- 2) Sozzi, G., Musso, K., Ratcliffe, C., Goldstraw, P., Pierotti, M. A. and Pastorino, U. : *Clin Cancer Res.*, **5**, 2689 (1999)
- 3) Silva, J. M., Dominguez, G., Garcia, J. M., Gonzalez, R., Villanueva, M. J., Navarro, F., Provencio, M., San, Martin, S., Espana, P. and Bonilla, F. : *Cancer Res.*, **59**, 3251 (1999)
- 4) Shao, Z. M., Wu, J., Shen, Z. Z. and Nguyen, M. : *Clin Cancer Res.*, **7**, 2222 (2001)

Allied products

Wako catalog No.	Description	Package Size
316-90025	TE (pH 8.0)	500 mL
200-14911	10 × TE (pH 8.0)	1 L × 10
318-90105	Distilled Water, Deionized, Sterile	500 mL
318-03231	Gene <i>Taq</i> NT (<i>Taq</i> DNA polymerase)	250 units
318-02871	Gene <i>Taq</i> (<i>Taq</i> DNA polymerase)	250 units
312-03491	P53 Primer Exon 2, 3	100 reactions
315-03501	P53 Primer Exon 4	100 reactions
312-03511	P53 Primer Exon 5	100 reactions
319-03521	P53 Primer Exon 6	100 reactions
316-03531	P53 Primer Exon 7	100 reactions
313-03541	P53 Primer Exon 8, 9	100 reactions
310-03551	P53 Primer Exon 10	100 reactions
317-03561	P53 Primer Exon 11	100 reactions

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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