〈For Research Use Only〉 Code No. 296-64303 (100 reactions)

for Genetic Research

## *S. pombe* Direct Transformation Kit *Wako*

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The *S. pombe* Direct Transformation Kit *Wako* is specially designed for rapid and easy transformation of the fission yeast, *Schizosaccharomyces pombe*.

The *S. pombe* Direct Transformation Protocol allows for successful transformation simply by mixing plasmid DNA and the kit reagents with overnight yeast cell cultures without complicated competent cell preparation.

The traditional yeast transformation methods require repeated centrifugation or washing steps, thus are problematic for high-throughput transformation. The *S. pombe* Direct Transformation Kit *Wako* is well suited for high-throughput transformation of a large number of yeast strains grown in 96-well plates. You may choose either a "96-well plate protocol" or a "tube protocol" for your purpose.

### 【Features】

- Direct transformation : One-step procedure by mixing your plasmid and the kit reagents with yeast culture.
- Simply : No need to prepare competent cells.
- High-throughput: Suitable for transformation of a large number of yeast strains using a 96-well plate.

### 【Transformation Efficiency】

≥ 250 cfu/μg

Transformation efficiency is calculated as the number of colonies per microgram of plasmid DNA.

### 【Kit contents】

(1) <i>Sp</i> Transformation Reagent	11.5 mL × 1
(2) Carrier DNA (5 μg/μL)	0.5 mL × 1

### 【Storage】

Store at -20°C

※The number of times of freeze-thaw should not exceed 10 times. If you do freeze-thaw more than 10 times, dispense the reagents and store at -20°C.

### 【Protocol】

#### I. Protocol 1 (96-well plate protocol)

##### I-1. Preparation of yeast culture

##### < Required Reagents and Equipment >

- Incubator (28-30°C)
- Sterile 96-well microtiter plate (Nunc Cat# : 267245)
- Microplate sealing tape
- Sterile 96-pin replicator (GENETIX, UK. Cat# : X5052)
- Square plate (14 × 10 cm) (Eikenkizai Co., Ltd. Cat# : AW2000)

- Sterile water
- YPD agar medium

[Preparation of YPD agar medium, 1 L]

Dissolve the following in distilled water to 1 L.

Yeast extract	10 g
Peptone	10 g
D(+)-Glucose	20 g

Add 20 g agar for plates.

Autoclave at 121°C for 20 minutes.

- YES medium

[Preparation of YES medium, 1 L]

Dissolve the following in distilled water to 1 L.

Yeast extract	10 g
D(+)-Glucose	30 g
L-Histidine	250 mg
L-Leucine	250 mg
L-Lysine	250 mg
Adenine	250 mg
Uracil	250 mg

Autoclave at 121°C for 20 minutes.

#### < Before Starting >

Inoculate a square YPD plate<sup>1)</sup> with *S. pombe* strains from glycerol stock using a sterile 96-pin replicator. Incubate the plate at 28-30°C for 2-7 days<sup>2)</sup>.

- 1) We recommend a YPD plate, because a YES plate shows a decreasing trend of transformation efficiency.
- 2) Cultured YPD plate can be stored at 4-10°C for 1 month without decreasing transformation efficiency.

#### < Yeast culture Protocol >

1. Aliquot 25 μL of YES medium into each well of a 96-well microtiter plate.
2. Pick *S. pombe* strains from the YPD plate culture using a 96-pin replicator and inoculate YES medium in a 96-well plate.
3. Seal the plate with a microplate sealing tape and incubate at 28-30°C for 18-24 hours without shaking.

#### I-2. Transformation

##### < Required Reagents and Equipment >

- Incubator (28-46°C)
- Multi channel pipettor
- Sterile 25 mL reserver
- Appropriate selective plates and plasmid
- Sterile water

#### < Before Starting >

Thaw the *Sp* Transformation Reagent and Carrier DNA completely at room temperature.

### < Transformation Protocol >

1. Add 1  $\mu\text{g}$  of plasmid DNA<sup>1)</sup> and 4  $\mu\text{L}$  of Carrier DNA to 90  $\mu\text{L}$  of *Sp* Transformation Reagent<sup>2)</sup>, and then add sterile water to give a final volume of 100  $\mu\text{L}$  per well. In the case for preparing multiple wells, the pre-mix can be prepared as follows.

	1 well	96 wells	X wells
<i>Sp</i> Transformation Reagent <sup>1)</sup>	90 $\mu\text{L}$	10.8 mL	$90 \times 1.25 \times X \mu\text{L}$
plasmid DNA <sup>2)</sup>	1 $\mu\text{g}$	120 $\mu\text{g}$	$1 \times 1.25 \times X \mu\text{g}$
Carrier DNA (5 $\mu\text{g}/\mu\text{L}$ )	4 $\mu\text{L}$	480 $\mu\text{L}$	$4 \times 1.25 \times X \mu\text{L}$

Sterile water	up to 100 $\mu\text{L}$	12 mL	$100 \times 1.25 \times X \mu\text{L}$
---------------	-------------------------	-------	--

2. Suspend the precipitated cells in the cultured 96-well microtiter plate using a plate mixer<sup>3)</sup>. Immediately add 100  $\mu\text{L}$  of the DNA mixture (Step 1) to each well of yeast culture using a multi-channel pipetter. Seal the plate and mix using a plate mixer<sup>3)</sup>.
3. Incubate at 46°C for 2 hours<sup>4)</sup>.
4. Spot 10  $\mu\text{L}$  of the mixture on the appropriate selection plate.
5. Incubate the selection plate at 28-30°C for 5-7 days.

- 1) The quality of DNA used for this method plays an important role in transformation efficiency. We recommend the use of the kits that have high purification performance for preparation of plasmid DNA.
- 2) Slow pipette operation is required because of the viscosity of this reagent.
- 3) If you do not have a plate mixer, a multi-channel pipetter can be used.
- 4) It is possible to incubate at 42°C for 4 hours or 37°C for 6 hours. You can choose the proper condition for your purpose.

## II. Protocol 2 (Tube protocol)

### II-1. Preparation of Yeast culture

#### < Required Reagents and Equipment >

- Shaking incubator (28-30°C)
- Sterile petri dish
- Sterile two position type tube (14 mL)
- YPD agar medium (See Section I for preparation (p. 2))
- YES medium (See Section I for preparation (p. 2))

#### < Before Starting >

Streak a YPD plate<sup>1)</sup> with *S. pombe* strain and incubate the plate at 28-30°C for 2-7 days<sup>2)</sup>.

- 1) We recommend a YPD plate, because a YES plate shows a decreasing trend of transformation efficiency.
- 2) Cultured YPD plate can be stored at 4-10°C for 1 month without decreasing transformation efficiency.

### < Yeast Culture Protocol >

1. Inoculate the colonies of your *S. pombe* strain into 2 mL of YES medium. Measure its OD at 600 nm. After adjusting the OD<sub>600</sub> to 0.2 with medium, transfer the 1 mL of its suspension to the 14 mL two position type tube. Cover the tube with the lid loosely.

Note : One mL of culture can be used for approximately 30 transformations (25  $\mu\text{L}$  per transformation). If you want to prepare over 30 transformations, increase the number of culture tubes but not the volume, because the volume of culture significantly affects the aeration and growth conditions.

- We highly recommend to perform yeast culture with a reciprocal or rotary shaking incubator at 150 rpm for 24 hours. However, yeast cell growth may vary from strain to strain, the value of OD<sub>600</sub> should be between 6.0 and 7.0 at the end of culture. Cell concentration plays an important role in transformation efficiency.
- The OD<sub>600</sub> should be measured by 10-fold to 50-fold dilution of the culture.

2. Incubate at 28-30°C for ~24 hours with vigorous shaking. The culture should reach a cell density to an OD<sub>600</sub> of 6.0-7.0.

## I-2. Transformation

### < Required Reagents and Equipment >

- Incubator or water bath (37°C)
- Incubator (28-30°C)
- Sterile microcentrifuge tube (1.5 mL)
- Appropriate selective plates and plasmid
- Sterile water

#### < Before Starting >

Thaw the *Sp* Transformation Reagent and Carrier DNA completely at room temperature.

### < Transformation Protocol >

1. Add 1  $\mu\text{g}$  of plasmid DNA<sup>1)</sup> and 4  $\mu\text{L}$  of Carrier DNA to 90  $\mu\text{L}$  of *Sp* Transformation Reagent<sup>2)</sup> in a microtube, and then add sterile water to give a final volume of 100  $\mu\text{L}$ . Mix by vortexing.
2. Add 25  $\mu\text{L}$  of yeast culture<sup>3)</sup> to the mixture and mix by vortexing.
3. Incubate at 37°C for 2 hours.
4. Spread on an appropriate selection plate using a sterile spreader. If you want to obtain single colonies, dilute the mixture to 2-fold to 10-fold with sterile water and spread 100  $\mu\text{L}$  of the diluted mixture on an appropriate selection plate.
5. Incubate the plate at 28-30°C for 5-7 days.

- 1) The quality of DNA used for this method plays an important role in transformation efficiency. We recommend the use of the kits that have high purification performance for preparation of plasmid DNA.
- 2) Slow pipette operation is required because of the viscosity of this reagent.
- 3) When yeast culture is concentrated 5-fold by centrifugation, transformation efficiency is increased to approximately 5-fold.

## 【Appendix (Flow Chart of the Experiments)】

### I. 96-well plate protocol

96 microtiter well  
 ↓  
 ← 25  $\mu$ L of YES medium  
 ← Inoculate yeast strains into each well  
 Incubate at 28-30°C for 18-24 hr. without shaking  
 Suspend precipitated cells using plate mixer  
 ← 100  $\mu$ L of plasmid /reagent solution<sup>1)</sup>  
 Mix with plate mixer  
 ↓  
 Incubate at 46°C for 2 hr.<sup>2)</sup>  
 Spotting 10  $\mu$ L  
 ↓  
 Incubate at 28-30°C for 5-7 days  
 Transformed yeast

### II. Tube protocol

2 mL of YES medium  
 ↓  
 ← Inoculate yeast colonies  
 Adjust OD<sub>600</sub> : 0.2 by medium  
 ↓  
 Transfer the 1 mL of the suspension to a 14 mL two position type tube  
 ↓  
 Incubate at 28-30°C with shaking until OD<sub>600</sub> : 6.0-7.0  
 Yeast culture...①

1.5 mL of microcentrifuge tube  
 ↓  
 ← 100  $\mu$ L of plasmid/reagent solution<sup>1)</sup>  
 Mix by vortexing  
 ← 25  $\mu$ L of yeast culture (①)  
 ↓  
 Incubate at 37°C for 2 hr.  
 Plating  
 ↓  
 Incubate at 28-30°C for 5-7 days  
 Transformed yeast

- 1) plasmid/reagent solution (per well, per tube)  
*Sp* Transformation Reagent 90  $\mu$ L  
 plasmid DNA 1  $\mu$ g  
 Carrier DNA 4  $\mu$ L  
 Sterile water up to 100  $\mu$ L  
 2) It is possible to incubate at 42°C for 4 hours or 37°C for 6 hours.

## 【Application】

### Higher transformation efficiency method (centrifugal concentration method)

2 mL of YES medium  
 ↓  
 ← Inoculate yeast colonies  
 Adjust OD<sub>600</sub> : 0.2 by medium  
 ↓  
 Transfer the 1 mL of the suspension to a 14 mL two position type tube  
 Incubate at 28-30°C with shaking until OD<sub>600</sub> : 6.0-7.0

500 × G, 5 min., RT  
 Remove 800  $\mu$ L of supernatant  
 Suspend by vortex mixer  
 centrifugal concentration (×5)  
 ↓  
 Yeast culture...①

1.5 mL of microcentrifuge tube  
 ↓  
 ← 100  $\mu$ L of plasmid/reagent solution  
 Mix by vortexing  
 ← 25  $\mu$ L of yeast culture (①)  
 ↓  
 Incubate at 37°C for 2 hr.  
 Plating  
 ↓  
 Incubate at 28-30°C for 5-7 days  
 Transformed yeast

## Allied products

Wako catalog No.	Description	Package Size
049-31165	D (+) -Glucose	500 g
084-00702	L-Histidine Hydrochloride Monohydrate	25 g
124-00852	L-Leucine	25 g
010-19612	Adenine Sulfate	25 g
212-00062	Uracil	25 g
010-15815	Ager, Powder	500 g
318-90105	Distilled Water, Deionized, Sterile	500 mL

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
 Telephone : + 81-6-6203-3741  
 Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
 Richmond, VA 23237  
 U.S.A.  
 Telephone : + 1-804-271-7677  
 Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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900mL の水に溶解  
Yeast Extract      10g  
Peptone              10g  
D (+)-Glucose      20g

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