**DNA Extractor<sup>®</sup> TIS Kit****〔Introduction〕**

The DNA Extractor<sup>®</sup> TIS Kit is specially designed for the extraction of DNA from human and animal parenchyma tissue. The base of this kit is a procedure using Sodium Iodide (NaI) as a chaotropic agent, which allows for a safe extraction without hazardous phenol/chloroform. The procedure using Sodium Iodide is an extraction procedure which causes less oxidation of DNA. Additionally, the DNA Extractor<sup>®</sup> TIS Kit is able to inhibit the oxidation of DNA by using an Oxidation Inhibitor. The DNA Extractor<sup>®</sup> TIS Kit is a most useful DNA extraction kit for the assay of 8-OHdG (8-Hydroxy-2'-deoxyguanosine) which is a marker of oxidative stress.

**〔Features〕**

- 1) Useful for the assay of markers of oxidative stress.
- 2) Based on the procedure using Sodium Iodide which causes less oxidation of DNA.
- 3) Inhibits the further oxidation of DNA by using an Oxidation Inhibitor.
- 4) Suitable for a wide range of applications which require the extraction of DNA from human and animal parenchyma tissue.

**〔Kit contents〕**

- |                               |                 |
|-------------------------------|-----------------|
| 1) Lysis Solution             | 75 mL × 2       |
| 2) Enzyme Reaction Solution   | 15 mL × 1       |
| 3) RNase Solution             | 50 $\mu$ L × 1  |
| 4) Protein Digestion Solution | 750 $\mu$ L × 1 |
| 5) Oxidation Inhibitor        | 350 $\mu$ L × 1 |
| 6) Sodium Iodide Solution     | 15 mL × 1       |
| 7) Alcohol Solution           | 30 mL × 1       |
| 8) PEG Solution               | 20 mL × 1       |

**Note :** In case precipitates appear in Enzyme Reaction Solution, solubilize them completely by warming the solution at 37°C or leaving the solution at room temperature.

**〔Materials needed〕****Reagents :**

- 1) 70% Ethanol : Code No. 059-07895
- 2) Deionized distilled water (DDW) : Code No. 316-90101
- 3) TE buffer (pH 8.0) : Code No. 314-90021

**Equipments :**

- 1) Microcentrifuge tubes (2 mL)
- 2) Microcentrifuge tubes (1.5 mL)
- 3) Homogenizer (Teflon homogenizer etc.)
- 4) High speed refrigerated microcentrifuge
- 5) Vortex mixer
- 6) Incubator (37°C)
- 7) Micropipette

**〔Storage〕**

Store at 2~10°C.

**〔Package size〕**

50 reactions

**〔Method〕**

The protocol contains vigorous mixing steps by a vortex mixer in order to finish the incubation and extraction steps quickly and inhibit the oxidation of DNA. If extraction of intact and long genomic DNA is desired, the mixing steps with the vortex mixer should be done more gently.

When DNA obtained by using Protocol 1 contains too much RNA, or sodium iodide, it is possible to get purer DNA using Protocol 2.

**〔Protocol 1〕**

- (1) Cut up the appropriate amount of tissue into small pieces and add 1 mL of Lysis Solution. Homogenize it in an ice water bath.

The following table shows the maximum amount of tissue that can be treated. Do not exceed this amount.

| Tissue          | Yield of DNA<br>(mg/100 mg of<br>tissue) | Maximum amount<br>of the tissue that<br>can be treated. | Predicted<br>yield of DNA |
|-----------------|------------------------------------------|---------------------------------------------------------|---------------------------|
| spleen          | 1-1.5 mg                                 | 30 mg                                                   | 300-450 $\mu$ g           |
| kidney          | 400-500 $\mu$ g                          | 100 mg                                                  | 400-500 $\mu$ g           |
| lung            | 400-500 $\mu$ g                          | 100 mg                                                  | 400-500 $\mu$ g           |
| pancreas        | 400-500 $\mu$ g                          | 100 mg                                                  | 400-500 $\mu$ g           |
| small intestine | 300-400 $\mu$ g                          | 100 mg                                                  | 300-400 $\mu$ g           |
| liver           | 250-300 $\mu$ g                          | 150 mg                                                  | approx. 400 $\mu$ g       |
| stomach         | 250-300 $\mu$ g                          | 150 mg                                                  | approx. 400 $\mu$ g       |
| testis          | approx. 200 $\mu$ g                      | 150 mg                                                  | approx. 300 $\mu$ g       |
| brain           | 100-150 $\mu$ g                          | 180 mg                                                  | 200-300 $\mu$ g           |
| heart           | approx. 100 $\mu$ g                      | 200 mg                                                  | approx. 200 $\mu$ g       |

The amount is based on the data for mouse tissue.

The maximum amount of the tissue which can be treated is 150~200 mg, and the maximum amount of DNA is 400~500  $\mu$ g. It becomes difficult to purify DNA when the amount of the tissue or DNA exceeds the maximum. We recommend the use the smallest amount of tissue in order to get good results when a large amount of DNA is not required.

- (2) Transfer the homogenate to a 2 mL microcentrifuge tube and centrifuge at  $600 \times G$  for 10 minutes at  $4^{\circ}C$ .

Centrifuge at  $5,000 \times G$  for 10 minutes at  $4^{\circ}C$  when the sample is brain tissue.  
(The condition of the centrifugation at Step 5 is  $600 \times G$ . This is the same as protocol.)

- (3) Decant and remove the supernatant from the tube then invert and place on paper towels to remove any excess liquid.  
(4) Add 1 mL of Lysis Solution to the pellet and vortex to detach and pull the pellet apart.  
(5) Centrifuge at  $600 \times G$  for 10 minutes at  $4^{\circ}C$ .  
(6) Repeat Step 3~5.  
(7) Decant and remove the supernatant from the tube then invert and place on paper towels to remove any excess liquid.  
(8) Add  $300 \mu L$  of Enzyme Reaction Solution,  $1 \mu L$  of RNase Solution and  $3 \mu L$  of Oxidation Inhibitor to the pellet, and vortex to detach and pull the pellet apart.  
(9) Incubate at  $37^{\circ}C$  for 10 minutes.  
(10) Add  $15 \mu L$  of Protein Digestion Solution and vortex to detach and pull the pellet apart.  
(11) Incubate at  $37^{\circ}C$  for 1 hour vortexing vigorously about every 20 minutes.

Frequent vortexing is more effective.

When the tissue has not dissolved completely, extend the incubation by 1 hour, and vortex occasionally in order to dissolve more quickly. We recommend vortexing frequently to dissolve the tissue quickly, because a long incubation time progresses the oxidation of DNA.

- (12) Centrifuge at  $10,000 \times G$  for 5 minutes at room temperature after the tissue has been dissolved completely.  
(13) Transfer the supernatant to a 1.5 mL microcentrifuge tube using a micropipette.  
Be careful not to take up the pellet.

The insoluble pellet is visible on the side which is subject to centrifugal force in the case of angle centrifugation. The pellet can be avoided by taking up the supernatant with a micropipette from the other side.

- (14) Add  $300 \mu L$  of Sodium Iodide Solution to the supernatant, mix, and flash centrifuge briefly. After that, add  $600 \mu L$  of Alcohol Solution and vortex.  
(15) Centrifuge at  $10,000 \times G$  for 10 minutes at room temperature.  
(16) Remove the supernatant from the tube then invert and place on paper towels for approx. 1 minute.  
Be careful not to drop the pellet while removing the supernatant.

A very slight amount of sodium iodide may remain in the DNA sample, and the DNA sample may show an UV (230 nm) absorption. The UV absorption of sodium iodide can be decreased by removing the supernatant carefully. Removing the sodium iodide solution which remains on the inner wall of the tube with tissue paper is more effective.

- (17) Add 1 mL of 70% Ethanol to the pellet and vortex to detach the pellet from the tube.  
(18) Centrifuge at  $10,000 \times G$  for 5 minutes at room temperature.  
(19) Repeat Step 16~18.  
(20) Remove the supernatant from the tube then invert and place on paper towels for approx. 5 minutes.  
Be careful not to drop the pellet while removing the supernatant.  
(21) Press the tube to paper towels in order to remove as much of the supernatant as possible, and leave the tube for approx. 10 minutes with the cap of the tube opened.  
Proceed to the next step while the pellet is still wet, because it is difficult to dissolve the pellet which is completely dry.  
(22) Add the appropriate amount of DDW or TE buffer, and vortex to detach the pellet from the tube.

The amount of DNA obtained varies between organs. Adjust the amount of DDW or TE buffer to prepare the appropriate concentration of DNA by referencing the table in Step 1. When assaying 8-OHdG using 「High Sensitive 8-OHdG Check (Code No.307-07921)」,  $1.5 \sim 1.8$  mg/mL is a benchmark.

- (23) Incubate at  $37^{\circ}C$  with occasionally vortexing to dissolve the pellet completely.

It may take longer to dissolve genomic DNA evenly. DNA may be dissolved by leaving it in the refrigerator over night, or with vigorous vortexing at high temperatures ( $50^{\circ}C \sim 60^{\circ}C$ ). Check the uniformity of the solution by visually checking whether the pellet (DNA) has been dissolved completely.

- (24) Measure the absorbance of OD 260 nm, and calculate the concentration of DNA.

The sample DNA should be stored at  $-80^{\circ}C$  because the DNA will continue to oxidize.

DNA obtained from spleen tissue may be brown in color. In the case of 8-OHdG assay, the color is removed from the DNA solution at the step of the centrifugal ultrafiltration after the treatment reaction with nuclease  $P_1$  and does not have an effect on 8-OHdG assay.

## 【Protocol 2】

- (1) Cut up the appropriate amount of tissue into small pieces and add 1 mL of Lysis Solution. Homogenize it in an ice water bath.

The table in Step 1 of Protocol 1 shows the maximum amount of tissue that can be treated. Do not exceed this amount.

The maximum amount of the tissue which can be treated is 150~200 mg, and the maximum amount of DNA is 400~500  $\mu$ g. It becomes difficult to purify DNA when the amount of the tissue or DNA exceeds the maximum. We recommend the use of the smallest amount of tissue in order to get good results when a large amount of DNA is not required.

- (2) Transfer the homogenate to a 2 mL microcentrifuge tube and centrifuge at  $600 \times G$  for 10 minutes at  $4^{\circ}C$ .

Centrifuge at  $5,000 \times G$  for 10 minutes at  $4^{\circ}C$  when the sample is brain tissue.  
(The condition of the centrifugation at Step 5 is  $600 \times G$ . This is the same as protocol.)

- (3) Decant and remove the supernatant from the tube then invert and place on paper towels.  
(4) Add 1 mL of Lysis Solution to the pellet and vortex to detach and pull the pellet apart.  
(5) Centrifuge at  $600 \times G$  for 10 minutes at  $4^{\circ}C$ .  
(6) Repeat Step 3~5.  
(7) Decant and remove the supernatant from the tube then invert and place on paper towels.  
(8) Add 300  $\mu$ L of Enzyme Reaction Solution, 15  $\mu$ L of Protein Digestion Solution and 3  $\mu$ L of Oxidation Inhibitor to the pellet, and vortex to detach and pull the pellet apart.  
(9) Incubate at  $37^{\circ}C$  for 1 hour vortexing vigorously about every 20 minutes.

Frequent vortexing is more effective.

When the tissue has not dissolved completely, extend the incubation by 1 hour, and vortex occasionally in order to dissolve more quickly. We recommend vortexing frequently to dissolve the tissue quickly, because a long incubation time progresses the oxidation of DNA.

- (10) Centrifuge at  $10,000 \times G$  for 5 minutes at room temperature after the tissue has been completely dissolved.  
(11) Transfer the supernatant to a 1.5 mL microcentrifuge tube using a micropipette.

Be careful not to take up the pellet.

The insoluble pellet is visible on the side which is subject to centrifugal force in the case of angle centrifugation. The pellet can be avoided by taking up the supernatant with a micropipette from the other side.

- (12) Add 300  $\mu$ L of Sodium Iodide Solution to the supernatant, mix, and flash centrifuge briefly. After that, add 600  $\mu$ L of Alcohol Solution and vortex.

- (13) Centrifuge at  $10,000 \times G$  for 10 minutes at room temperature.

- (14) Remove the supernatant from the tube then invert and place on paper towels for approx. 1 minute.

Be careful not to drop the pellet while removing the supernatant.

A very slight amount of sodium iodide may remain in the DNA sample, and the DNA sample may show an UV (230 nm) absorption. The UV absorption of sodium iodide can be decreased by removing the supernatant carefully. Removing the sodium iodide solution which remains on the inner wall of the tube with tissue paper is more effective.

- (15) Add 1 mL of 70% Ethanol to the pellet and vortex to detach the pellet from the tube.

- (16) Centrifuge at  $10,000 \times G$  for 5 minutes at room temperature.

- (17) Repeat Step 14~16.

- (18) Remove the supernatant from the tube then invert and place on paper towels for approx. 5 minutes.

Be careful not to drop the pellet while removing the supernatant.

- (19) Press the tube to paper towels in order to remove as much of the supernatant as possible.

- (20) Add 400  $\mu$ L of TE buffer(pH 8.0), 1  $\mu$ L of RNase Solution and 4  $\mu$ L of Oxidation Inhibitor, and vortex to detach the pellet from the tube.

- (21) Incubate at  $37^{\circ}C$  for 30 minutes vigorously vortexing every 10 minutes.

It is necessary to dissolve the pellet completely in order to finish the treatment of RNase. If the pellet is not dissolved completely after the incubation for 30 minutes, extend the incubation time by 5-10 minutes from when the pellet has been dissolved in the buffer uniformly.

- (22) After flash centrifuge, add 400  $\mu$ L of PEG Solution and vortex.

- (23) Centrifuge at  $10,000 \times G$  for 10 minutes at  $4^{\circ}C$ .

- (24) Remove the supernatant from the tube then invert and place on paper towels.

Be careful not to drop the pellet while removing the supernatant. Proceed to the next step even if some of the supernatant remains.

- (25) Add 1 mL of 70% Ethanol to the pellet and vortex to detach the pellet from the tube.

- (26) Centrifuge at  $10,000 \times G$  for 5 minutes at  $4^{\circ}C$ .

- (27) Remove the supernatant from the tube then invert and place on paper towels for approx. 5 minutes.  
Be careful not to drop the pellet while removing the supernatant.
- (28) Press the tube to paper towels in order to remove as much of the supernatant as possible, and leave the tube for approx. 10 minutes with the cap of the tube opened.  
Proceed to the next step while the pellet is still wet because it is difficult to dissolve the pellet which is completely dry.
- (29) Add the appropriate amount of DDW or TE buffer, and vortex to detach the pellet from the tube.

The amount of DNA obtained varies between organs. Adjust the amount of DDW or TE buffer to prepare the appropriate concentration of DNA by referencing the table in Step 1 of Protocol 1. When assaying 8-OHdG using 「High Sensitive 8-OHdG Check (Code No.307-07921)」, 1.5~1.8 mg/mL is a benchmark.

- (30) Incubate at 37 °C occasionally vortexing to dissolve the pellet completely.

It may take longer to dissolve genomic DNA evenly. DNA may be dissolved by leaving it in the refrigerator over night, or with vigorous vortexing at high temperatures (50 °C ~60 °C). Check the uniformity of the solution by visually checking whether the pellet (DNA) has been dissolved completely.


- (31) Measure the absorbance of OD260nm, and calculate the concentration of DNA.

The sample DNA should be stored at - 80 °C because the DNA will continue to oxidize.

DNA obtained from spleen tissue may be brown in color. In the case of 8-OHdG assay, the color is removed from the DNA solution at the step of the centrifugal ultrafiltration after the treatment reaction with nuclease P<sub>1</sub> and does not have an effect on 8-OHdG assay.

## Flow chart

### 【Protocol 1】

## 【Protocol 2】

Approx. 150 mg tissue  
|  
←1mL **Lysis Solution**  
|  
Homogenize  
|  
Homogenate/2mL microcentrifuge tube  
|  
600 × G, 4°C, 10 min.  
|  
Pellet  
|  
←1mL **Lysis Solution**  
|  
Vortex (Pull the pellet apart)  
|  
600 × G, 4°C, 10 min.  
|  
Pellet  
|  
←300 μL **Enzyme Reaction Solution**  
←15 μL **Protein Digestion Solution**  
←3 μL **Oxidation Inhibitor**  
|  
Vortex (Pull the pellet apart.)  
|  
37°C, 1hr (Vortex/20 min.)  
|  
10,000 × G, RT, 5 min.  
|  
Supernatant/1.5mL microcentrifuge tube  
|  
←300 μL **Sodium Iodide Solution**  
←600 μL **Alcohol Solution**  
|  
10,000 × G, RT, 10 min.  
|  
Pellet  
|  
←1mL 70% EtOH  
|  
10,000 × G, RT, 5 min.  
|  
Pellet (Nucleic acid)  
|  
←400 μL TE buffer  
←1 μL **RNase Solution**  
←4 μL **Oxidation Inhibitor**  
|  
37°C, 30 min. (Vortex/10 min.)  
|  
←400 μL **PEG Solution**  
|  
10,000 × G, 4°C, 10 min.  
|  
Pellet  
|  
←1mL 70% EtOH  
|  
10,000 × G, 4°C, 5 min.  
|  
Pellet (DNA)  
|  
←DDW or TE buffer  
DNA solution

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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