Tau ELISA kit *Wako*

[Introduction]

Tau is a microtubule-associated protein that is mainly expressed in neurons in the central nervous system and regulates stability of microtubules. In the brain in patients with Alzheimer's disease, aggregates of phosphorylated Tau or neurofibrillary tangles are formed, and extent of the formation has been reported to be correlated to severity of dementia. Tau, therefore, has been studied to investigate causes of Alzheimer's disease and to development drugs for treating the disease. On the other hand, the concentrations of Total Tau and phosphorylated Tau in cerebrospinal fluid are higher in Alzheimer's disease patients than in non-dementia subjects. This product is an ELISA kit that readily allows measurement of Tau. This kit recognizes total Tau whether it is phosphorylated or not.

[Kit Performance]

Standard curve range	4.1-1000 pg/mL
Reactive Tau	Total Tau
Sample	Human cerebrospinal fluid (CSF)
Sample volume	10 μ L *
Assay time	3 hours
Detection method	Luminescent system

*It is recommended that CSF sample volume is 50 μ L, considering accuracy of dilution.

[Kit Content]

Component	Condition	Volume or quantity
Antibody-coated Plate	Use as is	1 plate 96 wells (8 \times 12)
Tau Standard	Use after reconstitution	1 tube
Buffer	Use as is	60 mL/1 tube
Biotin-conjugated Antibody Solution	Use after preparation	100 μ L/1 tube
Peroxidase-conjugated Streptavidin Solution	Use after preparation	100 μ L/1 tube
Luminescent Reagent 1	Use as is	6 mL/1 tube
Luminescent Reagent 2	Use as is	6 mL/1 tube
Wash Solution (10 \times)	Use after preparation	100 mL/1 tube
Sample Buffer 1	Use as is	0.6 mL/1 tube
Sample Buffer 2	Use as is	30 mL/1 tube
Plate Seal	Use as is	3 sheets
Instruction manual		1 copy

[Measurement Principle]

The microplate is coated with anti Tau monoclonal antibody (epitope, C-terminal region of Tau corresponding to positions 412-432 a.a.). In each well, the standard solution or sample and biotin-conjugated anti Tau monoclonal antibody (epitope, middle region of Tau corresponding to positions 211-231 a.a.) are incubated to proceed with the antigen-antibody reaction. Furthermore, peroxidase-conjugated streptavidin is added to

proceed with the biotin-streptavidin binding. Finally, peroxidase activity in each well is measured to determine total Tau concentration in the sample.

Equipment or materials required but not supplied

- Purified water (distilled water)
- Tubes for dilution of a standard solution or sample
- Glass utensils for dilution of Wash Solution (graduated cylinder and beaker)
- Pipettes with disposable tips (one capable of pipetting 10 μ L of liquid accurately and one capable of pipetting 100 to 500 μ L)
- Repetitive dispenser pipette, capable of continuously dispensing 50 μ L and 100 μ L
- Water absorbing material such as paper towel (used to remove liquid remaining in the plate after rinsing)
- Stirrer (Vortex type)
- Microplate shaker (approximately 600 to 800 rpm)
- 96-well plate washer (preferable if available) or washing bottle
- 96-well plate reader (for luminescent measurement)
- Data processing software

[Preparation Methods of Reagents]

(1) Preparation of standard solution

Reconstitute the standard with the volume of purified water* stated in a separate sheet (standard stock solution, 10 ng/mL), and then mix with the Buffer included in the kit to prepare standard solutions as shown below.

*Because the volume of purified water to be added differs depending on the lot, see the specified volume in the annex

Concentration of the standard solution (pg/mL)	Volume of the preceding standard solution	Buffer
1000	Stock solution (10 ng/mL) : 50 μ L	450 μ L
400	Standard solution at 1000 pg/mL : 200 μ L	300 μ L
160	Standard solution at 400 pg/mL : 200 μ L	300 μ L
64.0	Standard solution at 160 pg/mL : 200 μ L	300 μ L
25.6	Standard solution at 64.0 pg/mL : 200 μ L	300 μ L
10.2	Standard solution at 25.6 pg/mL : 200 μ L	300 μ L
4.10	Standard solution at 10.2 pg/mL : 200 μ L	300 μ L
0.00	—	300 μ L

(2) Biotin-conjugated Antibody Solution

Dilute this solution 100 folds with the Buffer.

(3) Peroxidase-conjugated Streptavidin Solution

Dilute this solution 100 folds with the Buffer.

(4) Luminescent Reagent 1 and 2

Mix Luminescent Reagent 1 and 2 in equal volumes (1:1, vol./vol.) 15 to 30 minutes before use.

Store the mixture protected from light.

Example : Mix 6 mL each of Luminescent Reagent 1 and 2 (in a case where all of 96 wells are to be used)

(5) Wash Solution (10 \times)

Dilute this solution 10 fold with purified water (distilled water).

Example : Mix 100 mL of Wash Solution (10 \times) and 900 mL of purified water (distilled water) (in a case where all of 96 wells are to be used)

○ Use other reagents as they are.

[Stability and Storage Methods of Reagents]

(1) Antibody-coated Plate

Store unused antibody-coated strips (kept refrigerated with a seal covered) at 2°C to 10°C in a reclosable pack included in the kit. This component remains stable until the expiration date.

(2) Tau Standard

Reconstitute the standard with the volume of purified water stated in a separate sheet* (standard stock solution, 10 ng/mL). Proceed with the subsequent preparation using the Buffer included in the kit, which has been brought to room temperature. For the volume of purified water to be added, see the annex. (*because the volume of purified water to be added differs depending on the lot, see the specified volume in the annex) Use each standard solution immediately after the dilution, and do not store them.

(3) Buffer

If a part of the Buffer is used, transfer a volume slightly larger than that required to another container, close the lid tightly soon after use, keep the remaining Buffer cool without bringing it back to room temperature, and store it at 2°C to 10°C. This component remains stable until the expiration date.

(4) Biotin-conjugated Antibody Solution and (5) Peroxidase-conjugated Streptavidin Solution

If a divided part of the kit is used, dilute the stock solution immediately after taking out of the refrigerator, close the lid tightly soon after use, keep the remaining stock solution cool without bringing it back to room temperature, and store it at 2°C to 10°C. This component remains stable until the expiration date. Discard the remaining diluted solution after use.

(5) Luminescent Reagent 1 and 2

If a divided part of the kit is used, proceed with preparation immediately after taking the Luminescent Reagents out of the refrigerator, close the lid tightly soon after use, keep the remaining reagents cool without bringing them back to room temperature, and store them at 2°C to 10°C. This component remains stable until the expiration date.

(6) Wash Solution (10 ×)

Store the Wash Solution (10 ×) with the lid tightly closed at 2°C to 10°C. This component remains stable until the expiration date. Discard the remaining diluted wash solution after use.

(7) Sample Buffer 1 and 2

If a divided part of the kit is used, proceed with preparation immediately after taking the Sample Buffers out of the refrigerator, close the lid tightly soon after use, keep the remaining buffers cool without bringing them back to room temperature, and store them at 2°C to 10°C. This component remains stable until the expiration date.

[Sample Preparation Method]


1. To human cerebrospinal fluid (CSF) collected by the established procedure, add its 1/10 volume of Sample Buffer 1 that has been brought to room temperature, stir, and allow to stand at room temperature for 30 minutes. (preferably, to 50 μ L of CSF, add 5 μ L of Sample Buffer 1)
2. Then, further add its 4-fold volume of Sample Buffer 2 that has been brought to room temperature, stir, allow to stand at room temperature for 15 to 20 minutes, and use this mixture as the prepared sample (preferably, make 220 μ L). Note that Sample Buffer 2 has a blue color.

[Measuring Operations]

1. Fill each well of the antibody-coated plate with the previously diluted wash solution, and discard the solution to wash the well. Repeat the above washing operation 4 times in total. Then, tap the plate upside down gently against a paper towel to remove the remaining solution from wells.
2. Dispense 50 μ L of the diluted biotin-conjugated antibody solution to each of wells for standards and ones for samples.
3. Add 50 μ L of the standard solution at each concentration to each of wells for that standard solution.
4. Add 50 μ L of the sample prepared with sample buffers (prepared samples) to each of wells for that samples.
5. Shake the plate using a microplate shaker.
6. Attach a plate seal, and allow to stand at room temperature (20°C to 25°C) for 2.5 hours.
7. After end of the reaction, discard the reaction mixture. Fill each well with the diluted wash solution, and discard the solution to wash the well. Repeat the above washing operation 4 times in total. Then, tap the plate upside down gently against a paper towel to remove the remaining solution from wells.
8. Dispense 100 μ L of the diluted Peroxidase-conjugated streptavidin solution to each well. Shake the plate using a microplate shaker.
9. Attach a plate seal, and allow to stand at room temperature (20°C to 25°C) for 30 minutes.
10. After end of the reaction, discard the reaction mixture. Fill each well with the diluted wash solution, and discard the solution to wash the well. Repeat the above washing operation 4 times in total. Then, tap the plate upside down gently against a paper towel to remove the remaining solution from wells.
11. Dispense 100 μ L of the mixed luminescent reagent to each well. Shake the plate using a microplate shaker for 1 minute.
12. After stirring, measure the luminescent intensity with a 96-well plate reader (for luminescent measurement). It is recommended to perform the measurement within 10 to 20 minutes after addition of the combined luminescent reagent.
13. Prepare a standard curve by plotting the standard solution concentration (pg/mL) on the X-axis and the luminescent intensity on the Y-axis. Read the concentration (pg/mL) corresponding to the luminescent intensity of the diluted sample. Multiply the read concentration by a dilution factor for the sample (5.5) to obtain the measured value.

*For software-based computation, it is recommended to use cubic polynomial with 4 or 5 parameters.

Standard curve (example)

[Precautions for Use]


1. Please use polypropylene container. It is recommended to store samples at -35°C or lower if store for an extended period. Do not repeat freeze-thawing. Thaw the frozen sample just before use, and stir the thawed sample sufficiently. Prepare the samples before use.
2. Remove the turbidity or insoluble matters from a sample by centrifugation, etc., if any, before measurement.
3. If a sample is suspected of containing any interfering substance, measure the concerned samples at different dilution factors to check the dilution linearity. Dilute a prepared sample with the diluted wash solution ($1 \times$) if applicable.
4. This kit shall be used by an operator who has completed training for ELISA method or under appropriate supervision.
5. This kit shall be used by an operator who has a stably reproducible pipetting skill if the measurement is manually operated.
6. Wear gloves, glasses, and protective garment during preparation and operation with this kit.
7. Avoid skin contact with any reagent. In case where any reagent in this kit comes in contact with eyes, mouth, wound, or skin by mistake, provide emergency care, for instance, rinsing immediately with plenty of water, and seek medical advice if necessary.
8. Do not drink, eat, or smoke in a place where this kit is used.
9. Handle samples carefully as ones with a risk of infection. Note that this kit contains animal-derived substances.
10. Soak any used sample or consumables in a solution containing 1% formalin and 2% glutaraldehyde or $\geq 0.1\%$ sodium hypochlorite solution for at least 1 hour. Or autoclave them before disposition. Dispose used consumables and unused reagents in accordance with rules of the facility and regional regulations.
11. Do not mix reagents with different lot numbers for use. To prevent drying of the well surface, contamination with foreign matters, biased temperature distribution, or evaporation of dispensed reagents, affix the plate seal during static incubation in each step.
12. Note that any ELISA method is affected by measurement environment. Ensure that places for measurement operations and static incubation are kept at room temperature of 20°C to 25°C (on the bench or in an incubator). In addition, avoid the measurement in an environment with air flow (including one from an air-conditioning system) and low humidity.

[Example of measurement]

Spike-recovery test


Human CSF [1]	Unspiked	Spiked with Tau at 20 pg/mL		Spiked with Tau at 60 pg/mL		Spiked with Tau at 150 pg/mL	
	Measured value	Measured value	Recovery	Measured value	Recovery	Measured value	Recovery
	4.09 pg/mL	24.8 pg/mL	104%	58.1 pg/mL	90%	175 pg/mL	114%
Human CSF [2]	Unspiked	Spiked with Tau at 40 pg/mL		Spiked with Tau at 135 pg/mL		Spiked with Tau at 380 pg/mL	
	Measured value	Measured value	Recovery	Measured value	Recovery	Measured value	Recovery
	24.8 pg/mL	70.5 pg/mL	114%	160 pg/mL	100%	390 pg/mL	96.1%

Dilution linearity

Measurement of CSF samples

CSF samples from non-dementia subjects (non-dementia), patients with mild cognitive impairment (MCI), and those with Alzheimer's disease (AD) were measured with this ELISA kit.

Non-dementia1	430	MCI1	295	AD1	458
Non-dementia2	157	MCI2	294	AD2	231
Non-dementia3	246	MCI3	369	AD3	260
Non-dementia4	52.7	MCI4	142	AD4	826
Non-dementia5	14.1			AD5	230
Non-dementia6	397			AD6	670
Non-dementia7	53			AD7	1000
Non-dementia8	112			AD8	201
Non-dementia9	178			AD9	1011
Non-dementia10	203			AD10	886
Mean for non-dementia	184	Mean for MCI	275	Mean for AD	577

→ A difference was observed in measured value between non-dementia subjects and patients with Alzheimer's disease.

[Outline of Measurement Procedures]

- Bring the plate and reagents back to room temperature (20°C to 25°C).
- Dilution of Wash Solution concentrate : Dilute 10 folds with purified water that has been brought to room temperature.
- Dilution of the standard solution (example) : Reconstitute the standard with the volume of purified water stated in a separate sheet (standard stock solution, 10 ng/mL), and then mix with the Buffer of the kit to prepare standard solutions as shown below.

Example of dilutions	Concentration (pg/mL)	1000	400	160	64.0	25.6	10.2	4.10	0
Standard solution (μL)	Stock	50	200*	200*	200*	200*	200*	200*	0
Buffer (μL)		450	300	300	300	300	300	300	300

*: Preceding standard solution at a one-level higher concentration

- Preparation of diluted biotin-conjugated antibody solution (dilute the Biotin-conjugated Antibody Solution 100 folds with Buffer that has been brought to room temperature.)
- Preparation of working (prepared) samples
To a sample, add its 1/10 volume of Sample Buffer 1, stir, and allow to stand at room temperature for 30 minutes. Then, further add its 4-fold volume of Sample Buffer 2, stir, and allow to stand at room temperature for 15 to 20 minutes. (Preferably : 50 μL sample, 5 μL sample buffer 1, 220 μL sample buffer 2)

• Antibody-coated Plate

- ↓ Repeat washing operation 4 times (*[1])
- Diluted biotin-conjugated antibody solution 50 μL/well
- Standard solution or prepared sample 50 μL/well
- ↓ Stirring (*[2]) and then static incubation for reaction at room temperature (20°C to 25°C) for 2.5 hours (*[3])

↓

- *Preparation of diluted peroxidase-conjugated streptavidin solution (dilute Peroxidase-conjugated Streptavidin Solution 100 folds with Buffer that has been brought to room temperature)

- ↓ Repeat washing operation 4 times (*[1])

- Diluted peroxidase-conjugated streptavidin solution 100 μL/well

- ↓ Stirring (*[2]) and then static incubation for reaction at room temperature (20°C to 25°C) for 30 minutes (*[3])

↓

- *Preparation of combined luminescent reagent (mix Luminescent Reagents 1 and 2 in equal volumes [1:1, vol./vol.])

- ↓ Repeat washing operation 4 times (*[1])

- Mixed luminescent reagent 100 μL/well

- ↓ Stirring at room temperature (20 to 25°C) for 1 minute

- Measurement of luminescent intensity (within 10 to 20 minutes after addition)

(*[1]) In each washing operation, dispense the diluted wash solution to wells, shake the plate gently on the palm of your hand for approximately 10 seconds, and discard the solution. After 4 continuous washing operations, tap the plate upside down against a paper towel to remove the remaining solution completely from wells. After removing the diluted wash solution, dispense the next solution immediately with care so as not to dry the well surface. If a pipette is used to dispense the diluted wash solution, set the volume at 300 μL/well.

(*[2]) Stir at 600 to 800 rpm for 10 seconds, and repeat 3 times.

(*[3]) After stirring, affix a plate cover seal, and allow to stand. Peel the liner from a plate seal, and attach the seal with the adhesive side on the plate. Do not reuse any plate seal.

Storage Condition : Store at 2-10°C

Expiration date : Indicated on the label.

Package : For 96 tests

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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