SLP-HS Single Reagent Set II

《Research Use only, Not for use in diagnostic procedures》

【Preface】

The hemolymph of the silkworm *Bombyx mori* contains a self-defense mechanism termed the “prophenoloxidase (proPO) activating system” or “proPO cascade”. Upon invasion of microbe such as bacteria and fungi, the system participates in melanin formation observed in the body fluid of insects to protect them from the invader’s attack. The system is triggered by peptidoglycan (PG) from bacteria and (1 → 3)- β -D-glucan (β -glucan) from fungi and yeast, and consequently, proPO in the system is activated^{1) 2)}.

The system is thought to be a cascade reaction involving activation of multiple zymogens of proteinases.

The SLP (Silkworm Larvae Plasma) reagent is a lyophilized product containing factors of the proPO cascade, it is prepared from silkworm larvae body fluid and processed without melanin formation. The SLP reagent is activated by PG and/or β -glucan, then L-3,4-dihydroxyphenylalanine (DOPA) is oxidized and melanin is formed. Since PG is found on the cell wall of most bacteria and β -glucan is a component found in cell walls of many fungi, it is possible to detect various microorganisms by using melanin pigment formation of SLP reagent as an index³⁾.

【Features】

1. High sensitive detection of PG and β -glucan
2. Accurate and sensitive quantification of PG and β -glucan with Toxinometer

【Assay Principle】

The activation mechanism of SLP is considered as shown in Figure 1. Either PG or β -glucan binds to the respective recognition proteins termed PGRP or β GRP. The proPO cascade system is initiated by these reactions, and consequently, proPO is converted to phenoloxidase (PO)⁴⁾. The PO catalyzes oxidation of DOPA and is followed by the formation of melanin in the mixture. In the kinetic colorimetric assay with the Toxinometer, the amount of melanin pigment by SLP reaction is monitored as a change in the light transmission rate, and the reaction time until the light transmission rate reaches a preset threshold value is defined as activation time (Ta). SLP activating substances (PG and β -glucan) are quantified with Ta as an index.


Fig 1. Activation mechanism of SLP system⁴⁾
(modified from Ref.4)

【Kit Contents】

1. SLP-HS Reagent II (Lyophilized reagent containing Silkworm Larvae Plasma and DOPA) for 0.1 mL × 20 vials
Sensitivity : 10 pg/mL for PG, 1 pg/mL for β-glucan in 120 minutes.
2. SLP-Diluent 5 mL × 2 vials
3. Standard (Digested Peptidoglycan from *Staphylococcus aureus*) 0.5 mL × 1 vial

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1. Materials required but not supplied

- (1) Pipettes
- (2) Test tubes with aluminum caps for dilution of standard or sample
- (3) Vortex mixer
- (4) Toxinometer
- (5) Test water for dilution of standard or sample

NOTE : A test tube and the pipette should be dry heat sterilized by 250°C and more than two hours. Test water should be confirmed that they are no contamination and no interference for the measurement. The plastic products such as tips for micropipettes should be confirmed that they are no contamination and no interference for the measurement, or products shown in "Related Products" should be used.

2. Preparation of reagents

- (1) Dilution of peptidoglycan standard and sample
Vortex vigorously for approximately 120 seconds before opening the standard. Vigorously vortex each dilution step for about 20 seconds before preparing the next dilution.

**Table 1 Example of preparing dilution series of standard
(When the labeled concentration is 7.2 $\mu\text{g/mL}$)**

Standard solution (μL)	Test water (μL)	Prepared standard concentration (ng/mL)
100 (7.2 $\mu\text{g/mL}$)	620	1,000
100 (1,000 ng/mL)	900	100
100 (100 ng/mL)	900	10
100 (10 ng/mL)	900	1
100 (1 ng/mL)	900	0.1
100 (0.1 ng/mL)	900	0.01

(2) Preparation of SLP-HS Reagent II solution

Take SLP-HS Reagent II vials out of the package and put the vials in a rack. When powder of SLP-HS Reagent II is attached to the rubber stopper, knock on several times of the bottom of the vial. A little reagent attached to a rubber stopper does not influence the result of the measurement. Each rubber stopper is loosened by releasing the vacuum and slowly pulling up the rubber stopper slightly.

Dispense 0.1 mL of the SLP-Diluent into SLP-HS Reagent II vial and put the rubber stopper on the vial. Mix gently to dissolve the reagent without making bubbles. Please proceed to the step of "3. Measurement procedure" promptly after dissolving.

3. Measurement procedure

(1) Add 0.1 mL of dilution series of Standard, test water, or samples to SLP-HS Reagent II solution.

(2) Mix immediately so as not to make bubbles, and start to measure with the following conditions using Toxinometer at 30°C.

< Recommended measurement conditions for Toxinometer >

Threshold 90.0%
Count 3
Wait time 5 minutes
Measurement Time 120 minutes

(3) Calculate the concentration of SLP reactive substance in the sample from the activation time (T_a) of the standard and T_a of the sample as the concentration equivalent to peptidoglycan standard. A good fitting can be obtained from the logarithm (Log) of the peptidoglycan standard concentration on the x-axis and the double log (LogLog) of T_a on the y-axis with the quadratic regression of x as a calibration curve.


Fig. 2 Example of standard curve graph

【Usage notes】

1. Since this product reacts very sensitively to peptidoglycan and β -glucan, please pay attention to contamination from pipette, other instruments, test water etc. Please use equipment to be sterilized by dry heat at 250°C for 2 hours or more, or as described in "Related Products".
2. Do not use the reagent if it is discolored or insoluble materials are observed when dissolving.
3. Since the toxicity of this product is not confirmed well, please pay attention to handling so as not to inhale it.
4. The sensitivity of SLP reagent is suppressed by samples with a high ionic strength.
5. This kit is for research use only and not for use in diagnostic procedures.
6. Carefully open the vial so as not to get injured by the aluminium cap.

【Storage】 Keep at 2 ~ 10°C

【Expiration】 Indicate on the label (outer box)

【Package】 20 tests

【References】

- 1) Ashida, M. *Invertebrate Biological Defense* ed. by Natori, S. *et al.* Society Publishing Center, 111-126 (1992)
- 2) Ashida, M. and Yamazaki, H. I. *Molting and Metamorphosis* ed. by Ohnishi, E. and Ishizaki, H., Japan Sci. Soc. Press, Tokyo/Springer-Verlag, Berlin, 239-265 (1990)
- 3) Tsuchiya, M., Asahi, N., *et al.* *FEMS Immunol. Med. Microbiol.*, **15**, 129-134 (1996)
- 4) Ochiai, M. *SANSHI-KONCHU BIOTEC* **84** (3), 195-204 (2015)

【Related Products】

293-26551	Limulus Test Tube-S (12 × 75 mm, Endotoxin Free)	10 vials × 10
293-28251	Aluminium Cap-S (14.7 × 18 mm, Endotoxin Free)	10 pieces × 10
298-35031	Bio Clean Tip Wako® 1000 II	100 tips

291-35021	Bio Clean Tip Wako® 200 II	100 tips
294-35011	Bio Clean Tip Wako® Extend S II	100 tips

Manufactured by
FUJIFILM Healthcare Manufacturing Corporation

2613-2, Ogohara, Komono-Cho, Mie-Gun, Mie 510-1222, Japan

Distributed by
FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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