## Mature BDNF ELISA Kit *Wako*

### [1. Introduction]

Brain-derived neurotrophic factor (BDNF) is a member of NGF family of neurotrophic factors. It is known that BDNF is involved in neurogenesis and synaptogenesis. Therefore BDNF is expected as biomarker for diseases of nervous system. BDNF also has relation to heart failure and heart disease. As above, BDNF is a subject of study in various researches.

BDNF has precursor termed proBDNF, which is converted to mature BDNF (mBDNF) through the proteolytic removal of the N-terminal fragment by specific protease. It is reported that there are different functions between mBDNF and proBDNF.

Mature BDNF ELISA Kit *Wako* can specifically detect mBDNF.

### [2. Performance]

Principle of the assay	Sandwich method
Standard curve range	4.1 ~ 1,000 pg/mL
Specificity <sup>*1</sup>	Mature BDNF
Sample	Serum, Plasma (Human)
Sample volume	10 $\mu$ L
Assay time	4 hours
Detection method	Colorimetric method

<sup>\*1</sup> This kit has approximately 10% cross-reactivity to recombinant human proBDNF.

### [3. Material supplied]

Components	State	Volume or quantity
(A) Antibody-coated Plate	Use after washing	1 plate 96 wells (8 $\times$ 12)
(B) Mature BDNF Standard	Use after reconstitution	1 vial
(C) Buffer	Ready-to-use	60 mL/1 vial
(D) Biotin-conjugated Antibody Solution	Concentrated Use after dilution	100 $\mu$ L/1 vial
(E) Peroxidase-conjugated Streptavidin Solution	Concentrated Use after dilution	100 $\mu$ L/1 vial
(F) TMB Solution	Ready-to-use	6 mL/1 vial
(G) Stop Solution	Ready-to-use	6 mL/1 vial
(H) Wash Solution (10 $\times$ )	Concentrated Use after dilution	100 mL/1 vial
(I) Plate Seal	Ready-to-use	4 sheets

### [4. Measurement Principle]

The microplate is coated with anti mature BDNF monoclonal antibody. In each well, the standard solution or sample and biotin-conjugated anti BDNF polyclonal antibody are incubated to proceed with the antigen-antibody reaction. Furthermore, peroxidase-conjugated streptavidin is added to proceed with the biotin-streptavidin binding. Finally, peroxidase activity in each well is

measured to determine mature BDNF in the sample.

**[5. Equipment or Materials Required but not Supplied]**

- ☐ Purified water (distilled water)
- ☐ Test tube for dilution of a standard solution and samples
- ☐ Glass utensils for dilution of Wash Solution (e.g. graduated cylinders and beakers)
- ☐ Pipettes with disposable tips (one capable of pipetting 50  $\mu\text{L}$  of liquid accurately and one capable of pipetting 100 to 1,000  $\mu\text{L}$ )
- ☐ Dispenser capable of dispensing
- ☐ Water-absorbable material such as paper towel (to remove liquid remaining on a plate after washing)
- ☐ Mixer (Vortex type)
- ☐ Shaker for 96-well plate (approximately 600 to 1,200 rpm)
- ☐ Automatic washer for 96-well plate (If available) or washing bottle
- ☐ Microplate reader capable of measuring at  $450 \pm 10$  nm, with the correction wavelength set at 600-650 nm
- ☐ Software for data analysis

**[6. Reagent Preparation]**

Bring the reagents to room temperature prior to use. Reagents described as "Use after dilution" and "Use after reconstitution" in [3. Material Supplied] shall be prepared as follows. Prepare only as much reagent as needed on the day of the experiment.

6.1 Standard Solution

Reconstitute the standard with volume of purified water\*<sup>2</sup> described in a separate sheet (standard stock solution 10 ng/mL), and then mix with Buffer to prepare standard solutions as shown below.

\*<sup>2</sup> Because the volume of purified water to be added differs depending on the lot, see the specified volume in the separate sheet.

Concentration of the standard solution (pg/mL)	Volume of the preceding standard solution	Buffer
1000	Stock Solution (10 ng/mL) : 30 $\mu\text{L}$	270 $\mu\text{L}$
400	Standard Solution at 1000 pg/mL : 100 $\mu\text{L}$	150 $\mu\text{L}$
160	Standard Solution at 400 pg/mL : 100 $\mu\text{L}$	150 $\mu\text{L}$
64	Standard Solution at 160 pg/mL : 100 $\mu\text{L}$	150 $\mu\text{L}$
25.6	Standard Solution at 64.0 pg/mL : 100 $\mu\text{L}$	150 $\mu\text{L}$
10.2	Standard Solution at 25.6 pg/mL : 100 $\mu\text{L}$	150 $\mu\text{L}$
4.1	Standard Solution at 10.2 pg/mL : 100 $\mu\text{L}$	150 $\mu\text{L}$
0	—	150 $\mu\text{L}$

6.2 Biotin-conjugated Antibody Solution

Biotin-conjugated Antibody Solution is diluted 100-fold with Buffer.

6.3 Peroxidase-conjugated Streptavidin Solution

Peroxidase-conjugated Streptavidin Solution is diluted 100-fold with Buffer.

6.4 Wash Solution (1 $\times$ )

Wash Solution (10 $\times$ ) is diluted 10-fold with purified water (distilled water). (e.g.) When solution for 96-well reaction is prepared

Add 100mL of Wash Solution (10×) to 900 mL of purified water (distilled water) to prepare 1,000 mL of Wash Solution (1×).

#### **[7. Stability and Storage Method of Each Reagent]**

##### (A) Antibody-coated Plate

Unused antibody-immobilized strips (kept refrigerated and sealed) should be returned into the zip-seal bag provided in the kit and stored at 2-10°C. These are stable until the expiration date.

##### (B) Mature BDNF Standard

Reconstituted standard solution, which is standard stock solution 10 ng/mL should be stored at 2-10°C and should be used up within 2 weeks. Discard the remaining diluted standard solution after use.

##### (C) Buffer

If a part of Buffer is used, transfer the volume slightly greater than needed to another container, immediately close cap lid tightly without bringing the remaining Buffer to room temperature, and store at 2-10°C. It is stable until the expiration date.

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If the kit is divided for multiple assays, prepare these solutions by dilution just after taking out of the refrigerator, immediately close the cap tightly without bringing the remained stock solution to room temperature, and store at 2-10°C. It is stable until the expiration date. Discard the remaining diluted standard solution after use.

##### (E) Peroxidase-conjugated Streptavidin Solution

If the kit is divided for multiple assays, prepare these solutions by dilution just after taking out of the refrigerator, immediately close the cap tightly without bringing the remained stock solution to room temperature, and store at 2-10°C. It is stable until the expiration date. Discard the remaining diluted standard solution after use.

##### (F) TMB Solution

If a part of TMB solution is used, transfer the volume slightly greater than needed to another container, immediately close the lid tightly without bringing the remaining TMB solution to room temperature, and store at 2-10°C. It is stable until the expiration date.

##### (G) Stop Solution

Just after get the Stop Solution out of the refrigerator, dispense the required volume of the solution. The remained solution should be stored at 2-10°C with the cap tightened. It is stable until the expiration date.

##### (H) Wash Solution (10×)

Store the Wash Solution (10×) with the cap tightly closed at 2-10°C if applicable. It is stable until the expiration date. Discard the remaining diluted washing solution.

#### **[8. Sample Preparation]**

Serum : Dilute 20 to 100-fold with Buffer.

Plasma (EDTA) : Dilute 10 to 50 fold with Buffer.

### **[9. Assay Procedure]**

1. Discard the solution filled the Antibody-coated Plate with. Then wash 4 times with Wash Solution (1×). Invert the plate and gently blot it against clean paper towels.
  2. Add 50  $\mu$ L of diluted standard solution to each well.
  3. Add 50  $\mu$ L of diluted sample to each well.
  4. Agitate the plate on a microplate shaker.
  5. Cover with a Plate Seal. Incubate for 2 hours at room temperature (20-25°C).
  6. Discard the solution and wash 4 times with Wash Solution (1×). Invert the plate and gently blot it against clean paper towels.
  7. Add 50  $\mu$ L of Biotin-conjugated Antibody Solution to each well.
  8. Cover with a new Plate Seal. Incubate for 1 hour at room temperature (20-25°C).
  9. Repeat the wash as in step 6.
  10. Add 50  $\mu$ L of Peroxidase-conjugated Streptavidin Solution to each well. Then shake using microplate shaker.
  11. Cover with a new Plate Seal. Incubate for 30 minutes at room temperature (20-25°C) .
  12. Repeat the wash as in step 6.
  13. Add 50  $\mu$ L of TMB Solution to each well. Then shake using microplate shaker.
  14. Cover with a new Plate Seal. Incubate for 30 minutes at room temperature (20-25°C).
  15. Add 50  $\mu$ L of Stop Solution to each well.
  16. Measure the absorbance at 450 nm and 620 nm as reference wavelength (600-650 nm)<sup>\*3</sup> using a microplate reader after stirring microplate with microplate shaker.
- \*3 The range of 600-650 nm can be used as reference wavelength.

### **[10. Calculation of Results]**


Prepare a standard curve by plotting the absorbance readings (Y-axis) against values of standard concentration (X-axis). Determine the unknown sample concentration from the Standard Curve<sup>\*4</sup> and multiply the value by the dilution factor.

\*4 The use of a 3rd order regression curve for log-log plot or 4 or 5 parameters method for log normal plot in computer calculation is recommended.

Sample values obtained with this kit could be converted to the NIBSC/WHO BDNF ; WHO Reference Reagent Brain-derived neurotrophic factor (BDNF) code 96/534 Units, use multiplier coefficient<sup>\*5</sup> described in the separate sheet.

\*5 Multiplier coefficient differs depending on the lot, see the specified value in the separate sheet.

#### [11. Typical Data]

※If absorbance of highest standard concentration (1,000 pg/mL) is more than 3.0, the standard curves should be created except for the absorbance.

#### [12. Precautions]

- Wear protective gloves, clothing, eye, and face protection.
- Avoid contact of skin and mucous membranes with kit reagents or specimens. If any reagents come in contact with eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.
- Do not drink, eat, or smoke in a place where this kit is used.
- Handle the sample with proper care, being aware that the sample may have an infection risk. This kit contains animal-derived ingredients.
- Do not use reagents with different lot numbers together.
- TMB Solution is slightly yellow and clear liquid. Protect from light.
- It is recommended to store the sample frozen at -35°C or lower if long-term storage is intended. Avoid repeated freezing and thawing. Thaw the frozen sample just before assay and thoroughly agitate. Prepare the sample before use.
- Do not use a sample with hemolysis or containing high lipid.
- Centrifuge the sample to remove turbidity or insoluble matters where necessary before use for the assay.
- If presence of interfering substances is suspected for the sample, dilute the concerned one at multiple dilution ratios to check the dilution linearity. Dilute the prepared sample with the Buffer.
- When allowing the plate to stand in each step, always affix a plate sealer to protect the wells from drying, contamination with foreign matters, uneven temperature, and evaporation of a dispensed reagent.
- ELISA can be affected by assay environment. Ensure that room temperature at places for assay operation and incubation is strictly controlled at 20°C to 25°C (on a bench or in an incubator). Avoid performing assay under air flow (including air flow from an air-conditioning equipment) or in a low-humidity environment.

#### [13. Assay Procedure Summary]

- ☐ Bring the plate and reagents to room temperature (20-25°C) before use.
- ☐ Dilution of Wash Solution (10×) : Dilute Wash Solution (10×) 10-folds with purified water.

- ☐ Dilution of the standard solution: Reconstitute it with volume of purified water described in a separate sheet (standard stock solution 10 ng/mL), prepare the standard solutions by mixing with the Buffer to prepare as shown in a table below.

Example of dilution

Concentration (pg/mL)	1000	400	160	64.0	25.6	10.2	4.10	0
Standard Solution (μL)	30	100*	100*	100*	100*	100*	100	-
Diluent Buffer (μL)	270	150	150	150	150	150	150	150

\* : Standard solution at a 1-step higher concentration

- ☐ Antibody-coated Plate
- ☐ ↓ Washing 4 times (\*①)
- ☐ Sample or standard solution 50 μL/well
- ☐ ↓ Agitate (\*②) and incubate at room temperature (20-25°C) for 2 hour (\*③)
- ☐ \* Preparation of biotin-conjugated antibody solution (dilute 100- folds with Buffer)
- ☐ ↓ Washing 4 times (\*①)
- ☐ Biotin-conjugated antibody solution 50 μL/well
- ☐ ↓ Agitate (\*②) and incubate at room temperature (20-25°C) for 1 hour (\*③)
- ☐ \* Preparation of peroxidase-conjugated streptavidin solution (dilute 100- folds with Buffer)
- ☐ ↓ Washing 4 times (\*①)
- ☐ Peroxidase-conjugated streptavidin solution 50 μL/well
- ☐ ↓ Agitate (\*②) and incubate at room temperature (20-25°C) for 30 minutes (\*③)
- ☐ ↓ Washing 4 times (\*①)
- ☐ TMB Solution 50 μL/well
- ☐ ↓ Agitate (\*②) and incubate at room temperature (20-25°C) for 30 minutes (\*③)
- ☐ Stop Solution 50 μL/well
- ☐ ↓ Agitate (\*②)
- ☐ Measure absorbance at 450 nm (Reference wavelength 620 nm : 600 ~ 650 nm)

(\*①) In each washing operation, dispense the wash solution (1×) into wells, gently agitate the filled plate on the palm for about 10 seconds, and then empty the wells. After washing the wells 4 times consecutively, reverse the plate, and tap it against paper towel to remove the washing solution completely. After removal of the washing solution, immediately dispense the next solution with care not to dry the wells. Use of a pipet set at the liquid volume of 300 μL may be appropriate for dispensing the washing solution into each well.

(\*②) Three repeats of agitation at 600 to 1,200 rpm for 10 seconds may be appropriate.

(\*③) After agitation, cover with a plate sealer. Remove the liner from the plate sealer, and apply its adhesive side to the plate for affixation. Do not re-use any plate sealer.

**[Storage]** Store at 2-10°C

**[Expiration date]** Indicated on the label

**[Package]** For 96 assays

---

## **FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://flwk.fujifilm.co.jp>

### **FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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