**(R)-BiAC Derivatization Reagents Set**

This product is a set of reagents necessary for the chiral amino acids derivatization using (R)-BiAC (Code No. 025-19761).

**[Content of the reagent set]**


Code No.	Product Name	Volume	Amount
293-86011	Buffer for (R)-BiAC	3 mL	1
290-86021	Reaction Terminator	10 mL	1

**[Complementary products]**

Code No.	Product Name	Grade
012-19851	Acetonitrile	for LC/MS
025-19761	(R)-BiAC	for Amino Acid Analysis
235-64051	Wakopak® Ultra APDS TAG® φ 2.1 mm × 100 mm (D)	—

**[Principle]**


When (R)-BiAC solution is added and heated to a sample containing amino acids under basic conditions, a derivatized amino acid will be formed which is (R)-BiAC-tag bound to the amino group. The derivatized amino acid is separated into D- and L-amino acids by reversed-phase liquid chromatography due to the axially chiral biaryl skeleton and hydrophobicity of (R)-BiAC-tag, and detected by *m/z* in a mass spectrometer.

**[Operation procedures]****1. Preparation of reagent****1-1. (R)-BiAC solution**


Rinse the entire bottle of (R)-BiAC (Code No. 025-19761) with 1-2 mL of Acetonitrile (Code No. 012-19851) and transfer the solution into another container. Add Acetonitrile until the total volume becomes 5 mL and mix until it completely dissolved.

After dividing the prepared solution into the required number, store it at -20°C, avoid repeated freezing and thawing, and use it as soon as possible.

It can be stored at -20°C for about a week.

**1-2. Buffer for (R)-BiAC-acetonitrile solution**

Add 3 mL of Acetonitrile (Code No. 012-19851) to the Buffer for (R)-BiAC container and mix.

**2. Preparation of sample**

After centrifuging the sample solution, use the supernatant.

**3. Derivatization**

3-1. Add 30  $\mu$ L of Buffer for (R)-BiAC-acetonitrile solution and 10  $\mu$ L of (R)-BiAC solution into 15  $\mu$ L of the sample solution. Pipette 5 times, stir, and heat at 55°C for 10 min.

3-2. Cool the reaction solution up to room temperature, then mix it with 100  $\mu$ L of the Reaction Terminator.

**4. LC/MS Measurement**

4-1. Measure the solution prepared in procedure 3-2, using LC/MS or LC/MS/MS.

<Reference information>

Analysis examples are introduced on our website.

Details of LC/MS analysis conditions, etc. can be found below.

TOP > Analytical > Amino Acid, Peptides & Proteins > (R)-BiAC Method (Chiral Amino Acid LC/MS Analysis)

<https://labchem-wako.fujifilm.com/us/category/03157.html>

**[Instruction for use and handling]****1. Hazardous prevention**

1-1. Wear protective gloves, protective goggles or face shield, etc. Use respirator if necessary.

1-2. In case of contact, consult a physician.

1-3. In case of contact with eyes, skin and clothes etc., wash with plenty of water.

1-4. Safety Data Sheet is available for the detailed information.

**2. Precautions**

2-1. Follow the instruction for storing reagent and do not use chemicals past its expiration date.

2-2. Do not use the container and accessories in the kit for other purpose.

2-3. Use Acetonitrile, (R)-BiAC, and HPLC column sold separately.

**3. Disposal consideration**

Any disposal practice must be in compliance with country, local, state, and federal laws and regulations (contact country, local or state environmental agency for specific rules).

**[Storage / expiration date]**

Please check the label attached to the product.

**[References]**

Harada, M. and Karakawa, S. *et al.* : *J. Chromatogr. A*, **1593**, 91 (2019).

Harada, M. and Karakawa, S. *et. al.* : *Wako Junyaku Jiho*, **87** (1), 5 (2019).

Harada, M. and Karakawa, S. *et. al.* : *Symmetry*, **12**, 913 (2020).

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://ffwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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