LBIS™ Mouse Insulin ELISA Kit (High Sensitivity)

Please, read this instruction carefully before use.

1. Intended use

LBIS™ Mouse Insulin ELISA Kit (High Sensitivity) is a highly sensitive sandwich ELISA system for quantitative measurement of mouse insulin. This kit is intended for research use only.

2. Introduction

Insulin is a peptide hormone secreted from B cells of islet of Langerhans in the pancreas with a molecular weight of about 5800 and pI 5.4. It is consisted of 2 chains, A and B. It has 3 disulfide bonds formed between A6 and A11, A7 and B7, and A20 and B19. Insulin exists as a dimer molecule in acidic to neutral solution without Zn ion, and as a hexamer including two Zn ions in neutral solution if Zn ions are present. Main targets of insulin are liver, muscle, and adipose tissue. Insulin actions in these targets are as follows. In the liver, it promotes glycogenesis, protein synthesis, fatty acid synthesis, carbohydrate utilization, and inhibition of gluconeogenesis. In the muscle, it promotes membrane permeability for carbohydrates, amino acids and K ion, glycogenesis, protein synthesis, while inhibits protein degradation. In the adipose tissue, it promotes membrane permeability for glucose and fatty acid synthesis. A precursor of insulin, called proinsulin with a single polypeptide chain, is first synthesized in the cell, then sulfide bonds are formed, and finally by enzymatic cutting at two sites, active insulin and c-peptide (connecting peptide) are formed. Potency of an insulin preparation was originally determined by bioassay. However, whole body bioassay inevitably shows poor precision owing to individual variation.

WHO issued 1st International Standard for human insulin in 1986 which has the potency of 26 IU/mg (0.038 mg/IU). In the same year, 1st International Standard of bovine insulin, the potency of which is 25.7 IU/mg, and Porcine insulin 1st International Standard, 26 IU/mg, were provided. Before these standards, in 1974, 1st International Reference Preparation of human insulin for immunoassay was provided as 3 IU/ampoule. Based on the above data, if the biological activity of insulin per molecule is the same among various animal species, potencies of animal insulin might be calculated from their molecular weights. But, so far, we do not have experimental proof about this. As the molecular weights of insulin of various animals are nearly the same, and the differences are within 1%, there may be no critical fault if we think that the general potency of insulin is 26 IU/mg. Rat and mouse have two molecular species of insulin, type 1 and type 2. Amino acid sequences of these molecular species are same between rat and mouse. But as their ratios are different between these two animal species, it is recommended to use standard preparation derived from each animal.

3. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Mouse Insulin ELISA Kit (High Sensitivity), biotinylated anti-insulin antibody, and standard or sample are incubated in monoclonal anti-insulin coated wells to capture insulin bound with biotinylated anti-insulin antibody. After 2 hours' incubation and washing, HRP (horse radish peroxidase)-conjugated streptavidin is added, and incubated for 30 minutes. After washing, HRP-conjugated streptavidin remaining in wells are reacted with a TMB Solution for 30 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is proportional to insulin concentration. The standard curve is prepared by plotting absorbance against standard insulin concentrations. Insulin concentrations in unknown samples are determined using this standard curve.

4. Performance Characteristics

• Assay range

The assay range of the kit is 39 pg/mL - 2500 pg/mL.

If some samples show absorbance more than that of 2500 pg/mL standard, please repeat the assay after proper dilution of samples.

• Specificity

The antibodies used in this kit are specific to insulin. Cross-reactivity of the kit is shown below.

Substances	Cross-reactivity	Remarks
Mouse C-peptide	—	At 2.5 ng/mL
Mouse proinsulin	+	At 2.5 ng/mL
Rat insulin	+	At 2.5 ng/mL
Human insulin	+	At 2.5 ng/mL

- Precision of assay
Within assay variation (3 samples, 5 replicates assay) Mean CV was within 10%.
- Reproducibility
Between assay variation (3 samples, 3 days, assayed in triplicates) Mean CV was within 10%.
- Recovery test
Standard insulin was added in 3 concentrations to 2 serum samples and were assayed.
The recoveries were 95.8% - 109%.
- Dilution test
2 serum samples were serially diluted by 5 steps.
The dilution curves showed excellent linearity with $R^2 = 0.9994$ and 0.9991 .
- Correlation between U-type and ordinary type (TMB) : $R = 0.9983$.

5. Reference Assay Data

Serum insulin levels of fasted mice.

Mean assay value : 104 pg/mL, SD : 18.2 pg/mL

Mouse strain: BALB/c, 6 w males, starvation : 48 hours, n = 8

These data should be considered as guidance only. Each laboratory should establish its own normal and pathological reference ranges for insulin levels independently.

6. Precautions

- For professional use only, beginners are advised to use this kit under the guidance of experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Use clean laboratory glassware.
- Avoid contact with the acidic Stop Solution and TMB Solution, containing hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye and clothing protection when handling these reagents.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the Stop Solution because it is sulfuric acid. The Stop Solution and the substrate solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Do not drink, eat or smoke in the areas where assays are carried out.
- In treating assay samples of animal origin, be careful for possible biohazards.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- The materials must not be pipetted by mouth.
- Unused samples and used tips should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a Plate Seal supplied, during incubation.
- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.

7. Reagents supplied

Components	Use Status	Amount
(A) Antibody-coated Plate	Use after washing.	96 wells/1 plate
(B) Insulin Standard	Concentrated. Use after dilution.	500 μ L/1 bottle
(C) Buffer	Ready for use.	60 mL/1 bottle
(D) Biotin-conjugated Antibody Solution	Concentrated. Use after dilution.	10 μ L/1 bottle
(E) Peroxidase-conjugated Streptavidin Solution	Concentrated. Use after dilution.	20 μ L/1 bottle
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 $\times$)	Concentrated. Use after dilution.	100 mL/1 bottle
(J) Plate Seal	—	3 sheets

(D) Biotin-conjugated Antibody Solution, and (E) Peroxidase-conjugated Streptavidin Solution : Vials contain more than volumes shown in the list. You can easily take out 10 μ L and 20 μ L, respectively, from vials.

[Storage and Stability]

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C.

[(B) Insulin Standard (5 ng/mL)]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

[(C) Buffer] & [(F) TMB Solution]

If not opened, store at 2°C - 10°C. Once opened, we recommend using them as soon as possible to avoid influence by environmental condition.

[(D) Biotin-conjugated Antibody Solution] & [(E) Peroxidase-conjugated Streptavidin Solution]

Unused working solution (already diluted) should be disposed.

[(H) Stop Solution]

Close the stopper tightly and store at 2°C - 10°C.

[(I) Wash Solution (10×)]

The rest of undiluted buffer : Stored tightly closed at 2°C - 10°C. Dispose any unused diluted buffer.

8. Equipments or supplies required but not supplied ☐ Use as a check box

☐ Deionized water (or Distilled water). ☐ Test tubes for preparation of standard solution series. ☐ Glassware for dilution of Wash Solution (10×) (a graduated cylinder, a bottle). ☐ Pipettes (disposable tip type). One should be able to deliver 5 µL precisely, and another for 45 µL and 100 µL. ☐ Syringe-type repeating dispenser like Eppendorf multipipette plus which can dispense 50 µL. ☐ Paper towel to remove washing buffer remaining in wells. ☐ A vortex-type mixer. ☐ A shaker for 96 well-plate (600 rpm - 1200 rpm). ☐ An automatic washer for 96 well-plate (if available), or a washing bottle with a jet nozzle. ☐ A 96 well-plate reader (450 nm ± 10 nm, 620 nm : 600 nm - 650 nm). ☐ Software for data analysis.

9. Preparation of Samples

This kit is intended to measure insulin in mouse serum, plasma (preferably obtained with heparin), culture medium and tissue/cell extracts. The necessary sample volume for the standard procedure is 5 µL. Samples should be immediately assayed or stored below -35°C for several days. Defrosted samples should be mixed thoroughly for best results.

Hemolytic and hyperlipemic serum samples are not suitable.

*** To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis, do not use them for assay. Abnormal value might be obtained with hemolysis above 40 mg/dL with this kit.**

If presence of interfering substance is suspected, examine by dilution test at more than 2 points. Dilution of a sample should be made in a test tube using Buffer prior to adding them to wells. Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.

Storage and stability

Insulin samples will be inactivated if stored at 2°C - 10°C. If it is necessary to store sample in refrigerator (2°C - 10°C), add aprotinin at final concentration of 100 KIU/mL - 500 KIU/mL. (KIU : kallikrein inhibitor unit.) If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35°C. Avoid repeated freeze-thaw cycles.

10. Preparation of Reagents

◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.

◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Insulin Standard (5 ng/mL)]

Make a serial dilution of master standard (5 ng/mL) solution to prepare each standard solution (39 pg/mL - 2500 pg/mL).

*** Unit reduction for µIU/mL is 26 IU/mg. (Refer to 2. Introduction.)**

Volume of standard solution	(C) Buffer	Concentration (pg/mL)	Concentration (µIU/mL)
Original solution : 100 µL	100 µL	2500	65
2500 pg/mL solution : 100 µL	100 µL	1250	32.5
1250 pg/mL solution : 100 µL	100 µL	625	16.3
625 pg/mL solution : 100 µL	100 µL	313	8.13
313 pg/mL solution : 100 µL	100 µL	156	4.06
156 pg/mL solution : 100 µL	100 µL	78	2.03
78 pg/mL solution : 100 µL	100 µL	39	1.02
0 (Blank)	100 µL	0	0

[(D) Biotin-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the (C) Buffer to **1 : 4000**.

10 mL of the diluted solution is enough for 96 wells.

[(E) Peroxidase-conjugated Streptavidin Solution]

Prepare working solution by dilution of (E) with the (C) Buffer to **1 : 500**.

10 mL of the diluted solution is enough for 96 wells.

[(I) Wash Solution (10×)]

Dilute 1 volume of the concentrated Wash Solution (10×) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example: 100 mL of concentrated Wash Solution (10×) and 900 mL of deionized water (or distilled water).

11. Assay Procedure

Remove the cover sheet of the Antibody-coated Plate after bringing up to room temperature.

- (1) Wash the (A) Antibody-coated Plate by filling the wells with 300 μ L of washing buffer and discard 4 times (*①), then strike the plate upside-down onto several sheets of paper towel to remove residual buffer in the wells.
- (2) Pipette 45 μ L of (D) Biotin-conjugated Antibody Solution to all wells. Shake the plate gently on a plate shaker (*②).
- (3) Pipette 5 μ L of sample to the designated sample wells.
- (4) Pipette 5 μ L of standard solution to the wells designated for standards.
- (5) Shake the plate gently on a plate shaker (*②).
- (6) Stick a Plate Seal (*③) on the plate and incubate for 2 hours at 20°C - 25°C.
- (7) Discard the reaction mixture and rinse wells as step (1).
- (8) Pipette 50 μ L of (E) Peroxidase-conjugated Streptavidin Solution to all wells, and shake as step (5).
- (9) Stick a Plate Seal (*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
- (10) Discard the reaction mixture, and then wash the plate as step (1).
- (11) Pipette 50 μ L of (F) TMB Solution to wells, and shake as step (5).
- (12) Stick a Plate Seal (*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
- (13) Add 50 μ L of the Stop Solution to all wells and shake as step (5).
- (14) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm) immediately using a plate reader.

*Refer to 15. Summary of Assay Procedure for notes of *①, *② and *③.


12. Technical Tips

- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kit with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2°C - 10°C.
- Time reaction from the pipetting of the reagent to the first well.
- Dilution of the assay sample must be carried out using the Buffer provided in the kit.
- The TMB Solution should be almost clear pale yellow before use. It turns blue during reaction, and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the Antibody-coated Plate is module type of 8 wells $\times$ 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a Plate Seal and place the well plate in an incubator or a styrofoam box in each step of incubation.

13. Calculations

- (1) Prepare a standard curve using semi-logarithmic or two-way logarithmic section paper by plotting absorbance* (Y-axis) against insulin concentration (pg/mL) on X-axis. * Absorbance at 450 nm minus absorbance at 620 nm.
- (2) Using the standard curve, read the insulin concentration of a sample at its absorbance*, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the Buffer.

*We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.

Mouse Insulin assay standard curve (an example)
Absorbance may change due to assay environment.

14. Trouble Shooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not be added.
- 2) Reagents necessary for coloration such as Biotin-conjugated Antibody Solution, Peroxidase-conjugated Streptavidin Solution, or TMB Solution might not be added.
- 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotin-conjugated Antibody Solution or Peroxidase-conjugated Streptavidin Solution.
- 4) Contamination of peroxidase enzyme inhibitor(s).
- 5) Influence of the temperature under which the kits had been stored.
- 6) Excessive hard washing of the well plate.
- 7) Addition of TMB Solution soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- Intense coloration in all wells including blank

Possible explanations:

- 1) Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution.)
- 2) Overdeveloping. Incubation time with TMB Solution should be decreased before addition of Stop Solution.
- 3) Too high incubation temperature. Adjust the temperature to 20°C - 25°C.

- High coefficient of variation (CV)

Possible explanations:

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.

- Q-2 : I found there contains liquid in 96 well-plate when I opened the box. What is it?

A-2 : When we manufacture 96 well-plate, we insert preservation stabilizer in wells.


15. Summary of Assay Procedure ☐ : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

☐ Bring the well-plate and all reagents back to **20°C - 25°C for 2 hours**.

☐ Wash Solution (10×) concentrate must be diluted to **10 times** by deionized water (or distilled water) that returned to 20°C - 25°C.

☐ Insulin Standard dilution example :

☐ (D) Biotin-conjugated Antibody Solution : Dilute to **4000 times** by using (C) Buffer and use.

- ☐ Antibody-coated Plate
- ☐ ↓ Washing 4 times (*①) *⑤
- ☐ Biotin-conjugated Antibody Solution 45 µL
- ☐ ↓ Shaking (*②)
- ☐ Samples/Standards 5 µL
- ☐ ↓ Shaking (*②), Incubation for 2 hours at 20°C - 25°C. (Standing (*③))

(E) Peroxidase-conjugated Streptavidin Solution

☐ Dilute to **500 times** by using (C) Buffer and use.

(Dilute reagents during the first reaction.)

- ☐ ↓ Washing 4 times (*①) *⑤
- ☐ Peroxidase-conjugated Streptavidin Solution 50 µL
- ☐ ↓ Shaking (*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (*③))
- ☐ ↓ Washing 4 times (*①) *⑤
- ☐ TMB Solution 50 µL
(After dispensing, the color turns to blue depending on the concentration.)
- ☐ ↓ Shaking (*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (*③))
- ☐ Stop Solution 50 µL
(After dispensing, the color turns to yellow depending on the concentration.)
- ☐ ↓ Shaking (*②) (Immediately shake.)
- ☐ Measurement of absorbance (450 nm, Ref 620 nm) immediately *④
(Ref. wave cancels the dirt in the back of plate.)

*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume : 300 µL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the background tends to be high. If so, change washing frequency from 4 times to 5 - 8 times at the constant stroke after the reaction with Peroxidase-conjugated Streptavidin Solution. Standard of plate-washing pressure : 5 - 25 mL/min. (Adjust it depending on the nozzle's diameter.)

*② Guideline of shaking : 600 rpm - 1200 rpm 3 times for 10 seconds.

*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the Plate Seal used once.


*④ 600 nm - 650 nm can be used as reference wavelength.

*⑤ After removal of wash buffer, immediately dispense the next reagent.

Worksheet Example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	2500 pg/mL	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
B	1250 pg/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
C	625 pg/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
D	313 pg/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
E	156 pg/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
F	78 pg/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
G	39 pg/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
H	0 (Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40

Assay Worksheet

16. Storage and Expiration

The complete kit is stored at 2°C - 10°C (Do not freeze). Opened reagents should be used as soon as possible to avoid loss in optimal assay performance caused by storage environment.

LBIS™ Mouse Insulin ELISA Kit (High Sensitivity)

[Storage] Store at 2°C - 10°C (Do not freeze).

[Expiration date] Indicated on the label

[Package] For 96 tests

[Cat #] 296-89801

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://flwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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