Code No. 297-73101 (for 250 assays)

Protein Assay BCA Kit

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This product is used to determine the concentration of solubilized proteins, based on the method of Bicinchoninic acid (BCA). This method relies on the reduction of Cu^{2+} to Cu^+ by protein in an alkaline solution. BCA binds Cu^+ in an alkaline solution, and the complex exhibits absorption at 562 nm. The absorbance is linear with the concentration of solubilized proteins.

This product is for laboratory research use only ; use in any such application is the responsibility of the user.

Components

- Reagent A (BCA Solution)
- Reagent B (Copper (II) Sulfate Solution)

Preparation

Working solution

Mix Reagent A and Reagent B at ratio of 50 : 1 (preparation before use)

Standard solution

Dissolve Bovine Serum Albumin into phosphate buffer containing 0.02w/v% Sodium Azide. Dilute Albumin Solution with the same solvent as the one used for sample solution.

Assay Procedure with Test tube

- 1) Apply 0.1 mL of standard solution or unknown sample into test tubes.
- 2) Add 2 mL Working solution to each tube and mix.
- 3) Incubate tubes under the following conditions.

	Reactive conditions	Working range
Standard protocol	37 °C for 30 minutes	20~2,000 $\mu\text{g/mL}$
Room temperature protocol	Room temperature for 2 hours	20~2,000 $\mu\text{g/mL}$
Enhanced protocol	60 °C for 30 minutes	5~250 $\mu\text{g/mL}$

- 4) Cool tubes to room temperature.
- 5) Measure absorbance at 562 nm.

Assay Procedure with 96 well microplate

- 1) Apply 25 μL of standard solution or unknown sample into each well.
- 2) Add 200 μL Working solution to each well and mix.
- 3) Incubate microplate at 37 °C for 30 minutes.
- 4) Cool microplate to room temperature.
- 5) Measure absorbance at 562 nm with microplate reader.

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Storage : room temperature

Package :

Code No.	297-73101	
Package	for 250 assays	
Components	Reagent A	500 mL × 1 bottle
	Reagent B	25 mL × 1 bottle

* : Trans is the Japanese translation.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

FUJIFILM Wako Chemicals Europe GmbH

Fuggerstrasse 12
D-41468 Neuss
Germany
Telephone : + 49-2131-31110
Facsimile : + 49-2131-311100
<http://www.wako-chemicals.de>

2306KA2