Code No. 297-74201

For genetic research
MagCapture™ microRNA Isolation Kit, Mouse Ago2

Table of Contents	Page
Introduction	p.1
Outline of kit	p.1
Features	p.1
Kit contents	p.2
Storage	p.2
Outline of kit procedures	p.2
Additional materials required	p.2
Precautions for use	p.3
Procedure	p.3
Example of use	p.7
Patent	p.7
Related products	p.7

[Introduction]

MicroRNAs are endogenous non-coding RNAs of approximately 22-nucleotides, which play important roles in post-transcriptional regulation of gene expression by base-pairing to their target mRNAs. There are 1,000 or more kinds of microRNAs in humans and mice, and many researchers are identifying unknown microRNAs and investigating the function of microRNAs related to various disorders worldwide. Recently, microRNAs are reported to exist in not only intracellular but also extracellular body fluid such as blood and are noted as a candidate molecule for the clinical markers of disorders such as cancer. After being transcribed and processed, mature microRNAs are incorporated into Argonaute subfamily proteins of the core component of the RNA-induced silencing complex (RISC) for targeting mRNAs based on sequence complementation in 3'UTRs. In mouse, the Argonaute subfamily consists of four members (Ago1, Ago2, Ago3 and Ago4). Among them, Ago2 is the most abundant in a variety of cells and plays a major role in the microRNA pathway.

[Outline of kit]

MagCapture™ microRNA Isolation Kit, Mouse Ago2 can easily isolate Ago2-binding microRNAs and target mRNAs from mouse cell or tissue samples, by immunoprecipitation using an anti-mouse Ago2 monoclonal antibody. This kit also enables isolation of Ago2-binding microRNAs from biological fluids. (e.g., plasma, serum, and conditioned cell culture medium.) This kit also can isolate Ago2-binding microRNAs and target mRNAs from rat cell or tissue samples.

[Features]

- 1) Almost all reagents required for immunoprecipitation and purification of Ago2-binding microRNAs and target mRNAs are contained in MagCapture™ microRNA Isolation Kit, Mouse Ago2.
- 2) Because anti-mouse Ago2 antibody immobilized magnetic beads are adopted, handling of immunoprecipitation is improved from the preexisting kit (microRNA Isolation Kit, Mouse Ago2, Code No. 292-67301).
- 3) Two protocols for the RNA purification step are selectable - a one-step method for obtaining high-yield microRNAs by a simple procedure or a spin-column method for obtaining high-purity microRNAs.
- 4) Optimized protocols are configured for isolating microRNAs from cell, tissue, and biological fluids sample.
- 5) Since isolated RNAs using this kit contain Ago2-binding microRNAs and target mRNAs, the RNAs can be used in a variety of downstream applications, including identification or expression profiling of microRNAs and target mRNAs.

- 1/16 -


[Kit contents] (10 Reactions)

	Content
(1) Anti-Mouse Ago2 Magnetic Beads Solution	600 μ L \times 1 tube
(2) Cell Lysis Solution	20 mL \times 1 tube
(3) Washing Solution for IP	40 mL \times 1 tube
(4) Elution Solution I for IP	500 μ L \times 1 tube
(5) Elution Solution II for IP	500 μ L \times 1 tube
(6) Binding Solution for Spin Column	2 mL \times 1 tube
(7) Binding Enhancer for Spin Column	100 μ L \times 1 tube
(8) Washing Solution I for Spin Column	3 mL \times 1 tube
(9) Washing Solution II for Spin Column	4 mL \times 1 tube
(10) Elution Solution for Spin Column	1 mL \times 1 tube
(11) Spin Column/Collection Tube	10 tubes

[Storage]

Store at 2~10 °C.

[Outline of kit procedures]

[Additional materials required]

1. Equipment
 - 1) High speed micro centrifuge (Max > 20,000 g)
 - 2) Vortex mixer
 - 3) Tabletop centrifuge
 - 4) Magnetic stand
 - 5) Tube rotator
 - 6) Tube mixer
 - 7) Block heater
 - 8) Micro centrifugation tube (1.5 mL)
 - 9) Micro pipette
 - 10) Pipette tip
2. Reagents
 - 1) Isopropanol (Code No. 166-21671)
 - 2) Ethanol (Code No. 052-07221)

- 2/16 -

[Precautions for use]

1. Equipment

PCR tubes, microcentrifugation tubes, and pipette tips must be sterilized by autoclave. There is no problem with using commercially available RNase free products. We recommend the use gloves and a mask to avoid contamination of RNase.

2. Reagents

Additional reagents excluding kit must be sterilized by autoclave or filtration.

3. Handling of biohazardous waste

Please deal with human serum, human plasma, and the human tissue samples as a susceptibility infectious sample. This includes waste liquid and equipment such as PCR tubes, microcentrifugation tubes, pipette tips and gloves after experiment. These biohazardous wastes are must be disposed of according to the guidelines of the institution which belongs.

[Procedure]

1. Preparation of sample

1-A In the case of adherent cells

- 1) Culture the cells ($5 \times 10^6 \sim 2 \times 10^7$ cells).^{※1}
- 2) Remove the medium from culture dishes and wash cultured cells twice with PBS (-).
- 3) Dissociate the cells using a cell scraper or trypsin-EDTA solution from dishes and transfer the cell suspension into a new 1.5 mL centrifugation tube. In the case of trypsin-EDTA solution, add the medium including serum to the cell suspension immediately for suppressing damage by trypsin.
- 4) Centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 5) Remove the supernatant.
- 6) After adding 1 mL of PBS (-) and suspending the cell pellet, transfer the cell suspension into a new 1.5 mL centrifugation tube.
- 7) Centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 8) Remove the supernatant.
- 9) After adding 1 mL of PBS (-) and resuspending the cell pellet, centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 10) Remove the supernatant.
- 11) Add 1 mL of Cell Lysis Solution to the cell pellet ($5 \times 10^6 \sim 2 \times 10^7$ cells) and suspend cells by pipetting.^{※2}
- 12) Incubate on ice for 10 minutes.
- 13) Centrifuge the cell lysate at $20,000 \times g$ for 20 minutes at 4 °C.^{※3}
- 14) Transfer the supernatant into a new 1.5 mL micro centrifugation tube. ... [Cell lysate]

1-B In the case of suspension cultured cells

- 1) Culture the cells ($5 \times 10^6 \sim 2 \times 10^7$ cells).^{※1}
- 2) Centrifuge the culture solution at $200 \times g$ for 5 minutes at 4 °C.
- 3) Remove the supernatant.
- 4) After adding 1 mL of PBS (-) and suspending the cell pellet, transfer the cell suspension into a new 1.5 mL centrifugation tube.
- 5) Centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 6) Remove the supernatant.
- 7) After adding 1 mL of PBS (-) and re-suspending the cell pellet, centrifuge the cell suspension at $200 \times g$ for 5 minutes at 4 °C.
- 8) Remove the supernatant.
- 9) Add of 1 mL of Cell Lysis Solution to the cell pellet ($5 \times 10^6 \sim 2 \times 10^7$ cells) and suspend cells by pipetting.^{※2}
- 10) Incubate on ice for 10 minutes.
- 11) Centrifuge the cell lysate at $20,000 \times g$ for 20 minutes at 4 °C.^{※3}
- 12) Transfer the supernatant into a new 1.5 mL micro centrifugation tube. ... [Cell lysate]

1-C In the case of tissues

- 1) After adding 1mL of Cell Lysis Solution to the piece of frozen tissue (25~50 mg), homogenize the tissue by Teflon-lined homogenizer etc.^{※4}
- 2) Transfer the homogenized suspension into a new 1.5 mL centrifugation tube, and incubate it on ice for 10 minutes.
- 3) Centrifuge the cell lysate at $20,000 \times g$ for 20 minutes at 4 °C.^{※3}
- 4) Transfer the supernatant into a new 1.5 mL micro centrifugation tube. ... [Cell lysate]

1-D In the case of plasma or serum

- 1) Transfer the human plasma or serum (about 300 μ L) into a new 1.5mL centrifugation tube.
- 2) Centrifuge at $20,000 \times g$ for 10 minutes at 4 °C.
- 3) Transfer the 200 μ L of supernatant into a new 1.5 mL micro centrifugation tube. ... [Plasma or Serum sample]

1-E In the case of conditioned culture medium

- 1) Culture the cells.^{※5}
- 2) Transfer the 1 mL of conditioned culture medium into a new 1.5 mL centrifugation tube and centrifuge at $200 \times g$ for 5 minutes at 4 °C.
- 3) Transfer the supernatant into a new 1.5 mL centrifugation tube and centrifuge at $1,200 \times g$ for 10 minutes at 4 °C.
- 4) Transfer the supernatant into a new 1.5 mL centrifugation tube and centrifuge at $20,000 \times g$ for 10 minutes at 4 °C.
- 5) Transfer the supernatant into a new 1.5 mL centrifugation tube.^{※6} ... [Conditioned culture medium sample]

- ※1 Please examine appropriate culture conditions according to cell lines.
- ※2 Maximum cell numbers are about $1 \times 10^7 \sim 2 \times 10^7$.
- ※3 If a white floating layer appears in supernatant after centrifugation, please filtrate it using a centrifuge tube filter (pore size : 0.45 μ m).
- ※4 Please examine appropriate homogenizing condition according to the tissue. Maximum tissue weight is about 50 mg.
- ※5 Please examine appropriate culture conditions according to the cell lines. In order to avoid the contamination of microRNA in bovine serum, we recommend cultivation by using a serum-free medium.
- ※6 If there are few amounts of microRNAs in a culture supernatant, we recommend to condense a culture supernatant using a centrifugal type ultrafiltration filter (molecular weight 30000 cut) etc.

2. Prewashing of Anti-Human Ago2 Magnetic Beads

- 1) Transfer 60 μ L of Anti-Human Ago2 Magnetic Beads Solution into a new 1.5 mL centrifugation tube.^{※7}
- 2) Place this tube on the magnetic separation stand for 30 seconds to 1 minute.
- 3) After confirming the separation of beads and solution, discard the supernatant using a micro pipette.
- 4) Add 1 mL of Washing Solution for IP to magnetic beads and suspend it by vortexing.
- 5) After spin down of the solution, place this tube on the magnetic separation stand for 30 seconds to 1 minute.
- 6) After confirming the separation of beads and solution, discard the supernatant using a micro pipette. ... [Prewashed magnetic beads]

- ※7 Please use Anti-Human Ago2 Magnetic Beads Solution after completely suspending.

3. Immunoprecipitation

3-A In the case of cells or tissues

- 1) Add 1 mL of [Cell lysate] to [Prewashed magnetic beads] and suspend by vortexing.
- 2) Incubate with rotation for 2 hour or more at 4 °C.
- 3) After spin down of the solution, place this tube on the magnetic separation stand for 30 seconds to 1 minute.
- 4) After confirming the separation of beads and solution, discard the supernatant using a micro pipette.
- 5) Add 1 mL of Washing Solution for IP to beads and suspend by vortexing.
- 6) After spin down of the solution, place this tube on the magnetic separation stand for 30 seconds to 1 minute.
- 7) After confirming the separation of beads and solution, discard the supernatant using a micro pipette.
- 8) Repeat 5)-7) steps twice. ... [Immunoprecipitated magnetic beads]

3-B In the case of plasma or serum or conditioned culture medium

- 1) Add 200 μ L of [Plasma or Serum sample] or [Conditioned culture medium sample] to [Prewashed magnetic beads] and suspend by vortex mixer.
- 2) Incubate with rotation by tube mixer for 2 hours or more at 4 °C.※8
- 3) After spin down of the solution, place this tube on the magnetic separation stand for 30 seconds to 1 minute.
- 4) After confirming the separation of beads and solution, discard the supernatant using a micro pipette.
- 5) Add 1 mL of Washing Solution for IP to beads and suspend by vortexing.
- 6) After spin down of the solution, place this tube on the magnetic separation stand for 30 seconds to 1 minute.
- 7) After confirming the separation of beads and solution, discard the supernatant using a micro pipette.
- 8) Repeat 5)-7) steps twice. ... [Immunoprecipitated magnetic beads]

※8 Please adjust the appropriate rotating speed of a tube mixer as checking that beads are well suspended. (Indication : 1,200~1,400 rpm) Moreover, although the tube rotator can also be used, the yield of microRNA decreases compared with tube mixer.

4. Purification of microRNA from immunoprecipitated magnetic beads

4-A Spin-column method

This is a protocol for purifying the high purity RNA from immunoprecipitated magnetic beads.

Please prepare the following reagents in advance.

● [Washing Solution I for Spin Column (+ Isopropanol)]

Add 4.5 mL of isopropanol to Washing Solution I for Spin Column and mix it well.※9

● [Washing Solution II for Spin Column (+ Ethanol)]

Add 12 mL of ethanol to Washing Solution II for Spin Column and mix it well.※9

- 1) Add 50 μ L of Elution Solution I for IP to [Immunoprecipitated magnetic beads] and suspend by vortex mixer.※10
- 2) After spin down of the solution, place this tube on the magnetic separation stand for 30 seconds to 1 minute.
- 3) After confirming the separation of beads and solution, transfer the supernatant into a new 1.5 mL micro centrifugation tube.
- 4) Add 200 μ L of Binding Solution for Spin Column and 10 μ L of Binding Enhancer for Spin Column to the supernatant and mix them by vortexing.
- 5) Add 350 μ L of isopropanol to it and mix them by vortexing.
- 6) Add 50 μ L of Elution Solution for Spin Column into Spin Column/Collection Tube and centrifuge it at 8,000 \times g for 1 minute at room temperature.※11... [Pretreated column]

- 7) After spin down of the mixture of step 5) for 2~3 seconds※12 and stir this solution by pipetting,※13 apply it all on [Pretreated column].
- 8) Centrifuge at 8,000 \times g for 1 minute at room temperature.
- 9) Remove the Spin Column and discard the flow-through from the Collection Tube.
- 10) Place the Spin Column on to the Collection Tube.
- 11) Add 500 μ L of [Washing Solution I for Spin Column (+ isopropanol)] to the Spin Column.
- 12) Centrifuge at 8,000 \times g for 1 minute at room temperature.
- 13) Remove the Spin Column and discard the flow-through from the Collection Tube.
- 14) Place the Spin Column onto the Collection Tube.
- 15) Add 700 μ L of [Washing Solution II for Spin Column (+ ethanol)] to the Spin Column.
- 16) Centrifuge at 8,000 \times g for 1 minute at room temperature.
- 17) Remove the Spin Column and discard the flow-through from the Collection Tube.
- 18) Place the Spin Column onto the Collection Tube.
- 19) Repeat 15)-18) steps.
- 20) Centrifuge at 8,000 \times g for 3 minutes at room temperature.
- 21) Place the Spin Column onto a new 1.5 mL micro centrifugation tube.
- 22) Add carefully 50 μ L of Elution Solution for Spin Column to the center of membrane.
- 23) Incubate for 5 minutes at room temperature.
- 24) Centrifuge at 8,000 \times g for 1 minute at room temperature.
- 25) Please use the eluate in the 1.5mL micro centrifugation tube as the purified RNA sample.※14, ※15

※9 Store the prepared solution at 2~8 °C.

※10 Please use the Elution Solution I for IP after confirming that the components are entirely dissolved, because components are crystallized at a cool temperature.

※11 Please be sure to perform this operation because the pre-treatment of Spin Column elevate recovery amount of RNA.

※12 Please avoid a prolonged spin down.

※13 Because there is a possibility that RNA is precipitated, please be sure to mix by pipetting.

※14 Store the RNA sample at -80 °C.

※15 Because the amount of microRNAs obtained from cell or tissue samples is trace, it cannot be measured by spectrophotometer (A260/A280) or bioanalyzer. We recommend the detection of microRNAs by silver staining after Urea-PAGE or by quantitative PCR. Moreover, because the amount of microRNAs obtained from plasma or serum or conditioned cell culture medium is more trace, we recommend the detection of microRNAs by quantitative PCR. Because there is the difference of the amount of microRNA isolated from each kind of sample, we recommend the confirmation of the amount of microRNA using above mentioned methods after isolation.

4-B One-step method

This is a protocol for obtaining a high-yield RNA from immunoprecipitated magnetic beads more easily.


- 1) Add 50 μ L Elution Solution II for IP to [Immunoprecipitated magnetic beads] and suspend by vortexing.
- 2) Incubate for 90 seconds at 90 °C by using block heater immediately before beads are precipitated.※16
- 3) Place the incubated tube on ice for 1 minute or more.
- 4) After spin down of the solution, place this tube on the magnetic separation stand for 30 seconds to 1 minute.
- 5) After confirming the separation of beads and solution, transfer the supernatant into a new 1.5 mL micro centrifugation tube as the purified RNA sample.※14, ※15, ※17

- ※16 In the case of many specimen, we recommend the usage of Thermo Mixer.
- ※17 Although protein components are contained in the purified RNA sample, we have checked that it was satisfactory for quantitative PCR analysis.

【Example of use】

Purification of microRNAs from mouse P388D1 cells

Purification of microRNAs from mouse P388D1 cells (2×10^7 cells) was conducted by using MagCapture™ microRNA Isolation Kit, Mouse Ago2. The microRNA fraction purified by spin-column method or one-step method was detected by silver staining after Urea-PAGE.

Lane M: Molecular weight marker
 Lane 1 : Single strand RNA (23 base) 1ng
 Lane 2 : Purified RNA by Spin-column method
 Lane 3 : Purified RNA by One-step method

Figure : Detection of purified microRNAs from mouse P388D1 cells by silver staining

【Patent】

One-step method of the RNA purification is patent pending.

【Related Products】

Code No.	Product Name	Package
014-22023	Anti Mouse Ago2, Monoclonal Antibody	50 μ L
166-21671	2-Propanol	100 mL
052-07221	Ethanol (99.5%)	100 mL
316-90101	Distilled Water, Deionized, Sterile	100 mL
202-16931	0.05w/v% Trypsin-0.53 mmol/L EDTA ·4Na Solution with Phenol Red	100 mL
209-16941	0.25w/v% Trypsin-1 mmol/L EDTA ·4Na Solution with Phenol Red	100 mL
318-90041	5 \times TBE	1000 mL
194-15881	SuperSep™ RNA, 15%, 17well	5 gels
182-02711	R-Mark™ Small RNA Marker (20-50nt)	50 μ L

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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