## Neuron Dissociation Solutions S

Neuron Dissociation Solutions S are appropriate for dissociating and isolating neurons from rat or mouse central nervous system tissue. This product consists of three kinds of solutions (Enzyme Solution, Dispersion Solution, and Isolation Solution). You can easily isolate neurons in high cell viability with this product.

### [Kit components]

Component	Size
Enzyme Solution	2.5 mL × 10
Dispersion Solution	2.5 mL × 10
Isolation Solution	2.5 mL × 10

### [Reagent volume]

Experiment condition	Component	Volume
In the case of dissociating many tissues ("many tissues" means about 1~3 cerebral hemispheres from fetal mouse or rat)	Enzyme Solution	5.0 mL
	Dispersion Solution	3.0 mL
	Isolation Solution	5.0 mL
In the case of dissociating a few tissues ("few tissues" means about 1~8 hippocampuses from fetal mouse or rat)	Enzyme Solution	2.5 mL
	Dispersion Solution	2.5 mL
	Isolation Solution	2.5 mL

### [Procedure]

#### 〈Thawing and Washing〉

In case of using dissecting tissue

1. Thaw Neuron Dissociation Solutions S at room temperature or in a 37°C water bath and return it to a room temperature. If you need to excise the tissues in a short time, you can keep this product in a 37°C water bath.
2. Cut the tissues finely after removing meninges and transfer it to a 15 mL centrifuge tube.  
(We recommend that you cut the tissues in a Leibovitz's L-15 Medium or HBSS (+).)  
("cut finely" means that the section size should be about 1 ~ 4 mm<sup>2</sup>, but it is not necessary to rigidly observe this size.)


In case of using freezing nerve-cells (refer to Related Reagents)

1. Thaw Neuron Dissociation Solutions S at room temperature or in a 37°C water bath and return it to a room temperature.
2. Remove the cryotube from the liquid nitrogen storage container.
3. Thaw at room temperature. If the thawing time is shorter than 30 min, neurons can't thaw enough.
4. Remove the cryopreservation medium from the cryotube about 0.8 mL.

5. Add 1 mL of Washing Solution (Leibovitz's L-15 Medium or HBSS (+)) to the cryotube, then transfer to a transparent 15 mL centrifuge tube with the thawed tissues.
6. Add extra 4 mL of Washing Solution to the centrifuge tube.
7. Leave at room temperature for 5 min, and centrifuge at 300 rpm at room temperature for 1 min.
8. Remove the supernatant, and add 5 mL of Washing Solution.
9. Leave at room temperature for 5 min, and centrifuge at 300 rpm at room temperature for 1 min.

### 〈Dispersion〉

1. Remove the supernatant after the tissues sinking, and add Enzyme Solution, and incubate at 37°C for 30 min in a water bath. (In case of freezing neurons, the incubation time is 20 min.)
2. Suspend the tissues by pipetting in all Enzyme Solution, and disperse the tissues. The appearance remains unchanged in 30 min, but the enzyme reaction proceeds well. Suspend the tissues by pipetting repeatedly without bubbles. Suspend the tissues by pipetting again in procedure 4 after this procedure.
3. Centrifuge at 900~1,000 rpm at room temperature for 4 min.
4. Remove the supernatant, and add Dispersion Solution. Then suspend the tissues.  
In the case of dissociating a few tissues : Dispersion Solution 2.5 mL  
In the case of dissociating many tissues : Dispersion Solution 3.0 mL
5. Add all Isolation Solution to the centrifuge tube bottom quietly. This makes two separate solution layers. Turbulence of liquid surface does not affect the dissociation. (upper layer : Dispersion Solution, lower layer : Isolation Solution) (refer to the figure on the right)
6. Spin at 800 ~ 900 rpm at room temperature for 5 min.
7. Centrifuge the supernatant, and suspend the tissues with Neuron Culture Medium (Code No. 148-09671), etc., and culture onto the Poly-Lysine coated culture substratum. Add Ara-C (Day 3-5 : 5-10 μmol/L) when grial cells grow well.

### [Storage]

- Store at -80°C.
- Avoid repeating freeze-thaw.
- After opening, use as soon as possible.
- Avoid storing in environments with carbon dioxide such as dry ice.
- After opening, it may be affected by carbon dioxide in the storage environment. It is recommended to block the outside air with a clip with excellent airtightness.

【Related Reagents】

Code No.	Product Name
291-78001	Neuron Dissociation Solutions
148-09671	Neuron Culture Medium
128-06075	Leibovitz's L-15 Medium with L-Glutamine, Phenol Red and Sodium Pyruvate
082-09365	HBSS(+) with Phenol Red
082-10291	Hippocampus, from Mouse (embryonic day 16)
033-24871	Cerebral Cortex, from Rat (embryonic day 17)
036-24861	Cerebral Striatum, from Rat (embryonic day 17)
085-10301	Hippocampus, from Rat (embryonic day 19)
030-24881	Cerebral Cortex, from Mouse (embryonic day 15)

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : +1-804-271-7677  
Facsimile : +1-804-271-7791  
<http://www.wakousa.com>

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