LBIS™ Rat KLH-IgM ELISA Kit (TDAR)

Please, read this instruction carefully before use.

1. Intended use

LBIS™ Rat KLH-IgM ELISA Kit (TDAR) is a sandwich ELISA system for measurement of IgM-type anti-KLH (Keyhole limpet hemocyanin) antibody after inoculation of KLH to rat with high sensitivity. This is helpful in examining TDAR (T-cell dependent antibody reaction) in combination with LBIS™ Rat KLH-IgG ELISA Kit (TDAR). This is intended for research use only.

2. Introduction

ICH (International Conference on Harmonisation of Technical Requirements of Registration of Pharmaceuticals for Human Use) issues a guideline "Immunotoxicity Studies for Human Pharmaceuticals S8" in 2005. In the guideline, TDAR (T cell Dependent Antibody Reaction) is recommended in case target of a pharmaceutical's immunotoxicity is not identified. It says "The TDAR should be performed using a recognized T-cell dependent antigen (e.g., sheep red blood cells (SRBC) or keyhole limpet hemocyanin (KLH) that results in a robust antibody response". In this study, the production of IgM-type antibody caused by the primary response to e.g., KLH, and IgG-type antibody production by "class-switch" following the secondary response are tested. This kit makes it possible to measure IgM-type anti-KLH antibody in rat blood samples, and most suitable for TDAR test when used in combination with LBIS™ Rat KLH-IgG ELISA kit (TDAR).

3. Assay Principle

In LBIS™ Rat KLH-IgM ELISA Kit (TDAR), standards or samples are incubated in KLH coated wells to capture anti-KLH antibody. After 1 hour incubation and washing, HRP (horse radish peroxidase)-conjugated anti-rat IgM antibody is added and incubated for 1 hour together with captured anti-KLH-IgM. After washing, HRP-complex remaining in wells is reacted with a chromogen (TMB) for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to anti-KLH-IgM concentration. The standard curve is prepared by plotting absorbance against the standard concentrations. Anti-KLH-IgM concentrations in unknown samples are determined using this standard curve.

4. Performance Characteristics

- Assay Range
The assay range of the kit is 3.13 ng/mL - 200 ng/mL.
- Specificity
The Peroxidase-conjugated Antibody of this kit is specific to anti-rat IgM.
The cross-reactivity with anti-rat IgG is less than the detection limit.
- Precision of Assay
Within assay variation (2 samples, 5 replicates assay), the mean CV was within 10%.
- Reproducibility
Between assay variation (3 samples, 4 days, 4 replicates assay), the mean CV was within 10%.
- Recovery Test
Standard anti-KLH rat IgM was added in 3 concentrations to 2 serum samples and were assayed. The recoveries were 94.8% - 102%.
- Dilution Test
2 serum samples were serially diluted by 3 steps.
The dilution curves showed linearity with $R^2 = 0.998 - 0.999$.

5. Precautions

- For professional use only, beginners are advised to use this kit under the guidance of experienced person.
- Wear gloves and laboratory coats when handling assay materials.
- Do not drink, eat or smoke in the areas where assays are carried out.
- In treating assay samples of animal origin, be careful for possible biohazards.

- Wear gloves and goggles and clothing-protection when handling the Stop Solution (sulfuric acid) and the chromogen (hydrogen peroxide and tetramethylbenzidine). Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucous membrane. The Stop Solution and TMB Solution may cause skin/eyes irritation. In case of contact with these, wash the place thoroughly with enough water and seek medical attention if necessary.
- The materials must not be pipetted by mouth.
- Residual samples and used tips should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- Use clean laboratory glassware.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a Plate Seal supplied, during incubation.
- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.

6. Reagents supplied

Components	Use Status	Amount
(A) Antigen-coated Plate	Use after washing	96 wells/1 plate
(B) Anti-KLH Rat IgM Standard	Concentrated. Use after dilution	200 μ L/1 bottle
(C) Buffer	Ready for use.	100 mL/1 bottle
(D) Peroxidase-conjugated Antibody Solution	Concentrated. Use after dilution.	100 μ L/1 bottle
(F) TMB Solution	Ready for use.	12 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle
(I) Wash Solution (10 \times)	Concentrated. Use after dilution.	100 mL/1 bottle
Plate Seal	—	3 sheets

[Storage and Stability]

[(A) Antigen-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C.

[(B) Anti-KLH Rat IgM Standard (1,000 ng/mL)]

Standard solutions prepared above should be used as soon as possible and should not be stored.

[(C) Buffer] and [(F) TMB Solution]

If not opened, store at 2°C - 10°C. Once opened, we recommend using them as soon as possible to avoid influence by environmental condition.

[(D) Peroxidase-conjugated Antibody Solution]

Unused working solution (already diluted) should be disposed. The rest of the undiluted solution stored tightly closed at 2°C - 10°C.

[(H) Stop Solution]

Close the stopper tightly and store at 2°C - 10°C.

[(I) Wash Solution (10 \times)]

The rest of undiluted solution (unused) : close the cap tightly and store at 2°C - 10°C. Dispose any remaining diluted buffer.

7. Equipments required but not supplied ☐ Use as a check box

- ☐ Deionized water (or Distilled water).
- ☐ Test tubes for preparation of standard solution series.
- ☐ Glassware for dilution of washing buffer (a graduated cylinder, a bottle)
- ☐ Pipettes (disposable tip type). One should be able to deliver 10 μ L - 20 μ L precisely, and another for 50 μ L - 500 μ L.
- ☐ Syringe-type repeating dispenser like Eppendorf multipipette plus which can dispense 50 μ L.
- ☐ Paper towel to remove washing buffer remaining in wells.
- ☐ A vortex-type mixer.
- ☐ A shaker for 96 well-plate (600 rpm - 1,200 rpm)
- ☐ An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle.
- ☐ A 96 well-plate reader (450 nm \pm 10 nm, 620 nm : 600 nm - 650 nm)
- ☐ Software for data analysis.

8. Preparation of Samples

This kit is intended to measure IgM-type anti-KLH (Keyhole limpet hemocyanin) antibody in rat serum or plasma. The necessary sample volume for the standard procedure is 10 μ L. Samples should be immediately assayed or stored below -35°C for several days. Defrosted samples should be mixed thoroughly for best results. Don't repeat freeze and thaw.

Hemolytic and hyperlipemic samples are not suitable.

***To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis, do not use them for assay. Abnormal value might be obtained with hemolysis above 160 mg/dL with this kit.**

If presence of interfering substance is suspected, examine by dilution test at more than 2 points. Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.

Make sure to dilute samples more than 1 : 200 to avoid any nonspecific reaction. Recommended dilution rate is 1 : 200 - 20,000 depending on the antibody titer. Dilution should be carried out with the buffer solution of the kit using small test tubes before assay so as to be within the standard curve range. Dilution rate should be different depending on the immune or sampling conditions.

Example of dilution:	Rate	(20 \times)	200 \times	2,000 \times	20,000 \times
	Sample (μ L)	10	20*	20*	20*
	Buffer (μ L)	190	180	180	180

*One rank lower diluted sample

Storage and Stability

Sample is stable at 2°C - 10°C within a week. If you have to store assay samples for a longer period, snap-freeze samples and keep them below -35°C. Avoid repeated freezing and thawing cycles.

9. Preparation of Reagents

- ◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.
- ◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Anti-KLH Rat IgM Standard (1,000 ng/mL)]

Make a serial dilution of original standard solution to prepare each standard solution. Example is shown below.

Volume of standard solution	(C) Buffer	Concentration (ng/mL)
Original solution : 100 μ L	400 μ L	200
200 ng/mL solution : 200 μ L	200 μ L	100
100 ng/mL solution : 200 μ L	200 μ L	50
50 ng/mL solution : 200 μ L	200 μ L	25
25 ng/mL solution : 200 μ L	200 μ L	12.5
12.5 ng/mL solution : 200 μ L	200 μ L	6.25
6.25 ng/mL solution : 200 μ L	200 μ L	3.13
0 (Blank)	200 μ L	0

[(D) Peroxidase-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the (C) Buffer to **1 : 100**.

[(I) Wash Solution (10 \times)]

Dilute 1 volume of the concentrated (I) Wash Solution (10 \times) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of Wash Solution (10 \times) and 900 mL of deionized water (or distilled water).

10. Assay Procedure

Remove the cover sheet of the 96 well-plate after bringing up to room temperature.

- (1) Wash the (A) Antigen-coated Plate by filling the well with washing buffer and discard 3 times (*①), then strike the plate upside-down onto several layers of paper towels to remove residual buffer in the wells.
- (2) Pipette 50 μ L of standards or diluted samples to the wells designated for each.
- (3) Shake the plate gently on a plate shaker (*②).
- (4) Stick a Plate Seal on the plate and incubate for 1 hour at 20°C - 25°C (*③).
- (5) Discard the reaction mixture and rinse wells as step (1) (*①).
- (6) Pipette 50 μ L of (D) Peroxidase-conjugated Antibody Solution to all wells, and shake as step (3) (*②).
- (7) Stick a Plate Seal on the plate and incubate the plate for 1 hour at 20°C - 25°C (*③).
- (8) Discard the reaction mixture. Rinse wells as step (1) (*①).
- (9) Pipette 50 μ L of (F) TMB Solution to wells, and shake as step (3) (*②).
- (10) Stick a Plate Seal on the plate and incubate the plate for 20 minutes at 20°C - 25°C (*③).
- (11) Add 50 μ L of the (H) Stop Solution to all wells and shake as step (3) (*②).
- (12) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm : 600 nm - 650 nm can be used as reference wavelength) immediately using a plate reader.

*Refer to 14. Summary of Assay Procedure for notes of *①, *② and *③.


11. Technical Tips

- In manual operation, proficiency in pipetting technique is recommended.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2°C - 10°C.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- The TMB Solution should be almost colorless before use. It turns blue during reaction, and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated KLH, do not let the plate go dry.
- As the Antigen-coated Plate is module type of 8 wells \times 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a Plate Seal and place the well plate in an incubator or a styrofoam box in each step of incubation.

12. Calculations

- (1) Prepare a standard curve using a two-way logarithmic section paper by plotting absorbance* (Y-axis) against Anti-KLH-IgM concentration (ng/mL) on X-axis.
- (2) Using the standard curve, read the Anti-KLH-IgM concentration of a sample at its absorbance*, and multiply the assay value by dilution factor. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the buffer solution.

*We recommend the use of 3rd order regression curve for log-log plot, or 4 parameters or 5 parameters method for log-normal plot in computer calculation.

Anti KLH Rat IgM assay standard curve (an example)

Absorbance may change due to assay situation.

13. Trouble Shooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not be added.
- 2) Reagents necessary for coloration such as Peroxidase-conjugated Antibody Solution or TMB Solution might not be added.
- 3) Wrong reagents related to coloration might have been added. Wrong dilution of Peroxidase-conjugated Antibody Solution.
- 4) Contamination of enzyme inhibitor(s).
- 5) Influence of the temperature under which the kits had been stored.
- 6) Excessive hard washing of the well plate.
- 7) Addition of TMB Solution soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- Blank OD was higher than that of the lowest standard concentration (3.13 ng/mL).

Possible explanations :

Improper or inadequate washing. (Change washing frequency from 3 times to 4 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Antibody Solution.)

- High coefficient of variation (CV)

Possible explanations :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon

- Q-2 : I found there contains liquid in 96 well-plate when I opened the box. What is it?

A-2 : When we manufacture 96 well-plate, we insert preservation stabilizer in wells.


14. Summary of Assay Procedure ☐ : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

☐ Bring the well-plate and all reagents to **20°C - 25°C for 2 hours**.

☐ (I) Wash Solution (10×) must be diluted to **10 times** by purified water.

☐ Standard solution dilution example :

- | | | |
|--------------------------|--|-------|
| <input type="checkbox"/> | Antigen-coated Plate | |
| <input type="checkbox"/> | ↓ Washing 3 times (After removal of wash buffer, immediately dispense the next reagent.) | *① |
| <input type="checkbox"/> | Diluted Samples / Standards | 50 µL |
| <input type="checkbox"/> | ↓ Shaking, Incubation for 1 hour at 20°C - 25°C (Standing) | *②③ |
| <input type="checkbox"/> | Meanwhile, Dilute (D) Peroxidase-conjugated Antibody Solution to 1 : 100 with (C) Buffer returned to 20°C - 25°C. | |
| <input type="checkbox"/> | ↓ Washing 3 times (After removal of wash buffer, immediately dispense the next reagent.) | *① |
| <input type="checkbox"/> | Peroxidase-conjugated Antibody Solution | 50 µL |
| <input type="checkbox"/> | ↓ Shaking, Incubation for 1 hour at 20°C - 25°C (Standing) | *②③ |
| <input type="checkbox"/> | ↓ Washing 3 times (After removal of wash buffer, immediately dispense the next reagent.) | *① |
| <input type="checkbox"/> | TMB Solution
(After dispense, the color turns to blue depending on the concentration.) | 50 µL |
| <input type="checkbox"/> | ↓ Shaking, Incubation for 20 min at 20°C - 25°C (Standing) | *②③ |
| <input type="checkbox"/> | Stop Solution
(After dispense, the color turns to yellow depending on the concentration.) | 50 µL |
| <input type="checkbox"/> | ↓ Shaking (Immediately shake) | *② |
| <input type="checkbox"/> | Measurement of absorbance (450 nm, Ref 620 nm : 600 nm - 650 nm)
immediately using a plate reader. (Ref. wave cancels the dirt in the back of plate.) | |

*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 sec and remove the buffer. Guideline of washing volume : 300 µL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 3 times to 4 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Antibody Solution.

Standard of plate-washing pressure : 5 mL/min - 25 mL/min. (Adjust it depending on the nozzle's diameter.)


*② Guideline of shaking : 600 rpm - 1,200 rpm 3 times for 10 seconds.

*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the Plate Seal used once.

Worksheet Example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	200 ng/mL	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
B	100 ng/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
C	50 ng/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
D	25 ng/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
E	12.5 ng/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
F	6.25 ng/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
G	3.13 ng/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
H	0 (Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40

Assay Worksheet

15. Storage and Expiration

The intact kit is stored at 2°C - 10°C. Reagents, once opened, should be used as soon as possible to avoid losing its optimal assay performance by storage environment.

LBIS™ Rat KLH-IgM ELISA Kit (TDAR)

[Storage] Store at 2°C - 10°C (Do not freeze).

[Expiration date] Indicated on the container

[Package] For 96 tests

[Cat #] 297-91301

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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