## LabAssay™ Ammonia

Please, read this instruction carefully before use.

### 1. Introduction

Ammonia is converted to urea in the urea cycle and excreted in the urine.

LabAssay™ Ammonia is the reagent kit for assay of ammonia based on a colorimetric method using the Berthelot reaction. This kit is used for the quantitative determination of ammonia in mouse, rat and human blood. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

### 2. Kit Contents

	Components	State	Amount
(1)	Deproteinizing Reagent	liquid	100 mL/1 bottle
(2)	Chromogen Reagent A	liquid	50 mL/1 bottle
(3)	Chromogen Reagent B	liquid	25 mL/1 bottle
(4)	Chromogen Reagent C	liquid	50 mL/1 bottle
(5)	Standard Solution	liquid	15 mL/1 bottle
(6)	Dilute Solution For Standard	liquid	20 mL/1 bottle


※ Assay in a microplate format ; 700 tests. Assay in a test tube format ; 60 tests.

### 3. Storage and Expiration

When the complete kit is stored at 2°C - 10°C, the kit is stable until the expiration date shown on the label on the container.

### 4. Assay Principle

The Deproteinizing Reagent deactivates the enzymes and removes the components that inhibit color reaction in a sample. After adding Deproteinizing Reagent to a sample and centrifuging, add Phenol, Sodium Pentacyanonitrosylferrate (III) Dihydrate, and Sodium Hypochlorite to the supernatant. The indophenol blue pigment is obtained by the following schematic reaction. The amount of ammonia contained in the sample is determined by measuring the absorbance of the blue color.

※ Sodium pentacyanonitrosylferrate (III) dihydrate acts as catalyst

### 5. Materials and Apparatuses to be prepared

#### (For Microplate method)

- 96 wells microplate (transparent type) · Micropipette · Incubator maintaining at 37°C \*
- Plate mixer \*
- Microplate reader with a 630 nm wavelength filter (\* if the microplate reader is not equipped.)

#### (For Test Tube method)

- Test tube or microtube · Pipette · Incubator maintaining at 37°C
- Spectrophotometer with a 630 nm wavelength filter

### 6. Preparation of Reagents to be used

Before the measurement, transfer the required reagents to the tubes and adjust the reagents to the following temperatures.

- Return Ammonia Standard Solution, Deproteinizing Reagent, and Dilute Solution for Standard to room temperature \*.
- Incubate Chromogen Reagent A, B, and C at 37°C for an hour.

\* Room temperature : 20°C - 25°C

(1) Standard Solution (Microplate method)

Standard Solution is prepared by dilution of the provided Standard.

No.	Standard	Dilute Solution For Standard	Concentration
1	50 $\mu$ L	150 $\mu$ L	100 $\mu$ g/dL
2	100 $\mu$ L	100 $\mu$ L	200 $\mu$ g/dL
3	150 $\mu$ L	50 $\mu$ L	300 $\mu$ g/dL
4	undiluted	—	400 $\mu$ g/dL

(2) Standard Solution (Test Tube method)

Standard Solution is prepared by dilution of the provided Standard.

No.	Standard	Dilute Solution For Standard	Concentration
1	200 $\mu$ L	600 $\mu$ L	100 $\mu$ g/dL
2	400 $\mu$ L	400 $\mu$ L	200 $\mu$ g/dL
3	600 $\mu$ L	200 $\mu$ L	300 $\mu$ g/dL
4	undiluted	—	400 $\mu$ g/dL

## 7. Preparation of Samples

- Ammonia in whole blood increases at room temperature. Therefore, collect the blood and transfer to test tubes or microtubes containing Deproteinizing Reagent and centrifuge the samples to separate supernatant as quickly as possible.
- Store the supernatant on ice or at 2°C - 10°C and measure the sample within an hour.
- Do not use anticoagulants such as heparin and double oxalate.
- Amino acid preparations can cause a slightly negative effect on the assay.
- Amino acids in the blood may affect the assay.

## 8. Procedure

Before the measurement, transfer the required reagents to the tubes and adjust the reagents to the following temperatures.

- Return Standard Solution, Deproteinizing Reagent and Dilute Solution For Standard to room temperature \*.
- Incubate Chromogen Reagent A, B and C at 37°C for an hour.

\* Room temperature : 20°C - 25°C

(1) Assay in a Microplate format

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Blank
An experiment can be conducted in a microtube			
Deproteinizing Reagent	Deproteinizing reagent : Sample = 4 : 1 (v/v)	Deproteinizing reagent : Standard solution = 4 : 1 (v/v)	70 $\mu$ L
Sample			—
	Mix well. Centrifuge at 5000 $\times$ g, 4°C, for 15 min. Subsequent operations should be performed in the wells	Mix well. Subsequent operations should be performed in the wells.	
	Supernatant 70 $\mu$ L	Mixed solution 70 $\mu$ L	—
Chromogen Reagent A	70 $\mu$ L	70 $\mu$ L	70 $\mu$ L
Chromogen Reagent B	35 $\mu$ L	35 $\mu$ L	35 $\mu$ L
Chromogen Reagent C	70 $\mu$ L	70 $\mu$ L	70 $\mu$ L

After adding each reagent to a microplate, mix well and incubate at 37°C, for 20 min.

Mix well, then stand for 30 min at room temperature \*.

Mix well, then measure the absorbance of the test sample and standard solution at 630 nm with the blank as the control.

\* Room temperature : 20°C - 25°C

## (2) Assay in a Test Tube


Perform the assay in the tubes according to the following table scheme.

	Test	Standard	Blank
Operations should be performed in the microtubes.			
Deproteinizing Reagent	Deproteinizing reagent : Sample = 4 : 1 (v/v)	Deproteinizing reagent : Standard solution = 4 : 1 (v/v)	800 $\mu$ L
Sample			—
	Mix well. Centrifuge at 5000 $\times$ g, 4°C, for 15 min.	Mix well.	
	Transfer the supernatant (800 $\mu$ L) to a new tube.	Transfer mixed solution (800 $\mu$ L) to a new tube.	—
Chromogen Reagent A	800 $\mu$ L	800 $\mu$ L	800 $\mu$ L
Chromogen Reagent B	400 $\mu$ L	400 $\mu$ L	400 $\mu$ L
Chromogen Reagent C	800 $\mu$ L	800 $\mu$ L	800 $\mu$ L
After adding each reagent to a tube, mix well and incubate at 37°C, for 20 min. Mix well, and place in water to cool (20°C - 25°C) for 30 min. Mix well, then measure the absorbance of the test sample and standard solution with the blank solution as the control. Spectrophotometer : 630 nm.			


※ Assay in a test tube ; 60 tests.

## 9. Standard Curve

[Assay in a microplate]

[Assay in a test tube]

## 10. Performance

### (1) Sensitivity [assay in a test tube method]

- The absorbance is 0.010 to 0.070, when measuring purified water as a sample.
- The absorbance is 0.200 to 0.370, when measuring 400  $\mu$ g/dL ammonia as a sample.

### (2) Specificity

The ammonia concentration is less than  $\pm 35\%$ , when measuring the known concentration of control blood as a sample.

## 11. Usage Notes

### (1) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept

under the specified conditions.

- Do not use the containers and other materials in the package for any purpose other than those described herein.
- Reagents which are kept open for a long time cause dispersion of the measured values by absorbing ammonia. The opened vials should be capped and store at 2°C - 10°C.
- Avoid exposing reagents to direct sunlight during operation.
- Do not use the reagents at other reaction temperatures and reaction times other than described herein.
- Wash the equipment with a weak alkaline solution. Wash them with tap water. Rinse them well with deionized water. Tap water and distilled water may contain a lot of volatile ammonium salts. Pay attention to washing of experimental instruments.
- When collecting the supernatant after centrifuging, pay attention not to mix with suspended matter. Suspended matter may cause dispersion of the measured values.
- This kit is for research use only. Not for use in diagnostic procedure.

(2) Safety Precautions

- If reagents come into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult with a physician if necessary.
- Deproteinizing Reagent is an acid solution which is less than pH 3.
- Chromogen Reagent B and C are alkaline solutions which are more than pH 11.
- A pipette with a safety pipette filter should be used.

(3) Waste

- The waste should be processed appropriately according to the laws and regulations.
- Chromogen Reagent A contains 0.015% Sodium Pentacyanonitrosylferrate (III) Dihydrate and 40 g/L Phenol.
- All the devices including, reagents and vials, that come into contact with the specimen should be considered potentially infectious.

LabAssay™ Ammonia

[Storage]	Store at 2°C - 10°C (Do not freeze)
[Expiration date]	Indicated on the label
[Package]	For 700 tests
[Cat #]	297-94601

---

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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