LumiMAT™ Pyrogen Detection Kit

[Product information]

| Code No. | Product name | Size | Storage |
|-----------|----------------------------------------------|----------|---------|
| 297-96801 | LumiMAT™ Pyrogen Detection Kit - Reagent Set | 96 Tests | -20°C |
| 298-36991 | LumiMAT™ Pyrogen Detection Kit - Cells | 96 Tests | -80°C |

[Introduction]

This kit is used to detect pyrogen contamination in injectable pharmaceuticals and medical devices. This product utilizes a reporter assay system instead of the conventional ELISA system. The kit's very simple operation consists of adding luciferase substrate directly to cells that have been incubated with the samples to be tested. This unique procedure can reduce the risk of human errors and variation between tests. In addition, the assay is serum-free and is able to achieve constant performance independent of serum quality. The incubation time of the samples is only three hours, allowing for results to be obtained quickly. LumiMAT™ Pyrogen Detection Kit can detect a variety of non-endotoxin pyrogens as well as endotoxin.

[Assay principle]

This kit uses a monocytic cell line NOMO-1 in which a luciferase reporter gene is introduced to express luciferase protein in response to NF-κB signals activated by exposure to pyrogens. Because the expressed luciferase generates luminescence by reacting with the substrate, luminescence can be detected by a luminescence microplate reader to identify pyrogens in the tested samples.

[Product components]

LumiMAT™ Pyrogen Detection Kit - Reagent Set

| Component name | Size |
|-------------------------|--------|
| Assay medium | 20 mL |
| Dilution medium | 100 mL |
| Luciferase assay buffer | 12 mL |
| Luciferase substrate | 240 μL |

LumiMAT™ Pyrogen Detection Kit - Cells

| Component name | Size |
|----------------------------------------|--------|
| LumiMAT™ Pyrogen Detection Kit - Cells | 350 μL |

For single use only.

[Other reagents and instruments required]

- Safety cabinet
- CO₂ incubator
- Water bath
- Vortex mixer
- Electronic pipette controller and 10 mL disposal pipette, sterile, pyrogen-free
- Micropipette
- 200 μL pyrogen-free tip [BioCleanTip Wako™ Extend S II (#294-35011)]
- 1000 μL pyrogen-free tip [BioCleanTip Wako™ 1000 II (#298-35031)]
- Electronic multi-dispenser pipette [Eppendorf Multipette® E3 or equivalent] and 5 mL disposal pipette, sterile, pyrogen-free
- Pyrogen-free glass tube [Depyrogenated Dilution Tube (#DL-13100)]
- Aluminium Cap [Aluminium Cap-S (#293-28251)]
- Pharmacopeial standard endotoxin
- Water for bacterial endotoxin test [Lysate Reagent Water, 30 mL (#LRW-2030)]
- Ethanol for disinfection
- 15 mL centrifuge tube, sterile, pyrogen-free [Sumitomo Bakelite (#MS-56150) or equivalent]
- 96 well cell culture white plate [Nunc™ 96 well flat bottom (#136101)]
 - *In case using a plate other than the recommended one, confirm that the plate has no influence to assay in advance.
- Centrifuge
- Luminescence microplate reader

[Precautions]

- Both LumiMAT™ Pyrogen Detection Kit - Cells and LumiMAT™ Pyrogen Detection Kit - Reagent Set are necessary for assay.
- LumiMAT™ Pyrogen Detection Kit - Cells is very temperature sensitive product. Temperature increase in frozen condition significantly damage assay sensitivity. Store it in a -80°C or lower freezer immediately after receiving the product. When moving the product, put into a box of crashed-dry ice to cover the cell container with dry ice and avoid exposure to room temperature.
- LumiMAT™ Pyrogen Detection Kit - Cells are human derived, and gene modified cells. Handle these following the applicable regional, national, and local laws as well as regulations.
- This assay responds very sensitively to pyrogens. Contamination with microorganisms or pyrogens from the environment or procedure may lead to incorrect results. Use care and attention to ensure no contamination.
- The reaction of LumiMAT™ may be influenced by samples to be tested. Check the influence of the test samples before testing if necessary.
- This product is for single use only. Use all the reagents once they are opened. The performance is not guaranteed if the cells are used after culture growth. Be sure to use the cells immediately after thawing.

- It is recommended to regularly calibrate the luminescence microplate reader following the manufacturer's manual.
- Do not administrate the products to any living things.
- Follow your institute's safety manual and use appropriate safety precautions.
- Use appropriate personal protective equipment including protective gloves and protective eyeglasses.
- Disposal should be in accordance with applicable regional, national, and local laws and regulations by referring to the Safety Data Sheet.

[Assay procedure]

1. Preparation before assay


- (1) Thaw Dilution medium in advance (4°C over night or 37°C 1 hour). Store the thawed Dilution medium at 4°C until use.
- (2) Warm water bath to 37°C.
- (3) Warm Assay medium in a water bath for about 15 minutes.

2. Preparation of reference standard endotoxin

- (1) Reconstitute reference standard endotoxin to make 2000 EU/mL of solution by following the supplier's instruction.
- (2) Dilute the reference standard endotoxin solution with Dilution medium. The following table is an example dilution series. Use pyrogen-free vials which endotoxin do not adsorb for standard dilution.

| Tube name | Solution | Endotoxin concentration |
|-----------|------------------------------------------------------------|-------------------------|
| HC1 | 20 μ L Reference Standard +180 μ L Dilution medium | 200 EU/mL |
| HC2 | 100 μ L HC1 + 900 μ L Dilution medium | 20 EU/mL |
| HC3 | 100 μ L HC2 + 900 μ L Dilution medium | 2 EU/mL |
| PC1 | 800 μ L HC3 + 1200 μ L Dilution medium | 0.8 EU/mL |
| PC2 | 1000 μ L PC1 + 1000 μ L Dilution medium | 0.4 EU/mL |
| PC3 | 1000 μ L PC2 + 1000 μ L Dilution medium | 0.2 EU/mL |
| PC4 | 1000 μ L PC3 + 1000 μ L Dilution medium | 0.1 EU/mL |
| PC5 | 1000 μ L PC4 + 1000 μ L Dilution medium | 0.05 EU/mL |
| PC6 | 1000 μ L PC5 + 1000 μ L Dilution medium | 0.025 EU/mL |
| PC7 | 1000 μ L PC6 + 1000 μ L Dilution medium | 0.0125 EU/mL |

HC : High Concentration, PC : Positive Control

3. Sample preparation

Dilute samples not exceeding the maximum dilution factor (MVD) as calculated by the following expression.

$$\text{MVD} = (\text{CLC} \times \text{C}) / \text{Test Sensitivity}$$

- CLC : Contaminant limit concentration (EE/mg)

- C : Concentrations of test solutions (mg/mL)

- Test sensitivity : The lowest endotoxin reference standard concentration on the standard curve.

The CLC is the specified contaminant limit concentration for a product ; if unknown, it can be calculated using the following expression :

$$\text{CLC} = \text{K} / \text{M}$$

- K : Threshold pyrogenic dose per kilogram of body mass (EE/kg)

Administration method ① Intravenous : 5 EE/kg

② Intravenous (radioactive) : 2.5 EE/kg

③ Intrathecal : 0.2 EE/kg

- M : Maximum recommended bolus dose of product per kilogram of body mass

4. Transfer to plate


Add 50 μ L of each solution prepared in steps 2 and 3 to a 96-well white plate. An example of the plate layout is shown below.

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-------------------------|---|---|---|----------|---|---|---|-----------|----|----|----|
| A | Blank (Dilution medium) | | | | Sample 1 | | | | Sample 9 | | | |
| B | PC 7 | | | | Sample 2 | | | | Sample 10 | | | |
| C | PC 6 | | | | Sample 3 | | | | Sample 11 | | | |
| D | PC 5 | | | | Sample 4 | | | | Sample 12 | | | |
| E | PC 4 | | | | Sample 5 | | | | Sample 13 | | | |
| F | PC 3 | | | | Sample 6 | | | | Sample 14 | | | |
| G | PC 2 | | | | Sample 7 | | | | Sample 15 | | | |
| H | PC 1 | | | | Sample 8 | | | | Sample 16 | | | |

5. Reporter cell seeding

- (1) Add 9.5 mL of pre-warmed Assay medium to a 15 mL centrifuge tube.
- (2) Thaw LumiMAT™ Pyrogen Detection Kit - Cells in a 37°C water bath for 1 minute by gently swirling the vial. Take out the vial from the water bath when a small crump of ice remains.
- (3) Add 1 mL of Assay medium from step “(1)” to the cells, then mix by pipetting.

- (4) Transfer all the cell suspension to the 15 mL centrifuge tube prepared in step “(1)”.

- (5) Centrifuge the cells at 300 g for 5 minutes.
 (6) Remove the supernatant and suspend the cell pellet with 5.5 mL of Assay medium.
 (7) Add 50 μ L of cell suspension to each well in which standard/sample were dispensed.
 *Electronic multi-dispenser pipette is recommended over multi-channel pipette to minimize the variances among wells (Dispensing speed : 4 for Eppendorf Multipipette® E3).

6. Incubation

Incubate the plate for 3 hours in a 37°C, 5% CO₂ incubator.

7. Measurement

- (1) Add all of Luciferase substrate to Luciferase buffer and mix by inverting the vial. Bring Luciferase buffer to room temperature before use.
 (2) Add 100 μ L of the mixed solution to each well.
 *Electronic multi-dispenser pipette is recommended over multi-channel pipette to minimize the variances among wells (Dispensing speed : 4 for Eppendorf Multipipette® E3).
 (3) Measure the luminescent value with a luminescence microplate reader.

8. Calculation

Toximaster™ FQC1 Software for LumiMAT™ (sold separately) can be used for analysis. See the Toximaster™'s instruction manual for guidance on analysis by Toximaster™.

Calculate the pyrogen concentration as Endotoxin Equivalents (EE)/mL by comparing the luminescent value of the samples with the luminescent value of the reference standard endotoxin.

1. Plot the concentration of reference standard endotoxin on the X axis and the luminescent value on the Y axis and create a four parameters logistic (4PL) standard curve.

$$Y = D + \frac{A - D}{1 + (X/C)^B}$$

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2. Rearrange the expression above to the expression below, then substitute the luminescent values of the samples and calculate the “EE/mL” value.


$$X = C \times \{(A - Y)/(Y - D)\}^{(1/B)}$$

<calculation example>

A = 14732365, B = -1.3886, C = 0.37437, D = 61400, luminescent value of sample = 928030

$$X = 0.37437 \times \{(14732365 - 928030)/(928030 - 61400)\}^{(1/-1.3886)} = 0.0510 \text{ (EE/mL)}$$

<Standard curve example>

[Assay procedure video]

An assay procedure video is available on the product web page. Access the web page from the QR code below or search for product code No. (297-96801) on the Fujifilm Wako web site at <https://labchem-wako.fujifilm.com>.

<License related notice>

- * This product is patent pending.
 * CHO-Spica cells were used to express protein used in this kit.
 * NanoLuc® Technology is licensed from Promega Corporation. Licensed patents: U.S. Pat. No. 8557970 and U.S. Pat. No. 8669103, and all patents and patents pending which claim priority to the same priority application(s) as U.S. Pat. No. 8557970 and U.S. Pat. No. 8669103. NanoLuc® is a registered trademark of Promega Corporation.

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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※ NanoLuc® Technology is licensed from Promega Corporation. Licensed patents : U.S. Pat. No. 8557970 and U.S. Pat. No. 8669103, and all patents and patents pending which claim priority to the same priority application(s) as U.S. Pat. No. 8557970 and U.S. Pat. No. 8669103. NanoLuc® is a registered trademark of Promega Corporation.