## For Genetic Research RNA Quantification Kit for Purified EV

### [Introduction]

This kit enables highly sensitive quantification of trace amounts of total RNA contained in a sample by using RNA-binding fluorescent dyes. It can be used to quantify total RNA extracted from purified extracellular vesicles (EVs) derived from cell culture supernatants, serum, or plasma.

This kit can be used to quantify trace amounts of RNA derived from EVs extracted using the microRNA Extractor™ Kit for Purified EV (Code No. 294-84601).

This product is intended for laboratory use only.

### [Features]

- Quantification of trace amounts of RNA samples  
Sample volume required for measurement : 1 - 20  $\mu$ L
- High sensitivity (Quantification Range : 0.1 - 10 ng/ $\mu$ L)  
Ideal for quantifying trace amounts of RNA, such as RNA extracted from purified EV solutions
- Short measurement time  
Fluorescence-based measurement method with a reaction time of only 5 minutes

### [Kit contents (for 1,000 assays)]

This kit consists of 3 components.

| Item                           | Amount                      |
|--------------------------------|-----------------------------|
| Control RNA                    | 100 $\mu$ L $\times$ 1 vial |
| Fluorescent Reagent            | 100 $\mu$ L $\times$ 1 vial |
| Dilution Buffer (20 $\times$ ) | 2.5 mL $\times$ 1 vial      |

### [Other materials required]

- Nuclease-free water (Code No. 314-09291)
- 96-well plates
- Micropipettes
- Pipette tips
- Microplate shaker (if available)
- Microplate reader for 96-well plates (for fluorescence measurement)

### [Precautions]

#### 1. Tools

Microcentrifugation tubes, PCR tubes, pipette tips, and other experimental tools should be autoclaved, or use commercially available RNase-free products. Plastic gloves and masks should be worn during the experiments to avoid contamination of RNase.

#### 2. Reagents

Since the procedure primarily involves RNA quantification, handle all the reagents with great care, and strive to maintain the sterility of the working reagents during use.

### [Protocol]

Preparation of RNA samples derived from purified EVs and precautions

Prepare a purified EV solution in advance using the MagCapture™ Exosome Isolation Kit PS Ver.2 (Code No. 290-84103) or by some other means. This solution is used as the sample for extraction. Next, prepare a quantitative RNA sample using a MicroRNA Extractor™ Kit for purified EV (Code No. 294-84601) or by some other means.

#### Preparation of reagents

##### 1. 1 $\times$ Dilution Buffer

Add nuclease-free water to Dilution Buffer (20 $\times$ ) at a ratio of 1 : 20 and mix well.

##### 2. Fluorescent Buffer

Add 1 $\times$  Dilution Buffer to the Fluorescent Reagent at a ratio of 1 : 2,000 and mix well\*1.

##### 3. Control RNA Solution

Dilute the Control RNA (100 ng/ $\mu$ L) stepwise as shown below to prepare Control RNA Solutions of each concentration\*2.

| Concentration of the Control RNA Solution (ng/ $\mu$ L) | Volume of the Control RNA Solutions       | Volume of 1 $\times$ Dilution Buffer |
|---|---|--------------------------------------|
| 100   | Stock Control RNA : 100 $\mu$ L           | —                                    |
| 1   | 100 ng/ $\mu$ L solution : 2 $\mu$ L      | 198 $\mu$ L                          |
| 0.5   | 1 ng/ $\mu$ L solution : 50 $\mu$ L       | 50 $\mu$ L                           |
| 0.25  | 0.5 ng/ $\mu$ L solution : 50 $\mu$ L     | 50 $\mu$ L                           |
| 0.125   | 0.25 ng/ $\mu$ L solution : 50 $\mu$ L    | 50 $\mu$ L                           |
| 0.0625  | 0.125 ng/ $\mu$ L solution : 50 $\mu$ L   | 50 $\mu$ L                           |
| 0.03125   | 0.0625 ng/ $\mu$ L solution : 50 $\mu$ L  | 50 $\mu$ L                           |
| 0.015625  | 0.03125 ng/ $\mu$ L solution : 50 $\mu$ L | 50 $\mu$ L                           |

\*1 Protect from light by wrapping in aluminum foil, etc.

\*2 The table shows an example of stepwise dilution. Adjust the dilution concentration of the Control RNA Solutions for each sample.

#### Assay Procedure

- Dispense 200  $\mu$ L of Fluorescent Buffer into each well.
- Dispense 10  $\mu$ L of Control RNA Solution into the wells designated for standards.  
Add 10  $\mu$ L of 1 $\times$  Dilution Buffer to the blank.
- Dispense 1 - 20  $\mu$ L of pretreated-RNA samples into the wells designated for samples.
- Agitate the mixture in each well using a microplate shaker or by pipetting.
- Incubate at room temperature for 5 minutes in the dark.
- Determine the fluorescent intensity using the 96-well microplate reader (Excitation wavelength : 492 nm, Emission wavelength : 540 nm).

#### Calculation Method

- (1) Create a standard curve by plotting the concentration of Control RNA amount (ng) on the X-axis against the fluorescent intensity value on the Y-axis.
- (2) Read the concentration (ng) corresponding to the fluorescent intensity of the diluted sample.
- (3) Multiply the read concentration by the sample dilution factor to obtain the measured concentration.


Figure 1. Example of dispensing Control RNA Solution to a multiwell plate


Figure 2. Standard curve (an example)

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