<for Research Use Only> Code No. 290-68603 (10 reactions)  
298-68604 (40 reactions)

## Immuno-enhancer

### [Introduction]

Immuno-enhancer is a reagent which accelerates antigen-antibody reactions of Western blotting, dot blotting and ELISA by optimizing the reaction conditions. The effect is more noticeable when using a weakly-reactive antibody. It is possible to get high S/N ratio.

Immuno-enhancer contains Reagent A for the reaction of the primary antibody and Reagent B for the reaction of the secondary antibody. The reagents are ready-to-use, so can be used without dilution.

### [Features]

- 1) Enhancing of signals
- 2) High S/N ratio
- 3) Ready-to-use

### [Reagent contents]

	10 reactions	40 reactions
Reagent A	50 mL × 1 vial	200 mL × 1 vial
Reagent B	50 mL × 1 vial	200 mL × 1 vial

\*Package size is based on the use when 5 mL of Reagent A and 5 mL of Reagent B are used for 1 reaction.

### [Storage]

2 ~ 10 °C

### [Protocols]

#### Western Blotting

##### < Using primary antibody and secondary antibody >

- 1) Perform the electrophoresis of the samples. Transfer to the membrane and block the membrane following each protocol.

- It is possible to use the following blocking solutions: skim milk, casein, BSA, synthetics and so on. We recommend using of non-fat skim milk when using skim milk.
- Use the suitable blocking solution for the antibody and the assay because some antibodies may increase non-specific signals depending on the kind of the blocking solution.
- Wash the membrane after the blocking step because the residual blocking solution may cause background signals.  
Washing conditions : TBS-T (rinse × 2, wash for 15 min × 1, for 5 min × 2).

- 2) Dilute the primary antibody with Reagent A.

- Adjust the volume of Reagent A to the membrane.  
Ex : 5-10 mL / 10 × 10 cm membrane.
- Refer to the concentration of the antibody recommended by the supplier.

- 3) Incubate the membrane in the primary antibody solution at room temperature for 1 hour with shaking.

- Change the incubation time as necessary. Refer to the incubation time recommended by the supplier of the antibody.

- 4) Wash the membrane in TBS-T for 10 minutes with shaking.(3 times)
- 5) Dilute the secondary antibody with Reagent B.

- Adjust the volume of Reagent A to the membrane.  
Ex : 5-10 mL / 10 × 10 cm membrane.
- Refer to the concentration of the antibody recommended by the supplier.

- 6) Incubate the membrane in the secondary antibody solution at room temperature for 1 hour with shaking.

- Change the incubation time as necessary. Refer to the incubation time recommended by the supplier of the antibody.

- 7) Wash the membrane in TBS-T for 10 minutes with shaking.(3 times)
- 8) Proceed to the detection.

##### < Using labeled antibody, in case of direct detection >

- 1) Perform the electrophoresis of the samples. Transfer to the membrane and block the membrane following each protocol.

- It is possible to use the following blocking solutions: skim milk, casein, BSA, synthetics and so on. We recommend using of non-fat skim milk when using skim milk.
- Use the suitable blocking solution for the antibody and the assay because some antibodies may increase non-specific signals depending on the kind of the blocking solution.
- Wash the membrane after the blocking step because the residual blocking solution may cause background signals.  
Washing conditions : TBS-T (rinse × 2, wash for 15 min × 1, for 5 min × 2).

- 2) Dilute the antibody with Reagent B.

- Adjust the volume of Reagent B to the membrane.  
Ex : 5-10 mL / 10 × 10 cm membrane.
- Refer to the concentration of the antibody recommended by the supplier.
- When you can not get a good result with using Reagent B, please try to use Reagent A. Because some antibodies may be able to get better result with using Reagent A.

- 3) Incubate the membrane in the antibody solution at room temperature for 1 hour with shaking.
- 4) Wash the membrane in TBS-T for 10 minutes with shaking.(3 times)
- 5) Proceed to the detection.

### [Note]

Immuno-enhancer may not be able to be used for Western blotting of samples including the antibody for immunoprecipitation, because Immuno-enhancer accelerates antigen-antibody reactions therefore the secondary antibody reacts to the H-chain and L-chain of the antibody for immunoprecipitation.

### Dot Blotting

- 1) Spot the samples onto the membrane. Transfer to the membrane and block the membrane following each protocol.

- It is possible to use the following blocking solutions : skim milk, casein, BSA, synthetics and so on. We recommend using of non-fat skim milk when using skim milk.
- Use the suitable blocking solution for the antibody and the assay because some antibodies may increase non-specific signals depending on the kind of the blocking solution.
- Wash the membrane after the blocking step because the residual blocking solution may cause background signals.  
Washing conditions : TBS-T (rinse  $\times$  2, wash for 15 min  $\times$  1, for 5 min  $\times$  2).

- 2) Please follow the Western blotting protocol starting with Step 2).

### ELISA (Sandwich ELISA) in case of using 96-wells plate

#### < Using primary antibody and secondary antibody >

- 1) Coat the plate with antibody. Block the plate following ELISA protocol.

- It is possible to use the following blocking solutions : skim milk, casein, BSA, synthetics and so on. We recommend using of non-fat skim milk when using skim milk.
- Use the suitable blocking solution for the antibody and the assay because some antibodies may increase non-specific signals depending on the kind of the blocking solution.
- Wash the plate after the blocking step because the residual blocking solution may cause background signals.  
Washing conditions : PBS-T (250  $\mu$ L/well  $\times$  3).

- 2) Dilute the antigen in Reagent A.
- 3) Dilute the primary antibody in Reagent A.

- Refer to the concentration of the antibody recommended by the supplier.

- 4) Add 50  $\mu$ L of the antigen solution and 50  $\mu$ L of the primary antibody solution to each well and mix.
- 5) Seal the plate with a microplate sealing tape and incubate at 37°C for 1 hour.

- It is possible to incubate at 4°C for 16 hours.
- When performing the antigen reaction and the primary antibody reaction independently, wash the plate 3 times with 250  $\mu$ L of PBS-T per well after the antigen reaction.

- 6) Carefully remove the solution using a pipette.
- 7) Wash the plate with 250  $\mu$ L of PBS-T per well. (3 times )
- 8) Dilute the secondary antibody in Reagent B.

- Refer to the concentration of the antibody recommended by the supplier.

- 9) Add 100  $\mu$ L of the secondary antibody solution to each well.
- 10) Seal the plate with a microplate sealing tape and incubate at 37°C for 1 hour.

- It is possible to incubate at 4°C for 16 hours.

- 11) Carefully remove the solution using a pipette.

- 12) Wash the plate with 250  $\mu$ L of PBS-T per well. (3 times )
- 13) Proceed to the detection.

#### < Using labeled antibody, in case of direct detection >

- 1) Coat the plate with antibody. Block the plate following ELISA protocol.

- It is possible to use the following blocking solutions : skim milk, casein, BSA, synthetics and so on. We recommend using of non-fat skim milk when using skim milk.
- Use the suitable blocking solution for the antibody and the assay because some antibodies may increase non-specific signals depending on the kind of the blocking solution.
- Wash the plate after blocking step because the residual blocking solution may cause background signals.  
Washing conditions : PBS-T (250  $\mu$ L/well  $\times$  3).

- 2) Dilute the antigen in Reagent A.
- 3) Add 100  $\mu$ L of the antigen solution to each well.
- 4) Seal the plate with a microplate sealing tape and incubate at 37°C for 1 hour.

- It is possible to incubate at 4°C for 16 hours.

- 5) Carefully remove the solution using a pipette.
- 6) Wash the plate with 250  $\mu$ L of PBS-T per well. (3 times )
- 7) Dilute the labeled antibody in Reagent B.

- Refer to the concentration of the antibody recommended by the supplier.

- 8) Add 100  $\mu$ L of the secondary antibody solution to each well.
- 9) Seal the plate with a microplate sealing tape and incubate at 37°C for 1 hour.

- It is possible to incubate at 4°C for 16 hours.

- 10) Carefully remove the solution using a pipette.
- 11) Wash the plate with 250  $\mu$ L of PBS-T per well. (3 times )
- 12) Proceed to the detection.

## 【Troubleshooting】

Problem	Possible cause	Solution
Western blotting · Dot blotting		
Weak signal (Faint bands)	Insufficient protein	Increase the amount of total protein loaded on the gel. We recommend the examination of the optimum protein concentration by preparing a dilution series of protein.
	Insufficient antibody	Increase the antibody concentration. We recommend the examination of the optimum antibody concentration by preparing a dilution series of the antibody.
	Poor or incomplete transfer	Increase transfer time or voltage. High molecular weight proteins may require more time for transfer. Ensure that there is good contact between the membrane and gel. Altering transfer method from semi-dry to wet may improve the efficiency of transfer.
	Over transfer	Reduce transfer time or voltage for low molecular weight protein to avoid the protein going through the membrane. Altering the membrane from nitrocellulose to PVDF may get better results.
	Unsuitable blocking solution	Examine the suitable blocking solution for the antibody and the assay, because some antibodies may decrease the signals depending on the kind of blocking solution. Using skim milk including fat may cause weak signals.
Extra bands	Excessive antibody	Reduce the antibody concentration. We recommend the examination of the optimum antibody concentration by preparing a dilution series of the antibody.
	Excessive protein	Reduce the amount of protein loaded on the gel. We recommend the examination of the optimum protein concentration by preparing dilution series of protein.
	Insufficient blocking	Examine the suitable blocking solution for the antibody and the assay, because some antibodies may increase non-specific signals depending on the kind of the blocking solution. Do not use casein when using anti phosphor antibodies, because casein is a phosphorylated protein.
	Insufficient wash	Increase the number of washes.
High background	Excessive antibody	Reduce the antibody concentration. We recommend the examination of the optimum antibody concentration by preparing a dilution series of antibody.
	Excessive protein	Reduce the amount of protein loaded on the gel. We recommend the examination of the optimum protein concentration by preparing a dilution series of protein.
	Insufficient blocking	Examine the suitable blocking solution for the antibody and the assay, because some antibodies may increase non-specific signals depending on the kind of the blocking solution. Do not use casein when using anti phosphor antibodies, because casein is a phosphorylated protein.
	Insufficient wash	Increase the number of washes.
	Film overexposed	Reduce exposure time.

Problem	Possible cause	Solution
White band	Excessive protein	Reduce the amount of protein loaded on the gel. Excessive protein can cause extremely high levels of localized signal, and it results in rapid, complete consumption of substrate at the point. Therefore a part of the band cannot be observed when exposed to film.
	Excessive antibody	Reduce the antibody concentration. Excessive antibody can cause extremely high levels of localized signal, and it results in rapid, complete consumption of substrate at the point. Therefore a part of the band cannot be observed when exposed to film.
ELISA		
Weak signal	Insufficient antigen or antibody	Examine the optimum concentrations of antigen and antibody.
Excessive signal	Excessive antigen or antibody	Examine the optimum concentrations of antigen and antibody.
	Excessive reaction time	Reduce the antigen-antibody reaction time.
High background	Excessive antigen or antibody	Examine the optimum concentrations of antigen and antibody.
	Insufficient blocking	Examine the suitable blocking solution for the antibody and the assay, because some antibodies may increase non-specific signals depending on the kind of the blocking solution. Do not use casein when using anti phosphor antibodies, because casein is a phosphorylated protein.
	Insufficient wash	Increase the number of washes.
Data spread	ELISA plate	The absorption capacity may be very different between the plates. Additionally, some plates may be very different between lots. Be careful selecting the plates.

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
 Telephone : + 81-6-6203-3741  
 Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
 Richmond, VA 23237  
 U.S.A.  
 Telephone : + 1-804-271-7677  
 Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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