<For Research Use Only> Code No. 299-50601 (960 Tests)

for Cytotoxicity Assay

LDH-Cytotoxic Test *Wako*

■ Intended use

LDH-Cytotoxic Test *Wako* is used for the quantitative determination of cell death by measuring the LDH activity liberated from cells with cell membrane injuries.

■ Introduction

LDH-Cytotoxic Test *Wako* is a test kit for the determination of toxicity of various drugs, which utilizes the culture cells. LDH (lactate dehydrogenase) activity liberated from the cells with lysed cell membrane by drugs are determined with this kit. By using this kit, an objective measurement of toxicity can be made at a lower cost and with a simpler procedure compared to toxicity tests with experimental animals.

Furthermore, not only the natural cytotoxicity, following items concerning cell membrane injury can be measured as well.

1. Bioreactive substance (TNF and lymphotoxin, *etc.*) dependent cell injury.
2. Cell-mediated immunity (T-cells and macrophage, *etc.*) dependent cell injury.
3. Complement-dependent cell injury (Cytotoxin, *etc.*).
4. Antibody-dependent cell injury (Phagocytosis and K-cells, *etc.*).

This kit can also be used for the determination of survival rate of the cells.

■ Kit Contents


- | | | | |
|--------------|-----------------------------|-------------|------------------|
| 1. 296-50611 | Color Reagent | for 5 mL | 10 bottles |
| | · Nitrotetrazolium Blue | 3.7 mg/vial | |
| | · Diaphorase | | |
| | · NADH | | |
| 2. 293-50621 | Buffer Solution | 55 mL | 1 bottle |
| | · DL-Lithium lactate | 50 mg/mL | |
| 3. 290-50631 | Stop Solution | 55 mL | 1 bottle |
| | · 1 mol/L Hydrochloric acid | | |
| 4. 293-50641 | Microplates | | 10 plates |
| | | 96 wells | (Not sterilized) |

■ Characteristics

1. Measurement of alive cells not only dead cells is possible.
2. Highly sensitive determination of toxicity is possible by the direct measurement of dead cells.
3. Enzymatic measurement provides a quantitative determination of cell membrane injury with a high precision.
4. The detection sensitivity can be controlled by modifying the time for sample treatment or coloring reaction.
5. Because the culture supernatant is used as a specimen, there is no loss of the cells which enables prospective follow up of bioactivities.
6. Sterilized procedure is not necessary because of the short time procedure.

■ Principle of the method

LDH catalyzes the conversion of lactate to pyruvate in the presence of NAD. In the reaction lactate is oxidized to pyruvate and an equivalent amount of NAD is reduced to NADH. The NADH formed reduces nitrotetrazolium blue in the presence of diaphorase to produce purple-blue diformazan. The absorbance of this dye is measured at 560 nm (± 10 nm).

■ Detailed procedure

1. Materials required but not supplied

Medium
Cells
Sample
Phosphate buffered saline (PBS)
Sterilized 96 well microplate (for incubation):
for preparation of a Cell Plate and a Sample Plate
Tween 20
Centrifuge (with a rotor for microplate)
Dispenser (multichannel pipettes)
Microplate reader capable of reading absorbance at 560 $\pm$ 10 nm

2. Preparation of a Cell Plate

First, 96 well microplate with cells for sample treatment is prepared.

<Suspended cells>

- (1) Take cells into a centrifuge tube and centrifuge for 5 minutes at 1,000 rpm.
- (2) Remove the supernatant, suspend the cells in PBS, and centrifuge for 5 minutes at 1,000 rpm.
- (3) Remove the supernatant, dilute the cells to get a concentration of 1×10^5 cells/mL, and dispense 50 μ L of suspension into each well of the microplate (5,000 cells/well).
* Take care that a homogeneous suspension is dispensed into the well.
- (4) Use the plate prepared here as a Cell Plate for suspended cells.

<Adherent cells>

- (1) Take cells into a centrifuge tube and centrifuge for 5 minutes at 1,000 rpm.
- (2) Remove the supernatant, suspend the cells in the media, and centrifuge for 5 minutes at 1,000 rpm.
- (3) Remove the supernatant, dilute the cells to get a concentration of 1×10^5 cells/mL, and dispense 50 μ L of suspension into each well of the microplate (5,000 cells/well).
* Take care that a homogeneous suspension is dispensed into the well.
- (4) Incubate the plate with cells in 5% CO₂ incubator overnight at 37°C.
- (5) Remove the medium, wash twice with 100 μ L of PBS, and dispense 50 μ L of PBS into each well.
- (6) Use the plate prepared here as a Cell Plate for adherent cells.

3. Preparation of a Sample Plate

A 96 well microplate as a Sample Plate is prepared here before Sample Treatment.

The following procedure helps efficient dispensing of the sample into the Cell Plate using multichannel pipettes.

<Preparation of samples>

(1) Sample (S)

Prepare the samples of various concentrations diluted with PBS. Because samples are diluted by 2 (1 + 1) when sample treatment is made, prepare the samples with double concentration of the final concentration to be checked.

(2) Negative Control (NC)

Use PBS as a specimen for Negative Control.

(3) Positive Control (PC)

Use 0.2% Tween 20 dissolved in the PBS as a specimen for positive Control.

<Example of sample location>

All the samples with different concentration, Negative Control, and Positive Control should be located in the same positions as used for sample treatment.

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B		S	S	S	S	S	S	S	NC	NC	NC	
C		S	S	S	S	S	S	S	PC	PC	PC	
D												

4. Sample Treatment

(1) Dispensation of the samples

Transfer 50 μ L of each sample (including NC and PC) from Sample Plate to the respective well in the Cell Plate.

(2) Sample treatment

After dispensing the sample in (1), allow the plate to stand for 15 minutes at room temperature.

* Use the most suitable treatment condition for the experiment (e.g., usage of CO₂ incubator at 37°C, etc.). Detection sensitivity can be variable by changing the time for sample treatment.

(3) Mix the plate moderately.

(4) Centrifuge the plate for 3 minutes at 1,000 rpm.

(5) Take 50 μ L of supernatant from each well of the Sample Plate, dispense it into the plate included in the kit (Reaction Plate).

* Take care not to suck the cells when taking the supernatant.

5. Preparation of Coloring Solution

Dissolve the whole contents of Color Reagent with 5.0 mL of Buffer Solution and use this solution as Coloring Solution.

* One vial of Color Reagent is for use with one 96 well plate.

6. Coloring reaction

Dispense 50 μ L of Coloring Solution into the wells which contain supernatant from Sample Plate and allow it to stand for 45 minutes at room temperature.

* Detection sensitivity can be variable by changing the reaction time in this step.

7. Dilution of Stop Solution

Dilute the Stop Solution with the same volume of distilled or deionized water (1 + 1) and use this solution as Diluted Stop Solution.

8. Termination of coloring reaction

Dispense 100 μ L of diluted (1 + 1, with distilled or deionized water) Stop Solution into each well.

9. Measurement

Measure the absorbance of reaction mixture by a microplate reader at 560 $\pm$ 10 nm within 90 minutes after termination of coloring reaction.

* In the case that the absorbance for Positive Control exceeds 1.4, shorten the time for coloring reaction or dilute the supernatant obtained in step (5) of "4. Sample treatment" so that the absorbance of Positive Control becomes below 1.4.

10. Data analysis

Cell injury rate is calculated by using the following equation.


$$\text{Cell injury rate (\%)} = \frac{(S - N)}{(P - N)} \times 100 \quad \dots\dots\dots (1)$$

S: Absorbance of the sample

N: Absorbance of Negative Control

P: Absorbance of Positive Control

Prepare a graph similar to the one shown below by plotting the cell injury rate (ordinate) against sample concentration (abscissa) and determine LD₅₀ value. The LD₅₀ is defined as the sample concentration where the cell injury rate is 50%.

11. Correction of the results

LDH activity cannot accurately be measured when the sample contains some factor which influences the LDH activity is necessary.

When the sample inhibits or enhances the LDH activity, a correction should be made in the measurement of LDH activity (Coloring reaction). LDH solution prepared from the cell by the means like freezing and thawing of the cell is used to correct the LDH activity. Following absorbances are to be measured in the coloring reaction.

Cell injury rate is calculated by using the following equation.

$$\text{Cell injury rate (\%)} = \frac{(S'-N)}{(P-N)} \times 100 \dots\dots\dots(1)'$$

A: Absorbance obtained by the coloring reaction using LDH solution and sample with respective concentration as specimen.

B: Absorbance obtained by the coloring reaction using LDH solution and PBS.

In the case the sample is influentially colored, absorbance (C) for each sample obtained by the coloring reaction using PBS instead of color reagent should be deducted from S in equation (1) or (2).

■ Limitations of the procedure

- (1) Keep all the reagents under the specified conditions and do not use the reagent past the expiration date stated on each reagent container label.
- (2) Take care that the sample is not adhered to the side wall of the microplate and mix it well with Coloring Solution.
- (3) Colored reaction mixture after the addition of Reaction Terminator is stable for 3 hours under a light protection.
- (4) The Stop Solution contains 1N HCl and should be handled with a due caution.

■ Storage condition and expiry date

- (1) Storage: Keep at 2~10℃ in a dark place.
- (2) Expiry: Indicated on the label.

■ References

1. Korzeniewski, C. and Callewaert, D.M. An Enzyme-release Assay for Natural Cytotoxicity., *J. Immunol. Methods*, **64**: 313-320 (1983).
2. Decker, T. and Lohmann-Matthes, M.-L. A quick and simple method for the quantitation of lactate dehydrogenase release in measurements of cellular cytotoxicity and tumor necrosis factor (TNF) activity., *J. Immunol. Methods*, **15**: 61-69 (1988).

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
 Telephone : + 81-6-6203-3741
 Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
 Richmond, VA 23237
 U.S.A.
 Telephone : + 1-804-271-7677
 Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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