For Tissue Optical Clearing Reagent SCALEVIEW-S Trial Kit

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Dr. Atsushi Miyawaki et al. developed a sorbitol-based optical clearing method, *Scale S*.

SCALEVIEW-S is a high-performance clearing method that blessed with effective clearing and tissue/signal preservation. It is also noteworthy that this method is able to preserve highly fluorescence signal of fluorescent proteins and ultrastructure of cell membranes.

SCALEVIEW-S is a simple and reproducible method of accurate visualization of biological tissue.

【Kit contents】

This kit includes 6 components.

- | | |
|-------------------|-------------------|
| (1) SCALEVIEW-S0 | 100 mL × 1 bottle |
| (2) SCALEVIEW-S1 | 100 mL × 1 bottle |
| (3) SCALEVIEW-S2 | 100 mL × 1 bottle |
| (4) SCALEVIEW-S3 | 100 mL × 1 bottle |
| (5) SCALEVIEW-S4 | 100 mL × 1 bottle |
| (6) SCALEVIEW-SMt | 100 mL × 1 bottle |

Note : This kit needs an optional deScale Solution (Code No. 041-34425).

【Storage】

Store in the dark at 2-10℃

【Additional required materials】

Reagents

- 1) Paraformaldehyde (Code No. 160-16061)
- 2) 1 × PBS (-) (Code No. 164-25511)
- 3) Agarose
- 4) deScale Solution (Code No. 041-34425)

Equipments :

- 1) 5 mL/15 mL/30 mL/50 mL conical centrifuge tube
- 2) Culture dish (100 mm, 60 mm, and 35 mm diameter) (for each sample)
- 3) Seesaw shaker
- 4) Fluorescence microscope and Recommend Objective lenses for SCALEVIEW-S treated samples.
 1. Multiphoton microscope EVIDENT, model : FVMPE-RS
 - Multiphoton dedicated objectives :
 - EVIDENT, model : XLPLN10XSVMP : 10x, NA=0.6, WD=8 mm
 - EVIDENT, model : XLSLPLN25XGMP : 25x, NA=1.0, WD=8 mm
 2. Confocal microscope : EVIDENT, model : FV3000
 - Objectives :
 - EVIDENT, model : UPLSAPO10X2 : 10x, NA=0.4, WD=3.1 mm
 - EVIDENT, model : UCPLFLN20X : 20x, NA=0.7, WD=0.8-1.8 mm

【Procedure】

(Fixation)

A case for making mouse brain transparent according to reference (1).

- 1) Transcardially perfuse an anesthetised mouse with 4% paraformaldehyde (PFA)/PBS (pH 7.6 ~ 7.8).
- 2) Remove the whole brain and subject it to post-fixation in 4% PFA/PBS at 4℃ for 72 hours.
- 3) Wash the samples in PBS.
- 4) Prepare vibratome sections (optional, 0.2-3 mm thick)

(Clearing) (1-2 mm brain samples)

- 5) Transfer the sample into 5 mL of SCALEVIEW-S0 solution and incubate it with gentle shaking (10 rpm) at 35-40℃ for 30 min.
- 6) Transfer the sample into 5 mL of SCALEVIEW-S1 solution and incubate it with gentle shaking (10 rpm) at 35-40℃ for 30 min.
- 7) Transfer the sample into 5 mL of SCALEVIEW-S2 solution and incubate it with gentle shaking (10 rpm) at 35-40℃ for 30 min.
- 8) Transfer the sample into 5 mL of SCALEVIEW-S3 solution and incubate it with gentle shaking (10 rpm) at 35-40℃ for 30 min.
- 9) Transfer the sample into 5 mL of deScale Solution and incubate it with gentle shaking (10 rpm) at 4℃ for 3 hours × 2 times.
- 10) Transfer the sample into 5 mL of SCALEVIEW-S4 solution and incubate it with gentle shaking (10 rpm) at 35-40℃ for 12-24 hours.
- 11) The transparency can be evaluated by eye at this stage. If the sample is successfully cleared, the brain sample should look amber under a light source.
- 12) Transfer the sample into 5 mL of SCALEVIEW-SMt solution and incubate it with gentle shaking (10 rpm) at 35-40℃ for 1 hour.


Figure 1. Transmittance images of 1 mm brain slice before and after clearing with SCALEVIEW-S.

(Imaging)

- 13) Mount onto an imaging chamber. Use appropriate amount of SCALEVIEW-SMt solution to mount the sample. Observe the SCALEVIEW-S cleared samples using fluorescence microscopy.


Figure 2. Confocal imaging of adult hemisphere brain (Thy1-YFP-H mouse).

Note :

- For thicker brain slices (>2 mm) or hemisphere mouse brain whole mouse brain, incubation time in each step can become longer and SCALEVIEW-S solutions in each step can increase amount (e.g., hemisphere brain 25 mL, whole brain 40 mL each).
- SCALEVIEW-SMt mounted samples will not be solidified. If necessary for each sample, 1.5% melting agarose can be added to minimize movement artifacts during imaging.
- You can minimize movement artifacts even more by fixing the agarose gel and the observation container in place with super glue.

[Related products]

Code No.	Description	Size
193-18455	SCALEVIEW-A2	500 mL
196-18521	SCALEVIEW-S0	250 mL
194-18561	SCALEVIEW-S4	250 mL
041-34425	deScale Solution	500 mL
293-97001	AbScale-G Solution Kit	1 kit
294-35631	See Through Chamber, 0.3 mm thick	10 set
291-35641	See Through Chamber, 0.5 mm thick	10 set
295-35661	See Through Chamber, 1.0 mm thick	10 set
292-35671	See Through Chamber, 2.0 mm thick	10 set
299-35681	See Through Chamber, 3.0 mm thick	10 set

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