LBIS™ Rat Albumin ELISA Kit

Please, read this instruction carefully before use.

1. Intended use

LBIS™ Rat Albumin ELISA Kit is a sandwich ELISA system for quantitative measurement of rat albumin. This kit is intended for research use only.

2. Introduction

Albumin is mostly a simple hydrophilic protein present in cells and body fluids. Albumin is synthesized in the liver, and serum albumin (M.W. 69,000, pI 4.9) occupies 56% - 60% of total serum proteins. Because of its large population, albumin is very important in maintaining plasma osmotic pressure. Albumin can bind hydrophobic physiological substances e.g. fatty acids, bilirubin, and thyroxine and contributes the transfer of these substances. Concentration of albumin in serum is lowered in liver cirrhosis, malnutrition, and pyrexial diseases due to decreased biosynthesis or increased consumption of albumin in, and also by secretion into urine in renal damage. In normal human subjects, excretion of albumin into urine is very little, about 30 mg/day, but increased in glomerulonephritis, nephritic syndrome, and diabetic nephropathy. Urinary albumin sometimes increases in pyrexia, hypertension, congestive heart failure (CHF), and urinary tract infection (UTI). Even in healthy human, a transient upraise of urinary albumin is observed after hard exercise, muscle work, bathing in high temperature water, mental excitement, stress, intake of much protein, and before menstruation. These are called physiological or functional or sports proteinuria (albuminuria). Orthostatic proteinuria (albuminuria) is observed in teenagers. Rare hereditary analbuminemia in human is really albumin deficiency. A model analbuminemia rat has been established by Dr. Sumi Nagase which has been derived from Sprague-Dawley strain, and is called Nagase analbuminemia rat (NAR). Routine measurement of serum albumin is conveniently made with FUJIFILM Wako Pure Chemical Corporation's TIA assay kits for automatic analyzers. FUJIFILM Wako Pure Chemical Corporation's Albumin ELISA Kits can measure urinary albumin with high sensitivity, and are also applied to *in vitro* albumin biosynthesis system, checking of contamination with albumin in biological active substance preparations obtained by culture system, and liver transplantation experiments using NAR.

3. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Rat Albumin ELISA Kit, standards or diluted samples are incubated in monoclonal anti-albumin antibody-coated wells to capture albumin. After 1 hour incubation and washing, HRP (horse radish peroxidase)-labeled anti-albumin antibody is added, and incubated for 1 hour. After washing, HRP-complex remaining in wells is reacted with a TMB Solution for 20 minutes, and reaction is stopped by addition of acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to albumin concentration. The standard curve is prepared by plotting absorbance against standard albumin concentrations. Albumin concentrations in unknown samples are determined using this standard curve.

4. Performance Characteristics

• Assay range

The assay range of the kit is 50 ng/mL - 1000 ng/mL.

• Specificity

The antibodies used in this kit are specific to albumin.

Sample	Cross reaction	Sample	Cross reaction
Mouse albumin	Less than 5%	Bovine albumin	N.D.
Human albumin	N.D.	10% FCS	N.D.
Porcine albumin	N.D.		

*Cross reaction at 10000 ng/mL (Except FCS)

• Precision of assay

Within assay variation (3 samples, 5 replicates assay), Mean CV was within 10%

• Reproducibility

Between assay variation (3 samples, 3 days, 3 replicates assay), Mean CV was within 10%

• Recovery test

Standard albumin was added in 3 concentrations to serum sample and assayed.

The recoveries were 93.7% - 98.9%

- Dilution test
Serum sample was serially diluted by 3 steps.
The dilution curves showed linearity with $R^2 = 0.9951$ and 0.9984 .

5. Reference Assay Data

Rat albumin mean assay value : 36.0 mg/mL, SD : 4.45 mg/mL

Strain : CD(SD), 4 males, 10 week-old, serum

These data should be considered as guidance only. Each laboratory should establish its own normal and pathological reference ranges for albumin levels independently.

6. Precautions

- For professional use only. Beginners are advised to use this kit under the guidance of experienced person.
- Do not drink, eat or smoke in the areas where assays are carried out.
- This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. Especially be careful for the stop solution because it is sulfuric acid. The stop solution and the substrate solution may cause skin/eyes irritation. In case of contact with these wash skin/eyes thoroughly with water and seek medical attention, when necessary.
- Avoid contact with the acidic stop solution and TMB Solution containing hydrogen peroxide and tetramethylbenzidine. Wear gloves and eye and clothing protection when handling these reagents.
- The materials must not be pipetted by mouth.
- Unused samples and used tips should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by an autoclave before disposal.
- Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- Use clean laboratory glassware.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with a plate seal supplied, during incubation.
- ELISA can be easily affected by your laboratory environment. Room temperature should be at 20°C - 25°C strictly. Avoid airstream velocity over 0.4 m/sec. (including wind from air conditioner), and humidity less than 30%.

7. Reagents supplied

Components	Use Status	Amount	
		96 tests	96 tests × 2
(A) Antibody-coated Plate	Use after washing	96 wells/1 plate	96 wells/2 plates
(B) Rat Albumin Standard	Concentrated. Use after dilution.	150 μ L/1 bottle	150 μ L/2 bottles
(C) Buffer	Ready for use.	60 mL/1 bottle	120 mL/1 bottle
(D) Peroxidase-conjugated Antibody Solution	Concentrated. Use after dilution.	100 μ L/1 bottle	200 μ L/1 bottle
(F) TMB Solution	Ready for use.	12 mL/1 bottle	24 mL/1 bottle
(H) Stop Solution	Ready for use.	12 mL/1 bottle	24 mL/1 bottle
(I) Wash Solution (10×)	Concentrated. Use after dilution.	100 mL/1 bottle	100 mL/2 bottles
(J) Plate Seal	—	3 sheets	6 sheets

[Storage and Stability]

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip-seal originally used for well-plate container and store at 2°C - 10°C.

[(B) Rat Albumin Standard]

Standard solutions prepared above should be used as soon as possible, and should not be stored.

Dispose remaining prepared solution.

[(C) Buffer] and [(F) TMB Solution]

Use only volume you need for your assay. Remaining reagents should be stored at 2°C - 10°C fastening the cap tightly. It maintains stability until expiration date. Once opened, we recommend using as soon as possible to avoid influence by

environmental condition.

[(D) Peroxidase-conjugated Antibody Solution]

Unused working solution (already diluted) should be disposed.

[(H) Stop Solution]

Close the stopper tightly and store at 2°C - 10°C.

[(I) Wash Solution (10×)]

If not opened, store at 2°C - 10°C. Dispose any unused diluted buffer.

8. Equipments or supplies required but not supplied ☐ Use as a check box

☐ Purified water (distilled water). ☐ Test tubes for preparation of standard solution series.

☐ Glassware for dilution of Wash Solution (10×) (a graduated cylinder, a bottle).

☐ Pipettes (disposable tip type). One should be able to deliver 5 μ L precisely, and another for 50 μ L - 500 μ L. ☐ Syringe-type repeating dispenser like Eppendorf multipette plus which can dispense 50 μ L. ☐ Paper towel to remove washing buffer remaining in wells. ☐ A vortex-type mixer. ☐ A shaker for 96 well-plate (600 rpm - 1200 rpm). ☐ An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. ☐ A 96 well-plate reader (450 nm \pm 10 nm, 620 nm : 600 nm - 650 nm). ☐ Software for data analysis.

9. Preparation of Samples

This kit is intended to measure albumin in rat serum, plasma (heparin is recommended) or urine. The necessary sample volume for the standard procedure is 5 μ L. Dilute your samples with the buffer so as to be within the assay range, 50 ng/mL - 1000 ng/mL. The recommended dilution is 10000× - 50000× for serum or plasma, and 50× - 100× for urinary samples.

Samples should be immediately assayed or stored below -35°C until assay. Before starting assay, shake thawed samples sufficiently. Do not repeat freeze-and-thaw cycles.

Hemolytic and hyperlipemic serum samples are not suitable.

***To avoid influence of blood (high lipid or hemolysis, etc.), if your original samples have heavy chyle or hemolysis as the pictures below, do not use them for assay. Abnormal value might be obtained with hemolysis above 80 mg/dL with this kit.**

If presence of interfering substance is suspected, examine by dilution test at more than 2 points. Dilution of a sample should be made in a test tube using buffer solution prior to adding them to wells.

Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particulate matter.

10. Preparation of Reagents

◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.

◆ Prepare reagent solutions in appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Rat Albumin Standard]

Below is an example of preparing each standard solution.

Volume of standard solution	(C) Buffer	Concentration (ng/mL)
Original solution 50 μ L	450 μ L	1000
1000 ng/mL solution 400 μ L	100 μ L	800
800 ng/mL solution 300 μ L	100 μ L	600
600 ng/mL solution 200 μ L	100 μ L	400
400 ng/mL solution 100 μ L	100 μ L	200
200 ng/mL solution 100 μ L	100 μ L	100
100 ng/mL solution 100 μ L	100 μ L	50
Blank	100 μ L	0

[(D) Peroxidase-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the (C) Buffer to **1 : 100**.

[(I) Wash Solution (10×)]

Dilute 1 volume of the concentrated Wash Solution (10×) to **10 volume** with deionized water (or distilled water) to prepare working solution. Example : 100 mL of concentrated Wash Solution (10×) and 900 mL of deionized water (or distilled water).

11. Assay Procedure

Remove the cover sheet of the Antibody-coated Plate after bringing up to room temperature.

- (1) Wash the Antibody-coated Plate (A) by filling the wells with 300 μ L of washing buffer and discard 3 times (*①), then strike the plate upside-down onto folded several sheets of paper towel to remove residual buffer in the wells(*⑤).

- (2) Pipette 50 μ L of (C) Buffer to each well and shake the plate gently on a plate shaker (*②).
 - (3) Pipette 5 μ L of standard solution to the wells designated for standards, and 5 μ L of properly diluted samples to the designated sample wells.
 - (4) Shake the plate gently on a plate shaker (*②).
 - (5) Stick a Plate Seal (*③) on the plate and incubate for 1 hour at 20°C - 25°C.
 - (6) Discard the reaction mixture and rinse wells as step (1).
 - (7) Pipette 50 μ L of Peroxidase-conjugated Antibody Solution to all wells, and shake as step (4).
 - (8) Stick a Plate Seal (*③) on the plate and incubate the plate for 1 hour at 20°C - 25°C.
 - (9) Discard the reaction mixture and rinse wells as step (1).
 - (10) Pipette 50 μ L of TMB Solution to wells, and shake as step (4).
 - (11) Stick a Plate Seal (*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
 - (12) Add 50 μ L of the Stop Solution to all wells to stop the coloration.
 - (13) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (*④)) immediately using a plate reader.
- *Refer 15. Summary of Assay Procedure for notes of *①, *②, *③, *④, and *⑤.


12. Technical Tips

- In manual operation, proficiency in pipetting technique is recommended.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot number is the same, it is best not to mix the reagents with those that have been preserved for some period.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- Optimally, the reagent solutions of the kit should be used immediately after reconstitution. Otherwise, store them in a dark place at 2°C - 10°C.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- Dilution of the assay sample must be carried out using the buffer solution provided in the kit.
- The TMB Solution should be almost clear pale blue before use. It turns blue during reaction, and gives yellowish color after addition of Stop Solution. Greenish color means incomplete mixing.
- To avoid denaturation of the coated antibody, do not let the plate go dry.
- As the Antibody-coated Plate is module type of 8 wells \times 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- When ELISA has to be done under the airstream velocity over 0.4 m/sec. and the humidity less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a styrofoam box in each step of incubation.

13. Calculations

- (1) Prepare a standard curve by plotting standard concentration on X-log and absorbance on Y-axis.
- (2) Using the standard curve, read the albumin concentration of a sample at its absorbance*, and multiply the assay value by dilution factor if the sample has been diluted. Though the assay range is wide enough, in case the absorbance of some samples is higher than that of the highest standard, please repeat the assay after proper dilution of samples with the buffer solution.

*We recommend the use of 3rd order regression curve for log-log plot, or 4 or 5 parameters method for log-normal plot in computer calculation.

Rat Albumin assay standard curve (an example)

Absorbance may change due to assay environment.

14. Trouble Shooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not be added.
 - 2) Reagents necessary for coloration such as Peroxidase-conjugated Antibody Solution or TMB Solution might not be added.
 - 3) Wrong reagents related to coloration might have been added. Wrong dilution of Peroxidase-conjugated Antibody Solution.
 - 4) Contamination of enzyme inhibitor (s).
 - 5) Influence of the temperature under which the kits had been stored.
 - 6) Excessive hard washing of the well plate.
 - 7) Addition of TMB Solution soon after taking out from a refrigerator might cause poor coloration owing to low temperature.
- Blank OD was higher than that of the lowest standard concentration (50 ng/mL).

Possible explanations :

Improper or inadequate washing. (Change washing frequency from 3 times to 4 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Antibody Solution.)

- High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the other testing?

A-1 : Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator soon.

- Q-2 : I found there contains liquid in 96 well-plate when I opened the box. What is it?

A-2 : When we manufacture 96 well-plate, we insert preservation stabilizer in wells.


15. Summary of Assay Procedure ☐ : Use as a check box

*First, read this instruction manual carefully and start your assay after confirmation of details.

☐ Bring the well-plate and all reagents back to **20°C - 25°C for 2 hours**.

☐ Wash Solution (10×) concentrate must be diluted to **10 times** by deionized water (or distilled water) that returned to 20°C - 25°C.

☐ Standard solution dilution example :

- ☐ Antibody-coated Plate
- ☐ ↓ Washing 3 times (*①), (*⑤)
- ☐ Buffer 50 μ L
- ☐ ↓ Shaking (*②)
- ☐ (Diluted) Samples, or Standards 5 μ L
- ☐ ↓ Shaking (*②), Incubation for 1 hour at 20°C - 25°C. (Standing (*③))
- ☐ Dilute (D) Peroxidase-conjugated Antibody Solution to **1 : 100** with (C) Buffer returned to 20°C - 25°C.
- ☐ ↓ Washing 3 times (*①), (*⑤)
- ☐ Peroxidase-conjugated Antibody Solution (Diluted) 50 μ L
- ☐ ↓ Shaking (*②), Incubation for 1 hour at 20°C - 25°C. (Standing (*③))
- ☐ ↓ Washing 3 times (*①), (*⑤)
- ☐ TMB Solution 50 μ L
(After dispense, the color turns to blue depending on the concentration.)
- ☐ ↓ Shaking (*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (*③))
- ☐ Stop Solution 50 μ L
(After dispense, the color turns to yellow depending on the concentration.)
- ☐ ↓ Shaking (*②) (Immediately shake.)
- ☐ Measurement of absorbance (450 nm, Ref 620 nm (*④)) immediately
(Ref. wave cancels the dirt in the back of plate.)

*① After dispensing wash buffer to wells, lightly shake the plate on your palm for 10 seconds and remove the buffer. Guideline of washing volume : 300 μ L/well for an automatic washer and for a pipette if the washing buffer is added by pipette. In case of washing by using 8 channel pipette, sometimes the back ground tends to be high. If so, change washing frequency from 3 times to 4 - 6 times at the constant stroke after the reaction with Peroxidase-conjugated Antibody Solution.

Standard of plate-washing pressure : 5 mL/min - 25 mL/min. (Adjust it depending on the nozzle's diameter.)

*② Guideline of shaking : 600 rpm - 1200 rpm 3 times for 10 seconds.

*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal used once.


*④ 600 nm - 650 nm can be used as reference wavelength.

*⑤ After removal of wash buffer, immediately dispense the next reagent.

Worksheet Example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	1000 ng/mL	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
B	800 ng/mL	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
C	600 ng/mL	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
D	400 ng/mL	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
E	200 ng/mL	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
F	100 ng/mL	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
G	50 ng/mL	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
H	0 (Blank)	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40

Assay Worksheet

16. Storage and Expiration

When the complete kit is stored at 2°C - 10°C. Opened reagents should be used as soon as possible to avoid loss in optimal assay performance caused by storage environment.

LBIS™ Rat Albumin ELISA Kit

[Storage]	Store at 2 - 10°C
[Expiration date]	Indicated on the label
[Package]	For 96 tests / 96 tests × 2
[Cat #]	293-92501 (96 tests) 299-92503 (96 tests × 2)

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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