LabAssay™ NEFA (FFA)

Please, read this instruction carefully before use.

1. Introduction

Non-esterified fatty acid (NEFA) in serum binding albumin, is used as an important energy source of peripheral tissues. LabAssay™ NEFA is the reagent kit for assay of NEFA based on an enzymatic method using 3-methyl-*N*-ethyl-*N*-(β -hydroxyethyl)-aniline (MEHA) as a violet color agent. The kit is used for the quantitative determination of NEFA in mouse serum. It is a simultaneous multi-sample assay format using a microplate, but measurements can also be made using a test tube.

2. Kit Contents

	Components	State	Amount
(1)	Chromogen Reagent A	Freeze dry	For 10 mL/4 bottles
(2)	Solvent A (for Chromogen Reagent A)	Liquid	45 mL/1 bottle
(3)	Chromogen Reagent B	Freeze dry	For 20 mL/4 bottles
(4)	Solvent B (for Chromogen Reagent B)	Liquid	90 mL/1 bottle
(5)	Standard Solution	Liquid	7 mL/1 bottle


※ Assay in a microplate ; 500 tests. Assay in a test tube ; 40 tests.

3. Storage and Expiration

When the complete kit is stored at 2°C - 10°C, the kit is stable until the expiration date shown on the label on the container.

4. Assay Principle

Non-esterified fatty acid (NEFA) in the sample is converted to Acyl-CoA, AMP and Pyrophosphoric acid (PPi) by the action of Acyl-CoA synthetase (ACS), under coexistence with Coenzyme A (CoA) and Adenosine 5'-triphosphate disodium salt (ATP). The obtained Acyl-CoA is oxidized and yields 2,3-*trans*-enoyl-CoA and hydrogen peroxide by the action of Acyl-CoA oxidase (ACOD). In the presence of peroxidase (POD), the hydrogen peroxide yields a blue purple pigment by quantitative oxidation condensation with 3-methyl-*N*-ethyl-*N*-(β -hydroxyethyl)-aniline (MEHA) and 4-aminoantipyrine. NEFA concentration is obtained by measuring absorbance of the blue purple pigment.

5. Materials and Apparatuses to be Prepared

(For Microplate Method)

- 96 wells microplate (transparent type) · Micropipette · Incubator maintaining at 37°C
- Plate mixer
- Microplate reader with a 550 nm wavelength filter

(For Test Tube Method)

- Test tube or microtube · Pipette · Incubator maintaining at 37°C
- Spectrophotometer or colorimeter with a 550 nm wavelength filter

6. Preparation of Reagents to be used

① Color Reagent A :

Prepare Reagent A by mixing 1 vial of each Chromogen Reagent A (for 10 mL) and Solvent A (for Chromogen Reagent A) (10 mL). After reconstitution, the solution should be stored below 25°C and used the same day, or stored at 2°C - 10°C

and used within 5 days.

② Color Reagent B :

Prepare Reagent B by mixing 1 vial of each Chromogen Reagent B (for 20 mL) and Solvent B (for Chromogen Reagent B) (20 mL). After reconstitution, the solution should be stored below 25°C and used the same day, or stored at 2°C - 10°C and used within 5 days.

③ Standard solution (Microplate method)

Standard solution is prepared by dilution of the provided Standard Solution.

No.	Standard Solution	Distilled or deionized water	Sample volume	Concentration
1	80 μ L	120 μ L	4 μ L	0.40 mEq/L
2	120 μ L	80 μ L	4 μ L	0.60 mEq/L
3	160 μ L	40 μ L	4 μ L	0.80 mEq/L
4	Undiluted	—	4 μ L	1.00 mEq/L
5	Undiluted	—	8 μ L	1.97 mEq/L ^{*1}

*1 The test sample volume is usually 4 μ L, but 8 μ L is taken in this case. The NEFA concentration must be corrected accordingly as indicated in the table above.

7. Preparation of Samples

- Samples, especially mouse serum, should be used immediately after collecting, because enzymes such as lipoprotein lipase, phospholipase *etc.* hydrolyze lipids and form fatty acids. Freeze sample, when a serum is stored.
- Hemolysis gives a higher value on the assay.
- Bilirubin gives a lower value on the assay.
- Ascorbic acid does not have significant effects on the assay.

8. Procedure

(1) Assay in a Microplate

Perform the assay in the wells according to the following table scheme.

	Test	Standard	Reagent Blank	Sample Blank ^{*2}
Color Reagent A	80 μ L	80 μ L	80 μ L	80 μ L
Specimen	Serum 4 μ L	Standard 4 μ L	Distilled or Deionized water 4 μ L	—
Mix well ^{*3} and incubate at 37°C for 10 min.				
Color Reagent B	160 μ L	160 μ L	160 μ L	160 μ L
Specimen	—	—	—	Serum 4 μ L
Mix well ^{*3} and incubate at 37°C for 10min.				
After cooling the reaction solution to room temperature, measure the absorbance at 550 nm of the test sample and standard solution with the blank solution as the control within the next 30 min.				

*2 Sample blank is not required for common samples, but required for chyle or hemolyzed sample.

*3 Mix well with a plate mixer at 600 - 1,000 rpm for 1 minute. Be careful not to spill from the well.


(2) Assay in a Test Tube

Perform the assay in the test tubes of the microtubes according to the following table scheme.

	Test	Standard	Reagent Blank	Sample Blank ^{*2}
Color Reagent A	1.0 mL	1.0 mL	1.0 mL	1.0 mL
Specimen	Serum 0.05 mL	Standard 0.05 mL	Distilled or Deionized water 0.05 mL	—
Mix well and incubate at 37°C for 10 min.				
Color Reagent B	2.0 mL	2.0 mL	2.0 mL	2.0 mL
Specimen	—	—	—	Serum 0.05 mL
Mix well and incubate at 37°C for 10 min.				
After cooling the reaction solution to room temperature, measure the absorbance of the test sample and standard solution with the blank solution as the control within the next 30 min.				
Colorimeter with 550 nm filter				
Spectrophotometer : 550 nm				

*2 Sample blank is not required for common samples, but required for chyle or hemolyzed sample.

9. Standard Curve

10. Performance

(1) Sensitivity

- The absorbance is below 0.07, when measuring purified water as a specimen.
- The absorbance is 0.10 to 0.37, when measuring 1 mEq/L NEFA as a specimen.

(2) Specificity

The NEFA concentration is less than $\pm 15\%$, when measuring the known concentration of control serum as a sample.

11. Usage Notes

(1) Notes on the assay

- Do not use the reagents described above in any procedures other than those described herein. Performance cannot be guaranteed if the reagents are used in other procedures or for other purposes.
- Operate the measurement apparatuses according to operator's manuals under appropriate conditions. Consult the apparatus manufacturer for details.
- Keep the provided reagents under the indicated conditions before the expiration.
- Do not use reagents, which were frozen in error. Such reagents may give false results.
- After opening the reagents, it is recommended to use them immediately. The opened vials should be capped and kept under the specified conditions.
- Do not use the containers and other materials in the package for any purpose other than those described herein.
- The vial is stoppered at reduced pressure. Slowly remove the stopper in order not to blow the powder in the vial.
- Do not use the reagents at other reaction temperatures and reaction times than described herein.
- The reaction induced by Color Reagent A is completed in about 6 min at 37°C. Do not incubate more than 15 min.
- The reaction induced by Color Reagent B is completed in about 5 min at 37°C. When the incubation at 37°C is continued, absorbance decreases incrementally over time. Therefore, return it to room temperature immediately after the incubation for 10 min at 37°C.
- Avoid the direct sunlight during operation.
- This kit is for research use only. Not for diagnostic use.

(2) Safety Precautions

- If a reagent comes into contact with the mouth, eyes, or skin, immediately wash with a lot of water. Consult a physician if necessary.
- A pipette with a safety pipette filler should be used.
- Be careful not to cut yourself with the aluminum cap when removing it from the vial.

(3) Waste

- The waste should be processed appropriately according to regulations.
- All the devices including reagents and vials that come into contact with the specimen should be considered potentially infectious.
- Chromogen Reagent A and Solvent A (for Chromogen Reagent A) contain sodium azide as preservative (0.07% in Color A when reconstituted). Standard Solution contains 0.05% sodium azide. Sodium azide may react with lead or

copper plumbing to form explosive compounds. Even though the reagent contains minute quantity of sodium azide, drains should be flushed well with a large amount of water, when discarding the reagents.

LabAssayTM NEFA (FFA)

[Storage]

Store at 2°C - 10°C

[Expiration date]

Indicated on the container

[Package]

500 tests

[Cat #]

299-94301

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://ffwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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