LBIS™ Mouse IL-6 ELISA Kit

Please, read this instruction carefully before use.

1. Intended use

Mouse IL-6 is a secreted glycoprotein of 187 amino acids and it was discovered to be a cytokine that differentiates B cells into antibody-producing cells. It is also known to play an important role in various physiological processes such as immune response, inflammatory response, hematopoiesis, and the proliferation and differentiation of nervous system cells. After being produced by T cells, monocytes, fibroblasts, endothelial cells, keratinocytes, etc., it binds to the IL-6 receptor, which consists of two chains of IL-6R α and gp130, and signals are transduced. There is also a study that suggested that a possible correlation exists between IL-6 levels and disease activity in rheumatoid arthritis, and it is also attracting attention in the field of autoimmune diseases such as rheumatoid arthritis and inflammatory diseases.

This kit is a sandwich ELISA system for quantitative measurement of mouse IL-6. It can measure IL-6 in mouse serum (plasma) and culture supernatant, specifically and with high sensitivity. This product is intended for research use only. Not for use in diagnostics procedures.

2. Assay Principle

In FUJIFILM Wako Pure Chemical Corporation's LBIS™ Mouse IL-6 ELISA Kit, standards or samples are incubated in monoclonal antibody-coated wells to capture mouse IL-6. After 2 hours' incubation and washing, biotin-conjugated antibody is added and incubated for a further 1 hour to bind with captured mouse IL-6. After washing, peroxidase-conjugated streptavidin is added, and incubated for 30 minutes. After washing, bound peroxidase-conjugated streptavidin is reacted with a chromogen (TMB) for 20 minutes, and the reaction is stopped by addition of an acidic solution, and absorbance of yellow product is measured spectrophotometrically at 450 nm. The absorbance is nearly proportional to mouse IL-6 concentration. The standard curve is prepared by plotting absorbance against standard mouse IL-6 concentrations. Mouse IL-6 concentrations in unknown samples are determined using this standard curve.

3. Performance Characteristics

• Assay Range

2.05 pg/mL - 500 pg/mL

• Precision of assay (Within assay variation)

2 serum samples, 5 replicates assay

| n/ID | Sample1 (pg/mL) | Sample2 (pg/mL) |
|--------|--------------------|--------------------|
| 1 | 281 | 18.7 |
| 2 | 285 | 19.1 |
| 3 | 281 | 18.7 |
| 4 | 282 | 18.7 |
| 5 | 291 | 19.6 |
| mean | 284 | 19.0 |
| SD | 4.18 | 0.390 |
| CV (%) | 1.47 | 2.06 |

• Reproducibility (Between assay variation)

3 serum samples, 4 days, 3 replicates assay

| Day/ID | Sample1 (pg/mL) | Sample2 (pg/mL) | Sample3 (pg/mL) |
|--------|--------------------|--------------------|--------------------|
| 0 | 246 | 84.0 | 15.1 |
| 1 | 230 | 78.8 | 13.9 |
| 2 | 254 | 86.2 | 15.3 |
| 3 | 244 | 84.9 | 14.4 |
| mean | 244 | 83.5 | 14.7 |
| SD | 9.90 | 3.23 | 0.658 |
| CV (%) | 4.07 | 3.87 | 4.49 |

• Recovery test (Standard mouse IL-6 was added in 4 concentrations to a serum or plasma sample and assayed. The recoveries were 88.4% - 96.9%).


Serum sample

| Added (pg/mL) | Found (pg/mL) | Recovered (pg/mL) | Recovery (%) |
|------------------|------------------|----------------------|-----------------|
| 0 | 4.8 | — | — |
| 22.8 | 26.9 | 22.1 | 96.9 |
| 46.0 | 49.0 | 44.2 | 96.0 |
| 96.4 | 94.9 | 90.1 | 93.4 |
| 197 | 193 | 188 | 95.5 |

Plasma sample (EDTA)

| Added (pg/mL) | Found (pg/mL) | Recovered (pg/mL) | Recovery (%) |
|------------------|------------------|----------------------|-----------------|
| 0 | 4.4 | — | — |
| 28.1 | 30.6 | 26.2 | 93.0 |
| 56.0 | 55.0 | 50.5 | 90.2 |
| 116 | 107 | 102 | 88.4 |
| 227 | 206 | 202 | 88.8 |

- Dilution test (Serum sample and plasma sample underwent 3-step serial dilution.)

4. Technical tips/Precautions

- ELISA can be affected by the assay environment. Ensure that the room temperature at the places where the assay and incubation are performed is strictly controlled at 20°C - 25°C. Avoid performing assay under an airstream where the velocity exceeds 0.4 m/sec. and the humidity is less than 30%, seal the well plate with a plate seal and place the well plate in an incubator or a Styrofoam box in each incubation step.
- To avoid denaturation of the coated antibody, do not allow the plate dry out.
- In order to avoid dryness of wells, contamination by foreign substances and evaporation of dispensed reagents, never forget to cover the well plate with the supplied plate seal, during incubation.
- Be careful to avoid any contamination of assay samples and reagents. We recommend the use of disposal pipette tips, and 1 tip for 1 well.
- As the antibody-coated plate is a module type consisting of 8 wells × 12 strips, each strip can be separated by cutting the cover sheet with a knife and used independently.
- The chromogen (TMB) should be almost clear to pale yellow before use. It turns blue during the reaction, and gives a yellowish color after addition of the stop solution. A greenish color indicates incomplete mixing.
- Time the reaction from the pipetting of the reagent to the first well.
- Prepare a standard curve for each assay.
- For professional use only. Novices are advised to use this kit under the guidance of an experienced person. In manual operation, proficiency in pipetting technique is recommended.
- Wear protective gloves as well as protective clothing, eyewear, and face protection.
- Avoid kit reagents or specimens from coming into contact with the skin and mucous membranes. If any reagents come into contact with the eyes, skin, or mucous membranes, wash with copious amounts of water and contact a physician.
- Do not drink, eat, or smoke in the places where this kit is being used.
- The materials must not be pipetted by mouth.
- Do not combine reagents with different lot numbers.
- The reagents are prepared to give accurate results only when used in combination within the same box. Therefore, do not combine the reagents from kits with different lot numbers. Even if the lot numbers are the same, it is preferable not to mix reagents with those that have been retained for some period.
- Handle the sample carefully, in recognition of the fact that the sample may be associated with a risk of infection. This kit contains animal-derived ingredients.
- Residual samples and used tips should be sterilized before disposal.
- Dispose of consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.

5. Reagents supplied

| Components | Use Status | Amount |
|---|---|----------------------|
| (A) Antibody-coated Plate | Use after washing. | 96 wells/1 plate |
| (B) Mouse IL-6 Standard | Freeze-dried. Use after reconstitution. | 2 bottles |
| (C) Buffer | Ready for use. | 30 mL/1 bottle |
| (D) Biotin-conjugated Antibody Solution | Concentrated. Use after dilution. | 100 μ L/1 bottle |
| (E) Peroxidase-conjugated Streptavidin Solution | Concentrated. Use after dilution. | 100 μ L/1 bottle |
| (F) TMB Solution | Ready for use. | 12 mL/1 bottle |
| (G) Buffer for Standard and Sample | Ready for use. | 30 mL/1 bottle |
| (H) Stop Solution Be careful! | Ready for use. | 12 mL/1 bottle |
| (I) Wash Solution (10 \times) | Concentrated. Use after dilution. | 100 mL/1 bottle |
| (J) Plate seal | — | 4 sheets |

[Storage and stability]

[(A) Antibody-coated Plate]

If seal is not removed, put the strip back in a plastic bag with zip lock originally used as the well-plate container and store at 2°C - 10°C.

[(B) Mouse IL-6 Standard (Freeze-dried)]

The reconstituted standard solution (original standard solution) should be used within 2 weeks if stored in a refrigerator. Dispose of the remaining diluted standard solutions after use.

[(C) Buffer] & [(F) TMB Solution] & [(G) Buffer for Standard and Sample]

Use only the volume required to conduct the assay. Remaining reagents should be stored at 2°C - 10°C in a tightly capped container. Once opened, we recommend using as soon as possible to avoid being influenced by environmental conditions.

[(D) Biotin-conjugated Antibody Solution] & [(E) Peroxidase-conjugated Streptavidin Solution]

Remaining working solution (already diluted) should be disposed of. The rest of the undiluted solution (unused) : close the cap tightly and store at 2°C - 10°C. Once opened, we recommend using as soon as possible to avoid being influenced by environmental conditions.

[(H) Stop Solution]

Close the cap tightly and store at 2°C - 10°C.

[(I) Wash Solution (10 \times)]

The remaining undiluted solution (unused) : close the cap tightly and store at 2°C -10°C. Dispose of any remaining diluted buffer.

6. Equipment or supplies required but not supplied ☐ Use as a check box

☐ Deionized water (or distilled water) ☐ Test tubes for preparation of standard solution series. ☐ Glassware for dilution of (I) Wash Solution (10 \times) (a graduated cylinder, a bottle) ☐ Pipettes (disposable tip type). One pipette should be able to deliver 10 μ L - 100 μ L precisely, and another should be able to deliver 100 μ L - 1000 μ L. ☐ Syringe-type repeating dispenser like Eppendorf Multipette Plus which can dispense 100 μ L. ☐ Paper towel to remove washing buffer remaining in wells. ☐ A vortex-type mixer. ☐ A shaker for 96 well-plate (500 rpm -1200 rpm) ☐ An automatic washer for 96 well-plate (if available), or a wash bottle with a jet nozzle. ☐ A 96 well-plate reader (450 nm \pm 10 nm, 620 nm : 600 nm - 650 nm) ☐ Software for data analysis.

7. Preparation of Samples

- This kit is intended to measure IL-6 in mouse serum or plasma.
- Use EDTA as anticoagulant.
- Samples should be immediately assayed or stored below -35°C until assay. Before starting the assay, thoroughly shake thawed samples. Do not repeat freeze-and-thaw cycles.
- Please check in advance when using serum separation promoters, etc.
- Do not use hemolyzed samples or high-lipid samples.
- Centrifuge the sample to remove turbidity or insoluble matter where necessary before use in the assay.
- If presence of interfering substance is suspected, examine by dilution test at more than 2 points.
- Dilution of a sample should be made in a PP or PE test tube using (G) Buffer for Standard and Sample prior to adding them to wells.

8. Preparation of Reagents

- ◆ Bring all reagents of the kit to room temperature (20°C - 25°C) before use.
- ◆ Prepare reagent solutions in an appropriate volume for your assay. Do not store the diluted reagents.

[Concentrated reagents]

[(B) Mouse IL-6 Standard (Freeze-dried)]


Reconstitute the designated amount of purified water described in "Reconstitution of standard " to (B) Mouse IL-6 Standard, mix, and dissolve for 30 minutes to prepare the standard stock solution (5 ng/mL). Be sure to allow the solution to stand for 30 minutes, as accurate values may not be obtained.

Then prepare using the kit-attached (G) Buffer for Standard and Sample that has been allowed to warm up to room temperature.

*Please confirm "Reconstitution of standard" from this product page. Since the amount of purified water added varies by lot, please be sure to check each lot.

Below is an example of preparing each standard solution.

| Volume of standard solution | (G) Buffer for Standard and Sample | Concentration | Standard Curve |
|---------------------------------------|------------------------------------|---------------|----------------|
| Original standard solution 50 μ L | 450 μ L | 500 pg/mL | |
| 500 pg/mL solution 200 μ L | 300 μ L | 200 pg/mL | |
| 200 pg/mL solution 200 μ L | 300 μ L | 80 pg/mL | |
| 80 pg/mL solution 200 μ L | 300 μ L | 32 pg/mL | |
| 32 pg/mL solution 200 μ L | 300 μ L | 12.8 pg/mL | |
| 12.8 pg/mL solution 200 μ L | 300 μ L | 5.12 pg/mL | |
| 5.12 pg/mL solution 200 μ L | 300 μ L | 2.05 pg/mL | |
| Blank | 350 μ L | 0 | |

[(D) Biotin-conjugated Antibody Solution]

Prepare working solution by dilution of (D) with the (C) Buffer to **1 : 100**. 10 mL of the diluted solution is enough for 96 wells.

[(E) Peroxidase-conjugated Streptavidin Solution]

Prepare working solution by dilution of (E) with the (C) Buffer to **1 : 100**. 10 mL of the diluted solution is enough for 96 wells.

[(I) Wash Solution (10 \times)]

Dilute 1 volume of the (I) Wash Solution (10 \times) to **10% solution by volume** with deionized water (or distilled water) to prepare the washing solution.

Example : 100 mL of the (I) Wash Solution (10 \times) and 900 mL of deionized water (or distilled water).

Caution : Be sure to use the (I) Wash Solution (10 \times) included in this kit. Do not use the (I) Wash Solution (10 \times) included in another kit.


9. Assay Procedure

Remove the cover sheet of the (A) Antibody-coated Plate after bringing up to room temperature.

- (1) Wash the (A) Antibody-coated Plate by filling the wells with washing solution and discard 4 times (*①), then tap the plate upside-down onto several sheets of folded paper towels to remove residual washing solution in the wells (*⑤).
 - (2) Pipette **50 μ L** of standard solution to the wells designated for standards.
 - (3) Pipette **25 μ L** of (G) Buffer for Standard and Sample to the sample wells, and then add **25 μ L** of sample to the wells.
 - (4) Shake the plate gently on a plate shaker (*②).
 - (5) Affix a plate seal (*③) on the plate and incubate for 2 hours at 20°C - 25°C.
 - (6) Discard the reaction mixture and rinse the wells as described in step (1).
 - (7) Pipette **50 μ L** of biotin-conjugated antibody solution to all wells, and shake as described in step (4).
 - (8) Affix a plate seal (*③) on the plate and incubate the plate for 1 hour at 20°C - 25°C.
 - (9) Discard the reaction mixture and rinse wells as described in step (1).
 - (10) Pipette **50 μ L** of peroxidase-conjugated streptavidin solution to all wells, and shake as described in step (4).
 - (11) Affix a plate seal (*③) on the plate and incubate the plate for 30 minutes at 20°C - 25°C.
 - (12) Discard the reaction mixture and rinse wells as described in step (1).
 - (13) Pipette **50 μ L** of (F) TMB Solution to wells, and shake as described in step (4).
 - (14) Affix a plate seal (*③) on the plate and incubate the plate for 20 minutes at 20°C - 25°C.
 - (15) Add **50 μ L** of the (H) Stop Solution to all wells and shake as described in step (4).
 - (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm (*④)) immediately using a plate reader.
- *Refer to 13. Summary of Assay Procedure for notes on *①, *②, *③, *④, and *⑤.

10. Calculations

- (1) Create a standard curve by reducing the data using computer software capable of generating a 3rd order regression curve, or 4 or 5 parameters. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph. The data may be linearized by plotting the log of the mouse IL-6 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data.
- (2) Using the standard curve, read the mouse IL-6 concentration of samples at its absorbance, and multiply the assay value by the dilution factor if the sample has been diluted. The standard operation method is 5-fold dilution.

11. Calibration

The measured value obtained by this kit can be converted to the unit concentration based on NIBSC/WHO standard mouse IL-6 (Code : 93/730) by multiplying with the conversion factor specified in a separate sheet according to the following calculation formula.

NIBSC/WHO (93/730) unit concentration (U/mL) = conversion factor (*1) × measured mouse IL-6 value (pg/mL) (*2)

(*1) The exchange rate is different depending on the lot. Please see the specified conversion factor in the annex document.

(*2) Use the assay value that is multiplied by the dilution rate.

12. Troubleshooting

- Low absorbance in all wells

Possible explanations :

- 1) The standard or samples might not have been added.
- 2) Reagents necessary for coloration such as biotinylated anti-Mouse IL-6 antibody, HRP-conjugated streptavidin, or chromogen (TMB) might not have been added.
- 3) Wrong reagents related to coloration might have been added. Wrong dilution of biotinylated anti-Mouse IL-6 antibody or HRP-conjugated streptavidin.
- 4) Contamination of enzyme inhibitor.
- 5) Influence of the temperature under which the kits had been stored.
- 6) Excessively hard washing of the well plate.
- 7) Addition of chromogenic substrate reagent soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- Blank OD was higher than that of the lowest standard concentration (2.05 pg/mL).

Possible explanations :

Improper or inadequate washing. (Change washing frequency from 4 times to 5 - 8 times after the reaction with HRP-conjugated streptavidin.)

- High coefficient of variation (CV)

Possible explanation :

- 1) Improper or inadequate washing.
- 2) Improper mixing of standard or samples.
- 3) Pipetting at irregular intervals.

- Q-1 : Can I divide the plate to use it for the additional testing?

A-1 : Yes, cut off the clear seal on the plate with a cutter along the strip. Move the residual plate, which still has the seal on it, promptly to a refrigerator.

- Q-2 : I found that the 96 well-plate was empty when I opened the box.

A-2 : As this kit is the dried type, no preservation stabilizer is added.

- Q-3 : When I thawed the sample, there was some undissolved matter. Will this affect the measurement?

A-3 : There may be an impact. The measured value may be low or below the lower measurement limit.

- Q-4 : Can it be measured using plasma?

A-4 : It can also be measured in plasma. Use EDTA as an anticoagulant when collecting plasma.

13. Summary of Assay Procedure ☐ : Use as a check box

*First, read this instruction manual carefully and start your assay after confirming the details.

- ☐ Return the (A) Antibody-coated Plate and all reagents back to room temperature (20°C - 25°C).

This takes about **2 hours**.

- ☐ (I) Wash Solution (10×) must be diluted to **10 times** by deionized water (or distilled water) at room temperature (20°C - 25°C).

- ☐ Standard solution dilution example

Reconstitute the designated amount of purified water described in "Reconstitution of standard" * to (B) Mouse IL-6 Standard, mix, and dissolve for 30 minutes to prepare the standard stock solution (5 ng/mL). Be sure to allow the solution to stand for 30 minutes, as accurate values may not be obtained.

Then prepare using the kit-attached (G) Buffer for Standard and Sample that has been allowed to warm up to room temperature.

* Please confirm "Reconstitution of standard" from this product page. Since the amount of purified water added varies by lot, please be sure to check each lot.

Below is an example of preparing each standard solution.

| Conc.(pg/mL) | 500 | 200 | 80 | 32 | 12.8 | 5.12 | 2.05 | 0 |
|-------------------------------------|-----|------|------|------|------|------|------|-----|
| Std. sol.(μL) | 50 | 200* | 200* | 200* | 200* | 200* | 200 | — |
| Buffer for Standard and Sample (μL) | 450 | 300 | 300 | 300 | 300 | 300 | 300 | 350 |

*One rank higher standard.

- ☐ (A) Antibody-coated Plate
- ☐ ↓ Washing 4 times (*①), (*⑤)
- ☐ Samples (diluted sample : (G) Buffer for Standard and Sample **25 μL** + Sample **25 μL**), or Standards 50 μL
- ☐ ↓ Shaking (*②), Incubation for 2 hours at 20°C - 25°C. (Standing (*③))
- ☐ Dilute (D) Biotin-conjugated Antibody Solution to **1 : 100** with the (C) Buffer at room temperature (20°C - 25°C). (This should be prepared during incubation.)
- ☐ ↓ Washing 4 times (*①), (*⑤)
- ☐ Biotin-conjugated antibody solution (Diluted) 50 μL
- ☐ ↓ Shaking (*②), Incubation for 1 hour at 20°C - 25°C. (Standing (*③))
- ☐ Dilute (E) Peroxidase-conjugated Streptavidin Solution to **1 : 100** with the (C) Buffer at room temperature (20°C - 25°C). (This should be prepared during incubation.)
- ☐ ↓ Washing 4 times (*①), (*⑤)
- ☐ Peroxidase-conjugated streptavidin solution (Diluted) 50 μL
- ☐ ↓ Shaking (*②), Incubation for 30 minutes at 20°C - 25°C. (Standing (*③))
- ☐ ↓ Washing 4 times (*①), (*⑤)
- ☐ (F) TMB Solution 50 μL
- ☐ After dispensing, the color turns to blue depending on the concentration.
- ☐ ↓ Shaking (*②), Incubation for 20 minutes at 20°C - 25°C. (Standing (*③))
- ☐ (H) Stop Solution 50 μL
- ☐ After dispensing, the color turns to yellow depending on the concentration.
- ☐ ↓ Shaking (*②) Immediately shake.
- ☐ Measurement of absorbance (450 nm, Ref 620 nm (*④)) immediately.
Ref. wave cancels the dirt in the back of plate.

*① After dispensing wash buffer into the wells, lightly shake the plate on your palm for 10 sec and remove the buffer. Guideline for washing volume : 300 μL/well for an automatic washer and for a pipette if the washing buffer is added by pipette. Where washing is conducted using an 8 channel pipette, sometimes the background tends to be high. If so, change washing frequency from 4 times to 5 - 8 times after the reaction with peroxidase-conjugated streptavidin.

Standard of plate-washing pressure : 5 mL/min - 25 mL/min. (Adjust it depending on the nozzle's diameter.)


*② Guideline for shaking : 600 rpm - 1,200 rpm for 10 seconds × 3 times.

*③ Seal the plate during the reaction after shaking. Peel off the protective paper from the seal and affix the seal on the plate. Do not reuse a plate seal that has already been used once.

*④ 600 nm - 650 nm can be used as reference wavelength.

*⑤ After removal of wash buffer, immediately dispense the next reagent.

Assay Worksheet

14. Storage and Expiration

When the complete kit is stored at 2°C - 10°C (Do not freeze), the kit is stable until the expiration date shown on the label on the container. Opened reagents should be used as soon as possible to avoid less than optimal assay performance caused by storage environment.

LBIS™ Mouse IL-6 ELISA Kit

[Storage] Store at 2°C - 10°C (Do not freeze).

[Expiration date] Indicated on the label.

[Package] For 96 tests

[Cat #] 299-96001

FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan
Telephone : + 81-6-6203-3741
Facsimile : + 81-6-6201-5964
<http://fwk.fujifilm.co.jp>

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road
Richmond, VA 23237
U.S.A.
Telephone : + 1-804-271-7677
Facsimile : + 1-804-271-7791
<http://www.wakousa.com>

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