**FUJIFILM****Wako**

Code No. 299-97101 (1 kit)

**FM Fix (L)****【Introduction】**

Tissue fixation is an important procedure to protect biological specimens from degradation due to autolysis and decay, ensuring stable preservation of proteins and other components within the tissues. Formalin fixation is widely used in research and pathology, including microscopy specimen preparation, however, it can affect the steric structure and antigenicity of some proteins, leading to deterioration in staining quality and other properties.

FM Fix (L) is a fixative that has superior tissue staining properties while retaining tissue morphology and three-dimensional structure. It helps to maintain the three-dimensional structure of proteins even after tissue fixation, enabling the stable detection of antigens such as complements. In addition, it is characterized by minimal reduction of fluorescent protein signals such as GFP, making it suitable for tissue transparency techniques.

**【Features】**

- Superior tissue morphology maintenance and staining
- Fluorescent immunostaining of the complement provides staining images equivalent to or better than the acetone-fixation method after freezing
- Fluorescent protein signals such as GFP are less likely to weaken
- Applicable to SeeDB2, a tissue transparency method
- Allows detection of mRNA transcripts by in situ hybridization

**【Kit contents】**

This kit consists of two components.

Solution A (L) : 6 mL × 5 bottles

Solution B (L) : 42 mL × 5 bottles

**【Storage】**

Store at -20℃

**【Procedure】**

Note : Solution B (L) may turn slightly yellow, but this does not affect performance.

**1. Reagent Preparation**

Prepare the reagent just before use, and use the mixed solution on the same day.

- (1) Thaw each solution immediately before use.
- (2) After thawing, shake or stir gently to make sure there are no undissolved residues.
- (3) Add all of solution A (6 mL × 1 bottle) to solution B (42 mL × 1 bottle).
- (4) Close the container lid and shake gently to mix Solution A and Solution B.

— 1/4 —

**2. Fixation Procedure**

- (1) Immerse a tissue block (approx. 10 mL of mixed fixative solution for 1 g of tissue) in the prepared fixative solution for 24 to 48 hours.
  - \* Adjust the fixation time according to the size and type of tissue.
- (2) Follow the usual section preparation method according to the purpose of the experiment.

**【Precautions】**

When handling, wear a gas mask or hose mask, protective gloves, and glasses as necessary.

(For details, please refer to the handling precautions on the label.)

**FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://fwk.fujifilm.co.jp>

**FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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