



For rheumatology and tissue fibrosis research

Anti CTGF, Monoclonal Antibody

CTGF is a secreted protein of about 38 kDa produced from umbilical veins and vascular endothelial cells . CTGF is consist of four modules, and regulates chondrocyte proliferation / differentiation and intercellular adhesion. Some reports showed that CTGF is involved in various diseases such as cancer, renal failure, liver failure and rheumatism. This items are mouse monoclonal antibodies established by DNA immunization, specifically react with native form CTGF. it has neutralizing activity.

	Module1	Module2	Module3	Module4	
	CTGF				
Rheuma	tism Kid	lney failure	Liver fa	ilure Ca	ancer


Micamatism		Ridley failure =				
Epitope	Clone No.	Cross- reactivity	Subclass Applications		Concentration (1st lot)	
Module1	30D2	Human Mouse	Mouse IgG1	in vitro/ in vivo neutralization ELISA	0.8mg/mL	
Module2	2-3	Human, Mouse	Mouse IgG1	In vitro neutralization	0.9mg/mL	
Module3	3-54	Human, Mouse	Mouse IgG1	In vitro neutralization	1mg/mL	
Module4	4-69	Human, Mouse	Mouse IgG1	In vitro neutralization	0.7mg/mL	


Features


■ CTGF Neutralizing Activity Suppression of arthritis model animals (mice)


■ Each 4 clone recognizes each module 1, 2, 3 and 4, respectively.

in vitro neutralization assay - inhibition of osteoclastognenesis-










- ■Condition
 ①PBMCs are
 collected from whole
 human blood and
 CD14+ cell is isolated.
- ②CD14⁺ cells are cultured for 7days (37°C, 5%CO₂) with α-MEM (add 10% FBS) ,M-CSF (25ng/mL), sRANKL (40ng/mL) and anti-CTGF Abs.
- ③Osteoclasts (TRAP positive cells) were counted by optical microscopy.

Enhancement of osteoclastogenesis by CTGF were repressed all CTGF antibodies.

in vivo neutralization assay -suppression of arthritis model animals (mice)-

■ Anti CTGF Module1(30D2)



■ Condition

DBA/1J mice (Female, 7-8 weeks of age) were immunized intradermally at the base of the tail with 100 µL of emulsion (in which emulsified in equal amounts of bovine type II collagen and Freund's complete adjuvant (CFA)). The anti-CTGF antibody was intraperitoneally administered 1 week before the start of immunization, and administered once a week for a total of 8 times. Evaluation of limb arthritis was performed once a week until immunization for 8 weeks. (Refer to the table below for evaluation criteria)

Arthritis point criteria

point	Observation
О	No change
1	Regional reddening of the limbs, or swelling and reddening of one finger.
2	Swelling of limbs' joints, or swelling of two or more fingers
3	Medium – prominent swelling of hands and feet.
4	Remarkable swelling and reddening of limbs, or hardening.

* The above point is the total point of the limbs.

Mice treated with anti-CTGF module1 Ab significantly suppressed arthritis compared to the untreated mice with collagen induced arthritis.

This data was provided by Dr. Nozawa and Dr. Sekigawa of Institute for Environmental and Gender Specific Medicine, Juntendo University in Japan

8

7

Reference

15

Arthritis score (Max; 16)

9

0

3

- 1) Ikawa, Y., et al.: J. Cell Physiol., 216, 3, 680(2008).
- 2) Nozawa, K., et al.: Arthritis Rheum., 65, 6, 1477 (2013).

5 6 Week after immunization

- 3) Miyashita, T., et al.: Autoimmunity., 49, 2,109 (2016).
- 4) Miyazaki O, et al.: Ann. Clin. Biochem., 47, 205(2010).

Product Name	Wako Cat. No	Pkg. Size	Grade	Storage Condition
Anti CTGF Module1, Monoclonal Antibody(30D2)	012-27421 018-27423	10µL 50µL	For Immunochemistry	Keep at -20°C
Anti CTGF Module2, Monoclonal Antibody(2-3)	019-27431 015-27433	10µL 50µL	For Immunochemistry	Keep at -20°c
Anti CTGF Module3, Monoclonal Antibody(3-54)	016-27441 012-27443	10μL 50μL	For Immunochemistry	Keep at -20°C
Anti CTGF Module4, Monoclonal Antibody(4-69)	013-27451 019-27453	10µL 50µL	For Immunochemistry	Keep at -20°C

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: https://labchem-wako.fujifilm.com/ / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site https://labchem-wako.fujifilm.com



FUJIFILM Wako Pure Chemical Corporation 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Tel: +81 6 6203 3741 Fax: +81 6 6203 1999

ffwk-cservise@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road, Richmond, VA 23237, U.S.A. Toll-Free (U.S. only): +1 877 714 1920 Tel: +1 804 271 7677 Fax: +1 804 271 7791 wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH Fuggerstr 12, 41468 Neuss, Germany Tel: +49 2131 311 0 Fax: +49 2131 311 100 labchem_wkeu@fujifilm.com