(109 × 210mm Size)

**FUJIFILM**

**Wako**

Code No. 298-80101

## High Molecular Amyloid $\beta$ Oligomer ELISA Kit *Wako*

### [Introduction]

Amyloid  $\beta$  ( $A\beta$ ) is a peptide consisting of about 40 amino acids and accumulates in the brain of Alzheimer disease patients forming a structure known as senile plaques. Therefore,  $A\beta$  is considered as a causal factor of Alzheimer disease. In recent years, an  $A\beta$  oligomer hypothesis has been focused. The hypothesis is that soluble aggregates called  $A\beta$  oligomers resulting from polymerization of  $A\beta$  monomers, causes synaptic failure in nerve cells, not insolubilized and accumulated  $A\beta$  fibrils. There are several types of  $A\beta$  oligomers due to the difference in molecular weight, and low molecular  $A\beta$  oligomers (such as low-n oligomer) and high molecular  $A\beta$  oligomers (such as ADDL,  $A\beta$  \*56, proto-fibril) have been reported. This product is an ELISA kit that can specifically measure high molecular  $A\beta$  oligomers.

### [Performance]

Standard curve range	0.41 - 100 pM *Calculated on the basis of 16-mer MAP peptide
Reactive $A\beta$ oligomers	$\geq$ 9-mer $A\beta$ oligomer
Sample	Human cerebrospinal fluid in vitro $A\beta$ oligomer
Sample volume	25 $\mu$ L (4 times dilution)
Assay time	4.5 hours
Detection method	Luminescence method

### [Reagents supplied]


Components	Condition	Volume
(1) Antibody (BAN50)-coated Plate	Use as it is.	1 plate, 96 wells (8 × 12)
(2) Standard Solution (10 nM)	Use after preparation.	100 $\mu$ L/1 vial
(3) Buffer	Use as it is.	60 mL/1 bottle
(4) Biotin-conjugated Antibody (BAN50) Solution	Use after preparation.	100 $\mu$ L/1 vial
(5) Peroxidase-conjugated Streptavidin Solution	Use after preparation.	100 $\mu$ L/1 vial
(6) Luminescent Reagent 1	Use as it is.	6 mL/1 bottle
(7) Luminescent Reagent 2	Use as it is.	6 mL/1 bottle
(8) Wash Solution (10 × )	Use after preparation.	100 mL/1 bottle
(9) Plate Seal	Use as it is.	3 sheets
Instruction manual		1 copy

### [Measurement principle]

The microplate is coated with anti  $A\beta$  antibody (BAN50). The standard solution or a sample is added in a well to bind  $A\beta$  monomer and  $A\beta$  oligomer to the antibody. Then, the biotin-conjugated anti  $A\beta$  monoclonal antibody (BAN50) is allowed to react to bind only to the high molecular  $A\beta$  oligomer in the well.

In addition, the peroxidase-conjugated streptavidin is reacted. Finally, the peroxidase activity in the well is measured to determine the concentration of high molecular A $\beta$  oligomer in the sample. As a standard, a 16-mer MAP peptide with 16 epitopes of anti-A $\beta$  antibody (BAN50) is used and the measured value is the value calculated on the basis of 16-mer MAP peptide.

○ Standard (16-mer MAP peptide)

**[Equipment or materials required but not supplied]**

- ☐ Purified water (distilled water)
- ☐ Tubes for dilution of the standard solution/sample
- ☐ Glass apparatuses (volumetric cylinders, beakers) for dilution of the wash solution
- ☐ Tip replaceable pipets (pipets with a disposable tip that can accurately pipet 10  $\mu$ L, and in the range between 200 and 500  $\mu$ L)
- ☐ Continuous dispensing pipet that can dispense 100  $\mu$ L continuously
- ☐ Water absorbent paper towel (to eliminate liquid remaining in plates after washing)
- ☐ Agitator (vortex type)
- ☐ Microplate shaker (about 600 - 800 rpm)
- ☐ 96-well plate washer (preferred if available) or washing bottle
- ☐ 96-well plate reader (for luminescence measurement)
- ☐ Data processing software

**[Preparation of reagents]**

(1) Standard Solution (10 nM)

Mix the Standard Solution (10 nM) and the Buffer as described below to prepare the standard solutions.

Concentration (pM) of standard solution	Volume of the standard solution	Buffer
100	Standard Solution (10 nM) : 10 $\mu$ L	990 $\mu$ L
40.0	100 pM standard solution : 200 $\mu$ L	300 $\mu$ L
16.0	40.0 pM standard solution : 200 $\mu$ L	300 $\mu$ L
6.40	16.0 pM standard solution : 200 $\mu$ L	300 $\mu$ L
2.56	6.40 pM standard solution : 200 $\mu$ L	300 $\mu$ L
1.02	2.56 pM standard solution : 200 $\mu$ L	300 $\mu$ L
0.41	1.02 pM standard solution : 200 $\mu$ L	300 $\mu$ L
0.00	—	300 $\mu$ L

(2) Biotin-conjugated Antibody (BAN50) Solution

Dilute the solution 100 times with the Buffer.

- (3) **Peroxidase-conjugated Streptavidin Solution**  
Dilute the solution 100 times with the Buffer.
- (4) **Luminescent Reagent 1 and (5) Luminescent Reagent 2**  
Mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1 : 1 (vol./vol.) 15 to 30 minutes before use. Store the reagent solution away from light until use.  
Example : Mix 6 mL of Luminescent Reagent 1 and 6 mL of Luminescent Reagent 2 (when all of 96 wells are used)
- (6) **Washing Solution (10 ×)**  
Dilute the solution 10 times with purified water (distilled water).  
Example : Dilute 100 mL of the concentrated wash solution (10 ×) with 900 mL of purified water (distilled water) (when all of 96 wells are used)

☐ Use other reagents as they are.

**[Storage and stability of reagents]**


- (1) **Antibody (BAN50)-coated Plate**  
Return unused (under a sealed and refrigerated condition) antibody-coated strips in the zip-sealed pack in the container and store at 2 to 10 °C. The stability is maintained within the expiration date.
- (2) **Standard Solution (10 nM)**  
When the kit is divided and used, take out the standard solution from the refrigerator and dilute it immediately before use and keep the remaining stock solution with the cap tightly closed at 2 to 10 °C without returning to room temperature. The stability is maintained within the expiration date. Use the respective diluted standard solutions immediately and do not store them after use.
- (3) **Buffer**  
If a part of the solution is used, transfer a volume slightly greater than the necessary volume in another container. Promptly close the cap of the container containing the remaining solution tightly and store at 2 to 10 °C without returning to the room temperature. The stability is maintained within the expiration date.
- (4) **Biotin-conjugated Antibody (BAN50) Solution and (5) Peroxidase-conjugated Streptavidin Solution**  
When the kit is divided and used, take out the solutions from the refrigerator and dilute them immediately and keep the remaining stock solutions with the caps tightly closed at 2 to 10 °C without returning to room temperature. The stability is maintained within the expiration date. Discard the diluted solutions remaining after use.
- (6) **Luminescent Reagent 1 and (7) Luminescent Reagent 2**  
When the kit is divided and used, take out the reagents from the refrigerator and prepare them. Keep the remaining luminescent reagents with the caps tightly closed at 2 to 10 °C without returning to room temperature. The stability is maintained within the expiration date.
- (8) **Washing Solution (10 ×)**  
Store the Wash Solution (10 ×) at 2 to 10 °C with its cap tightly closed. The stability is maintained within the expiration date. Discard the diluted wash solution remaining after use.

#### [Measurement procedure]

1. Wash the Antibody (BAN50)-coated Plate 4 times by filling the wells with the prepared wash solution. Then, invert the plate on a paper towel and pat lightly to eliminate the solution remaining in the wells.
2. Dispense 100  $\mu$ L each of the standard solutions at the respective concentrations in the reference standard measurement wells.
3. Dispense 100  $\mu$ L each of the sample diluted with the Buffer in the sample measurement wells (2~4-time dilution in the standard operating procedure).
4. Shake the plate using a microplate shaker.
5. Attach the Plate Seal and leave to stand at room temperature (20 to 25  $^{\circ}$ C) for 2 hours.
6. After completion of the reaction, discard the reaction solution, fill the wells with the wash solution, and wash the plate 4 times. Then, invert the plate on a paper towel and pat lightly to eliminate the solution remaining in the wells.
7. Dispense 100  $\mu$ L each of the biotin-conjugated anti A $\beta$  (BAN50) antibody solution to the wells. Shake the plate using the microplate shaker.
8. Attach the Plate Seal and leave to stand at room temperature (20 to 25  $^{\circ}$ C) for 2 hours.
9. After completion of the reaction, discard the reaction solution, fill the wells with the wash solution, and wash the plate 4 times. Then, invert the plate on a paper towel and pat lightly to eliminate the solution remaining in the wells.
10. Dispense 100  $\mu$ L each of the peroxidase-conjugated streptavidin solution to the wells. Shake the plate using the microplate shaker.
11. Attach the Plate Seal and leave to stand at room temperature (20 to 25  $^{\circ}$ C) for 30 minutes.
12. After completion of the reaction, discard the reaction solution, fill the wells with the wash solution, and wash the plate 4 times. Then, invert the plate on a paper towel and pat lightly to eliminate the solution remaining in the wells.
13. Dispense 100  $\mu$ L each of the mixed luminescent reagents to the wells. Shake the plate for 1 minute using the microplate shaker.
14. After agitation, determine the luminescent intensity using the 96-well microplate reader (for luminescence measurement). It is recommended to perform the measurement 10 to 20 minutes after addition of the luminescent reagent.
15. Create a standard curve by plotting the concentrations (pM) of the standard solutions along the X axis and the luminescence intensity along the Y axis. Read the concentration (pM) corresponding to the luminescence intensity of the diluted sample. Multiply the concentration read with the sample dilution factor to obtain a measured value. Convert the measured value to the concentration based on the 16-mer MAP peptide.

\*Use of a tertiary polynomial equation with 4 or 5 parameters is recommended for the arithmetic processing using computer software.

#### Calibration curve

\*The concentration is a value converted on the basis of 16-mer MAP peptide.


### [Precautions for use]



1. Please use container made of polypropylene. If other containers are used, A $\beta$  may adsorb to the container and the measured value may decrease.
2. Storage in a freezer at  $-35^{\circ}\text{C}$  or lower is recommended when a sample is stored for a long period of time. Do not repeat freezing and thawing. Thaw the frozen sample and thoroughly shake just before measurement. Dilute the sample just before use.
3. Do not use a sample with hemolysis or containing high lipid.
4. Centrifuge a sample with opacity and insoluble matter before measurement.
5. With a sample suspected of the influence of an interfering substance, check the linearity of dilution of the same sample using 2 or more different dilution factors.
6. When a sample is diluted, dilute it with the Buffer in advance using test tubes (PP, PE) and dispense in the measurement wells. The dilution factor is 2~4-time dilution in the standard operating procedure.
7. This kit should be used by a person who has completed the training for the ELISA or under his/her supervision.
8. For measurement by the manual operation, the kit should be used by a person who is well-trained for stable reproducibility of pipetting operation.
9. Wear gloves, goggles and protective gown during preparation and operation of the kit.
10. Take care not to touch the reagents to the skin. If any of the reagents of this kit has adhered to the eye, mouth, wound or skin by error, promptly take emergent measures such as washing thoroughly with tap water and visit a physician if necessary.
11. Do not eat, drink or smoke at a place where this kit is used.
12. Handle samples with thorough care as a substance with a risk of infection. This kit contains animal-derived ingredients.
13. Immerse used samples and used consumables in 1% formalin, 2% glutaraldehyde or sodium hypochlorite solution at 0.1% or higher for at least 1 hour. Or, autoclave these before disposal. Discard the used consumables and unused drugs in accordance with the rules of the testing site and local regulations.
14. Do not mix the reagent of different lot numbers. Be sure to affix the Plate seal during standing for reaction in each step to prevent well drying, contamination of foreign materials, temperature bias and evaporation of the dispensed reagents.
15. The ELISA is affected by the measurement environment. Strictly observe measurement procedure and room temperature of 20 to 25  $^{\circ}\text{C}$  (on the experimental table or incubator) at the place for reaction on standing. Never perform the measurement in an environment with air blowing (including circulation from the air-conditioner) and low humidity.

### [Example of measurement]

#### Reactivity with in vitro A $\beta$ oligomer

After overnight incubation of A $\beta$  1-42 peptides in a refrigerator, separate the fractions containing peptides not greater than 30 kDa (8-mer and smaller) and that containing peptides equal to and greater than 30 kDa (9-mer and larger) and the respective fractions were measured using this ELISA kit.


→No signal was detected in the fraction containing 8-mer and smaller (< 30 kDa), whereas a strong signal was exhibited in the fraction containing 9-mer and greater (> 30 kDa).

#### Measurement of a cerebrospinal fluid sample

Cerebrospinal fluids of healthy person (healthy), patients with mild cognitive impairment (MCI) and patients with Alzheimer disease (AD) were measured with this ELISA kit.

	Aβ oligomer(pM)	MMSE score	ADAS score
Healthy①	0.98	-	-
Healthy②	1.99	-	-
Healthy③	1.69	-	-
Healthy④	2.67	-	-
Healthy⑤	4.59	-	-
Healthy⑥	1.73	-	-
MCI	4.59	26	12
AD①	4.4	20	16
AD②	4.11	20	34
AD③	4.95	26	23
AD④	5.84	24	22
AD⑤	5.74	15	N/A

→There was a difference between healthy person and Alzheimer's disease patients.

#### [Outline of the measurement procedure]

- ☐ Return the plate and reagents thoroughly to room temperature (20 - 25 °C).
- ☐ Dilution of the concentrated wash solution : Dilute 10 times with purified water returned to room temperature.
- ☐ Dilution of the standard solution (example) : Dilute with the Buffer returned to room temperature.

Concentration (pM)		100	40.0	16.0	6.40	2.56	1.02	0.41	0
Example of dilution	Standard solution (μL)	Stock solution: 10							
	Buffer (μL)	990	300	300	300	300	300	300	300

\* : Standard solution of one higher concentration

- ☐ Antibody (BAN50)-coated Plate
- ☐ ↓ Washing 4 times (\*(i))
- ☐ Sample or standard solution 100 μL/well
- ☐ ↓ Shake (\*(ii)), 2-hour reaction at room temperature (20-25 °C), standing (\*(iii))
- ☐ ↓
- ☐ \*Preparation of the biotin-conjugated antibody (BAN50) solution (dilute 100 times with the Buffer returned to room temperature)

- ☐ ↓ Washing 4 times (\*i))
- ☐ Biotin-conjugated antibody (BAN50) solution 100 µL/well
- ☐ ↓ Shake (\*ii)), 2-hour reaction at room temperature (20 - 25 °C), standing (\*iii))
- ↓
- ☐ \*Preparation of the peroxidase-conjugated streptavidin solution (dilute 100 times with the Buffer returned to room temperature)
- ☐ ↓ Washing 4 times (\*i))
- ☐ Peroxidase-conjugated streptavidin solution 100 µL/well
- ☐ ↓ Shake (\*ii)), 30-minute reaction at room temperature (20 - 25 °C), standing (\*iii))
- ↓
- ☐ \*Preparation of the luminescent reagent (mix Luminescent Reagent 1 and Luminescent Reagent 2 in 1 : 1 (vol./vol.))
- ☐ ↓ Washing 4 times (\*i))
- ☐ Luminescent reagent 100 µL/well
- ☐ ↓ Shake for 1 minute at room temperature (20 - 25 °C)
- ☐ Measurement of luminescence intensity (to be measured between 10 and 20 minutes)

(\*i)) In each washing, dispense the wash solution in the wells, shake lightly for 10 seconds on the palm and discard. After 4 consecutive washings, invert the plate on a paper towel and pat to eliminate the wash solution completely. After elimination of the wash solution, promptly dispense the subsequent solution taking care not to dry the wells. The volume of the wash solution to be pipetted to the plate is about 300 µL/well.

(\*ii)) Shake 3 times each at about 600 - 800 rpm for 10 seconds.

(\*iii)) After completion of agitation, attach a Plate Seal and leave to stand. Peel the protective paper from the Plate Seal and attach the adhesive side to the plate. Never re-use the Plate Seal used once.

#### **[Storage condition]**

Store at 2 to 10 °C.

#### **[Expiration date]**

Indicated on the label

#### **[Packaging]**

For 96 tests

#### **[References]**

- 1) Fukumoto, H., *et al.* : *FASEB J.*, **24**, 2716 (2010).
- 2) Kasai, T., *et al.* : *Biochem. Biophys. Res. Commun.*, **422**, 3, 375 (2012).
- 3) Kasai, T., *et al.* : *Neurosci. Lett.*, **551**, 17 (2013).
- 4) Chambers, J. K., *et al.* : *Acta. Neuropathol. Commun.*, **3**, 78 (2015).
- 5) Kasai, T., *et al.* : *PLoS One*, **12**, 4, e0174630 (2017).

---

## **FUJIFILM Wako Pure Chemical Corporation**

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
Telephone : + 81-6-6203-3741  
Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

### **FUJIFILM Wako Chemicals U.S.A. Corporation**

1600 Bellwood Road  
Richmond, VA 23237  
U.S.A.  
Telephone : + 1-804-271-7677  
Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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