





Acetyl CoA plays a central role in energy metabolism. When energy is obtained from the breakdown of carbohydrates, amino acids, and lipids, it is mostly oxidized as acetyl CoA in the TCA circuit. RH-NH₂ has a tetramethylrhodamine backbone, which is known to localize to mitochondria, and acetylation of the amino group increases the fluorescence The fluorescence intensity increases when the amino group is acetylated (fluorescence quantum yield in PBS buffer (1% DMSO): RH-NH₂ 4.9 × 10⁻⁴, RH-NHAc 0.12).

RH-NH₂



When HeLa S3 cells are treated with $RH-NH_2$ (2 μ M), they emit weak fluorescence induced by acetyl CoA localized in mitochondria, as shown below, and when the accelerator tributylphosphine (PBu₃) is added, the fluorescence intensity increases within 1 minute.



Feature

- Acetyl CoA in cells can be visualized.
- 1 mg package for easy use in assays

<Reference>

Komatsu, H. et al. : Chem. Commun., 49, 2876 (2013).

Product List

Code No.	Product Name	Structure	CAS RN [®]	Package Size
			Storage conditions	
181-03521	RH-NH ₂	$H_3C_{N} \leftarrow CF_3COO^-$	1427006-06-1	• 1 mg
			Keep at -20 degrees C.	

Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: https://labchem-wako.fujifilm.com/ / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us.

FUJIFILM Wako Laboratory Chemicals site https://labchem-wako.fujifilm.com

FUJIFILM Wako Pure Chemical Corporation 1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Tel: +81 6 6203 3741 Fax: +81 6 6203 1999 ffwk-cservise@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation 1600 Bellwood Road, Richmond, VA 23237, U.S.A. Toll-Free (U.S. only): +1 877 714 1920 Tel: +1 804 271 7677 Fax: +1 804 271 7791 wkuslabchem@fujifilm.com

Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road, Tsuen Wan, N.T., Hong Kong Tel: +852-2799-9019 Fax: +852-2799-9808 wkhk.info@fujifilm.com



FUJIFILM Wako Chemicals Europe GmbH Fuggerstr 12, 41468 Neuss, Germany Tel: +49 2131 311 0 Fax: +49 2131 311 100 labchem_wkeu@fujifilm.com

FUJIFILM Wako Chemicals (Hong Kong) Limited FUJIFILM Wako (Guangzhou) Trading Corporation Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road, Guangzhou, 510095, China Tel: +86-20-8732-6381(Guangzhou) Tel: +86-21-6288-4751(Shanghai) Tel: +86-10-6413-6388(Beijing) wkgz.info@fujifilm.com