




CEPT Cocktail is highly cytoprotective in single cell cloning

CEPT is a cocktail of small molecules that exhibit cytoprotective effects when added to cell culture. It particularly excels under high-stress conditions such as single cell cloning. Our research focused on the efficacy of the CEPT Cocktail in single cell cloning of human iPS cells, using an automated cell dispensing machine.





Product information

CEPT is a cocktail of small molecules developed at the National Institutes of Health (NIH). It contains four components, and CEPT is an acronym for them: Chroman 1, Emricasan, Polyamines, and Trans-ISRIB.

Features

- · Protects human ES/iPS cells from stresses including DNA damage, helping to maintain cell structure and function
- Filter-sterilized, ready-to-use cocktail solution
- One of the CultureSur™ series products that are tested for endotoxin contamination and are mycoplasma negative





Please check our website for more details.

Wako CEPT Search

Wako Cat. No.	Product Name	Grade	Package Size
033-26071	CultureSure™ CEPT Cocktail (1,000×)	for Cell Culture	300 μL

Related Products

Wako Cat. No.	Product Name	Grade	Package Size
197-17571	StemSure® hPSC Medium Δ	for Cell Culture	100 mL
193-17573			100 mL×4
064-05381	Fibroblast Growth Factor (basic) (FGF-basic / bFGF / FGF2), Human, recombinant, Animal-derived-free(154aa) (powder)	for Cell biology	50 μg
068-05384			100 μg
060-05383			1 mg

We also have many other reagents for regenerative medicine research. For details, please visit our website.

- Medium
- Cytokine / Growth Factor
- **Small Molecules**
- Stem Cell Detection / Removal

Wako stem cell culture Search



Listed products are intended for laboratory research use only, and not to be used for drug, food or human use. / Please visit FUJIFILM Wako Laboratory Chemicals site: https://labchem-wako.fujifilm.com/ / This leaflet may contain products that cannot be exported to your country due to regulations. / Bulk quote requests for some products are welcomed. Please contact us

FUJIFILM Wako Laboratory Chemicals site https://labchem-wako.fujifilm.com

1-2, Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Tel: +81 6 6203 3741 Fax: +81 6 6203 1999 ffwk-cservise@fujifilm.com

FUJIFILM Wako Chemicals U.S.A. Corporation 1600 Bellwood Road, Richmond, VA 23237, U.S.A. тои- венwooa коаd, кксптопа, VA 23237, U Toll-Free (U.S. only): +1 877 714 1920 Tel: +1 804 271 7677 Fax: +1 804 271 7791 wkuslabchem@fujifilm.com

FUJIFILM Wako Chemicals (Hong Kong) Limited FUJIFILM Wako (Guangzhou) Trading Corporation Room 1111, 11/F, International Trade Centre, 11-19 Sha Tsui Road, Tsuen Wan, N.T., Hong Kong Tel: +852-2799-9019 Fax: +852-2799-9808 wkhk.info@fujifilm.com

FUJIFILM Wako Chemicals Europe GmbH Fuggerstr 12, 41468 Neuss, Germany Tel: +49 2131 311 0 Fax: +49 2131 311 100 labchem_wkeu@fujifilm.com

Room 3003, 30/F., Dong Shan Plaza 69, Xian Lie Zhong Road, Guangzhou, 510095, China Tel: +86-20-8732-6381(Guangzhou) Tel: +86-21-6288-4751(Shanghai) Tel: +86-10-6413-6388(Beijing) wkgz.info@fujifilm.com

^{*}This product is manufactured and sold under license from NIH.