

A novel affinity-based method for the isolation of highly purified extracellular vesicles

Masayuki Yamane⁴, Wataru Nakai¹, Takeshi Yoshida^{1,2}, Diego Diez³, Yuji Miyatake^{1,2}, Takahiro Nishibu⁴, Ryo Ukekawa⁴, Naoko Imawaka⁴, Ken Naruse⁴, Yoshifusa Sadamura⁴ & Rikinari Hanayama^{1,2,5}

1 Laboratory of Immune Network, WPI Immunology Frontier Research Center (IFReC), Osaka University, Japan

2 Department of Immunology, Kanazawa University Graduate School of Medical Sciences, Japan

3 Quantitative Immunology Research Unit, WPI Immunology Frontier Research Center (IFReC), Osaka University, Japan

4 Wako Pure Chemical Industries Ltd, Japan

5 PRESTO, Japan Science and Technology Agency (JST), Japan




Abstract

Extracellular vesicles (EVs) such as exosomes and microvesicles serve as messengers of intercellular network, allowing exchange of cellular components between cells. EVs carry lipids, proteins, and nucleic acids derived from their producing cells, and have potential as biomarkers specific to cell types and even cellular states. However, conventional methods (such as ultracentrifugation or polymeric precipitation) for isolating EVs have disadvantages regarding purity and feasibility. Here, we have developed a novel method for EV purification by using Tim4 protein, which specifically binds the phosphatidylserine displayed on the surface of EVs. Because the binding is Ca²⁺-dependent, intact EVs can be easily released from Tim4 by adding Ca²⁺ chelators. Tim4 purification, which we have applied to cell conditioned media and biofluids, is capable of yielding EVs of a higher purity than those obtained using conventional methods. Tim4 protein can also be used as a powerful tool for quantification of EVs in both ELISA and flow cytometry formats. Therefore, the affinity of Tim4 for EVs will find abundant applications in EV studies.

MagCapture™ Exosome Isolation Kit PS

Tim4-affiity purification method



Method	Tim4-affinity method (MagCapture™ Exosome Isolation Kit PS)	Ultracentrifugation	Antibody-based affinity purification
EV's Purity	■■■	■■	■■■
State of vesicles	Intact	Intact	Not Intact
Operability	Easy and Stable	Easy	Easy and Stable
Recovery amount	■■■	■■	■■

MagCpature™ Exosome Isolation Kit PS





Sample Type: cell culture supernatant, serum, plasma, urine, etc.

◆ MagCpature™ Exosome Isolation Kit PS can purify any EVs which expose phosphatidylserine on the outer surface of their lipid bilayer. It has been confirmed that human, mouse, and bovine EVs can be purified by this isolation kit.



Code No.	Product Name	Package Size
299-77603	MagCapture™ Exosome Isolation Kit PS	2 tests
293-77601	MagCapture™ Exosome Isolation Kit PS	10 tests

Particle analysis of sEVs purifeid by the Tim4-affinity methods

Tim4-affiity purification



Ultracentrifugation




The appearance of sEVs isolated by the Tim4-affinity purification method matched the typical saucer-like shape as previous reported*, and **almost no contaminants could be observed**. In contrast, sEVs isolated by UC were accompanied by a large number of small precipitates probably derived from supplements added in advance into the medium.

*Raposp, G. *et al.* (1996)

sEVs from 10K sup of COLO201 cells isolated by each method were examined by nanoparticle tracking analysis (NTA) using NanoSight and transmission electron microscope (TEM)

PS Capture™ Exosome ELISA Kit


Tim4-based ELISA method




Code No.	Product Name	Package Size
297-79201	PS Capture™ Exosome ELISA Kit (Anti-mouse IgG POD)	96 Reactions

The yield and purity of sEVs

Tim4-based ELISA




BCA protein assay



sEVs from 10K sup of COLO201 cells isolated by each method were examined by Tim4-based ELISA analysis and BCA protein assay


The Tim4-affinity purification method could isolate **about twice as many sEVs** as the Ultracentrifugation method. In contrast, equal amount of proteins were detected in sEV fractions using BCA protein assay.

The Uptake of sEV




sEVs from COLO201 cells were labeled by the fluorescence dye and the uptake of these sEVs by HeLa cells was examined.

sEVs purified by the Tim4-affinity purification method were **more efficiently incorporated** into the HeLa cells than that by the ultracentrifugation.



Microarray analysis of miRNA



miRNA microarray analysis of sEVs from COLO201 cells was performed by the 3D-Gene (TORAY)

miRNA microarray analysis of sEVs from COLO201 cells revealed **high correlation of miRNA profiles between Tim4-affinity purification method and ultracentrifugation method**.

Conclusions

- Tim4-affinity purification method could be used for the **efficient isolation** of sEVs
- Tim4-affinity purification method can purify **intact** extracellular vesicles.
- sEVs purified by the Tim4-affinity purification method were **efficiently taken up** by the recipient cells.
- High correlation of microRNA profiles between Tim4-affinity purification method and ultracentrifugation was revealed.