For Research Use Only

GLP-1 ELISA Kit *Wako*,  
High Sensitive

Code No.299-75501

**Please read these instructions carefully before use.**

## Table of Contents

1. Intended use.....	3
2. Storage and expiration.....	3
3. Principle of the assay.....	3
4. Precautions.....	3
5. Reagents supplied.....	4
6. Equipment or materials required but not supplied.....	4
7. Preparation of reagents.....	4
8. Technical tips.....	6
9. Preparation of samples.....	6
10. Assay procedure.....	6
11. Calculations.....	7
12. Performance characteristics.....	7
13. Troubleshooting.....	8
14. Summary of assay procedures and a check list.....	9
15. Assay worksheet.....	10

### **1. Intended use**

GLP-1 ELISA Kit *Wako*, High Sensitive, is a sandwich ELISA system for quantitative measurement of mouse/rat GLP-1. For research use only. Not for human or animal therapeutic or diagnostic use.

### **2. Storage and expiration**

Store the kit at 2-10°C (do not freeze). The kit is stable at 2-10°C until the expiration date shown on the label on the box. Once opened, the reagents should be used as soon as possible to avoid lower optimal assay performance caused by storage environment.

### **3. Principle of the assay**

Standards of known GLP-1 content or samples are incubated in monoclonal anti GLP-1 antibody-coated wells to capture GLP-1. After 2 hours incubation, the wells are washed and biotin-labeled anti GLP-1 antibody is added. The wells are incubated further for 1 hour to bind with captured GLP-1. After washing, HRP (horse radish peroxidase)-labeled avidin is added and the wells are incubated for 30 minutes. After washing, HRP-complex remaining in wells is reacted with a chromogenic substrate (TMB). The reaction is stopped by addition of an acidic solution and absorbance of yellow product is measured spectrophotometrically at 450 nm (reference wavelength, 620 nm). The absorbance is proportional to the GLP-1 concentration. The standard curve is prepared by plotting absorbance against standard GLP-1 concentrations. GLP-1 concentrations in unknown samples are determined using this standard curve.

### **4. Precautions**

- For professional use only. Beginners are advised to use this kit under the guidance of an experienced person.
- Wear gloves and laboratory coats, and use clean laboratory glassware when handling the assay materials.
- Be careful not to allow the reagent solutions of the kit to touch the skin, eyes and mucus membranes. In case of contact with these, wash skin/eyes thoroughly with water and seek medical attention when necessary.
- Do not drink, eat or smoke in the areas where assays are carried out.
- In treating assay samples of human or animal origin, be careful for possible biohazards. This kit contains components of animal origin. These materials should be handled as potentially infectious.
- Used samples and tips should be rinsed in 1% formalin, 2% glutaraldehyde, or more than 0.1% sodium hypochlorite solution for more than 1 hour, or be treated by autoclaving before disposal. Dispose consumable materials and unused contents in accordance with applicable regional/national regulatory requirements.
- The materials must not be pipetted by mouth.
- Do not mix the reagents from various lots.
- In order to avoid dryness of wells, contamination of foreign substances and evaporation of dispensed reagents during incubation, never forget to cover the well plate with a plate seal supplied.
- ELISA can be easily affected by the laboratory environment. Room temperature should be strictly maintained at 20-25°C. Avoid airstream velocity over 0.4 m/sec including wind from an air conditioner, and humidity less than 30%.

## 5. Reagents supplied

Code No.	Components	State	Amount
296-75511	Antibody-coated Plate	Use after washing	96 (8×12) wells/1 plate
293-75521	GLP-1 Standard	Lyophilized Use after reconstitution	1 bottle
290-75531	Buffer Solution	Ready for use	60 ml/1 bottle
297-75541	Biotinylated Anti GLP-1 Antibody Solution	Concentrated Use after dilution	100 µl/1 vial
294-75551	Peroxidase-conjugated Streptavidin Solution	Concentrated Use after dilution	100 µl/1 vial
291-75561	TMB Solution	Ready for use	12 ml/1 bottle
298-75571	Stop Solution	Ready for use	12 ml/1 bottle
295-75581	Wash Solution (x10)	Concentrated Use after dilution	100 ml/1 bottle
292-75591	Plate Seal	—	4 sheets
-	Instruction Manual	—	1 copy

## 6. Equipment or materials required but not supplied Use as a check box

- Purified water (distilled water).
- Test tubes for preparation of standard solution series.
- Glassware for dilution of Wash Solution (a graduated cylinder and a bottle)
- Pipettes (disposable tip type). One should be able to deliver 50 µl precisely, and another for 100-500 µl.
- Syringe-type repeating dispenser like Eppendorf multipipette plus which can dispense 50 µl.
- Paper towel to remove washing buffer remaining in wells.
- A vortex-type mixer.
- A shaker for 96 well-plate (600-1,200 rpm).
- An automatic washer for 96 well-plate (if available), or a washing bottle with a jet nozzle.
- A 96 well-plate reader (450 nm ±10 nm, 620 nm ( 600-650 nm)).
- Software for data analysis, if available.

## 7. Preparation of reagents

- ◆ Bring all reagents of the kit to room temperature (20-25°C) before use.
- ◆ Prepare reagent solutions in appropriate volumes for the assay just before use. Do not store the diluted reagents.

### **[Concentrated reagents]**

#### ●GLP-1 Standard

Reconstitute the standard with the volume of purified water\* stated in a separate sheet (GLP-1 standard solution, 300 pM), and then mix with the Buffer in the kit to prepare standard solutions as shown below. \*Because the volume of purified water to be added differs depending on the lot, see the specified volume in the annex.

Volume of standard solution	Buffer Solution	Concentration (pM)	Concentration (pg/ml)
GLP-1 standard solution (300 pM): 50 µl	450 µl	30	100
30 pM solution: 120 µl	120 µl	15	50
15 pM solution: 120 µl	120 µl	7.5	25
7.5 pM solution: 120 µl	120 µl	3.75	12.5
3.75 pM solution: 120 µl	120 µl	1.88	6.2
1.88 pM solution: 120 µl	120 µl	0.94	3.1
0 pM (Blank)	120 µl	0.0	0.0

#### ●Biotinylated Anti GLP-1 Antibody Solution

Prepare working solution by dilution of **Biotinylated Anti GLP-1 Antibody Solution** with the **Buffer Solution** to 1:100.

#### ●Peroxidase-conjugated Streptavidin Solution

Prepare working solution by dilution of **Peroxidase-conjugated Streptavidin Solution** with the **Buffer Solution** to 1:100.

#### ●Wash Solution (x10)

Dilute 1 volume of the **Wash Solution (x10)** to 10 volume with purified water (or deionized water) to prepare washing buffer.

Example: 100 ml of the concentrated Wash Solution (x10) and 900 ml of purified water (or deionized water) in case all wells of a 96-well plate are used.

### **[Storage and stability]**

#### ●Antibody-coated Plate

Put the unused strips with seals back in a plastic bag with zip-seal which is originally used for a well-plate container and store at 2-10°C. The strips will be stable until the expiration date, but note that the reagents in the kit need to be used as soon as possible.

#### ●GLP-1 Standard

When dividing and using the kit, do not return the GLP-1 standard solution(300pM) to room temperature, immediately close the lid tightly and store at 2-10 °C. It remains stable until the expiration date. Please use each diluted standard solution immediately, don't keep it.

#### ●Buffer Solution and TMB Solution

If not opened, store at 2-10°C. It remains stable until the expiration date. Once opened, use them as soon as possible to avoid influence by environmental conditions.

#### ●Biotinylated Anti GLP-1 Antibody Solution and Peroxidase-conjugated Sreptavidin Solution

The rest Biotinylated Anti GLP-1 Antibody Solution or Peroxidase-conjugated Sreptavidin Solution is stored at 2-10°C in a container with the stopper closed tightly. It remains stable until the expiration date. Residual working solution, which is already diluted, should be disposed.

#### ●Stop Solution

In case the rest of Stop Solution is kept, close the stopper tightly and store at 2-10 °C. It remains

stable until the expiration date.

● Wash Solution (x10)

In case the rest of Wash Solution is kept, close the stopper tightly and store at 2-10°C. It remains stable until the expiration date. Dispose any unused diluted Wash Solution.

### 8. Technical tips

- Be careful to avoid any contamination of assay samples and reagents. The use of disposal pipette tips, and 1 tip for 1 well are recommended.
- The **TMB Solution** should be almost colorless or clear pale blue before use. Protect it from the light.
- The **Stop Solution** should be almost colorless before use.

### 9. Preparation of samples

- Samples should be immediately assayed. If assay samples have to be stored for a longer period, snap-freeze the samples and keep them below –35°C. Avoid repeated freeze-thaw cycles. Thaw the frozen samples and mix them thoroughly before use. Dilute samples immediately before the assay.
- When a serum separating medium is used, check previously if it works.
- Hemolyzed or hyperlipidemic serum samples are not suitable.
- Turbid samples or those containing insoluble materials should be centrifuged before testing to remove any particular matter.
- If presence of an interfering substance is suspected, examine with a dilution test at more than 2 points.
- Before starting the assay, dilute samples 5x in test tubes (PP or PE). Multiply the assay value by this dilution rate (5x) in calculation.

#### ■ Worksheet example

	Strip 1&2	Strip 3&4	Strip 5&6	Strip 7&8	Strip 9&10	Strip 11&12
A	30 pM	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
B	15 pM	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
C	7.5 pM	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
D	3.75 pM	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
E	1.88 pM	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
F	0.94 pM	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
G	0 pM	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
H	Positive control	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40

### 10. Assay procedure

- All reagents should be brought to room temperature before use.
  - Allow the 96-well plate to reach room temperature. Remove the cover sheet of the plate.
- (1) Wash the **Antibody-coated Plate** by filling the wells with washing buffer and discard. Repeat the wash and discard steps 3 times (\*①), then turn the plate upside-down onto several layers of paper towels to remove residual buffer in the wells.
  - (2) Pipette 50 µl of standard solution to the wells designated for standards.
  - (3) Pipette 50 µl of diluted sample to the designated sample wells (5x dilution is recommended).

- (4) Shake the plate gently on a plate shaker (\*②).
- (5) Stick a plate seal (\*③) on the plate and incubate for 2 hours at 20-25°C.
- (6) Discard the reaction mixture and rinse wells as step (1).
- (7) Pipette 50 µl of Biotinylated Anti GLP-1 Antibody Solution to all wells, and shake as step (4).
- (8) Stick a plate seal (\*③) on the plate and incubate the plate for 1 hour at 20-25°C.
- (9) Discard the reaction mixture and rinse wells as step (1).
- (10) Pipette 50µl of Peroxidase-conjugated Streptavidin Solution to all wells, and shake the plate as step (4).
- (11) Stick a plate seal (\*③) on the plate and incubate for 30 minutes.
- (12) Discard the reaction mixture and rinse wells as step (1).
- (13) Pipette 50 µl of TMB Solution to wells, and shake as step (4).
- (14) Stick a plate seal (\*③) on the plate and incubate the plate for 30 minutes at 20-25°C.
- (15) Add 50 µl of Stop Solution to all wells and shake the plate as step (4).
- (16) Measure the absorbance of each well at 450 nm (reference wavelength, 620 nm) within 30 minutes using a plate reader.

\*Refer to the page 10 for notes of \*①, \*② and \*③.


### 11. Calculations

- (1) Prepare a standard curve using paper by plotting readings of absorbance\* (Y-axis) against values of GLP-1 concentration (pM) on X-axis.

\*Absorbance at 450 nm minus absorbance at 620 nm.

- (2) Using the standard curve, read the GLP-1 concentration of a sample at its absorbance\*, and multiply the assay value by dilution factor if the sample is diluted.

- In case the absorbance of a sample is higher than that of the highest standard, repeat the assay after proper dilution of samples with Buffer Solution.
- The use of a 3rd order regression curve for log-log plot or 4 or 5 parameters method for log-normal plot in computer calculation is recommended.

This standard curve is for the purpose of illustration only. The readings of absorbance may change due to the assay environment.

### 12. Performance characteristics

- Assay range  
The assay range of the kit is 0.94 ~ 30 pM (for 5x dilution of sample, 4.7 ~ 150 pM).  
If some samples show absorbance more than that of 30 pM standard, repeat the assay after proper dilution of samples.
- Precision of assay  
Within assay variation (2 samples, 5 replicates assay) mean CV is less than 10 %.
- Reproducibility  
Between assay variation (3 samples, 3 replicates assay within 4 days) mean CV is less than 10 %.
- Dilution test  
Two plasma samples were serially diluted by 3 steps and assay was performed.  
The dilution curves showed excellent linearity,  $R^2 = 0.999$ .

### 13. Troubleshooting

- Low absorbance in all wells

Possible explanations:

- 1) The standard or samples might not be added.
- 2) Reagents necessary for coloration such as Biotinylated Anti GLP-1 Antibody Solution, Peroxidase-conjugated Streptavidin Solution or TMB Solution might not be added.
- 3) Wrong reagents related to coloration might have been added. Wrong dilution of Biotinylated Anti GLP-1 Antibody Solution or Peroxidase-conjugated Streptavidin Solution.
- 4) Contamination of enzyme inhibitor(s).
- 5) Influence of the temperature under which the kits had been stored.
- 6) Excessive hard washing of the well plate.
- 7) Addition of TMB Solution soon after taking out from a refrigerator might cause poor coloration owing to low temperature.

- The absorbance of blank is higher than that of the lowest standard concentration (0.94 pM).

Possible explanations:

Improper or inadequate washing.

- High coefficient of variation (CV)

Possible explanation:

- 1) Improper or inadequate washing.
- 2) Improper mixing of standards, positive controls or samples. (Frozen sample should be agitated sufficiently.)
- 3) Pipetting at irregular intervals.

- Q-1: The plate can be divided to use for the other testing?

A-1: Yes, cut off the clear seal on the plate with cutter along strip. Put the residual plate, which is still the seal on, in a refrigerator immediately.

- Q-2: What is the liquid in 96 well-plates found when the box is opened?

A-2: When 96 well-plates are manufactured, a preservation stabilizer is added to the wells.

- Q-3: When samples are thawed, a cloud-like substance appeared. Does this influence the assay?

A-3: It may lower the assay results or be out of the detection limit. When adding aprotinin and heparin together, fibrin may appear in samples. Raise the heparin concentration or use EDTA-2Na.

**14. Summary of assay procedure and a check list**  : Use as a check box

**\*First, read this instruction manual carefully and start the assay after confirmation of details.**

- Bring the well-plate and all reagents to 20-25°C. It will take about 2 hours
- Wash Solution must be diluted 10 times by purified water of room temperature (20-25°C ) to prepare washing buffer.
- Make GLP-1 standard solution : Add the volume of purified water stated in a separate sheet, shake and dissolve it completely (GLP-1 standard solution, 300 pM).
- Preparation of diluted GLP-1 standard solution:

Concentration (pM)	30	15	7.5	3.75	1.88	0.94	0
GLP-1 standard solution (300 pM) (μl)	50	120*	120*	120*	120*	120*	0
Buffer solution (μl)	450	120	120	120	120	120	120

\*One rank higher standard.

		<b>Precautions &amp; related info</b>
<input type="checkbox"/>	<b>Antibody -coated Plate</b>	After removal of washing buffer, dispense the next reagent immediately.
<input type="checkbox"/>	↓ Wash 3 times (*①).	
<input type="checkbox"/>	<b>Diluted Samples/Standards</b> 50 μl	Dilute reagents during the first reaction. After removal of washing buffer, dispense the next reagent immediately.
<input type="checkbox"/>	↓ Shake (*②), incubate (leave to stand still (*③)) for 2 hours at 20-25°C.	
<input type="checkbox"/>	Dilute Biotinylated Anti GLP-1 Antibody Solution 100 times by using Buffer Solution	
<input type="checkbox"/>	↓ Wash 3 times (*①).	Dilute reagents during the second reaction. After removal of washing buffer, dispense the next reagent immediately.
<input type="checkbox"/>	<b>Biotinylated Anti GLP-1 Antibody Solution</b> 50 μl	
<input type="checkbox"/>	↓ Shake (*②), incubate (leave to stand still (*③)) for 1 hour at 20-25°C.	
<input type="checkbox"/>	Dilute Peroxidase-conjugated Streptavidin Solution 100 times by using Buffer Solution.	After removal of washing buffer, dispense the next reagent immediately.
<input type="checkbox"/>	↓ Wash 3 times (*①).	
<input type="checkbox"/>	<b>Peroxidase-conjugated Streptavidin Solution</b> 50 μl	
<input type="checkbox"/>	↓ Shake (*②), incubate (leave to stand still (*③)) for 30 minutes at 20-25°C.	After removal of washing buffer, dispense the next reagent immediately.
<input type="checkbox"/>	↓ Wash 3 times (*①).	
<input type="checkbox"/>	<b>TMB Solution</b> 50 μl	Shake the plate immediately
<input type="checkbox"/>	↓ Shake (*②), incubate (leave to stand still (*③)) for 30 minutes at 20-25°C.	
<input type="checkbox"/>	<b>Stop Solution</b> 50 μl	Measurement of absorbance (450 nm, Ref 620 nm)
<input type="checkbox"/>	↓ Shake (*②).	

- \*①After dispensing washing buffer to wells, shake the plate lightly on the palm for 10 sec and remove the washing buffer. Guideline of washing volume: 300  $\mu$ l/well washing buffer is added by pipette. In case of using an 8 channel pipette, the back ground sometimes tends to be high. In such a case, increase washing frequency from 3 times to 4-6 times at the constant stroke after the reaction with HRP-conjugated Streptavidin. Recommended flow rate for a plate washer: 5-25ml/min.
- \*②Guideline of shaking: 600-1,200 rpm for 10 seconds x 3 times.
- \*③The plate should be sealed during the reaction. Peel off the protective paper from the seal and stick the seal on the plate. Do not reuse the plate seal.

### 15. Assay worksheet

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

[Storage condition] Store the kit at 2-10°C (do not freeze).

[Term of validity] 6 months from production (expiration date is indicated on the container).

---

## FUJIFILM Wako Pure Chemical Corporation

1-2, Doshomachi 3-Chome, Chuo-Ku, Osaka 540-8605, Japan  
 Telephone : + 81-6-6203-3741  
 Facsimile : + 81-6-6201-5964  
<http://www.wako-chem.co.jp>

### FUJIFILM Wako Chemicals U.S.A. Corporation

1600 Bellwood Road  
 Richmond, VA 23237  
 U.S.A.  
 Telephone : + 1-804-271-7677  
 Facsimile : + 1-804-271-7791  
<http://www.wakousa.com>

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