

For life science research only
Not for use in diagnostic procedures



iMatrix-411

Product No. 892 041 350 µg
Product No. 892 042 1,050 µg

Version 001
Store at 2-15 °C
Protect from light

Background Information

Laminin-411 is well known to bind to the integrin α6β1 which is located on the cell surface. iMatrix-411 is recombinant Laminin411-E8 fragments which domain is known the binding region to integrin α6β1. iMatrix-411 increases the induction frequency to the blood vessel endothelial cells from pluripotent stem cells.

Content

Recombinant Human Laminin411-E8 Fragments

Amount

175 µg / tube (892 041: 2 tubes, 892 042: 6 tubes)

Form

Liquid solution

Product Information

iMatrix-411 is recombinant human Laminin411-E8 fragments expressed by CHO-S cell (Life Technologies).

Storage and Stability

The liquid solution is stable at +2 to +15 °C until the expiration date printed on the label.
Protect from light.

iMatrix-411 is stable at 4 °C
for 2 years from the manufacturing date.

Activity

The dissociation constant of the binding activity with integrin α6β1 is less than 10 nM.

Application

iMatrix-411 is able to be used as cell culture substrate for various cell types including ES/iPS cells.

Procedure

- 1) Dilute the solution with sterile PBS(-). Coat dishes with 0.5 µg/cm².
* For example, for one well of a 6-well plate (9.6 cm²/well), add 9.6 µL of iMatrix-411 (4.8 µg) in 2 mL of PBS(-). Add 2 mL of diluted iMatrix-411 solution to the well.
- 2) Incubate for 1 h at 37 °C, 3 h at room temperature, or over night at 4 °C.
- 3) Remove remaining fluid from the coated surface. No rinse is needed.
- 4) Immediately plate the cells at desired density.
* Don't allow the plate to dry.
* The optimum coating concentration depends on cell lines, from 0.1 to 1.5 µg/cm².

References

Laminin 411 and 511 promote the cholangiocyte differentiation of human induced pluripotent stem cells
Kazuo Takayama, et al.

Laminin-guided highly efficient endothelial commitment from human pluripotent stem cells
Ryo Ohta, et al.

Regulatory Disclaimer

For life science research only. Not for use in diagnostic procedures.

Contact and Support

E-mail: info@matrixome.co.jp

Designed by Matrixome, Inc.
3-2 Yamadaoka,
Suita-shi, Osaka 565-0871, Japan
Institute for Protein Research, Osaka University

Manufactured by Nippi, Incorporated
1-1-1 Senjumidori-cho, Adachi-ku, Tokyo 120-8601, JAPAN